

specTM

X-Ray Diffraction Software

USER MANUAL and TUTORIALS

REFERENCE MANUAL

STANDARD MACROS GUIDE

FOUR-CIRCLE REFERENCE

ADMINISTRATOR'S GUIDE

HARDWARE REFERENCE

© 1989,1990,1992,1993,1994,1995,1996,1997,1998,1999,2000
by Certified Scientific Software. All rights reserved.

Portions derived from appendices to the Doctoral thesis of Alan Braslau,
Department of Physics, Harvard University, 1988.

This is version 2.8 of the `spec` documentation, printed August 16, 2001,
describing features of release 3.03 to 4.05 of the software.

`spec` and `C-PLOT` are trademarks of Certified Scientific Software.
All other trademarks and registered trademarks are the property of their respective
owners.

The material in this manual is furnished for informational use only, is subject to
change without notice and should not be construed as a commitment by Certified Sci-
entific Software. Certified Scientific Software assumes no responsibility or liability
for any errors or inaccuracies that may appear in this manual. The software de-
scribed in this manual is furnished under license and may only be used or copied in
accordance with the terms of such license.

spec

X-Ray Diffraction Software

Certified Scientific Software

PO Box 390640 □ Cambridge, Massachusetts 02139 □ (617) 576-1610

FAX: (617) 497-4242 □ spec@certif.com

<http://www.certif.com>

MANUAL SUMMARY

USER MANUAL and TUTORIALS

A beginner's guide to diffractometer operation and a tutorial on the basic features of the `spec` user interface.

REFERENCE MANUAL

An overview of the internal structure of `spec` and a complete description of all the built-in keywords, operators, grammar rules, commands and functions.

STANDARD MACROS GUIDE

Tips on writing macros, examples of macros from the standard library and a detailed description of the structure of the standard scan macros.

FOUR-CIRCLE REFERENCE

Describes the special functions, variables and macros used to operate the standard four-circle diffractometer. The various modes are described, along with the use of sectors, cut points and frozen angles. Also contains an explanation of how to fit the orientation matrix when the lattice parameters are unknown.

ADMINISTRATOR'S GUIDE

How to install and update the `spec` package. Explains the hardware configuration and motor settings file and how to use security features to protect motors from being moved by unqualified users.

HARDWARE REFERENCE

Information on the specific hardware devices and interfaces supported by `spec`.

TABLE OF CONTENTS

USER MANUAL AND TUTORIALS	1
Introduction	3
Beginner's Guide to Diffractometer Operation	4
Starting Up	4
Using the Printer and Data Files	6
Setting Motor Positions and Moving Motors	7
Counting	11
Scans	12
Introduction To the spec User Interface	14
spec as a Calculator	14
Command Recall (History)	15
Controlling Output To the Printer and Data Files	17
Using Variables	18
Flow Control	20
Macro Facility	22
Command Files	25
Status and Help	26
UNIX Commands	27
Moving Motors	28
Diffractometer Geometry	31
Counting	32
CAMAC, GPIB and Serial	34
Using spec with C-PLOT and Other UNIX Utilities	35
Standard Data File Format	35
Scans.4	36
Contents	40
Showscans	41
REFERENCE MANUAL	43
Introduction	45
Internal Structure Of spec	45
Syntax Description	46
Comments	46
Identifiers	46

Arrays	47
Keywords	48
Numeric Constants	49
String Constants	49
String Patterns and Wild Cards	50
Command Recall (History)	50
Starting Up	51
Keyboard Interrupts	54
Cleanup Macros	54
Exiting	55
Variables	55
Operators	59
Flow Control	62
Grammar Rules	64
Built-In Functions and Commands	68
Utility Functions and Commands	69
Input and Output	81
Commands For Variables	89
Commands For Macros	90
Hardware Functions and Commands	94
Data Handling and Analysis Functions	106
String and Number Functions	118
 STANDARD MACRO GUIDE	 121
Introduction	123
Some Tips	124
Utility Macros	127
UNIX Commands	127
Basic Aliases	127
Basic Utility Macros	128
Reading From Command Files	130
Saving To Output Devices	132
Start-up Macros	133
Motor Macros	135
Counting Macros	140
Plotting Macros	142

Reciprocal Space Macros	143
Scan Macros	145
Scan Miscellany	145
Motor Scans	146
Basic Reciprocal Space Scans	147
Special Reciprocal Space Scans	147
Temperature Scans	148
Powder Mode	148
Customizing Scan Output	148
Temperature Control Macros	149
Printer Initialization Macros	153
The Scan Macros In Detail	154
Standard Data-File Format	162
FOUR-CIRCLE REFERENCE	163
Introduction	165
Diffractometer Alignment	166
Orientation Matrix	168
Four-Circle Modes	169
Freezing Angles	172
Sectors	173
Cut Points	174
Four-Circle Files	175
Four-Circle Variables	175
Four-Circle Functions	177
Four-Circle Macros	178
Zone Macros	179
Least-Squares Refinement of Lattice Parameters	180
ADMINISTRATOR'S GUIDE	185
Introduction	187
Quick Install	187
Steps For Installing spec	187
Extracting the Distribution	188
Installing the spec Program Files	189
Installing Drivers	193

Selecting the Hardware Configuration	194
Fine Tuning the Hardware Configuration	195
Adding Site-Dependent Help Files	195
Adding Site-Dependent C Code	195
Updating spec	197
Installed Files	198
File Hierarchy	198
The CAMAC Debugging Program	199
Enabling I/O Ports On PC Platforms (not <i>linux</i>)	200
Enabling I/O Ports On PC Platforms Running <i>linux</i>	200
Enabling E/ISA I/O Ports On HP 700 Systems	201
The Configuration Editor	202
The Settings File	203
The Config File	203
Security Issues	210
Extra Protection	210
 HARDWARE REFERENCE	 213
Introduction	215
Interface Controllers and General Input/Output	215
CAMAC Controllers	216
GPIB Controllers	219
VME Controllers	226
Serial (RS-232C) Ports	227
Generalized CAMAC I/O	227
PC Port Input/Output	228
Motor Controllers	229
Timers and Counters	251
Multichannel Data Acquisition Devices	260
 REFERENCES	 267
 INDEX	 268
 MANUAL PAGE	 279

Preface For the Interim Manual

We apologize for the work-in-progress nature of this manual. We hope to have a completely up-to-date version someday in either this or the next century. Until then, we hope this version will be of some use.

Presently, the *Reference Manual* and *Administrator's Guide* included in this interim manual are more or less current with respect to release 4 of spec. Also, portions of the *Four-Circle Reference* have been recently rewritten. In addition, a new *Hardware Reference* has been started. However, the *User Manual and Tutorials* could use a major rewrite, and the *Standard Macros Guide* does not reflect many revisions to the standard macros.

At any rate, this version of the manual should certainly give you the flavor of spec. For details on exactly how macros are constructed, you should probably consult the source code for the macros in your current distribution. The most recent on-line *changes* help files should be perused for the most up-to-date information on new features.

Also, please note that our World Wide Web site at <http://www.certif.com> contains both online and downloadable PostScript copies of this manual, and will surely be updated with more interim versions, before the next *official* version is available.

Finally, if there are questions or ambiguities that cannot otherwise be resolved, by all means contact us at CSS for the final word.

Thanks for your patience.

G.S. 6/13/99

USER MANUAL AND TUTORIALS

Introduction

`spec` is a UNIX-based software package for instrument control and data acquisition widely used for X-ray diffraction at synchrotrons around the world and in university, national and industrial laboratories. Developed in 1986 for X-ray diffraction experiments, `spec`'s portability, flexibility and power are winning it increasing application as general-purpose data-acquisition software. `spec` is available on a wide range of UNIX platforms and supports numerous hardware configurations. Features include:

- Built-in code to regulate motor controllers and detection electronics using VME, CAMAC, GPIB, RS-232, PC-board and ethernet-socket interfaces.
- Generalized access for VME, CAMAC, GPIB, RS-232, PC I/O ports and socket I/O to read from and write to user devices.
- Sophisticated user interface with command interpreter, complete with variables, looping and flow control, allowing for creative experiment automation.
- Language uses a familiar C-like syntax.
- A command-file facility allows complicated or commonly used command sequences to be immediately called up.
- An easy-to-use macro facility, with a large library of predefined macros. Macros can be readily modified with any text editor to suit experiments.
- Scans, data-file formats, etc. are not built into the compiled program but defined as easily modified macros.
- High-resolution, real-time data plots are available on X Windows and Sunview systems.
- Macro libraries and geometry-calculation routines support two-, four-, five- and six-circle diffractometers, kappa diffractometers, many liquid surface X-ray diffractometers and other configurations. The standard four-circle diffractometer supports many advanced modes and includes features such as least-squares refinement of the lattice parameters. New geometry configurations can be easily created.
- Hardware configuration employs a spread-sheet-style interface to select device names, addresses, CAMAC slot assignments, motor parameters, etc.
- Security features let site administrators restrict access to particular motors (such as those at a synchrotron beam-line front end).
- Available for most UNIX Systems, widely used on Linux PC platforms and UNIX workstations, including SUN (both SunOS 4.x and Solaris 2.x), HP 700 series and IBM RS/6000.

Beginner's Guide to Diffractometer Operation

Starting Up

In this introduction to the basics of `spec`, as in most of this manual, the standard four-circle diffractometer configuration is used in the examples. Other specialized diffractometer geometries are available, but all configurations function in a similar manner.

To start up the four-circle version of the `spec` package from a UNIX shell, type:

```
% fourc
```

(In this manual, input from the keyboard will be indicated in **bold-faced type**.) You will see output similar to the following:

```
Welcome to 'fourc' Release 3.01.01
Copyright (c) 1987,1988,1989,1990,1991,1992,1993,1994,1995
by Certified Scientific Software. All rights reserved.

(Portions derived from a program developed at Harvard University.)

Using '/usr/local/lib/spec.d/' for auxiliary file directory (SPEC.D).

Getting configuration parameters from 'SPEC.D/fourc/config'.
Using '/dev/ca00' for CAMAC device.
Using National Instruments GPIB.
Using 8 motors.

=
Spec Hot Line: (617) 576-1610.
Type "h changes" for info on latest changes.
See spec manual for tutorials.
=

Opening "SPEC.D/standard.mac" at input nest level 1.
Warning: No open data file. Using "/dev/null".
Type "startup" to initialize data file, etc.

Opening "SPEC.D/four.mac" at input nest level 1.
Warning: Using default lattice constants.
(UB recalculated from orientation reflections and lattice.)

1.FOURC>
```

The welcome message identifies the geometry configuration (`fourc`) and the release number of the program (3.01.01). The directory name that contains `spec`'s auxiliary

files is then identified (/usr/local/lib/spec.d). That name is assigned to spec's interval variable named SPECD.

A *configuration* file is then read to obtain the hardware configuration (hardware devices and types, stepper motor parameters, etc.). Various messages are printed as the specific hardware devices are initialized. The number of motors described in the configuration file is reported. During the start-up hardware configuration, spec reads the current diffractometer angle settings from the motor-controller hardware registers and verifies that they agree with the positions stored in a *settings* file associated with the diffractometer. If there is a discrepancy, you will see output similar to the following:

```
E500 at 0 steps (.1 user), spec at 24431 steps (12.3155 user)
on motor 1, slot 6, "Theta".  Modify the E500 registers?
```

The E500 is one of many different motor controllers available. Since the controller shows 0 steps, it probably has been powered down, and the program value is probably correct. Type *yes* or *y* to modify the controller registers. If you are uncertain what to do, the safest thing is to immediately terminate the program without updating the motor *settings* files by typing the quit control character (usually a ^v on IBM AIX platforms and a ^\ on most others), and then seek help.

A *news* file is displayed each time the program starts up. In this example, the Spec Hot Line message is from that file. The spec administrator can keep the news file up-to-date with messages for local users.

The first time you run spec, standard command files from the auxiliary file directory are automatically read (SPECD/four.mac and SPECD/standard.mac). These files contain the standard macro definitions used to operate the diffractometer. There are also some commands that assign default values to the variables used in the macros. The displayed *warning* message about no data file being open is produced by these standard macros, along with the message that suggests running the startup macro.

Finally, you are prompted for input. The prompt indicates the geometry configuration and includes a prepended command sequence number that is used with the command recall (or history) feature. You can exit the program by typing a ^D (control-D) at the prompt:

```
1.FOURC> ^D
Bill's state is stored for /dev/console.

%
```

The closing message confirms that your spec state is saved. The spec state consists of all your current macro definitions, variables, open output files and command history. Each user has a unique state associated with a particular terminal. Your saved state is automatically restored the next time you run spec from the same terminal.

(See page 52 in the *Reference Manual* to see how to start `spec` with a state from another user or terminal.)

`spec` is built around an interpreter that has a C-like syntax and recognizes over a hundred built-in commands and function names. However, you will typically be invoking the standard macros. These macros are written to do specific jobs using the built-in commands and functions and require minimum keyboard input. If you just want to move motors, count photons and do scans, you will only have to learn the few standard macros presented in this *Beginner's Guide*.

Using the Printer and Data Files

The standard macros in `spec` are designed to keep records of the experiment in progress on a printer and in a data file, although neither is required. The `startup` macro will prompt for a printer and a data file, along with asking for many other parameters and options. For now, enter information just for the printer and data file and accept the current values for the other parameters.

```
1.FOURC> startup
```

```
Enter <return> for no change in the displayed parameters.  
The names of start-up macros that can be invoked separately  
are shown in parenthesis above a set of queries.  
Type ^C to return to command level without finishing.  
(Interrupting one of the specialized start-up macros will  
likely undo any changes entered for its associated parameters.)
```

```
(newsample)
```

```
Title for scan headers (fourc)? cu 110
```

```
(newfile)
```

```
Data file (/dev/null)? cu110/94_01_31.a
```

```
Using "cu110/94_01_31.a". Next scan is number 1.
```

```
Last scan # (0)? <return>
```

```
Use a printer for scan output (NO)? y
```

```
Printer device (/dev/null)? /dev/lp
```

```
(And so on ...)
```

```
2.FOURC>
```

When prompting for input, `spec` generally displays the default or current response in parentheses. Simply hitting `<return>` makes that selection.

You can use the `newfile` macro directly to open (or reopen) a data file. Usage is `newfile [filename [scan_number]]`. (As is the convention in this manual, the square brackets indicate optional arguments, and the *Courier Oblique* typeface

denotes variable parameters you supply.) The optional argument *scan_number* is the number of the last scan and should be specified when appending to an existing data file.

The standard `spec` macros allow you to use a printer to record scan data and other status information. Not all users use a printer, though, as the information is also stored in data files.

When using a printer, `spec` generates output for a 132-column wide format. Most `spec` users use 8½" wide paper with their printer set to compressed mode. The `initfx` macro sends the correct programming sequence to put an Epson printer into compressed mode. The `initdw` macro does the same for a Decwriter. (Other macros are available for other printers – type `lsdef init*` from `spec` for a list.) You could also use printer switches to select compressed mode.

Also, when using the printer, you should set the top-of-form position correctly. That way, each scan will begin at the top of a new page, and it will be much easier to locate scans when thumbing through the data printout later.

Use the `comment` macro (also available as `com`) to insert arbitrary comments in the data file and on the printer. For example,

```
2.FOURC> com Absorber inserted in front of detector
Mon Feb 15 01:41:52 1994. Absorber inserted in front of detector.
3.FOURC>
```

Setting Motor Positions and Moving Motors

When `spec` is used to control an X-ray diffractometer, the `wh` (where) macro is available to show the positions of the most interesting angles and the diffractometer position in reciprocal space coordinates. With the four-circle diffractometer, the output is as follows:

```
3.FOURC> wh
H K L = 0 0 1
Alpha = 30 Beta = 30 Azimuth = 90
Omega = 0 Lambda = 1.54

Two Theta      Theta      Chi      Phi
60.0000      30.0000     -90.0000    0.0000

4.FOURC>
```

The incident and scattered angles for surface diffraction (ALPHA and BETA), the

AZIMUTH angle used in advanced modes (see the *Four-Circle Reference*) and the incident X-ray wavelength, `LAMBDA`, are also listed.

The angular positions listed above are the *user* angles. You set user angles during diffractometer alignment to satisfy the premises of the geometry calculations, such as the positions of the zeroes of the angles. *Dial* angles keep track of hardware limits and prevent complete loss of angles from alignment errors or computer failure. The dial angles are generally made to agree with a physical indicator on each motor, such as a dial. User angles are related to the dial angles through the equation:

$$user = sign \times dial + offset$$

Redefining a user angle changes the internal value of *offset*. Dial angles are directly proportional to the values contained in the hardware controller registers. The *sign* of motion is set in the configuration file by the `spec` administrator and normally isn't changed.

The `set_dial motor position` macro is used to set the dial position of a motor. The argument *motor* is the motor number or mnemonic. All motors have short mnemonics, such as `tth`, `th`, `chi`, and `phi`.

```
4.FOURC> set_dial tth 24.526
```

```
Mon Feb 15 01:42:10 1994. Two Theta dial reset from 0 to 24.526.
```

```
5.FOURC>
```

The `set motor position` macro is used to set the *user* position of a motor (i.e., to change *offset*). If you had a slit motor with mnemonic `ts1`, you might enter

```
5.FOURC> set ts1 .5
```

```
Mon Feb 15 01:43:31 1994. Top Slit1 reset from 0 to .5.
```

```
6.FOURC>
```

The `wa` (where all) macro lists both the user and dial positions of all configured motors.

```
6.FOURC> wa
```

```
Current Positions (user, dial)
```

Two Theta	Theta	Chi	Phi	Top Slit1	Bot Slit1
tth	th	chi	phi	ts1	bs1
24.6310	12.3155	90.0000	0.0000	0.5000	-0.5000
24.5260	12.2155	89.7865	0.0950	0.5000	-0.5000

```
7.FOURC>
```

spec also keeps track of software motor limits. These limits are always checked before any motors are moved. The `lm` macro lists these limits in both user and dial angles, as well as the current positions of the motors.

```
7.FOURC> lm
```

```
USER Limits (high, current, low):
```

Two Theta	Theta	Chi	Phi	Top Slit1	Bot Slit1
180.1050	90.1000	135.2135	179.9050	5.0000	0.0000
24.6310	12.3155	90.0000	0.0000	0.5000	-0.5000
-179.8950	-89.9000	-134.7865	-180.0950	0.0000	-5.0000

```
DIAL Limits (high, current, low):
```

Two Theta	Theta	Chi	Phi	Top Slit1	Bot Slit1
180.0000	90.0000	135.0000	180.0000	5.0000	0.0000
24.5260	12.2155	89.7865	0.0950	0.5000	-0.5000
-180.0000	-90.0000	-135.0000	-180.0000	0.0000	-5.0000

```
8.FOURC>
```

The macro `set_lm motor low high` changes the software limits for a single motor. The values for *low* and *high* are given in *user* angles (although they are stored internally in *dial* angles).

The `wm motor [motor ...]` macro lists complete information for up to six motors given as arguments.

```
8.FOURC> wm tth th
```

	Two Theta	Theta
	tth	th
User		
High	180.1050	90.1000
Current	24.6310	12.3155
Low	-179.8950	-89.9000
Dial		
High	180.0000	90.0000
Current	24.5260	12.2155
Low	-180.0000	-90.0000

```
9.FOURC>
```

Once the diffractometer has been aligned, you can move to any allowed reciprocal space position using the `br H K L (Bragg)` macro.

```

9.FOURC> br 2 0 0

10.FOURC> wh

H = 2   K = 0   L = 0
ALPHA = -25.251   BETA = 25.251   AZIMUTH = 90   LAMBDA = 1.54

Two Theta      Theta      Chi      Phi
   50.5030     25.2515     90.0000     0.0000

11.FOURC>

```

You can see where the motors would move for particular values of (H, K, L) using the `ca H K L (calculate)` macro.

```

11.FOURC> ca 2 1 1

Calculated Positions:

H = 2   K = 1   L = 1
ALPHA = -25.252   BETA = 25.252   AZIMUTH = -90   LAMBDA = 1.54

Two Theta      Theta      Chi      Phi
   62.9960     31.4980     54.7355    135.0000

12.FOURC>

```

Conversely, `spec` will display the (H, K, L) that corresponds to a particular set of motor positions using the `ci tth th chi phi (calculate inverse)` macro.

A single motor may be moved in real space using the `mv motor position` macro. For example,

```

12.FOURC> mv tth 50

13.FOURC>

```

will move the 2θ motor to 50° . You are prompted for more input immediately, even though the motors are still moving.

You can tell when the motor has stopped moving by using the `w` macro. The program will pause until the motor has stopped moving and then generate a beep on the terminal. Alternatively, you can have the motor position displayed on the screen as it is moving by invoking the `umv (updated-move)` macro instead of `mv`. To stop the motors before they have finished moving, type the interrupt character, usually a `^c`.

You can use the `mvr motor relative_position` macro to move a motor relative to its current position.

```
13.FOURC> mvr th 1
```

```
14.FOURC>
```

will move θ by one degree.

The *tw motor delta* (tweak) macro is useful when lining up the diffractometer or when searching for the beam.

```
14.FOURC> tw th .1
```

```
Indicate direction with + (or p) or - (or n) or enter  
new step size. Type something else (or ^C) to quit.
```

```
th = 26.2515, which way (+)? <return>  
th = 26.3515, which way (+)? <return>  
th = 26.4515, which way (+)? <return>  
th = 26.5515, which way (+)? <return>  
th = 26.6515, which way (+)? -  
th = 26.5515, which way (-)? <return>  
th = 26.4515, which way (-)? ^C
```

```
15.FOURC>
```

Each time you hit *<return>*, the motor moves *delta* in the plus or minus direction.

Counting

You count photons using the *ct* macro. Without arguments, this macro counts for the time set by the variable *COUNT*, which is typically one second. According to the convention used in the standard macros, positive count times indicate counting to seconds and negative count times indicate counting to monitor counts.

```
15.FOURC> ct 10
```

```
Mon Feb 15 01:45:12 1994
```

```
Seconds = 10  
Monitor = 389387 (38939/s)  
Detector = 192041 (19204/s)
```

```
16.FOURC> ct -40000
```

```
Mon Feb 15 01:45:28 1994
```

```
Seconds = 1.027  
Monitor = 40000 (38948/s)  
Detector = 19756 (19237/s)
```

```
17.FOURC>
```

Type a `^c` to abort counting. The macro `show_cnts` will display the current scaler contents. The `uct` macro will update the screen with the current scaler contents during the counting period.

Scans

Scans in `spec` are built of macros. Many different standard scans are available. Absolute-position motor scans such as `ascan`, `a2scan` and `a3scan` move one, two or three motors at a time. Relative-position motor scans are `lup` (or `dscan`), `d2scan` and `d3scan`. The relative-position scans all return the motors to their starting positions after the last point. Two motors can be scanned over a grid of points using the `mesh scan`.

Simple reciprocal space scans are `hscan`, `kscan` and `lscan`. The `hklscan` macro moves the diffractometer along an arbitrary straight line in reciprocal space. Scans such as `hkcircle` or `hkradial` describe other trajectories. The `hklmesh` scan measures intensities over a grid of reciprocal-space points.

If you do not know the arguments for a scan or how a scan is used, you can call up its usage by typing its name with no arguments.

```
17.FOURC> ascan
Usage:  ascan motor start finish intervals time

18.FOURC> hscan
Usage:  hscan start finish intervals time

19.FOURC>
```

When the program does a scan such as `hscan`, the following happens: the program waits for motors to stop moving, calculates (H, K, L) for the current position and then scans H , holding K and L fixed for a reciprocal space scan along the H direction.

```
19.FOURC> hscan .9 1.1 20 1
Mon Feb 15 01:50:33 1994

Scan   1   Mon Feb 15 01:50:33 1994   file = cull10/90_01_31.a
hklscan 0.9 1.1 0 0 0 0 20 1

#       H       K       L Detector   Monitor   Seconds
0       0.9     0       0       2604    38939    1
1       0.91    0       0       3822    38820    1
2       0.92    0       0       5295    39034    1
3       0.93    0       0       7259    38789    1
4       0.94    0       0       9298    38804    1
5       0.95    0       0      11505    38909    1
6       0.96    0       0      13907    38821    1
7       0.97    0       0      16022    39110    1
```

8	0.98	0	0	17603	38839	1
9	0.99	0	0	18834	38950	1
10	1	0	0	19103	38917	1
11	1.01	0	0	18701	39013	1
12	1.02	0	0	17652	39135	1
13	1.03	0	0	16011	38836	1
14	1.04	0	0	13848	38901	1
15	1.05	0	0	11585	38933	1
16	1.06	0	0	9302	39022	1
17	1.07	0	0	7237	39205	1
18	1.08	0	0	5324	38957	1
19	1.09	0	0	3780	38801	1
20	1.1	0	0	2580	38975	1

```
Peak at 1 is 19103   FWHM at 1 is 0.05   COM is 1
Sum = 231272   Ave.Mon./Time = 38921   Ave.Temp. = 0C
28 second
```

```
20.FOURC>
```

The output shown is what would generally appear on the screen. More detailed output is sent to the printer. Also, a complete scan header and the data points are stored in the data file. A rudimentary plot can be produced on the printer at the end of the scan. Typing `splot` will produce a plot of the data on the screen. Typing `pts` will list the data on the screen.

The `setplot` macro configures how the data will be displayed during and at the conclusion of scans.

```
20.FOURC> setplot
```

```
Do real-time screen plots during scans (NO)? y
Do screen plot after scan (YES)? <return>
Do printer plot after scan (NO)? y
```

```
21.FOURC>
```

Scans can be aborted by typing `^c`. Typing `scan_on` restarts an aborted scan at the current point.

Introduction To the spec User Interface

spec as a Calculator

In some respects, the `spec` user interface behaves like a BASIC language interpreter that uses the C language syntax. For example, you can easily print strings and the results of arithmetic expressions:

```
1.FOURC> p 2+2, sqrt(3), "2^16 =", 1<<16
4 1.73205 2^16 = 65536

2.FOURC>
```

(The `p` macro is defined as `print`, a built-in command.) You do not need to search for your calculator, as all the standard operators and functions are available.

The arithmetic operators (`=`, `*`, `/`, `%`, `+`, `-`, `++`, `--`, `+=`, `-=`, `*=`, `/=`, `%=`), the relational operators (`>`, `<`, `<=`, `>=`, `==`, `!=`), the boolean operators (`!`, `&&`, `||`), the bitwise operators (`>>`, `<<`, `~`, `&`, `^`, `|`, `>>=`, `<<=`, `&=`, `^=`, `|=`) and the ternary operator (`? :`) are all available. Parentheses can be used for grouping within expressions. See the *Reference Manual* for a description of all the operators and their rules of precedence.

The most useful standard C math functions are included, such as `sin()`, `cos()`, `tan()`, `asin()`, `acos()`, `atan()`, `exp()`, `log()`, `log10()`, `pow()`, `sqrt()`, and `fabs()`. Conversions functions such as `deg()` and `rad()` convert between degrees and radians, while `bcd()` and `dcb()` convert between decimal and binary-coded decimal. A `rand()` function to return random numbers is also provided.

Numbers can be entered in decimal, octal or hexadecimal notation, just as in C.

```
2.FOURC> p 100, 0100, 0x100
100 64 256

3.FOURC>
```

Special string functions also exist. The `date()` function provides the current date and time as a string:

```
3.FOURC> p date()
Mon Feb 15 02:13:13 1994

4.FOURC>
```

The `date()` function can also take an argument that is the number of seconds from the UNIX epoch.


```
4.FOURC> p date(1e9)
Sat Sep 8 21:46:40 2001
```

```
5.FOURC> p date(0)
Wed Dec 31 19:00:00 1969
```

```
6.FOURC> p int(time()), date(time())
729760917 Mon Feb 15 02:21:57 1994
```

```
7.FOURC>
```

The second example shows the (Eastern Standard Time) moment of the UNIX epoch. The function `time()` returns the number of seconds since that moment, including a fractional part with a resolution determined by the system clock. The difference of subsequent calls to `time()` can, for example, give a reasonable elapsed time for each point in a scan.

The function `input()` reads a string from the keyboard. An optional argument will be printed first. For example, a macro or command file might prompt you for information:

```
7.FOURC> TITLE = input("Please enter a title: ")
Please enter a title: Au (001) Sample #1
```

```
8.FOURC>
```

Other string functions such as `index()`, `substr()`, `length(s)` and `sprintf(format, [args])` are also available. See the *Reference Manual* for details.

Command Recall (History)

A command recall (or *history*) feature lets you recall previously typed commands. `spec`'s command recall implements a subset of the features of the standard `cs` history mechanism. When using command recall, note that only keyboard input is saved, command recall cannot be used in command files, and the command recall characters must occur at the beginning of a new line.

When you run `spec` interactively, a command sequence number is always prepended to the prompt. The `history` command lists by number the commands that can be recalled. (At present, only the most recent 1000 commands are available for recall.)

```

8.FOURC> history
  1 p 2+2, sqrt(3), "2^16 =", 1<<16
  2 p 100, 0100, 0x100
  3 p date()
  4 p date(1e9)
  5 p date(0)
  6 p time(), date(time())
  7 TITLE = input("Please enter a title: ")
  8 history

```

```
9.FOURC>
```

To use command recall, type `!!` or `!-1` to recall the previous command. Typing `!-2` will recall the second previous command. Type `!2` to recall command number 2. Also, `!TI` will recall the last command beginning with the string `TI`.

```

9.FOURC> !2
p 100, 0100, 0x100
100 64 256

```

```
10.FOURC>
```

Notice that the recalled command is first printed and then executed.

Recalled commands can be modified by appending text.

```

10.FOURC> !p , "= 100, 0100 and 0x100."
p 100, 0100, 0x100 "= 100, 0100 and 0x100."
100 64 256 = 100, 0100 and 0x100.

```

```
11.FOURC>
```

Arbitrary substitutions to recalled commands are allowed using the `:s/left/right/modifier`, as in

```

11.FOURC> !-1:s/./,respectively./
p 100, 0100, 0x100 "= 100, 0100 and 0x100, respectively."
100 64 256 = 100, 0100 and 0x100, respectively.

```

```
12.FOURC>
```

You can also use a circumflex `^` to make a substitution on the most recent command, just as with the standard UNIX `csh`.

```

12.FOURC> ^=^are^
p 100, 0100, 0x100 "are 100, 0100 and 0x100, respectively."
100 64 256 are 100, 0100 and 0x100, respectively.

```

```
13.FOURC>
```

Controlling Output To the Printer and Data Files

`spec`'s output facility is unusual. Output files and devices, including the screen, are turned on or off for output. The output of each printing command, whether generated by a user command or internally, is sent to all the turned-on devices.

`open("filename")` opens a file or device to append output. The current contents of existing files are never erased. The `on("filename")` function turns on printing to the file or device and opens the file if `open()` was not previously called. The `off("filename")` function ends printing to that file or device, and `close("filename")` closes the file or device and removes the name from the program's table of file pointers. The name `"tty"` is special when used as an argument to these functions. It always refers to your current terminal.

Whenever there is an error or a `^c` interrupt, all files (except *log* files) are turned off, and output to the terminal is turned on. A log file is used for debugging purposes and is any file that begins with `log`. Output to all *on* files and devices is automatically copied to a log file.

To get the status of all open files, type:

```
13.FOURC> on()
'tty' has output ON.
'/usr/alan/default.dat' has output OFF.
'/dev/null' has output OFF.

14.FOURC>
```

If you change `spec`'s current directory, you can reference open files either by the name with which the files were opened or by the correct path name relative to the new directory.

The standard macros use three output devices: the screen, a printer and a data file. The `ont`, `offt`, `onp`, `offp`, `ond` and `offd` macros are usually used to simplify controlling output to these devices, where `ont` is defined as `on("tty")`, etc. Typical usage is

```
ond; offt; printf("#S %d\n", ++SCAN_N); offd; ont
```

Often, printing commands are placed between `onp` and `offp` to direct the output both to the screen and the printer. For instance,

```
1.FOURC> onp; p "This is also being printed on the printer."; offp
This is also being printed on the printer.

2.FOURC>
```

Formatted printing is available using the `printf()` and `fprintf()` functions. The

format specifications are the same as for the C-language routine and can be found in the *printf()* write-up in any C reference manual.

```
2.FOURC> printf("The square root of two is %.12g.\n", sqrt(2))
The square root of two is 1.41421356237.

3.FOURC>
```

Using Variables

spec's variables can be used as both strings and as double-precision floating-point numbers. Variables are not declared, but come into existence through usage. Some variables are built-in, though, and of these, some have preassigned values. The variable `PI` is an example.

```
3.FOURC> {
4.more> k = 2 * PI / 1.54
5.more> print k
6.more> }
4.07999

7.FOURC>
```

Curly brackets (`{` and `}`) are used to delimit a block to be interpreted together, since variables are local to interpreted blocks. Notice that the prompt indicates the program is expecting further input before interpreting and taking action.

```
7.FOURC> print k
0

8.FOURC>
```

The value of `k` disappeared because `k` was local to the previous statement block. New variables start off with a value of zero.

A variable may be declared `global` to hold its value outside an interpreted block:

```
8.FOURC> global Lambda CuKa

9.FOURC> Lambda = 1.54

10.FOURC> CuKa = "Copper K-alpha"

11.FOURC> print CuKa, "=", Lambda
Copper K-alpha = 1.54

12.FOURC>
```

By convention, global variables in the standard macro package use capital letters or begin with an underscore. Variables can be made constant to protect them from accidental reassignment,

```
12.FOURC> constant Lambda 1.54
```

```
13.FOURC> Lambda = 1.7
Trying to assign to a constant 'Lambda'.
```

```
14.FOURC>
```

Variables defined as constant are automatically global.

Most built-in variables with preassigned values are of the *immutable* type and cannot be changed at all:

```
14.FOURC> PI = 1
Trying to assign to an immutable 'PI'.
```

```
15.FOURC>
```

Some built-in variables, such as `DEBUG`, can be changed by the user. Another variable, the `A[]` array, may be filled by the program with the current motor positions or can be set to target motor positions. For a list of all current symbols, type:

```
15.FOURC> syms
( Built-In/Global/Local Array Number String Constant/Immutable )
8256 A (BA...) 80 SLIT_W (G.NS.) 80 _f1 (G.NS.)
96 ADMIN (G..S.) 96 SPEC (B..SI) 80 _f2 (G.NS.)
80 BG (G.N..) 112 SPECD (B..SI) 80 _f3 (G.NS.)
96 COLS (B.N..) 80 TEMP_CS (G.NS.) 80 _fx (G.NS.)
80 COUNT (G.N..) 80 TEMP_SP (G.NS.) 80 _g1 (G.NS.)
96 COUNTERS (B.N.I) 96 TERAMP_MIN (G.N..) 80 _g2 (G.NS.)
96 COUNT_TIME (G.NS.) 96 TERM (B..S.) 96 _hkl_col (G.N.C)
128 CP_FILTER (G..S.) 80 TIME (G.NS.) 80 _m (G.NS.)
112 CWD (B..SI) 96 TIME_END (G.NS.) 80 _m1 (G.NS.)
128 DATAFILE (G..S.) 96 TITLE (G..S.) 80 _m2 (G.NS.)
112 DATA_DIR (G..S.) 80 T_AV (G.NS.) 80 _m3 (G.NS.)
80 DATE (G.NS.) 80 T_HI_SP (G.N..) 80 _n1 (G.NS.)
80 DEBUG (B.N..) 80 T_L (G.NS.) 80 _n2 (G.NS.)
80 DEGC (G.NS.) 80 T_LO_SP (G.N..) 80 _nm (G.NS.)
80 DEGC_SP (G.NS.) 2736 U (BA...) 80 _numgeo (G.N..)
80 DET (G.N..) 1056 UB (BA...) 80 _pmot (G.NS.)
80 DOFILE (G.NS.) 592 UNITS (GA...) 96 _pmotflag (G.NS.)
96 DO_DIR (G..S.) 80 UPDATE (G.N..) 96 _pre_chk (G.N..)
80 EPOCH (G.N.C) 96 USER (B..SI) 80 _pwid (G.NS.)
80 FPRNT (G.NS.) 96 USER_CHK_ACQ (G.NS.) 80 _reg_f (G.NS.)
1296 G (BA...) 96 USER_CHK_COUNT (G.NS.) 80 _reg_i (G.NS.)
96 GS_file (G..S.) 96 USER_CHK_MOVE (G.NS.) 80 _reg_n (G.NS.)
80 GS_scan (G.N..) 80 VFMT (G.NS.) 80 _reg_s (G.NS.)
80 GS_xcol (G.N..) 80 VPRNT (G.NS.) 112 _reg_scan (G..S.)
80 GS_ycol (G.N..) 80 X_L (G.NS.) 80 _reg_t (G.NS.)
```

96	GTERM	(B..S.)	80	Y_L	(G.NS.)	80	_s	(G.NS.)
80	HEADING	(G.NS.)	816	Z	(BA...)	80	_s1	(G.NS.)
112	HOME	(B..SI)	80	_1	(G.N..)	80	_s2	(G.NS.)
96	MAIL	(G..S.)	80	_2	(G.NS.)	80	_s3	(G.NS.)
80	MODES	(G.N.C)	80	_3	(G.NS.)	80	_sleep	(G.NS.)
80	MON	(G.N..)	80	_4	(G.NS.)	80	_stime	(G.NS.)
96	MON_RATE	(G.NS.)	80	_5	(G.NS.)	80	_stypc	(G.NS.)
80	MOTORS	(B.N.I)	80	_6	(G.NS.)	80	_sx	(G.NS.)
80	MT_AV	(G.NS.)	80	_7	(G.NS.)	96	_upd_flg	(G.NS.)
80	NPTS	(G.NS.)	80	_8	(G.NS.)	80	bg_m	(G.NS.)
80	PFMT	(G.NS.)	80	_9	(G.NS.)	80	bg_pts	(G.N..)
80	PI	(B.N.I)	80	_LAMBDA	(G.NS.)	80	bg_yI	(G.NS.)
96	PLOT_MODE	(G.N..)	96	_bad_lim	(G.NS.)	80	chi	(B.N.I)
80	PL_G	(G.N..)	80	_c1	(G.NS.)	96	chk_thresh	(G.NS.)
80	PL_G1	(G.N..)	80	_c2	(G.NS.)	80	det	(B.N.I)
80	PL_X	(G.N..)	80	_c3	(G.NS.)	1184	gmodes	(GA...)
80	PL_Y	(G.N..)	80	_c4	(G.NS.)	656	mA	(GA...)
80	PPRNT	(G.NS.)	80	_cols	(G.NS.)	80	mon	(B.N.I)
112	PRINTER	(G..S.)	80	_const	(G.NS.)	80	phi	(B.N.I)
1856	Q	(BA...)	80	_cp	(G.NS.)	128	rplot_col	(G..S.)
80	REFLEX	(G.NS.)	80	_ctime	(G.NS.)	80	sec	(B.N.I)
80	ROWS	(B.N..)	80	_d	(G.NS.)	128	splot_col	(G..S.)
10576	S	(BA...)	80	_d1	(G.NS.)	80	th	(B.N.I)
80	SCAN_N	(G.NS.)	80	_d2	(G.NS.)	80	tth	(B.N.I)
80	SLIT_H	(G.NS.)	80	_d3	(G.NS.)			
80	SLIT_N	(G.N.C)	80	_f	(G.NS.)			

Memory usage is 41088 bytes.

16.FOURC>

The number preceding each name is the number of bytes of memory the variable consumes. All the global variables in the list above come from the standard start-up macro files. Those variables that begin with an underscore are internal to the standard macro package.

Flow Control

Flow control allows you to construct complex scans and other macros to control experiments and take data. The syntax of the flow control is very similar to standard C. For example, to list all the motor positions, you can use the following loop:

```

16.FOURC> for (i = 0; i < MOTORS; i++) {
17.more>     printf("Motor %d %-10s = %g\n", i, motor_name(i), A[i])
18.more> }
Motor 0  Two Theta = 3
Motor 1      Theta = 1.5
Motor 2      Chi = 0
Motor 3      Phi = 0

19.FOURC>

```

As in C, the `for` statement contains within parentheses three optional expressions separated by semicolons. The first expression is executed before entering the loop. The second is the test done before each pass of the loop — if it evaluates false, the loop is terminated. The third expression is executed at the end of each loop.

The conditional statements

```
if ( condition ) statement
```

and

```
if ( condition ) statement else statement
```

are also available. For example, to test whether a variable has been assigned a value, you could examine the return value of the built-in `whatis()` function (described on page 75 in the *Reference Manual*).

```

19.FOURC> if (whatis("DATAFILE")>>16&0x8000)
20.more> print "Warning, Data file is uninitialized!"
21.more> ;
Warning, Data file is uninitialized!

22.FOURC>

```

When there is a solitary `if` statement, a semicolon, extra newline or some other command must be read before the `if` statement will be executed, as it is not clear to the parser whether an `else` statement will follow. If there is an `else` statement, it must follow the `if` portion of the statement on the next line.

```

22.FOURC> if (whatis("DATAFILE")>>16&0x8000)
23.more> print "Warning, Data file is uninitialized!"
24.more> else
25.more> print "Data is being stored in", DATAFILE
Data is being stored in /usr/alan/default.dat

26.FOURC>

```

The `while` construction is also available. Usage is

```

26.FOURC> while (wait(0x22)) {
27.more>     getangles
28.more>     printf("%10.4f\r", A[tth])
29.more> }

30.FOURC>

```

As in C, `continue` and `break` statements may be used within loops. The statement `exit` can be used anywhere and always causes a jump back to the command level.

Macro Facility

One of `spec`'s most powerful features is its provision for defining macros. Through macros, you can simplify use of the diffractometer as well as determine the style and format of the data output. Through the macro facility, you can customize the environment to include any enhancements or specialized requirements for your experiment. Standard macro sets included in the `spec` package support conventional two-circle, four-circle and z-axis diffractometers along with some specialized liquid surface diffractometers. These macros control the measurement and the recording of experimental data and establish a standard format for ASCII data files.

An example of a simple macro that can be used to record a comment on the printer is

```

30.FOURC> def com '
31.quot>     on(PRINTER)
32.quot>     printf("$*\n")
33.quot>     off(PRINTER)
34.quot> '

35.FOURC>

```

Notice the prompt shows the program is expecting the quote to be closed. The variable `PRINTER` contains a string naming the printer device used to document the diffractometer operation. To use the above macro, type:

```

35.FOURC> com This is a comment.
This is a comment.

36.FOURC>

```

The text `This is a comment` is substituted in the `printf()` function for the symbol `$*` and is printed on both the screen and the printer.

Each argument following a macro call is available to the macro using `$1`, `$2`, ..., where `$1` refers to the first argument, and so on. Up to 25 arguments may be used. An argument is a string of characters separated by *white space* (spaces and tabs) or enclosed in single or double quotes, `$*` represents all the arguments, and `$#` is the number of arguments.

When a macro definition contains argument substitution, and you invoke that macro with more arguments than needed, the extra arguments you typed up to the next `;`, `}` or newline disappear. However, if the macro does not use argument substitution in its definition, text typed following the macro invocation is not thrown away.

To see what a macro contains, use the command `prdef` to print out the macro definition.

```
36.FOURC> prdef com
def com '
    on(PRINTER)
    printf("$*\n")
    off(PRINTER)
'
37.FOURC>
```

Notice that the form of the definition, if written to a file, would be suitable for reading back in as a macro definition.

The standard macro library is read automatically the first time you run `spec` or when you start the program with the `-f` flag. You can get a listing of all the currently defined macros with the command `lsdef`.

```
37.FOURC> lsdef
ALPHA      (4)  bug      (275)  hklscan  (1639)  qdo      (10)
AZIMUTH    (4)  ca       (182)  hkradial (334)   rplot    (10)
BETA       (4)  calcA    (7)    hlcircle (339)   rplot_res (112)
CEN        (10)  calcE    (7)    hlradial (333)   save     (432)
Escan     (1448) calcG    (7)    hscan    (139)   savegeo  (1121)
F_ALPHA    (4)  calcHKL  (7)    initdw   (42)    saveslits (99)
F_AZIMUTH  (5)  calcL    (8)    initfx   (39)    saveusr   (0)
F_BETA     (4)  calcM    (7)    initnec  (44)    savmac   (118)
F_OMEGA    (5)  calcZ    (7)    initoki  (39)    scan_head (5)
F_PHI      (5)  cat      (14)   inittemp (33)    scan_loop (5)
Fheader    (0)  cd       (11)   klcircle (339)   scan_move (5)
Flabel     (2)  ci       (177)  klradial (333)   scan_on   (192)
Fout       (2)  cl       (22)   kscan    (139)   scan_plot (0)
Ftail      (0)  com      (12)   l         (16)    scan_tail (5)
H          (4)  comment  (184)  less     (15)    set       (344)
K          (4)  config   (109)  lm       (539)   set_E     (314)
L          (4)  count    (6)    lp_plot  (674)   set_dial  (649)
LAMBDA     (4)  ct       (47)   ls       (13)    set_lm    (349)
OMEGA      (4)  cuts     (764)  lscan    (139)   setaz     (448)
Pheader    (0)  cz       (177)  lup      (419)   setlat    (764)
Plabel     (2)  d        (12)   mail     (16)    setmode   (927)
Pout       (2)  d2scan   (564)  measuretemp (1)    setmono   (368)
RtoT_0     (162) d3scan   (688)  mesh     (1221)  setplot   (1119)
RtoT_1     (162) debug    (212)  mk       (175)   setpowder (867)
RtoT_2     (161) do       (9)    move_E   (208)   setscans  (67)
RtoT_3     (163) dscan    (95)   move_em  (8)     setsector (1341)
TtoR_0     (160) dtscan   (143)  mv       (175)   setslit   (464)
```

TtoR_1	(160)	dumbplot	(334)	mvd	(192)	setslits	(469)
TtoR_2	(159)	end_reflex	(131)	mvr	(191)	settemp	(217)
TtoR_3	(161)	freeze	(456)	mz	(268)	show_cnts	(327)
_check0	(240)	g_aa	(4)	ned	(14)	showslits	(114)
_chk_lim	(266)	g_al	(4)	newfile	(1165)	showtemp	(133)
_cleanup2	(0)	g_bb	(4)	newmac	(270)	splot	(8)
_cleanup3	(0)	g_be	(4)	offd	(13)	splot_res	(150)
_count	(158)	g_cc	(4)	offp	(12)	startgeo	(45)
_do	(479)	g_chi0	(5)	offsim	(126)	starttemp	(225)
_getcut	(8)	g_chil	(5)	offt	(10)	startup	(245)
_head	(1363)	g_frz	(4)	ond	(12)	te	(105)
_hkl_lim	(84)	g_ga	(4)	onp	(11)	teramp	(458)
_hklline	(1254)	g_h0	(5)	onsim	(127)	th2th	(141)
_hklmesh	(638)	g_h1	(5)	ont	(9)	tscan	(732)
_loop	(413)	g_haz	(5)	or0	(549)	tw	(742)
_mo_loop	(175)	g_k0	(5)	or1	(551)	u	(10)
_mot	(124)	g_k1	(5)	or_swap	(320)	uan	(23)
_move	(37)	g_kaz	(5)	p	(8)	ubr	(31)
_pcount	(132)	g_l0	(5)	pa	(1301)	uct	(364)
_plot_scale	(392)	g_l1	(5)	pl	(176)	umk	(31)
_pmove	(131)	g_laz	(5)	pl_CFWMH	(10)	umv	(19)
_scanabort	(102)	g_mo_d	(5)	pl_COM	(10)	umvr	(20)
_setcut	(8)	g_mo_s	(5)	pl_FWHM	(10)	unfreeze	(39)
_settemp	(1)	g_mode	(4)	pl_LHMX	(10)	upl	(24)
_tail	(83)	g_om0	(4)	pl_MAX	(10)	uwm	(694)
_update1	(204)	g_om1	(5)	pl_MAXX	(10)	vi	(13)
_update2	(255)	g_phi0	(5)	pl_MIN	(10)	vt52_rplot	(1190)
_update4	(370)	g_phil	(5)	pl_MINX	(10)	vt52plot	(914)
_var	(171)	g_sect	(4)	pl_SUM	(10)	w	(12)
a2scan	(1172)	get_E	(74)	pl_SUMSQ	(11)	wa	(231)
a3scan	(1439)	getvar	(65)	pl_UHMX	(10)	waitall	(7)
add_reflex	(408)	gpset	(123)	pl_xMAX	(11)	waitcount	(7)
an	(173)	gt101_rplot	(1190)	pl_xMIN	(11)	waitmove	(7)
ansi_rplot	(1196)	gt101plot	(915)	plot	(40)	wh	(48)
ansiplot	(933)	h	(4)	plot_res	(435)	whats	(744)
ascan	(865)	help	(24)	prcmd	(42)	wm	(1232)
beep	(12)	hi	(7)	pts	(87)	yesno	(162)
beg_reflex	(283)	hkcircle	(339)	pwd	(11)		
br	(175)	hklmesh	(1095)	qcomment	(176)		

38.FOURC>

Each macro is listed as well as the number of characters in its definition. Some macros have zero length — their definitions are assigned during the course of an experiment. In the standard library, macros that are only used within other macros (and not meant to be referenced directly by the user) begin with an underscore.

Another macro handling command allows you to remove a macro definition.

```
38.FOURC> undef com
```

```
39.FOURC> prdef com
com: undefined.
```

```
40.FOURC>
```

There are several special macro names. If a macro named `cleanup` is defined, it will be automatically invoked whenever there is an error or `^C` interrupt. This macro can be defined to print a message, update a file, return motors to a starting position, etc. For example, in the standard macro library, something like the following is defined for the duration of a scan:

```
def cleanup '
    comment "Scan aborted after %g points." NPTS
    undef cleanup
,
```

Similarly, a macro named `cleanup1` can be defined, which behaves the same way. However, if `cleanup` exists, it will be run first.

Also `begin_mac`, `end_mac` and `prompt_mac` have special meaning. (text forthcoming ...)

Command Files

Macros are generally defined and maintained using the command file facility. In addition, sequences of experimental scans are often called up using command files. Command files are ASCII files of text, created with any of the UNIX text editors, and contain input just as it would be typed interactively. Command files are read line by line by `spec` when invoked with the functions `dofile()` or `qdofile()`. For example,

```
40.FOURC> dofile("spec.mac")
Opened command file 'spec.mac' at level 1.

FOURC.1> (Commands from file echoed as read ...)
```

The `.1` extension to the prompt indicates the level of nesting. Command files can be nested to five levels.

The function `qdofile()` is identical to the function `dofile()` except that the commands are not echoed as they are read.

```
40.FOURC> qdofile("spec.mac")
Opened command file 'spec.mac' at level 1.

41.FOURC>
```

If a file named *spec.mac* exists in your current directory, it is read as a command file each time you run *spec*. You can have private initialization code and macros in this file.

Two standard macros are defined to simplify reading command files. The macros *do* and *qdo* are normally used in place of the above functions. In addition to supplying the parentheses and quotation marks around the file name and recording the *do* command on the printer and in the data file, these macros also allow you to repeat the last command file when a dot is given as the argument:

```
41.FOURC> qdo .

qdo spec.mac
Opened command file 'spec.mac' at level 1.

42.FOURC>
```

Just as with keyboard input, comments can be included in a command file. Everything following a # up to a newline is ignored by *spec*.

Status and Help

Below is a summary of the diagnostic commands — most have been previously mentioned.

```
syms    lists built-in, global and local symbols.
lsdef   lists the names and sizes of macros.
prdef   prints out macro definitions.
lscmd   lists built-in keywords and functions.
```

All the above commands can take *pattern* arguments, employing the metacharacters *?* and ***. In a pattern argument, *?* stands for any single character, while *** stands for any string of characters, including the null string. For example, *lsdef ???* lists all the three-letter macros:

```
42.FOURC> lsdef ???

CEN (10)  com (12)  ned (14)  or0 (549)  qdo (10)  uct (364)
_do (479) lup (419)  ond (12)  or1 (551)  set (344) umk (31)
bug (275) mvd (192)  onp (11)  pts (87)  uan (23)  umv (19)
cat (14)  mvr (191)  ont (9)   pwd (11)  ubr (31)  upl (24)

43.FOURC>
```

Likewise, *lsdef *scan* shows all the macro names that end in *scan*.

```
43.FOURC> lsdef *scan
```

```
Escan (1448)  ascan (865)  dscan (95)   hscan (139)  tscan (732)
a2scan (1172) d2scan (564) dtscan (143) kscan (139)
a3scan (1439) d3scan (688) hklscan(1639) lscan (139)
```

```
44.FOURC>
```

A single * matches everything.

An on-line *help* facility exists to display files from *spec's help* directory:

```
44.FOURC> help
```

```
Help is available on the following subjects (type "h subject"):
```

```
386      config      files      geometry  macros     powder     sizes
ackno    counting    flow      gpib      news      print      spec
angles  debug      fourc    help     pdp      serial     syms
changes  do          functions history  plot      simulate   syntax
```

```
45.FOURC>
```

`help` (and simply `h`) are macros that use the built-in `gethelp()` function to print files contained in the `help` directory. The above listing is also produced by typing `gethelp("help")`. The command `gethelp("news")` is automatically executed each time the program starts up. New help files particular to a site may be added to the help directory.

UNIX Commands

The easiest way to write macro definitions is to use a standard text editor to create a command file, and the easiest way to get at the text editor is through the `unix()` function that spawns subshells.

```
45.FOURC> unix("vi macro.defs")
"macro.defs" 3 lines, 20 characters
```

```
46.FOURC> qdo macro.defs
Opened input file 'macro.defs' at level 1.
```

```
47.FOURC>
```

Any UNIX command may be spawned as in the above example. Because this is so useful, a macro has been written to simplify the syntax. You could type:

```

47.FOURC> u vi macro.defs
"macro.defs" 3 lines, 20 characters

48.FOURC> qdo macro.defs
Opened input file 'macro.defs' at level 1.

49.FOURC>

```

The `unix()` command (or the `u` macro) with no argument will spawn a subshell. You return to `spec` upon exiting the subshell. `spec` uses the shell environment variable `SHELL` or `shell`, if set, to select the type of UNIX shell. By default, `/bin/sh` is used. With arguments, `unix()` uses `/bin/sh` to execute the one-line command. For some common UNIX commands, macros such as the following are defined in the standard library.

```

def cat 'unix("cat $*")'
def ls  'unix("ls $*")'
def l   'unix("ls -l $*")'
def vi  'unix("vi $*")'

```

The working directory of `spec` can be changed as with the shell.

```

49.FOURC> cd data
Now in 'data'

50.FOURC>

```

The macro `cd` used above is defined using the built-in function `chdir()`. Only the working directory of the program `spec` is changed; the shell from which you started `spec` is not touched.

Moving Motors

A primary purpose of `spec` is to manipulate an X-ray diffractometer according to a calculated geometry. The automation of the angular settings is accomplished through the use of motor controllers interfaced to the computer. `spec` can be configured to control any number of motors.

As explained earlier, motor positions are referred to as dial positions and user positions. The diffractometer is operated in user positions. Dial positions are used to provide a stable point of reference. The two differ possibly by a sign and/or an offset. Dial positions should be set to always agree with the physical dials of the diffractometer motors. The user positions are then set in the line-up procedure of the

diffractometer. For example, they may be set to zero at the direct beam. The relations between the two positions are:

$$\begin{aligned} dial &= hardware_register / steps_per_unit \\ user &= sign \times dial + offset \end{aligned}$$

The *hardware_register* contains the value maintained by the stepper motor controller. The value of *steps_per_unit* is assigned in the hardware configuration file, as is *sign*. The latter must be chosen to agree with the conventions of the built-in geometry calculations.

The motor positions are often placed in the `A[]` array. The array is built-in and its elements can be used like any other variables. What makes it special, however, are the commands that use the array to convey the positions of the diffractometer motors. For example, the command `getangles` sets all of the elements of the `A[]` array to the current motor positions in user angles, while the `move_all` command sends the motors to the positions contained in `A[]` (in user angles). Typical usage is,

```
50.FOURC> waitmove      # Make sure no motors are active.
51.FOURC> getangles     # load A[] with user angles.
52.FOURC> A[0] = 3      # move motor #0 to 3.
53.FOURC> move_all      # start the move.
54.FOURC>
```

(The `#` symbols introduce comments.) It is important to first wait for any previous motions to complete. Then `getangles` is used to load the angle array with the current positions. Only the values for the motors to be moved are reassigned before using `move_all` to set the motors in motion.

A macro that would list the user positions of all the configured motors might be:

```
54.FOURC> def wa '
55.quot>     getangles
56.quot>     for (i = 0; i < MOTORS; i++)
57.quot>         printf("%9.9s = %g\n", motor_name(i), A[i])
58.quot> '

59.FOURC> wa
Two Theta = 3
Theta = 1.5
Chi = 0
Phi = 0

60.FOURC>
```

The `motor_name()` function returns the motor name assigned in the configuration file.

The motor positions are stored in three locations: in program memory, on the computer's hard disk and in the hardware registers associated with the motor controller. The program manipulates the angles in its memory. The values on the hard disk are updated every time a motor is moved but are only read when the program is initialized, or after the `reconfig` command is invoked. The controller registers count the actual number of steps moved and should be the true positions of the motors (unless a motor was switched off). Before each motor is moved, the controller registers are compared with program memory. If there are discrepancies, you are notified and asked which value is correct. The `sync` command can also be used to synchronize the controller registers with program memory. The angles can get out-of-sync by moving the motors with manual controls, by turning off the power to motor controllers or perhaps by a computer crash.

Although the motor controllers work in steps, it is much more convenient to use real units such as degrees (or, for linear motion, millimeters). The user and dial angles are in these units, converted from steps by the step-size parameters that are read from the configuration file.

The `chg_dial(motor, dial_angle)` function sets the dial register to *dial_angle* for one motor. The `chg_offset(motor, user_angle)` function sets the offset used to convert from dial positions to user positions for one motor. Often during the line-up procedure you will want to zero a particular angle:

```
60.FOURC> chg_offset(th, 0) # set motor theta to zero
61.FOURC>
```

The `set` macro includes the above and documents the change on the printer and in the data file.

```
61.FOURC> set th 0
Wed Aug 19 11:53:33 1987. Theta reset from 1.5 to 0
62.FOURC>
```

Dial and user settings may also be set by the `spec` administrator using the program `edconf`. See page 202 in the *Administrator's Guide* for further details.

Usually, diffractometer motions have a limit of travel, beyond which physical damage may occur (or a hardware limit switch tripped). Software limits therefore exist to prevent you from accidentally moving a motor out of range. The lower and upper limits for each motor are contained in internal arrays accessed through the

`set_lim(motor_number, low, high)` and `get_lim(motor_number, flag)` functions. With the latter, the lower limit is returned if `flag` is less than zero, otherwise the upper limit is returned.

If a `move_all` command would take any motor outside of its limits, an error message is printed and no motors are moved. The limit values are stored in dial angles since they correspond to physical limitations to the motion. The limit values are therefore preserved as the user-angle offsets are changed. The `set_lm` macro can be used to set a single motor's limits:

```
62.FOURC> set_lm tth 0 360
```

```
Two theta limits set to 0 360 (dial units).
```

```
63.FOURC>
```

The angle arguments to the macro `set_lm` are given in user angles. (In this example, dial and user angles are the same.)

Diffraction Geometry

You can operate a two-circle diffractometer in terms of angles alone. However, for a four-circle diffractometer (and others such as the z-axis or liquid-surface diffractometers) it makes more sense to work in three-dimensional reciprocal space coordinates. It is therefore necessary to be able to calculate angles according to the diffractometer geometry.

`spec` is designed to accommodate a variety of diffractometer configurations. The particular calculations are contained in geometry code (the source for which is included in the standard `spec` package.) accessible through the `calc()` function. The arguments to `calc()` determine the particular code that is called. For example, a `calcA` macro is defined as `calc(1)`. Its purpose is to load the `A[]` array with the angles corresponding to the current reciprocal space coordinates. The four-circle configuration represents the three reciprocal space coordinates as the first elements of the built-in array `Q[]`. For convenience, the following definitions are made:

```
def H 'Q[0]'  
def K 'Q[1]'  
def L 'Q[2]'
```

Thus, to move to a position in reciprocal space such as the point [100], the appropriate commands would be

```
63.FOURC> H = 1; K = L = 0; waitmove; getangles; calcA; move_all
```

```
64.FOURC>
```

Whenever the `move_all` command is used, it is important that the `A[]` array contain the current motor positions for all motors except the ones to be moved. In the above example, the `getangles` command loads `A[]` with the current positions after the `waitmove` ensures all motors have stopped. The `calcA` changes the appropriate elements of the `A[]` array and the `move_all` starts the motors.

Often, you might change a single angle, or several angles, and then wonder where in reciprocal space the diffractometer is set. The `calcHKL` macro will take the positions in the `A[]` array and set the variables `H`, `K`, and `L` to the calculated coordinates. For example:

```
64.FOURC> waitmove; getangles; calcHKL; print H, K, L
1 0 0

65.FOURC>
```

The command `getangles` loads the `A[]` array with the current positions.

Counting

Another important function of the diffractometer program is to measure the scattered X-ray intensities. `spec` supports several types of timers, scalers and multichannel analyzers (MCAs). Timers control the count time. Scalers count detected photons. MCAs accumulate many channels of counts and are used with energy-dispersive detectors and positional-sensitive detectors.

To count the number of X rays incident on the detector per second, the counting hardware must be able to accumulate detector counts accurately within a fixed time period. The scaler hardware is gated by a clock that operates independently of the computer. Thus, the response time of the computer to interrupts (real-time events) does not affect the accuracy of the count. `spec` programs and starts the clock and senses when the clock period, and hence the counting, has ended. `spec` can then read the contents of the scalers and save the measurement in a data file.

Clearing the scalers and starting the clock is accomplished by the function `tcount(seconds)`. To count for one second, type:

```
65.FOURC> tcount(1)

66.FOURC>
```

The contents of the scalers are accessed through the built-in `S[]` array. The hardware scalers are read and their contents loaded into the scaler array by the `getcounts` command. A second, associated string array, `S_NA[]`, is defined in the standard macros and identifies each scaler:

```

66.FOURC> getcounts; printf("%s = %g\n%s = %g\n%s = %g\n",\
67.cont> S_NA[0], S[0]/1000, S_NA[1], S[1], S_NA[2], S[2])
seconds = 1
monitor = 347
detector = 35031

68.FOURC>

```

The first scaler, labeled `seconds`, is usually fed a 1 kHz signal, so it actually tracks milliseconds and is therefore divided by 1000. The number of scalers available depends on the particular hardware and the number of detectors and monitors used. The default scaler channel numbering for the first three scalers puts a 1 kHz time signal in scaler 0, monitor counts in scaler 1 and detector counts in scaler 2.

You can also count to a fixed number of monitor pulses, rather than to a fixed time period.

```

68.FOURC> mcount(1e4)

69.FOURC> getcounts; printf("%s = %.1f\n%s = %g\n%s = %g\n",\
70.cont> S_NA[0], S[0]/1000, S_NA[1], S[1], S_NA[2], S[2])
seconds = 28.8
monitor = 10000
detector = 1.00954e+6

71.FOURC>

```

Counting is asynchronous, i.e., the `tcoun`() and `mcoun`() functions return immediately after starting the clock. They do not wait until the counting period is over. Use the `wait`() function to determine when counting is finished.

A useful macro has been written to count and print the scaler contents:

```

71.FOURC> ct 5

Thu Aug 20 19:11:51 1987
Seconds = 5 Detector = 175103 (35020.6/s) Monitor = 1730 (346/s)

72.FOURC>

```

If the argument is omitted, a default count time (stored in the global variable `COUNT`) is used. A positive argument to `ct` signifies seconds; a negative argument signifies monitor counts.

```

72.FOURC> ct -10000

Thu Aug 20 19:13:42 1987
Seconds = 28.3 Detector = 1.0434e6 (36869.3/s) Monitor = 10000 (353.36/s)

73.FOURC>

```

CAMAC, GPIB and Serial

Besides the built-in hardware support for moving motors and counting using CAMAC (IEEE-583), GPIB (IEEE-488) and serial (RS-232C) interfaces, `spec` provides generalized input/output support over these interfaces from the command level.

The `ca_get(device, subaddress)` command will return the contents of the addressed module register ($F=0$ in the standard *FNA* CAMAC command code), while `ca_put(data, device, subaddress)` will write the 24-bit value `data` to the addressed module register ($F=16$). The `device` argument is the I/O module index from the configuration file and can be 0, 1 or 2. The CAMAC slot number of the module is set in the configuration file. The `subaddress` argument is the module's subaddress (the *A* in the *FNA*).

`spec` allows you to send a string of characters to a GPIB instrument at any GPIB address and to read a string of characters from any instrument. When reading characters, `spec` will accept a string terminated by either a newline or by carriage return—newline, or you can specify the number of bytes to be read. For example, to initialize a particular voltmeter having GPIB address 12, you would issue the command:

```
73.FOURC> gpib_put(12, "D0R0Z0B0T0K1M0G1X")
74.FOURC>
```

That instrument might then be read with:

```
74.FOURC> {k_ohms = gpib_get(12); print k_ohms}
100.024
75.FOURC>
```

The command

```
75.FOURC> x = gpib_get(12, 4)
76.FOURC>
```

would read 4 bytes from device 12 and not look for terminators.

When sending strings using `gpib_put()`, you cannot send null bytes. Usually a device that requires null lower order bits in a data byte will ignore the high order (parity) bit. In this case, you can usually set that highest bit to avoid sending a null byte.

The `ser_get(device, n)` and `ser_put(device, string)` functions access the serial interface, where `device` is the index from the configuration file and can be 0, 1 or 2. In `ser_get()`, `n` is the most number of bytes to read. The function will return after reading one line (terminated by a newline or carriage return) from the device, even if

the number of bytes is less than n . In `ser_put()`, *string* contains the characters to be written.

Using spec with C-PLOT and Other UNIX Utilities

Standard Data File Format

The format of the data files used by `spec` is set at the macro level. The files are ASCII, so they can be easily manipulated by other UNIX utilities such as `grep`, `sed`, `awk` or any of the editors. The format of files produced by the standard macros is described here.

When opened with the `newfile` macro, the following header is written to initialize the data file:

```
#F /tmp/data
#E 729994936
#D Wed Feb 17 19:22:16 1994
#C cu 110 User = bill
#O0 Two Theta Theta Chi Phi
```

Information or control lines begin with a `#` character, with the character in the second column indicating the type of information that follows. The first line of the data file header contains the name by which the file was opened. The next line is the number of seconds from the UNIX epoch as returned by the `time()` function. In the data that will follow, each scan point will include a field containing the number of seconds elapsed since that file creation time. Next in the header is a line containing the date as returned by the `date()` function, then a comment line. That is followed by a line containing all the motor names in use. Each motor name is separated from the other by two spaces. What will then follow will be various comment lines created by the `comment` macros, user defined entries and scan data.

Each scan has a header that looks like the following, always beginning with a blank line:

```
#S 1 hklscan 0.9 1.1 0 0 0 0 20 1
#D Wed Feb 17 19:25:55 1994
#T 1 (Seconds)
#G0 0 0 0 0 0 1 0 0 0 0 0
#G1 1 1 1 90 90 90 3 3 3 90 90 90 1 0 0 0 1 0 60 30 0 0 0 0 60 30 0 -90 0 0 0
#Q .9 0 0
#P0 29.745 29.745 90 0
#N 7
#L H K L Epoch Seconds Monitor Detector
```

The first line of the scan header contains the scan number followed by the basic scan

name and its arguments. (Scans invoked as `hscan` or `kscan` execute the basic `hklscan`, scans such as `lup` or `dscan` execute `ascan`, etc.) The next line is the date and time the scan was started. Following that, the `#T` control line indicates that the scan was counting to time and for what duration at each point. A `#M` would indicate counting to monitor.

The next lines give information that describe the diffractometer configuration at the start of the scan. Following `#G0` are the current values of the four-circle parameters, which are defined four-circle macro file, `macros/fourc.src`. Following `#G1` are the current parameters describing the crystal lattice and orientation matrix. The identification of these parameters is in the macro file `macros/ub.mac`. The `#Q` line gives the H, K, L coordinates at the start of the scan, while the `#P0` line gives the motor positions at the start of the scan with each column corresponding to the motor names in the `#00` line of the scan header. The `7` after `#N` indicates there will be seven columns of data in the scan, and the `#L` line gives the names for each column, each name separated from the other by two spaces. The detector counts are always placed in the last column, preceded by the monitor counts, if counting to time, or the elapsed time for that data point, if counting to monitor. The `EPOCH` column has the number of seconds elapsed from the time of the `#E` at the start of the file.

Following the scan header is the scan data. These are just lines of space-separated numbers that correspond to the column headers given with `#L`. Intervening `#C` comment lines may lie within the rows of data if, for example, the scan was aborted and then restarted with the `scan_on` macro. Otherwise, the data continues until the next non-comment control line or blank line.

You can develop your own programs and scripts to extract data from the `spec` data file, or you may want to use the `scans.4` user function that is part of the C-PLOT package, or the stand-alone `scans` program that is based on `scans.4`.

Scans.4

The C-PLOT user function `scans.4` reads in files of ASCII data according to a modest set of conventions. In particular, `scans.4` manipulates the X-ray scattering data from the `spec` data files, doing scan averaging, background subtraction, data normalization and error bar calculation. The C-language source code to `scans.4` is always available on your system, is liberally commented and should always be consulted if there is any question as to what operations are being done on your raw data points.

The `scans.4` user function can be invoked from C-PLOT either as,

```
PLOT-> fn scans.4
```

or

PLOT-> **fn scans.4 options scan-numbers**

Use the second form when running from command files. The possible options are:

.	Use same options as last time.
-i	Initialize, used to start up function and return.
-f <i>filename</i>	Select scan file name.
-p	Print scan file contents.
+e or -e	Calculate (or don't) error bars from statistics.
+s or -s	Sort and merge (or don't) data by x values.
+d or -d	Collect (or don't) 3 columns of data.
+r or -r	Rerange (or don't) plot axis for each new data set.
+S or -S	Retrieve scans by scan (or file) number.
+v or -v	Print (or don't) each line of scan file (verbose).
+n or -n	Normalize (or don't) data points.
-m	Turn on +n flag and normalize to monitor counts.
-t	Turn on +n flag and normalize to time.
x=#	Set column for x values.
y=#	Set column for y values.
z=#	Set column for z values and turn on +d flag.

The default options are:

-f data +esSn -rvd -m x=1 y=-1

Retrieving Scans By Scan Number or File Position Number

Scans can be retrieved by entering either the scan number (option +s, default) or the file position number (option -s). Scan numbers are determined by the #s lines in the file. The file position number is the sequence position of the scan in the file, irrespective of scan number. Normally, the scan number and the file position number are the same.

When selecting by scan numbers, if there is more than one scan with the same number, the last one is retrieved. Specify which instance of a repeated scan to retrieve by using the *scan.sub* syntax. For example, selecting scan 10.3 retrieves the third instance of scan number 10.

Negative numbers count back from the end of the file and are always considered to be file position numbers. For example,

```
fn . -1
```

will always return the last scan in the file. You can read in multiple scans by giving several scan numbers as arguments. You can read in a group of consecutive scans with

```
fn . 3-7 10-14
```

The above reads in scans 3 through 7 and 10 through 14.

Merging Scans and Background Subtraction

The default `+s` option causes the data points from all the scans read in to be sorted by x values and data points with the same x value averaged. If data is to be normalized and/or error bars calculated, the appropriate weight is given to the count time for each point.

Scan numbers that end with the letter `b` are used as background scans. The sort-and-merge flag should be in effect when using background scans. Entering

```
fn . 12b 13 14b 15b 16 17b
```

for example, or

```
fn . 13-15 16-27b
```

will merge the data from the background scans with the data scans, subtracting the background counts from the data counts at each x value. When doing merging and background subtraction, the x values must be identical for the data points to be merged.

File Conventions

scans.4 only uses some of the control lines in the standard spec data files described earlier.

The control conventions used by *scans.4* are:

#S <i>number</i>	Starts a new scan. <i>number</i> is the user's numbering scheme.
#M <i>number</i>	Indicates data was taken counting to <i>number</i> monitor counts.
#T <i>number</i>	Indicates data was taken counting for <i>number</i> seconds.
#N <i>number</i>	Indicates there are <i>number</i> columns of data.

The following control lines do nothing, although they will be printed to the screen while reading a scan.

#C <i>comment</i> ...	Conventionally a comment.
#D <i>date</i>	Conventionally the date.
#L <i>lab1 lab2</i> ...	Conventionally data column labels, with each label separated by two spaces.

Data Columns

By default, *x* values are taken from the first column, *y* values from the last column. Monitor counts are always taken from the column prior to the *y* column.

When entering column numbers, a negative number counts backwards from the last column. If the column for *x* is zero, the value put in for *x* is just the row number of the point with respect to the start of the data for the current scan.

More Details

After *scans.4* reads and indexes a data file, it remembers the file length. If you answer affirmative to the `Change modes?` query, *scans.4* will add to the index if the file has lengthened.

If you give a dot (.) as the command line argument or in response to `Scans/options` query, the previous argument or option string will be used. That is, the string is remembered, not the options chosen interactively using `Change modes?` For instance, if you enter a long sequence of scan numbers and read in the scans, then change something via `Change modes?`, you can simply enter a dot in response to `Scans/options` and recover the previous sequence of scan numbers.

When you do enter a string of flags and scan numbers, the modes set by the flags only apply to the scans that follow the flags and not the preceding scans.

The Index File

Indexing a long ASCII data file to find at what byte offset each scan begins takes time. Once *scans.4* has indexed a file, it saves the index information in a binary-format *index* file. The name of the index file is formed by appending *.I* to the original data-file name. As long as the index file is more recent than the data file, *scans.4* will take the index information from the index file.

Normalization and Error Bars

The values returned as error bars are those due to counting statistics (the square root of the number of counts). When the counts are derived from the algebraic combination of detector, background and monitor counts, the error bars are calculated using the appropriate *propagation of errors* formalism.¹

Contents

A *contents* program is included in the *spec* package. The program tries to print summary scan information from *spec*'s standard data files. Usage is

```
contents [ options ] file1 [ file2 ... ]
```

Current options are

-o <i>output</i>	Name of output file, otherwise standard output is used.
-s <i>start</i>	Starting number for first scan of first file.
-p <i>page</i>	Lines per page.
-d	Send control codes appropriate for DecWriter II.
-c	Print #C comment lines.

¹P.R. Bevington, *Data Reduction and Error Analysis for the Physical Sciences*, (McGraw-Hill, New York, 1969), p. 64.

Showscans

The *showscans* program is a shell script that serves as a front end for an *awk* script called *show.awk*, which is installed in the *spec* auxiliary file directory. The *awk* script prepares a C-PLOT command file that will be automatically run to make summary plots of scans from *spec* data files. Twelve scan plots are placed on each page. The text is very small, so a high resolution display device is recommended.

Usage is

```
showscans [options] [file_options] file [[file_options] file ...]
```

where the global *options* relate to C-PLOT commands and are:

-eb <i>error_bars</i>	Error bar mode (default is 0, off).
-ft <i>font</i>	Font code (default is font 2).
-sy <i>symbol</i>	Symbol code (default is L, a line).
-zi <i>filter</i>	Graphics filter (default is <i>psfilter</i>).
-w	Wait for keyboard <return> after each page
-x	Shorthand for -zi @x11 -rotate@ -w

The *file_options* apply only to the next file on the command line and are:

-f <i>from_scan_number</i>	Starting scan number for the following file.
-t <i>to_scan_number</i>	Ending scan number for the following file.

In no scan numbers are specified, all the scans within a file will be processed.

File names that end with *.I* are ignored as they are assumed to be *scans.4* index files, allowing you to use metacharacters to specify filenames a little more freely.

The *awk* script recognizes all the standard scan headings produced by the standard *spec* scan macros and uses that information to select the *x*-axis label and the column for the independent variable in the data file. The title of each of the small plots is the starting *Q* vector of the displayed scan, taken from the #*Q* line of the data files.

If the monitor counts are zero for any point in a scan, that scan is not plotted. The *show.awk* file can be edited to disable that feature.

REFERENCE MANUAL

Introduction

The material contained in this part of the documentation describes most of the built-in features of the `spec` program, that is, those parts of `spec` that are compiled into the code and cannot be changed at the installed site. These features include the user interface and the general hardware support, but do not include application-dependent features, such as geometry code and macro libraries.

Internal Structure Of `spec`

This section briefly explains `spec`'s internal structure to give an overview of how it is constructed.

First, consider how a user's input gets interpreted. The initial translation of characters typed at the keyboard (or read from a command file) is done by the *input preprocessor*, which keeps track of the input sources and handles command recall (or history) substitution.

The input text is then broken into *tokens* by the *lexical analyzer*. Tokens represent the different kinds of input objects, such as *keywords*, *operators*, *variables*, *number constants* and *string constants*. When the lexical analyzer recognizes a predefined *macro* name, its definition, possibly with argument substitution, is pushed back onto the input stream to become further input to the lexical analyzer.

The *parser* in `spec` repeatedly calls the *lexical analyzer* to obtain tokens. The *parser* contains a set of *grammar rules* that determine acceptable sequences of tokens. A *syntax error* occurs when input violates these rules. When enough input is read to satisfy a complete set of rules, the parser returns a parse tree (or mini-program) to the *executor*. The executor code then steps through the parse tree, calling the appropriate internal functions for each action.

The macro-definition command `def` is an exception to the above rules. As soon as the parser recognizes the `def` command sequence, the macro name and its definition are saved and made available to the lexical analyzer, even while the parser is still building the mini-program. A different command, `rdef`, defers storing the macro definition until the mini-program is executed. The `rdef` command is useful when some flow control logic needs to be run to decide what definition to assign to the macro.

Understanding the difference between the parse phase and execution phase of `spec` is important. Each time the command-level prompt is given, a new parse tree will be created. If several semicolon-separated commands are given on the same line, a separate parse tree will be created for each. However, curly brackets can be used to

group any number of lines together to form just one parse tree. A significant consequence of the parse tree mechanism is limitation of the scope of a nonglobal variable to the statement block in which it is referenced.

`spec` may detect error conditions during each of the phases described above. Most of these errors (and the interrupt character, usually a `^c`) reset `spec` to the command-level prompt.

`spec` manages a memory arena that is used as storage for parse trees, macro and variable symbol tables and command history. In earlier versions of `spec`, this memory arena was fixed in size. `spec` now allows the arena to grow up to the maximum data-space size allowed by the system. The `memstat` command displays `spec`'s memory usage.

`spec` also manages some auxiliary files. The user's state file stores the important parts of the memory arena. This file is used to let the user exit and restart `spec` without losing the program's state. Another file stores the current data points and data group configuration, so that they, too, will be available after exiting and restarting `spec`. These files are placed in the *userfiles* subdirectory of each configuration's auxiliary file area.

Syntax Description

Comments

A `#` introduces a comment. Everything following a `#` on an input line is ignored (unless the `#` is within a string). Comments are retained in macro definitions and are counted in the macro length.

Identifiers

An identifier is a name — it can be a variable name, a macro name or an array name. An identifier may begin with the letters `a-z`, `A-Z` or `_` (underscore). The remaining letters in a name may be those characters or the digits `0-9`. There is no limit to the number of characters in a variable name. In the syntax rules described later, such names are represented by the term *identifier*.

Arrays

Arrays are formed using square brackets, as in *identifier* [*expression*]. *spec* has two kinds of arrays: associative arrays and data arrays.

Associative Arrays

Associative arrays have been around since the earliest versions of *spec*. For associative arrays, the array index *expression* can be any numeric or string-valued constant or expression, although the internal code always uses the string value. Thus `A["12"]` refers to the same element as `A[12]`. Elements of associative arrays can contain numbers or strings. Currently, all built-in global arrays, such as those that hold motor positions and scaler counts, are associative arrays.

Associative arrays may also be two dimensional, as of release 4.01.11. Internally, the two array indices for each element are stored as a single string formed using the string value of the first index, followed by the character `\034`, followed by the string value of the second index. You can access such a 2D array element using a single index constructed according to the above recipe. That is, `arr["list"]["one"]` refers to the same item as `arr["list\034one"]`.

Data Arrays

The second kind of array is the *data* array introduced in *spec* release 4.00. The data arrays store and manipulate their contents far more efficiently than associative arrays, and also use a far more intuitive syntax than the data group facility (described later) for which they are an alternative. For data arrays, the array index must be a number, and array elements are numbered consecutively starting at zero. Data arrays must be explicitly declared using the `array` keyword. The data values of data-array elements are not currently saved in the user's state file, unlike associative array elements and the contents of data groups.

Keywords

The following names are reserved, being either grammar keywords or the names of built-in commands or functions. New reserved names may be added. The list can be obtained using the built-in command `lscmd`. Parentheses after a name indicate a function.

<code>acos()</code>	<code>data_get()</code>	<code>gethelp()</code>	<code>motor_num()</code>	<code>sleep()</code>
<code>array</code>	<code>data_grp()</code>	<code>getline()</code>	<code>motor_par()</code>	<code>sock_io()</code>
<code>array_dump()</code>	<code>data_info()</code>	<code>getval()</code>	<code>move_all</code>	<code>spec_par()</code>
<code>array_fit()</code>	<code>data_nput()</code>	<code>global</code>	<code>move_cnt</code>	<code>split()</code>
<code>array_op()</code>	<code>data_pipe()</code>	<code>gpib_cntl()</code>	<code>off()</code>	<code>sprintf()</code>
<code>array_pipe()</code>	<code>data_plot()</code>	<code>gpib_get()</code>	<code>on()</code>	<code>sqrt()</code>
<code>array_plot()</code>	<code>data_put()</code>	<code>gpib_poll()</code>	<code>open()</code>	<code>strand()</code>
<code>array_read()</code>	<code>data_read()</code>	<code>gpib_put()</code>	<code>plot_cntl()</code>	<code>sscanf()</code>
<code>asc()</code>	<code>data_uop()</code>	<code>history</code>	<code>plot_move()</code>	<code>stop()</code>
<code>asin()</code>	<code>date()</code>	<code>if</code>	<code>plot_range()</code>	<code>string</code>
<code>atan()</code>	<code>dcb()</code>	<code>image_get()</code>	<code>port_get()</code>	<code>substr()</code>
<code>atan2()</code>	<code>def</code>	<code>image_par()</code>	<code>port_getw()</code>	<code>syms</code>
<code>bcd()</code>	<code>deg()</code>	<code>image_put()</code>	<code>port_put()</code>	<code>sync</code>
<code>break</code>	<code>delete</code>	<code>in</code>	<code>port_putw()</code>	<code>tan()</code>
<code>byte</code>	<code>dial()</code>	<code>index()</code>	<code>pow()</code>	<code>tcount()</code>
<code>ca_cntl()</code>	<code>dofile()</code>	<code>input()</code>	<code>prdef</code>	<code>time()</code>
<code>ca_fna()</code>	<code>double</code>	<code>int()</code>	<code>print</code>	<code>tty_cntl()</code>
<code>ca_get()</code>	<code>else</code>	<code>length()</code>	<code>printf()</code>	<code>tty_fmt()</code>
<code>ca_put()</code>	<code>eprint</code>	<code>local</code>	<code>qdofile()</code>	<code>tty_move()</code>
<code>calc()</code>	<code>eprintf()</code>	<code>log()</code>	<code>quit</code>	<code>ubyte</code>
<code>cdef()</code>	<code>exit</code>	<code>log10()</code>	<code>rad()</code>	<code>ulong</code>
<code>chdir()</code>	<code>exp()</code>	<code>long</code>	<code>rand()</code>	<code>undef</code>
<code>chg_dial()</code>	<code>exp10()</code>	<code>lscmd</code>	<code>rdef</code>	<code>unglobal</code>
<code>chg_offset()</code>	<code>extern</code>	<code>lsdef</code>	<code>read_motors()</code>	<code>unix()</code>
<code>close()</code>	<code>fabs()</code>	<code>mca_get()</code>	<code>reconfig</code>	<code>user()</code>
<code>cnt_mne()</code>	<code>file_info()</code>	<code>mca_par()</code>	<code>return</code>	<code>ushort</code>
<code>cnt_name()</code>	<code>float</code>	<code>mca_put()</code>	<code>savstate</code>	<code>vme_get()</code>
<code>cnt_num()</code>	<code>fmt_close()</code>	<code>mca_sel()</code>	<code>ser_get()</code>	<code>vme_get32()</code>
<code>constant</code>	<code>fmt_read()</code>	<code>mca_sget()</code>	<code>ser_par()</code>	<code>vme_move()</code>
<code>continue</code>	<code>fmt_write()</code>	<code>mca_spar()</code>	<code>ser_put()</code>	<code>vme_put()</code>
<code>cos()</code>	<code>for</code>	<code>mca_sput()</code>	<code>set_lim()</code>	<code>vme_put32()</code>
<code>counter_par()</code>	<code>fprintf()</code>	<code>mcounct()</code>	<code>set_sim()</code>	<code>wait()</code>
<code>data_anal()</code>	<code>get_lim()</code>	<code>memstat</code>	<code>shared</code>	<code>whatis()</code>
<code>data_bop()</code>	<code>getcounts</code>	<code>motor_mne()</code>	<code>short</code>	<code>while</code>
<code>data_dump()</code>	<code>getenv()</code>	<code>motor_name()</code>	<code>sin()</code>	<code>yesno()</code>
<code>data_fit()</code>				

Numeric Constants

Numeric constants can be integers or floating point numbers. Integer constants are considered octal if the first digit is 0. Valid digits in the rest of the constant are 0 through 7. A hexadecimal constant begins with 0x or 0X and is followed by digits or the letters a or A through f or F, which have values 10 through 15. Otherwise, a sequence of digits is a decimal constant.

Floating-point constants have an integer part, a decimal point, a fraction part, an e or E and an optionally signed exponent. The integer part or the fraction part, but not both, may be missing. The decimal point or the e and exponent, but not both, may be missing.

Octal and hexadecimal constants can have values from $\pm 2^{31}$. As spec stores number values internally in double-precision format, the range of other integer and floating constants is determined by the range of double-precision numbers allowed by the host computer's compiler and architecture. To conserve resources, spec by default stores the numbers used in the data group facility in single-precision format. Double-precision can be selected at the spec administrator's option, though, when spec is installed.

The following are valid numeric constants.

```
65535  0177777  0xFFFF
+1066  1.066e3  1.066e+3
```

String Constants

Strings are delimited by pairs of single or double quotes. The following escape sequences, introduced by a backslash, can be included within strings to represent certain special characters:

<code>\a</code>	attention, audible alert (bell)
<code>\b</code>	back space
<code>\f</code>	form feed
<code>\n</code>	newline
<code>\r</code>	carriage return
<code>\t</code>	horizontal tab
<code>\\</code>	backslash
<code>\'</code>	single quote

<code>\"</code>	double quote
<code>\ooo</code>	octal code
<code>\[xx]</code>	tty control code

Tty control codes are only recognized when embedded in strings passed to `spec`'s built-in functions `tty_move()` and `tty_fmt()`. The recognized strings for `xx` are described in the description of the `tty_cntl()` function on page 87.

For any other character `x`, `\x` is just that character. The sequence `\ooo` represents one to three octal digits that have the ASCII value of a single character. For example, `\033` represents the *escape* character.

A character string can be continued over more than one line by using a `\` at the end of a line. On the other hand, new lines not preceded with a `\` are inserted literally into the string.

String Patterns and Wild Cards

For the commands `lscmd`, `lsdef`, `prdef` and `syms`, if the characters `?` or `*` appear in the arguments, the argument is taken as a pattern. Only information about those commands, macros or symbols that match the pattern is displayed. In forming the pattern, the character `?` matches any single character, while `*` matches any string of characters, including the null string.

Command Recall (History)

`spec` standard history feature lets you recall previous commands¹. Examples of the recognized syntax are:

¹`spec` is usually installed with the optional GNU *readline* history library, at the `spec` administrator's option. With the *readline* library, the history syntax is greatly expanded. In addition, features such as command-line editing, command completion and tilde expansion become available. See the on-line *readline* help file for detailed information.

!!	Redo the previous command.
!14	Redo command number 14.
!-2	Redo the second to previous command.
!asc	Redo the last command that began with <code>asc</code> .
!asc -10000	As above and append -10000 to the command.
history	List the last 1000 commands.
history <i>N</i>	List the last <i>N</i> commands.
history - <i>N</i>	List the last <i>N</i> commands in reverse order.

The command number is prepended to the `spec` prompt as an aid in using the history feature. Only commands typed at the keyboard are remembered for history, and no more than 1000 commands are retained. The history feature cannot be used in command files.

Command recall must occur at the beginning of a line, although initial white space is allowed. Text may follow the command-recall word to extend that command.

Appending `:s/left/right/` to a recalled command will modify the first occurrence of the string `left` in the recalled command to the characters `right`. The delimiter of the left and right strings may be any character. The final delimiter is optional. If `left` is empty, the last entered left string is used as the pattern to match.

In addition, `^left^right^` at the start of a line is shorthand for `!-1:s/left/right/`. In this case, the circumflex (`^`) must be used as the delimiter. The final delimiter is optional.

The history is saved in the state file when exiting `spec`. Restarting `spec` reads in the saved history from the state file.

Starting Up

When you run `spec`, you invoke it using a name such as *fourc*, *twoc*, *surf*, *spec*, etc. That name determines the kind of geometry code that will be available and which macro and configuration files in the auxiliary file directory will be used, as explained below.

The following command line options are recognized by `spec`:

- f** – Fresh start. All symbols are set to their default values and the standard macros are read to establish the default state.
- F** – Clean and fresh start. All symbols are set to their default values but no command files are read and no macros are defined. Only the built-in commands are available.

- s** – Simulation mode. No hardware commands are issued. Simulation mode cannot be turned off after entering the program.
- y** – Yes, change motor controller registers initially if they disagree with the *settings* file. Normally, *spec* requires you to confirm such a change. This flag would be useful if you know controller power had been turned off, and the controller's current positions should be updated with the software positions.
- o *option=value*** – Initializes the `spec_par()` *option* to *value*. Type `spec_par("?")` from the *spec* prompt to see the available options.
- p *fd pid*** – Indicates that *spec* input is coming from a pipe from another program. The argument *fd* is the file descriptor that *spec* should use for standard input. The argument *pid* is the process ID of the spawning process. If *fd* is zero, *spec* will not re-echo input from the file descriptor to *spec*'s standard output.
- q** – Indicates that *spec* should operate in quiet mode and allow output to all devices to be turned off. This option is only valid when used with the **-p** option.
- g *geometry*** – Force loading of macro files and activation of geometry calculations for the specified geometry, while using the configuration files taken from the name by which *spec* is invoked.
- N *my_name*** – Use *my_name* instead of the name by which *spec* was invoked to establish the command prompt and the name of the directory in *SPECD* in which the configuration-dependent files exist. This command also sets the geometry to *my_name*. Follow this option with the **-g** option to choose a different name for the geometry.
- u *user*** – Use *user*'s last saved state as the current user's initial state.
- t *tty*** – Use the current user (or *user*'s) last saved state from the terminal specified by *tty* as the initial state. The terminal can be specified as **-t /dev/tty01** or **-t *tty01***. Pseudo-tty names, such as */dev/ttyp0*, */dev/ttyp1*, etc., are saved as */dev/ttyp#*, since there is no special significance to the number.
- T *fake_tty*** – This option creates a user state associated with *fake_tty*, which may be any name. This option allows you to bypass the locking feature that prevents multiple instances of *spec* to be started by the same user from the same terminal.
- d *debug*** – Sets the initial value of the debugging variable `DEBUG` to *debug*, which may be either in decimal or hexadecimal (with a leading 0x) format. The available debugging categories are described on page 57.
- D *direc*** – Use *direc* as the auxiliary file directory, instead of the compiled-in name (usually */usr/local/lib/spec.d*) or the value of the *SPECD* environment variable.

In some installations, `spec` is installed as a set-user-id root process, to allow certain calls that allow privileged access to hardware device registers. The first thing `spec` does on start up is to set the effective user and group ids to that of the real user, so there is no danger of the user spawning subshells or creating files as root. The root effective id is only used for the duration of the calls that enable the privileged access.

`spec` then performs other initialization tasks, including obtaining values for its internal variables `DISPLAY`, `GTERM`, `TERM`, `HOME` and `PAGER` from variables of the same name in the process environment. It also obtains the value of the `SHELL` environment variable for use with the `unix()` function.

`spec` then reads the hardware configuration from the appropriate *config* file from the auxiliary file directory. The path name of that file is *SPECD/spec/config*, where *SPECD* is the auxiliary file directory, established when `spec` is installed (or by the `-D` invocation flag, or by the *SPECD* environment variable), and *spec* is the name by which `spec` is invoked.

The first time a user starts `spec`, up to seven macro files are automatically read. The path names of these files are

- SPECD/standard.mac*
- SPECD/geom.mac*
- SPECD/spec/geom.mac*
- SPECD/site_f.mac*
- SPECD/site.mac*
- SPECD/spec/conf.mac*
- ./spec.mac*

where *SPECD* is the auxiliary file directory, as described above, *geom* matches the first four letters of the name by which `spec` was invoked and *spec* is the complete name by which `spec` was invoked. The files are only read if they exist. The files *SPECD/standard.mac*, *SPECD/geom.mac*, *SPECD/spec/geom.mac* and *SPECD/site_f.mac* are only read if the user is starting `spec` for the first time or has invoked `spec` with the `-f` (fresh start) flag.

Each time `spec` starts up, if a macro named `begin_mac` is defined, that macro will be run after reading any startup command files, but before input is read from the keyboard.

After reading the start-up macro or command files and possibly running `begin_mac`, `spec` prompts the user for commands.

Keyboard Interrupts

On UNIX systems, two different asynchronous signals can be sent from the keyboard to programs. These signals are *interrupt* and *quit*. A `^C` is usually used to generate the *interrupt* signal. A `^\` generates a *quit* on many versions of UNIX, although `^v` is used on some others. The control key assignments are arbitrary and can be changed using the *stty* command from the UNIX shell. A common problem for new users is that their default UNIX environment has either no keys or archaic keys assigned to these signals. To display the current key assignments on BSD-type systems, type `stty everything` from the shell. On System V systems, type `stty -a`.

With *spec*, the *interrupt* key halts all activity, including asynchronous motor motion or counting, and closes all command files. All output files and devices (except *log*, *dlog* and *elog* files) are closed. On keyboard interrupts (and command and syntax errors), *cleanup* macros, as described below may be run.

Typing the *quit* character will asynchronously terminate *spec* without saving the user's state. However, if motors are moving, the program will wait for them to halt and then update the *settings* file.

Cleanup Macros

On keyboard interrupts (and command and syntax errors), if macros named `cleanup` or `cleanup1` have been defined, their definitions are read as input. Typical uses of the `cleanup` macro are to return motors to starting positions and/or to insert comments in data files after aborted scans. After running the `cleanup` and/or `cleanup1` macros, *spec* gives the standard prompt and waits for the next command from the keyboard. If there is another `^C` interrupt or error while the commands contained in the `cleanup` or `cleanup1` macros are being executed, the macro definitions are removed.

As of release 4.03.13, the additional `cleanup_once` and `cleanup_always` macro names are recognized. Each, if defined, will be pushed on to the input stream on `^C` interrupts and command and syntax errors, as above. The `cleanup_once` definition, though, is always removed before the next main prompt is issued. On the other hand, the `cleanup_always` definition will never be automatically removed. Definitions for these macros should be constructed using the `cdef()` (chained macro) function (described on page 91) in order to allow various clean-up actions to be added and removed during the flow of arbitrary statement blocks.

Exiting

A `spec` session is normally terminated by typing `^D` at the start of an input line. As of release 3.01.04, the `quit` command will also terminate `spec`, but only when typed directly from the keyboard or when read from a command file. The `quit` command can't be included in a macro.

As of release 4.02.04, if there is a macro named `end_mac` defined, it will be run after the `^D` or `quit` command.

When exiting, `spec` saves the user's state, and if any motors are moving, waits for them to halt and then updates the `settings` file.

If `spec` appears hung, typing the `quit` character (often `^\`) should terminate the program. If `spec` is waiting for unresponsive hardware to indicate it is finished with a move or a count, typing the `quit` or `interrupt` character again may be necessary.

Variables

A variable is brought into existence simply by using it. Variables assigned values at the top level (outside of any curly-bracketed statement block) are automatically made global. Otherwise, unless explicitly given *global* attributes, the scope of a variable lies only within the statement or statement block in which it occurs.

The possible attributes of a variable are as follows:

Local	A symbol with scope only in the statement block in which it appears.
Global	A symbol with scope that carries across separate statement blocks. All built-in symbols are global.
Built-in	A symbol that is compiled into <code>spec</code> and that cannot be removed.
Constant	A global symbol that cannot have its value altered by ordinary assignment. Such a symbol can only be changed using the <code>constant</code> command.
Immutable	Certain built-in symbols and motor and counter mnemonics from the hardware <code>config</code> file, which cannot have their values altered.

Variables can have string, number or array type and may have both string and number types simultaneously. The `print` command always prints the string representation of a variable. The formatted printing commands `printf`, `fprintf` and `sprintf`

print the string or number value, depending on the format specification. If the string value cannot be interpreted as a number, its number value is zero. All number values are maintained as double-precision floating-point quantities.

Uninitialized variables have a number value of zero and a string value of "" (the null string). Although associative array indices are internally derived from the string value of the index expression, if the index is an uninitialized variable, its value for purposes of indexing the array is the string "0".

Functions such as `input()`, `getline()`, `gpib_get()` and `ser_get()` return string values that possibly represent numbers. When the string is used in a number context, automatic string to number conversion takes place. The conversion rules require that there are no extraneous characters in the string. An initial `0x` or `0X` introduces a hexadecimal integer. An initial `0` introduces an octal constant, unless there is a decimal point or an `e` or `E` exponential string, in which case the entire string is interpreted as a decimal number, leading zero notwithstanding.

Built-In Variables

The following is a list of most of `spec`'s built-in variables. Some site-dependent code, along with most of the standard diffractometer geometry code, will create additional built-in variables. Also, motor and counter mnemonics entered in the hardware configuration file become built-in variables when the *config* file is read by `spec`.

- `A[]` — is an array dimensioned to the number of motors as obtained from the *config* file. The function `read_motors(0)` fills the array with user angles. The user can assign values to any of the elements. The commands `move_all` and `move_cnt` use the values in the array. Also, the various site-dependent, geometry-specific calculations, accessed through the `calc()` user-hook function, base their results on the values in this array or place new values in it.
- `COLS` — is a number-valued variable set to the number of text columns on the user terminal or workstation window. The value is used for formatting text-mode plots and on-line help files. `COLS` is generally automatically assigned a value from the system terminal data base when `spec` starts up, or, if available, by using the `TIOCGWINSZ` command in the `ioctl()` system call whenever a value for `COLS` is needed.
- `COUNTERS` — is the number of counters being used as determined from reading the *config* file.
- `CWD` — is a string-valued variable that contains the name of the user's current working directory. It is assigned a value when `spec` starts up, and is updated each time the `chdir()` function is executed.

`DEBUG` — is a user-assignable numeric variable that determines the level of debugging messages printed. The level is determined by the sum of the values given in this table:

Hex	Decimal	What is shown
0x1	1	Show input tokens during parsing.
0x2	2	Show node execution while running.
0x4	4	Show node allocate and free.
0x8	8	Show symbol table creation and lookup.
0x10	16	Show value get and set.
0x20	32	Show misc info.
0x40	64	Show hardware related info.
0x80	128	Show more hardware related info.
0x100	256	Show macro substitution.
0x200	512	Show memory allocate and free.
0x400	1024	Show input after macro expansion.
0x800	2048	Print warning messages.
0x1000	4096	Show certain low level hardware info.
0x2000	8192	Show data array allocation.
0x4000		Show signal blocking.
0x8000		Show sleeps and other low level hardware info.
0x10000		Show input file queuing.
0x20000		Show readable runtime debugging.
0x40000		Print input context on execution-time errors.

A change in the value of `DEBUG` only takes effect after the complete mini-program in which the new value is assigned has been executed. If a debugging log file is open (any file that begins with the characters `dlog` or ends with the characters `.dlog`), debugging messages are only written to that file, not to the screen or any other file or device.

`DISPLAY` — is a user-assignable, string-valued variable. Its value at the time an X-Windows graphics filter process is spawned with the `plot_cntl("open")` function (with `GTERM` set to `"x11"`) determines on which host and screen the plot window will be displayed. The initial value for `DISPLAY` is taken from the environment variable of the same name.

`FRESH` — is a built-in variable that has an initial nonzero value if `spec` was invoked with the `-f` (fresh start) flag or if a fresh start was forced by an incompatible state file version. The value is zero otherwise, and is set to zero in any case after all start-up command files and macros have been read and their commands executed. In the standard start-up macros, the value of `FRESH` is

checked to see if initial default parameter assignments should be made.

- GTERM — is a user-assignable, string-valued variable containing a value describing the display type to use for high-resolution graphics. Its value will be taken from an environment variable of the same name when `spec` starts up, if such a variable exists. Currently supported GTERM types are `vga`, `vpc`, `ega`, `epc`, `herc` and `cga` for various PC console graphics (not *linux*); `x11` for X Window System graphics; and `sun` for SunView graphics. If GTERM is not set in the environment or has not been assigned a value, it defaults to `x11`.
- HOME — is string valued and is initialized to the user's home directory as taken from the environment variable *HOME*. If not found in the environment, its value is set to the current directory.
- MOTORS — is the number of motors being used as determined from reading the *config* file.
- PAGER — can contain the name of a program to use for paging through `spec` help files. Such a program can be an alternative to `spec`'s built-in paging code, allowing, for example, searching and line-by-line perusal. On start up, `spec` examines the environment variables `SPEC_PAGER` and `PAGER`, in turn. If one exists, its value is assigned to `PAGER`. Also, if a version of the *less* utility greater than version 332 is installed on the `spec` platform, the standard macros will assign to `PAGER` a string that invokes *less* with special flags. Otherwise, `PAGER` is left unset.
- PI — is a number-valued symbol with the value 3.14159...
- ROWS — is a number-valued variable set to the number of text rows on the user terminal or workstation window. The value is used for formatting text-mode plots and on-line help files. `ROWS` is generally automatically assigned a value from the system terminal data base when `spec` starts up, or, if available, by using the `TIOCGWINSZ` command in the *ioctl()* system call whenever a value for `ROWS` is needed.
- S[] — is an array that will be filled with the hardware scaler contents when the command `getcounts` is executed.
- SPEC — is string valued and set to the name by which `spec` is invoked, such as *fourc*.
- SPECD — is string valued and set to `spec`'s auxiliary file directory. The default name is compiled in when `spec` is installed, but can be overridden by the `-D` invocation option or by the *SPECD* environment variable.
- TERM — is a user-assignable, string-valued variable. It is initialized to the user's terminal type as taken from the environment variable *TERM*. If not found in the environment, it is set to terminal-type *vt100*.

`USER` — is string valued and is set to the login name of the current user.

`VERSION` — is string valued and is set to the version number of `spec`, as in 3.03.11.

The values of `A[]`, `S[]`, `DEBUG`, `TERM`, `GTERM`, `ROWS` and `COLS` can be changed by the user. The other variables are immutable and cannot have their value changed.

Motor mnemonics obtained from the *config* file become built-in, immutable variables. User-added code, such as the X-ray diffractometer geometry code, typically creates other built-in variables, such as `G[]`, `Q[]`, `U[]` and `Z[]`.

Operators

The following tables summarize the operators available in `spec`. (Almost all these operators work the same as in the C language, so a C-language reference manual could be consulted to provide more detailed information on the use of unfamiliar operators.) Operators that require integral operands use the integer part of nonintegral operands. The precedence rules give the evaluation order when multiple operators appear in the same expression.

Unary Operators

Unary operators have equal precedence, are evaluated right to left and have higher precedence than the binary operators.

Operator	Name	Result
-	Unary minus	Negative of operand.
+	Unary plus	Operand.
!	Logical negation	1 if operand is zero, 0 otherwise.
++	Increment	Increment the operand by one.
--	Decrement	Decrement the operand by one.
~	Bitwise not	One's complement of operand. (A tilde.)

When used as prefix operators, `++` and `--` operate before the result is used in the expression. When used as postfix operators, `++` and `--` are done after the current value of the operand is used in the expression.

Indirection Operator

The `@` character is a special unary operator that provides indirection when used in front of a symbol name (available as of release 4.03.01). The behavior is similar to

the C language * indirection operator. In `spec`, the `@` operator allows reference to a variable whose name is the string value of another variable. For example:

```
1.FOURC> a = "b"; b = PI; print a, @a
b 3.14159

2.FOURC>
```

Binary Operators

All binary operators have lower precedence than the unary operators. Binary operator precedence is indicated in the grammar rules that are listed on page 64, where higher precedence operators are listed first, and operators with the same precedence are listed on the same line. Binary operators with the same precedence are evaluated left to right.

Operator	Name	Result (L is left operand, R is right)
*	Multiplication	$\mathbf{L} \times \mathbf{R}$
/	Division	\mathbf{L} / \mathbf{R}
%	Modulus	Remainder of \mathbf{L} / \mathbf{R} (operands are integers).
+	Addition	$\mathbf{L} + \mathbf{R}$
-	Subtraction	$\mathbf{L} - \mathbf{R}$
<	Less than	1 if $\mathbf{L} < \mathbf{R}$ (or 0 if not).
>	Greater than	1 if $\mathbf{L} > \mathbf{R}$ (or 0 if not).
<=	Less than or equal	1 if $\mathbf{L} \leq \mathbf{R}$ (or 0 if not).
>=	Greater than or equal	1 if $\mathbf{L} \geq \mathbf{R}$ (or 0 if not).
==	Logical equality	1 if \mathbf{L} is equal to \mathbf{R} (or 0 if not).
!=	Logical inequality	1 if \mathbf{L} is not equal to \mathbf{R} (or 0 if it is).
&&	Logical and	1 if both \mathbf{L} and \mathbf{R} are nonzero, otherwise 0.
	Logical or	1 if either \mathbf{L} or \mathbf{R} are nonzero, otherwise 0.
<<	Bitwise left shift	\mathbf{L} shifted left by \mathbf{R} bits (both integers).
>>	Bitwise right shift	\mathbf{L} shifted right by \mathbf{R} bits (both integers).
&	Bitwise and	“Bitwise and” of integers \mathbf{L} and \mathbf{R} .
^	Bitwise exclusive or	“Bitwise exclusive or” of integers \mathbf{L} and \mathbf{R} .
	Bitwise or	“Bitwise or” of integers \mathbf{L} and \mathbf{R} .
	Concatenation	\mathbf{LR}

If either of **L** or **R** are strings, the relational operators `<`, `>`, `<=`, and `>=` use the lexicographic comparison provided by the C subroutine `strcmp()`.

The concatenation operator comes into effect when expressions are combined separated only by space characters. The resulting expression is the concatenation of the string values of the constituent expressions. Concatenation is only allowed on the right side of an assignment operator, in the arguments to the `print` command and in the value assigned to a variable in a `constant` statement. Concatenation has lower precedence than the other operators. For example,

```
2.FOURC> print "ab" "cd" 1 2 + 3 4
abcd154

3.FOURC>
```

Assignment Operators

Operator	Name	Result
=	Equals	L = R
+=	Plus equals	L = L + R
-=	Minus equals	L = L - R
*=	Times equals	L = L × R
/=	Divide equals	L = L / R
%=	Mod equals	L = L % R
<<=	Left-shift equals	L = L << R
>>=	Right-shift equals	L = L >> R
&=	Bitwise-and equals	L = L & R
^=	Bitwise-exclusive-or equals	L = L ^ R
=	Bitwise-or equals	L = L R

Ternary Operator

`spec`, like the C language, has one ternary operator, which works in a manner similar to the if-else construction and uses the characters `?` and `:`. Its use is

$$expression_1 ? expression_2 : expression_3$$

where the result of this entire expression is $expression_2$ if $expression_1$ is nonzero, otherwise $expression_3$.

Flow Control

Conditional Statement

The forms of the conditional statement are

```
if ( expression ) statement
if ( expression ) statement else statement
```

The *expression* is evaluated in both cases. If nonzero, the first *statement* is executed. In the second form, if *expression* is zero, the second *statement* is executed. Each *else* is always matched with the last *else-less if* within the same statement block.

While Statement

The form for the while statement is

```
while ( expression ) statement
```

The *expression* is evaluated and the *statement* executed until *expression* is zero.

For Statement

The *for* statement comes in two forms. The first is

```
for ( expression1 ; expression2 ; expression3 ) statement
```

The flow expressed in this statement can be thought of as

```
expression1
while ( expression2 ) {
    statement
    expression3
}
```

Any of the expressions can be missing. A missing *expression*₂ is equivalent to

```
while ( 1 )
```


The second form of the `for` statement concerns associative arrays. A syntactical feature for these arrays (forms of which are available in the standard *awk* and *gawk* utilities) is available in `spec` as of release 4.01.05. The construction

```
for ( identifier in assoc-array ) statement
```

will run through each element of the associative array *assoc_array* assigning to *identifier* the string value of the index of each element. Note, though, the order in which the elements are assigned should be considered arbitrary.

For two-dimensional associative arrays, the construction

```
for ( identifier in assoc-array[expr] ) statement
```

will step through each element of *assoc-array* having *expr* as the first index.

Break Statement

The statement `break` terminates the smallest enclosing `while` or `for` loop.

Continue Statement

The statement `continue` passes control to the loop-test portion of the smallest enclosing `while` or `for` loop.

Exit Statement

The statement `exit` terminates execution of the current parse tree and jumps control back to command level just as if an error occurred or a `^C` was typed.

Grammar Rules

This syntax summary defines how all the built-in keywords, operators and commands of `spec` can be combined. These grammar rules are similar to those given in standard C-language manuals. Operators are listed in order of precedence, with the highest precedence operators listed first.

The following terms are used in the grammar rules:

- lvalue* - “Left value”, something on the left side of an equals sign.
- binop* - A binary operator (\neq , etc.).
- asgnop* - An assignment operator ($=$, $+=$, etc.).
- assoc-array* - An associative (original-style) array.
- assoc-elem-list* - A *space-* or *comma-*separated list of associative array elements.
- identifier* - A variable.
- identifier-list* - A *space-* or *comma-*separated list of identifiers.
- pattern* - An alphanumeric string possibly containing the metacharacters `?` or `*`.
- pattern-list-opt* - An optional *space-*separated list of patterns.
- expression-list* - A *comma-*separated list of expressions.
- expr-opt* - An optional expression.
- `[;]` - A semicolon or a newline. (A semicolon after a statement is optional if the statement is followed by a newline.)

Note, that in the following list, the entry

expression in *assoc-array*

is included in the rules of what constitutes an *expression*. This is a special expression that evaluates to nonzero (or **true**) if *assoc-array*[*expr*] is an existing element of the array, and zero (or **false**) otherwise. For a two-dimensional associative array,

*expr*₁ in *assoc-array*[*expr*₂]

is nonzero if *assoc-array*[*expr*₂][*expr*₁] is an element of the array.

These are the grammar rules:

expression:

lvalue
numeric-constant
string-constant
(expression)
function (expression-list)
- expression
! expression
~ expression
++ lvalue
-- lvalue
lvalue ++
lvalue --
expression ? expression : expression
expression binop expression
expression in assoc-array
lvalue asgnop expression
expression , expression
expression expression

lvalue:

identifier
identifier [expression]
identifier [expression] [expression]

binop:

** / %*
+ -
>> <<
> < <= >=
== !=
&
^
|
&&
||

asgnop:

*= += -= *= /= %= >>= <<= &= ^= |=*

data-array-type

string
byte
ubyte
short
ushort
long
ulong
float
double

data-array-declaration:

array *identifier*[*expression*]
data-array-type array *identifier*[*expression*]
array *identifier*[*expression*][*expression*]
data-array-type array *identifier*[*expression*][*expression*]

compound-statement:

{ *statement-list* }

statement-list:

statement
statement statement-list

statement:

compound-statement
expression [;]
if (*expression*) *statement*
if (*expression*) *statement* else *statement*
while (*expression*) *statement*
for (*expr-opt* ; *expr-opt* ; *expr-opt*) *statement*
for (*identifier* in *assoc-array*) *statement*
break [;]
continue [;]
exit [;]
history *expr-opt* [;]
lscmd *pattern-list-opt* [;]
print *expression-list* [;]
global *identifier-list* [;]
constant *identifier expression* [;]
unglobal *identifier-list* [;]
delete *assoc-elem-list* [;]

```
local identifier-list [;]
syms pattern-list-opt [;]
data-array-declaration [;]
shared data-array-declaration [;]
def identifier string-constant [;]
rdef identifier expression [;]
undef identifier-list [;]
prdef pattern-list-opt [;]
lsdef pattern-list-opt [;]
memstat [;]
savstate [;]
reconfig [;]
getcounts [;]
move_all [;]
move_cnt [;]
sync [;]
```

Built-In Functions and Commands

These built-in functions and commands are described in the following sections.

UTILITY FUNCTIONS AND COMMANDS							
<i>System Functions</i>				<i>Miscellaneous</i>			
chdir()	time()	gethelp()	file_info()	sleep()	calc()	lscmd	savstate
unix()	date()	getenv()	whatis()	spec_par()	history	memstat	sock_io()

KEYBOARD AND FILE INPUT, SCREEN AND FILE OUTPUT						
<i>Controlling Output Files</i>		<i>Reading From Files</i>	<i>Keyboard Input</i>	<i>Text Output</i>		<i>Binary Input/Output</i>
open()	on()	getline()	input()	print	fprintf()	fmt_read()
close()	off()	dofile()	yesno()	printf()	tty_cntl()	fmt_write()
		qdofile()	getval()	eprint	tty_move()	fmt_close()
				eprintf()	tty_fmt()	

VARIABLES			
global	constant	delete	syms
unglobal	local	array	

MACROS		
def	prdef	undef
rdef	lsdef	cdef()

HARDWARE FUNCTIONS AND COMMANDS						
<i>MCA</i>	<i>Images</i>	<i>CAMAC</i>	<i>GPIB</i>	<i>SERIAL</i>	<i>PORT I/O</i>	<i>VME</i>
mca_sel()	image_par()	ca_get()	gpib_get()	ser_get()	port_get()	vme_get()
mca_par()	image_get()	ca_put()	gpib_put()	ser_put()	port_getw()	vme_put()
mca_get()	image_put()	ca_fna()	gpib_poll()	ser_par()	port_put()	vme_get32()
mca_put()		ca_cntl()	gpib_cntl()		port_putw()	vme_put32()
mca_spar()						vme_move()
mca_sget()						
mca_sput()						

<i>Misc</i>	<i>Counting</i>		<i>Moving and Motors</i>			
reconfig	mcount()	cnt_mne()	move_all	motor_mne()	read_motors()	
set_sim()	tcount()	cnt_name()	move_cnt	motor_name()	dial()	user()
wait()	getcounts	cnt_num()	sync	motor_num()	chg_dial()	chg_offset()
stop()		counter_par()		motor_par()	get_lim()	set_lim()

DATA HANDLING AND PLOTTING FUNCTIONS					
array_dump()	array_op()	data_grp()	data_dump()	data_bop()	plot_cntl()
array_read()	array_fit()	data_info()	data_read()	data_uop()	plot_move()
array_pipe()		data_nput()	data_pipe()	data_anal()	plot_range()
array_plot()		data_get()	data_plot()	data_fit()	
		data_put()			

STRING AND NUMBER FUNCTIONS							
<i>Math Functions</i>				<i>String</i>		<i>Conversion</i>	
exp()	log()	sin()	cos()	tan()	index()	length()	asc() deg()
exp10()	log10()	asin()	acos()	atan()	split()	sprintf()	bcd() rad()
srand()	pow()			atan2()	substr()	sscanf()	dcb()
rand()	sqrt()	fabs()	int()				

Utility Functions and Commands

All functions return a number or string value that may be used in an expression. The return values of some functions are only of use in conditional expressions, so their return values are given as **true** or **false**. The corresponding number values are 1 and 0, respectively.

System Functions

`chdir()` – Changes `spec`'s current working directory to to the user's home directory as obtained from the user's environment variable *HOME*. Returns **true** or **false** according to whether the command was successful or not. The value of the built-in string variable `CWD` is updated to the current working directory.

`chdir(directory)` – As above, but changes to the directory *directory*, which must be a string constant or expression.

`unix()` – Spawns an interactive subshell using the program obtained from the user's environment variable *SHELL* (or *shell*). Uses */bin/sh* if the environment variable is not set. Returns the integer exit status of the shell.

`unix(command)` – As above, but uses */bin/sh* to execute the one-line command *command*, which must be a string constant or expression. Returns the integer exit status of the command.

`unix(command, str [, len])` – As above, but the argument *str* is the name of a variable in which to place the string output from the command in the first argument. The maximum length of the string is 4096 bytes (including a null byte). The optional third argument can be used to specify a larger size.

`time()` – Returns the current epoch in seconds. The UNIX epoch is the number of seconds from January 1, 1970, 00:00:00 GMT. The value returned includes a fractional part with the resolution dependent on the host platform. Millisecond resolution is standard, although on 80X86 systems only hundredth of a second resolution is returned.

`date()` – Returns a string containing the current date as

```
Thu Jul 23 21:12:39 EDT 1992 .
```

`date(seconds)` – As above, but the returned string represents the epoch given by *seconds*. See `time()` above.

Miscellaneous

`gethelp(topic)` – Formats and displays the help file with the name *topic* to fit the current screen size (as defined by the built-in variables `ROWS` and `COLS`). The argument *topic* is a string constant or expression. If *topic* is a path name containing the `/` character, that is the path name used for locating the file. Otherwise, the file must be in the subdirectory *help* of `spec`'s auxiliary file directory, given by the global variable `SPECD`. The function returns nonzero if the file could not be opened.

`spec` will also read a file named `SPECD/help/.links` that contains a list of help topics and file names. If a file corresponding to *topic* is listed in *.links*, that file will be displayed.

Normally, the help file is displayed one screenful at a time. At the end of each page, the user may type `<space>`, `<return>`, `y` or `Y` to see the next page, `b` to go back a page or anything else to return to the `spec` prompt.

If the variable `PAGER` is set, the command therein is executed with the formatted help file as input. Typically, `PAGER` is set to a pagination utility such as the *more* or *less* programs.

The help file contents are only written to the screen or display window and not to any other turned on output files or devices.

The syntax recognized by `spec`'s built-in formatting uses the conventions of the traditional UNIX text formatters (*nroff*, *troff*, *ditroff*, etc.), although only a small subset of the formatting directives are supported. Those are described below.

The `spec` distribution includes utilities for preparing both PostScript and HTML (hyper text mark-up language) versions of the help files using either the standard *ditroff* or the GNU *groff* packages.

Input text that has a dot on the first line is a formatting command. Unrecognized formatting commands are ignored. A command such as `.if t ...` can be used to insert commands to be executed only when the help files are being printed with *troff*, etc., as the `.if` is not recognized by `spec`.

The following commands are supported:

- `.\` – Begins a comment. The rest of the text on the line is ignored.
- `.bp` – Begins a new page.
- `.BP` – Same as above.
- `.ne N` – Indicates *N* lines are needed on the current page. If there isn't enough room on the current page, a new page is started.

- .sp *N* — spaces *N* lines. Without an argument, one line is skipped. The argument may contain a decimal point and a unit indicator for the benefit of *nroff/troff*, although with *spec*, only the integer part of the argument is used.
- .ta *N* — Start words every *N* letters, somewhat like using tabs.
- .br — Causes a line break. The current line is finished and a new line is started. Useful for forcing a newline in fill mode.
- .fi — Switch to fill mode. Words are collected from the input until a full line of output is available.
- .nf — Switch to no-fill mode. Each newline in the input causes a newline in the output.
- .TH *name category date SCCS_info* — Should be on the first line of a help file. The arguments are used in the page header (and footer, when printed) and are, respectively, the name of the help file, the category (such as macro, function, reference, etc.), the date the file was last revised and the source-code control information.
- .SH *subhead ...* — Does a subhead. There can't be more than six arguments, although you can use double quotes to put more than one word in an argument. Subheads are placed all the way to the left of the page.
- .HP *N* — Starts a hanging paragraph. The optional argument *N* specifies the indent level. The default is 0.
- .hp *N* — Can be used within a *.HP* block to increase or decrease the indent. The optional argument *N* specifies the indent level. This command will skip a line.
- .PP — Starts a paragraph.
- .LP — Same as above.
- .(L *F N I* — Starts a no-fill block. The optional argument *F* indicates the font to use and may be R, I, C, B or O. (See below for a description of the font code letters.) The optional argument *N* is the number of lines needed. If there is not enough room to print that many lines on the current page, a new page will be started. The optional argument *I* is the number of spaces of extra indentation to add.
- .)L — Ends the no-fill block started by *.(L*.
- .Q — Puts quotation marks around the arguments. If no arguments, puts quotation marks around text on the next input line.
- .NE — Only relevant for screen displays, this command indicates the screen should not be erased before displaying the first page.
- # — Starts a comment if it is the first character on a line. The rest of the text on the line is ignored. This comment notation is only for files that will not be run off using *nroff/troff*.

- .BL – Begin a list of items.
- .EL – Ends a list of items.
- .LI *item* – Adds *item* to a list of items. For HTML documents, one list item appears per line of output. For non-HTML documents, list items are separated by tabs, and will fill the row according to the current tab stops.
- .DL + or .DL – – For HTML output, indicates the start or end of a definition block.
- .DD *text* – For HTML output, indicates the start of a definition. For the other formats, the optional argument is simply copied verbatim.
- .HR *text* [C [URL]] – For HTML output, creates a hypertext link. If one or two arguments are present, names the link . If the second argument is C, outputs *text* in Courier font. If the third argument is present, uses that as the URL link. For non-HTML formats, the text is copied verbatim in the indicated font.

The help files use the following fonts for the purposes indicated. In the printed documentation, the New Century Schoolbook typeface is used for the Roman, italic and bold fonts. When used with `spec` and `C-PLOT`, the directives to switch to bold, italic or one of the Courier typefaces result in highlighted text on a video terminal.

Roman (R) – Used for normal text.

Italic (I) – Used for file and path names, including system commands, along with book titles, etc.

Bold (B) – Used for section heads. Occasionally used for emphasis. Sometimes used to display arguments to system commands.

Courier (C) – Used to represent literal text that might be output by the computer program or might be typed by a user at the keyboard.

Courier Bold (CB) – Occasionally used to distinguish text a user types from text output by the computer.

Courier Oblique (O) – Used to represent variable computer text, such as the argument to a function.

“Quoted Text” (Q) – Not really a font, but used to place double quotes around a bit of text.

All the fonts, except Roman, display as highlighted text on the video terminal.

The formatting commands:

```
.B .C .I .O
```

put their arguments in the associated font. If there are no arguments to the above commands, the next line of input is put in that font.

The formatting commands:

```
.RB .BR .IR .CR .OR .CB
.RI .BI .IB .CI .OC
.RC .IC .CO
.RO
```

place the first argument in the first font, the second argument in the second font and append an optional third argument in the regular text font, except for `.CB`, which uses Courier Bold and Roman for the two fonts. The space characters between the arguments are removed. Use double quotes to include space characters within arguments. A space character is added after the last argument. The formatting commands:

```
.rb .ir
.ri .ci .oc
.rc .ic .co
.ro
```

work as above, except a space characters is not appended after the last argument.

The `@` symbol is used in the help file sources to represent literal double quote characters, which makes it easier to insert literal double quote characters into formatting command arguments. To enter a literal `@` use both a backslash and the troff `.tr` directive, as in

```
.tr @@
a literal \@ appears
.tr @"
```

The following special character sequences are recognized:

<code>\"</code>	A double quote
<code>\&</code>	A zero-width space
<code>\0</code>	A digit-width space
<code>\ </code>	A sixth-em space (no space on screen)
<code>\(**</code>	An asterisk (*)
<code>\(em</code>	A long dash (—)
<code>\(hy</code>	A hyphen (-)
<code>\(mu</code>	The multiplication sign (×)
<code>\(space)</code>	An unpaddable space
<code>*(CP</code>	The string C-PLOT
<code>*(SP</code>	The string spec
<code>*(UN</code>	The string UNIX
<code>\-</code>	A minus sign (-)

<code>*-</code>	A long dash (—)
<code>\@</code>	An @ sign
<code>\\</code>	A backslash
<code>\e</code>	A backslash
<code>\f1</code>	Switch to font R
<code>\fC</code>	Switch to font C
<code>\fP</code>	Switch to previous font
<code>\f*(#B</code>	Switch to font B
<code>\f*(#C</code>	Switch to font C
<code>\f*(#I</code>	Switch to font I
<code>\f*(#O</code>	Switch to font O

`getenv(string)` — Returns the value of the environment variable represented by the string *string*. If the environment variable is unset, the null string is returned. Environment variables are exported to `spec` by the invoking shell program.

`file_info(filename [, cmd])` — Returns information on the file or device named *filename*. With a single *filename* argument, `file_info()` returns **true** if the file or device exists. If the argument *filename* is the string "?", the possible values for *cmd* are listed. If *filename* is the string ".", `spec` uses the information from the last `stat()` system call made using the previous argument for *filename*, avoiding the overhead associated with an additional system call.

Possible values for *cmd* and the information returned follow. Note that the first set of commands parallel the contents of the data structure returned by the `stat()` system call, while the second set of commands mimic the arguments to the `test` utility available in the shell.

- "dev" — The device number on which *filename* resides.
- "ino" — The inode number of *filename*.
- "mode" — A number coding the access modes and file attributes.
- "nlink" — The number of hard links for *filename*.
- "uid" — The user id of the owner.
- "gid" — The group id of the owner.
- "rdev" — The device ID if *filename* is a block or character device.
- "size" — The size in bytes of *filename*.
- "atime" — The time when *filename*'s data was last accessed.
- "mtime" — The time when *filename*'s data was last modified.
- "ctime" — The time when *filename*'s attributes were last modified.
- "isreg" or "-f" — Returns **true** if *filename* is a regular file.
- "isdir" or "-d" — Returns **true** if *filename* is a directory.

- "ischr" or "-c" — Returns **true** if *filename* is a character device.
- "isblk" or "-b" — Returns **true** if *filename* is a block device.
- "islnk" or "-h" or "-L" — Returns **true** if *filename* is a symbolic link.
- "isfifo" or "-p" — Returns **true** if *filename* is a *named pipe* (sometimes called a *fifo*).
- "issock" or "-S" — Returns **true** if *filename* is a socket.
- "-e" — Returns **true** if *filename* exists.
- "-s" — Returns **true** if the size of *filename* is greater than zero.
- "-r" — Returns **true** if *filename* is readable.
- "-w" — Returns **true** if *filename* is writable.
- "-x" — Returns **true** if *filename* is executable.
- "-o" — Returns **true** if *filename* is owned by you.
- "-G" — Returns **true** if *filename* is owned by your group.
- "-u" — Returns **true** if *filename* is setuid mode.
- "-g" — Returns **true** if *filename* is setgid mode.
- "-k" — Returns **true** if *filename* has its sticky bit set.

Finally, the usage `file_info(pid, "alive")` will return **true** if the process associated with the process ID *pid* exists.

`whatis(string)` — Returns a number that encodes what the identifier in *string* represents. The return value is a two-word (32-bit) integer, with the low word containing a code for the type of object and the high word containing more information for certain objects.

High Word	Low Word	Meaning
0	0	Not a command, macro or keyword.
0	1	Command or keyword.
Length	2	Macro name (length is in bytes).
0x0010	4	Number-valued symbol.
0x0020	4	String-valued symbol.
0x0040	4	Constant-valued symbol.
0x0100	4	Associative array name.
0x0200	4	Built-in symbol.
0x0400	4	Global symbol.
0x0800	4	Unset symbol.
0x2000	4	Immutable symbol.
0x4000	4	Local symbol.
0x8000	4	Associative array element.
0x0001	4	A data array or data array element.

Most type-4 symbols have more than one of the high-word bits set.

The following code uses the `whatis()` function to determine whether or not a value has been assigned to the `PRINTER` variable:

```
3.FOURC> if ((whatis("PRINTER")>>16)&0x0800) {
4.more>     printf("Warning:  No open printer.\n")
5.more>     open(PRINTER = "/dev/null")
6.more> }

7.FOURC>
```

`whatis(string, 1)` – With two arguments, returns a printable string explaining in words what kind of thing the thing named in *string* is to `spec`.

`sleep(t)` – Suspends execution for a minimum of *t* seconds, where *t* may be non-integral. Actual sleep times may vary depending on other activity on the system and the resolution of the system clock. Returns **true**. Can be interrupted with `^C`. Fractional-second sleeps are available on all currently supported systems, although on the System V/386 UNIX systems, such sleeps are available only if the special *nap* driver has been installed.

`spec_par(par [, value])` – Sets internal parameters. Typing `spec_par("?")` lists the available parameters. The currently available parameters are:

"`auto_file_close`" – The auto-file-close option is available to automatically close output files that haven't been accessed for some interval of time. The parameter units are hours, and the parameter can have nonintegral values. When the auto-close option is enabled, each time an `on()`, `off()`, `open()`, `close()` or `fprintf()` function is called, `spec` will check its list of opened output files. Any files which have not been written to for the length of time given by *value* hours will be closed. Enabling this option can help prevent errors when your macros or commands do not close files when appropriate, resulting in `spec` running out of resources to open additional files.

As files are opened automatically when sent output, auto-close mode operates transparently for the most part. However, if you change to a different working directory between the time the file is first opened and subsequently automatically closed, and if the file is not opened by an absolute path name, the next time you refer to the file, `spec` will reopen it using a relative path based on the current directory.

If *value* is zero, the mode is disabled. By default, the mode is initially disabled.

"`auto_hdw_poll`" – When automatic hardware polling is turned on, `spec` will automatically poll busy motor controllers, timers and acquisition devices to determine when they are finished. For interrupt-driven devices (certain motor controllers and timers), this mode is irrelevant. For some

polled devices, `spec` needs to perform an action, such as starting a motor backlash move, when the device is finished with its current business. Without automatic hardware polling, a call to the `wait()` function is necessary. The default is for this mode to be on. A reason to turn it off may be to reduce the amount of debugging output during hardware debugging.

- "`check_file_names`" – The check-file-names option can prevent you from (accidentally) creating files with names containing nonstandard characters. When enabled, if a file name passed to the `on()`, `open()` or `fprintf()` functions contains any of the characters `()[]{}|$\`*?;!&<>`, the space character, any control characters or any characters with the eighth bit set, and the file doesn't already exist, `spec` will print an error message and reset to command level. By default, this mode is on.
- "`eelog_timestamp`" – The time interval for the optional time stamps for the `eelog` error file capability is set using this option. The units of the "`eelog_timestamp`" parameter are minutes. The default value is five minutes. Note, time stamps are only added before a command or error message is logged, so that the interval between time stamps can be greater than that specified if no commands are being typed or errors generated.
- "`flush_interval`" – The flush-interval parameter controls how often `spec` flushes output to the hard disk or other output device. Traditionally, `spec` flushed output after each print command. However, as some users report that this frequent flushing introduces considerable delays when the output device is to an NFS-mounted file system, the flushing interval can now be controlled. The "`flush_interval`" parameter specifies how many seconds to allow between output buffer flushing. The default value is 0.5 seconds. If the interval is set to zero, the traditional frequent-flushing behavior will be restored. Output to the screen is still flushed immediately. Output is also flushed each time the main `spec` prompt is issued.
- "`HKL_rounding`" – Traditionally, `spec` rounded values for `H`, `K`, and `L` (and other geometry values derived from motor positions) to five significant digits for configurations using reciprocal space calculations. As of release 4.03.01, the rounding is now turned off by default. It can be turned on using the command `spec_par("HKL_rounding", 1e5)` where the argument indicates the magnitude of the rounding, i.e., one part in `1e5`, for example. Note, values with an absolute value less than `1e-10` are still rounded to zero whether or not the optional rounding is turned on.

- "`hdw_poll_interval`" – When the `wait()` function is called to wait for polled motors, timers or other acquisition devices to finish, `spec` sleeps for a small interval between each check of the hardware. Use this `spec_par()` option to change that interval. The units of the parameter are milliseconds, and the default sleep time is 10msec. A value of zero is allowed, though not recommended if the computer is being used for anything else.
- "`keep_going`" – Normally, when taking commands from a command file, `spec` resets to command level and the main interactive prompt when there are syntax errors in the file, certain floating point exceptions, references to unconfigured hardware, hardware access errors, along with a number of other errors. When the "`keep_going`" option is set, `spec` will keep reading and executing commands from a command file no matter what errors occur. When there is an error, the next line from the current command file will be read. Note, depending on where the error is in a file, reading subsequent lines may generate more errors, particularly if the error occurs inside a statement block.
- "`modify_step_size`" – Normally, `spec` doesn't allow users to modify the motor step-size parameter with the `motor_par()` function, as the consequences are generally undesirable. However, in the rare circumstance that it is necessary, this parameter allows you to enable such modifications.
- "`parse_units`" – When parsing of units is turned on, numbers typed as input to `spec`'s parser with one of the following suffixes appended will automatically be multiplied by the corresponding factor.

1r	57.2958	radian
1mr	0.0572958	milliradian
1d	1	degree
1md	0.001	millidegree
1mm	1	millimeter
1um	0.001	micrometer
1m	0.0166667	minute
1s	0.000277778	second

Note, however, suffixes on numbers converted from strings or entered using the `getval()` function are not parsed. The only known use for unit-suffix parsing is with the user-contributed macros in the file `macros/units.mac`. These macros require that unit suffixes be supplied for all motor position arguments in the standard `spec` macros. The default is for this mode to be off.

"show_prdef_files" — When this mode is on, the source file for each macro definition is displayed with the `prdef` command. The default is for this mode to be on.

"specwiz" — Allows `spec` administrators to gain access to motors locked out in the `config` file to ordinary users. This feature allows the administrator to position the motors without having first to go into the configuration editor to change access modes. Entering `spec_par("specwiz", 1)` causes `spec` to prompt for the "Wizard's password". If entered correctly, the characters `_WIZ` are appended to the `spec` prompt to remind the administrator of the special powers, and motors configured with protected status can be moved. Entering `spec_par("specwiz", 0)` disables the special mode.

`spec` looks for the encrypted password belonging to the `spec_wiz` (or `specwiz`) user in the files `SPECD/passwd`, `/etc/shadow`, and `/etc/passwd` in turn. If a shadow password file is used, ordinary users won't be able to read it, and the normal password file won't contain the encrypted password.

The `spec` distribution includes a `wiz_passwd` utility (as of release 4.03.07) that can be used to create a `passwd` file in the `spec` auxiliary file directory that contains just the entry for the `spec_wiz` user. The `SPECD/passwd` file should probably be owned and writable only by root or the `spec` administrator.

"use_sem_undo" — This flag relates to whether the `SEM_UNDO` flag is set when semaphores are used. It exists to get around a memory leak bug observed with some releases of the Solaris 2 operating system. The flag should be ignored unless you are instructed otherwise by CSS.

`calc(i)` — Calls a user-added function having code `i`. Codes are assigned in the distribution source file `u_hook.c`. Returns a user-supplied value or zero if there is no user-supplied value.

`calc(i, x)` — As above, but passes the argument `x` to the user-supplied function.

The geometry calculations that transform motor positions to reciprocal space coordinates and vice versa are implemented using calls to `calc()`. A description of the particular calls implemented for the four-circle diffractometer are in the *Four-Circle Reference*. See page 195 in the *Administrator's Guide* for information on how to include other user-added functions in the program.

`history` — This command lists the 1000 most recent commands typed at the keyboard. See the description of the command-recall (history) feature on page 50.

`history N` – As above, but only prints the N most recent commands. If N is negative, the commands are printed in reverse order.

`lscmd` – This command lists all the built-in commands, functions and keywords.

`lscmd pattern ...` – As above, except only names matching *pattern* are listed.

`memstat` – Shows the current memory usage.

`savstate` – Writes the current state to the state file.

The *state* file contains the variables, macro-file bitmap, output file names and command history unique to each user, terminal and diffractometer. The state file preserves the current situation when the user leaves the program, so that the same situation can be restored when the user later returns to the program.

`spec` can be invoked with a `-u user` flag and a `-t tty` flag. These flags instruct the program to initialize the current user's state from the files associated with the other user and/or terminal. Subsequent `savstate` commands access the user's natural state file.

The `savstate` command does not save the state files until the entire parsed mini-program in which the command occurs has been run.

Note, as of release 4.05.01, a new `show_state` utility is included in the `spec` distribution. This utility can display the contents of `spec` state files. Type `show_state` – to see the utility's usage message.

`sock_io(host:port, cmd [, arg])` – Connects and communicates with sockets created by other processes on the local or a remote host. The *host* can be specified by a symbolic name or by an IP address. The string *cmd* is one of the following:

"?" – Lists the command options available.

"show" – Lists the current open sockets and their status.

"connect" – Opens the socket to the specified host and port. Returns **true** for success and **false** for failure.

"listen" – Sets up a socket for listening, allowing another instance of `spec` or some other program to make a connection.

"close" – Closes the socket associated with the specified host and port.

"send" – Sends the string *arg* to the specified socket, opening it if not already connected.

"read" – Reads from the specified socket, opening it if not already connected. Currently, the message to be read can be no longer than 1,024 bytes. With no argument, reads until a newline or until the timeout expires. An optional string argument can specify an end-of-string (EOS) delimiter, in which case the socket will be read until characters matching the

EOS sequence are read or until the timeout expires. Trailing newlines and carriage returns will be removed from the returned string. Also, if the optional argument is a number, it specifies the number of bytes to be read. That number of bytes (or less, if the timeout expires) will be read and returned.

"timeout" — With no argument, returns the current read-timeout value in seconds. Otherwise, sets the read timeout to *arg*. A negative timeout value causes the socket read to block until it is satisfied. The default timeout value is five seconds.

Note, the "connect" command doesn't generally need to be used, as the connection will be made if needed with the other commands. However, the return value from the "connect" command may be used to test if a socket can be created. A connection remains open until the "close" argument is used. Simultaneous connections to multiple ports are allowed.

The following example connects to the *echo* service on port 7 of the local host.

```
7.FOURC> sock_io("localhost:7", "send", "This is a test.\n")
8.FOURC> print sock_io("localhost:7", "read")
This is a test.
9.FOURC>
```

Input and Output

Controlling Output Files

With this group of functions, the names "tty" and "/dev/tty", when used for *filename*, are special and refers to the user's terminal. The names "null" and "/dev/null" are also special and when used as an output device, result in no output. The name "pipe" is also special, but only when *spec* is invoked with the **-p** flag, where it refers to the special data stream from *spec* to a front-end program.

`open()` — Lists all open files, including their directories, and indicates which files are currently turned on for output. Returns zero.

`open(filename)` — Makes *filename*, which is a string constant or expression, available for output. Files are opened to append. Returns zero for success, -1 if the file can not be opened or if there are too many open files. If the `spec_par()` "check_file_name" option is on, and if *filename* contains any of the characters `()[]{}|$\`*?;!&<>`, *spec* will print an error message and reset to command level, unless the file already exists.

`close(filename)` – Closes *filename* and removes it from the table of files available for output. Returns zero for success, -1 if the file wasn't open. Files should be closed before modifying them with editors.

`on()` – Lists all open files and indicates which ones are currently turned on for output.

`on(filename)` – Turns on *filename* for output. All messages, except for some error and debugging messages, but including all `print` and `printf()` output, are sent to all turned-on devices. If *filename* has not been made available for output with the `open()` function, it will be opened. Returns zero for success, -1 if the file can't be opened or if there are too many open files.

`off(filename)` – Turns off output to *filename*, but keeps it in the list of files available for output. If this was the last turned-on file or device, `tty` is turned back on automatically. Returns zero for success, -1 if the file wasn't open.

`spec` remembers the directory the files are in when they are first opened. If the user changes `spec`'s current directory, open files may be referenced either by the name with which the files were opened or by the correct path name relative to the current directory. If an open file disappears from the file system (for example, if a user removes the file using a subshell), the next time the file is written to, `spec` prints a warning message and creates a new instance of the file.

Files should be closed before attempting to edit them outside of `spec`.

Errors during parsing or execution of commands, or typing a `^c` turns off all open files except *log* files (see next section).

Log Files

Special *log* files can be created for debugging or archiving purposes. A regular *log* file is a file with a name beginning with the characters `log` or ending with the characters `.log`. All output sent to any device is sent to a *log* file.

A special type of *log* file is a debugging *log* file. That is a file with a name that begins

with `dlog` or ends with `.dlog`. Debugging messages (generated when a value is assigned to the built-in variable `DEBUG`) are only written to such a file, not to the screen or any other file or device.

Another kind of *log* file is available (as of release 4.05.01) that records typed commands, error messages and optional time stamps. The file is intended to be useful to administrators trying to diagnose user problems. A file opened for output named `elog` or that ends in `.elog` is such a file. Commands entered at the `spec` prompt are logged prefixed by a `#C`. Error messages produced by the built-in C code, or generated by the `eprint` or `eprintf()` built-in keywords, are logged prefixed by a `#E`. If time stamps are enabled (via the `spec_par()` "elog_timestamp" option), the UNIX epoch and the corresponding date string are logged (at the time-stamp interval) prefixed by a `#T`.

Reading From Files

The first function below is for reading strings from a file one line at a time. The second and third functions cause `spec` to switch its source of command input from the keyboard to the specified files.

`getline(file [, arg])` – This function reads successive lines from the ASCII file *file* each time it is called and returns the string so obtained, including the trailing newline. If *arg* is the string "open", the function returns zero if the file can be opened for reading. Otherwise `-1` is returned. If *arg* is "close", the file is closed and zero is returned. If *arg* is zero, the first line of the file is returned. If only the first argument is present, the next line of the file is read and returned. At the end of the file, a `-1` is returned.

The previous file, if any, is closed and the new file is opened automatically when the filename argument changes (at least in this preliminary implementation).

`dofile(file [, line_num|search_pattern])` – Queues the file *file* for reading commands. *file* must be a string constant or expression. Returns nonzero if the file doesn't exist or permit read access. As of release 4.05.01, an optional second argument can specify a line number or a text pattern that will be used to locate the point in the file to begin reading. If the argument is an integer, the number specifies at which line to start reading the file. (Currently, only positive integers are allowed.) If the argument is anything else, it is considered a search string, and text is read from the file starting at the first line containing that search string. The metacharacters `*`, which matches any string, and `?`, which matches any single character, are allowed in the search string. Initial and trailing white space is ignored in the file.

`qdofile(file)` – As above, but does not show the contents of the file on the screen as the file is read.

Normally, most errors that occur while reading from a command file cause `spec` to reset to the level of the main interactive prompt with any open command files then closed. However, as of release 4.05.01, the `spec_par()` option "keep_going" is available to override that behavior. See `spec_par()` on page 78.

When a command file is opened within a statement block, the source of input isn't switched to the command file until all the commands in the statement block are executed. Thus it isn't possible to execute commands from a command file within a statement block. Note, though, the `getline()` (on page 83) function is available to scan strings from files.

When multiple command files are queued on a single command line, the input source can only change after the current line is exhausted, as the following example demonstrates:

```
9.FOURC> print "hi";dofile("file1");dofile("file2");print "bye"
hi
bye

FOURC.2> Reading "file2".
print "This is text from file2"
This is text from file2

FOURC.1> Reading "file1".
print "This is text from file1"
This is text from file1

10.FOURC>
```

Here *file1* contains the single line `print "This is text from file1"`, and *file2* contains `print "This is text from file2"`.

Keyboard Input

`input()` – Reads a line of input from the keyboard. Leading white space and the trailing newline are removed and the string is returned. Returns the null string "" if only white space was entered. Example:

```

10.FOURC> def change_it '{
11.quot>     local it
12.quot>     printf("Change it? ");
13.quot>     if ((it=input()) != "")
14.quot>         change_mac(it)
15.quot> }'

16.FOURC>

```

`input(prompt)` – As above, but prompts with the string *prompt*. Examples:

```

16.FOURC> input("Hit return when ready ... ")
Hit return when ready ... <return>

17.FOURC>

```

`input(n)` – This function behaves differently depending on whether the input source is the keyboard or a pipe from another program (where `spec` is invoked with the `-p fd pid` option, with nonzero *fd*.)

In the usual case, if *n* is less than or equal to zero, the tty state is set to “cbreak” mode and input echo is turned off. Then `input()` checks to see if the user has typed a character and immediately returns a null string if nothing has been typed. Otherwise, it returns a string containing the single (or first) character the user typed. If *n* is less than zero, the cbreak, no-echo mode remains in effect when `input()` returns. If *n* is greater than zero, the normal tty state is restored (as it is also if there is an error, if the user types `^C` or if the user enters the `exit` command). Also, no characters are read and the null string is returned. The normal state is also restored before the next main prompt is issued, whether due to an error, a `^C`, or through the normal flow of the program.

On the other hand, when `spec` is invoked with the `-p fd pid` option, with nonzero *fd*, `input()` reads nothing but does return the number of characters available to be read. If *n* is nonzero, `input()` simply reads and returns a line of text, as if it had been invoked with no argument.

`yesno(val)` – Reads a line of input from the keyboard. The function returns 1 if the user answers with a string beginning with `Y`, `y` or `1`. The value of *val* is returned if the user simply enters return. Otherwise the function returns zero.

`yesno(prompt, val)` – As above, but prompts the user with the string *prompt*. The characters “ (YES)? ” are appended to the prompt string if *val* is nonzero. Otherwise the characters “ (NO)? ” are added.

`getval(val)` – Reads a line of input from the keyboard. If the user enters a value, that value is returned. The value of *val* is returned if the user simply enters

return. The function works with both number and string values.

`getval(prompt, val)` – As above, but prompts the user with the string *prompt*. The string is printed followed by the current value of *val* in parenthesis, a question mark and a space. For example,

```
17.FOURC> DATAFILE = getval("Data file", DATAFILE)
Data file (pt100.133)? <return>

18.FOURC>
```

Text Output

`print a [, b ...]` – Prints the string value of each argument, adding a space between each string. If the argument is an associative array, each element of the array is printed in a list, as in:

```
18.FOURC> print mA
mA["0"] = 2
mA["1"] = 0
mA["2"] = 1
mA["3"] = 3

19.FOURC>
```

If the argument is a data array, the contents of the array are printed in a compressed format, as in:

```
19.FOURC> array data[64][64]; data[0] = 1; data[1] = 2
20.FOURC> print data
{{1 <64 repeats>}, {2 <64 repeats>}, {0 <64 repeats>} <62 repeats>}

21.FOURC>
```

`eprint a [, b ...]` – As above, except that if an error-log file is open, the generated string will also be written to that file prefixed by the #E characters.

`printf(format [, a ...])` – Does formatted printing on the turned-on output devices. *format* contains the format specifications for any following arguments. See the description of *printf()* in any C-language manual. Returns **true**.

`eprintf(format [, a ...])` – As above, except that if an error-log file is open, the generated string will also be written to that file prefixed by the #E characters.

`fprintf(file_name, format [, a ...])` – Does formatted printing on *file_name*. All other devices (except *log* files) are turned off while the string is printed. The specified file is opened, if necessary, and remains open until closed with the `close()` function.

`tty_cntl(cmd)` — Sends terminal-specific escape sequences to the display. The sequences are only written to the "tty" device and only if it is turned on for output. The sequences are obtained from the system terminal-capability data base using the value of the environmental variable `TERM`. The following values for `cmd` are recognized:

- "ho" — Move the cursor to the home position (upper left corner).
- "cl" — Clear the screen.
- "ce" — Clear to the end of the line.
- "cd" — Clear from current position to the end of the screen.
- "so" — Start text stand-out mode.
- "se" — End text stand-out mode.
- "md" — Start bold (intensified) mode.
- "me" — End bold mode.
- "us" — Start underline mode.
- "ue" — End underline mode.
- "mb" — Start blink mode. (Note, xterms don't blink.)
- "mh" — Start half-bright mode.
- "mr" — Start reverse video mode.
- "up" — Move cursor up one line.
- "do" — Move cursor down one line.
- "le" — Move cursor left one space.
- "nd" — Move cursor right one space (nondestructive).
- "resized?" — A special option that updates the `ROWS` and `COLS` variables in the event the window size has changed and returns a nonzero value if the window size has changed since the last call to `tty_cntl("resized?")`.

Returns **true** if `cmd` is recognized, otherwise returns **false**.

`tty_move(x, y [, string])` — Moves the cursor to column `x` and row `y` of the display, where column 0, row 0 is the upper left corner of the screen. If the third argument `string` is present, it is written as a label at the given position. The sequences and string are only written to the "tty" device and only if it is turned on for output. Special tty control sequences of the form `\[xx]`, where `xx` is one of the codes listed for the `tty_cntl()` function above, can be used with `string`. Negative `x` or `y` position the cursor relative to the left or bottom edges of the screen, respectively. Relative moves are possible by adding ± 1000 to `x` or `y` position arguments. Both coordinates must specify either relative or absolute moves. If one coordinate specifies a relative move, the absolute move in the other coordinate will be ignored. Please note, not all terminal types support relative moves. Returns **true**.

`tty_fmt(x, y, wid, string)` – Writes the string *string* to the screen starting at column *x* and row *y*, where column 0, row 0 is the upper left corner of the screen. The string is only written to the "tty" device and only if it is turned on for output. If *string* is longer than the width given by *wid*, the string is split at space characters such that no line is longer than *wid*. Newlines in the string are retained, however. The function will truncate words that are wider than *wid* and drop lines that would go off the bottom of the screen. Special tty control sequences of the form `\[xx]`, where *xx* is one of the codes listed for the `tty_cntl()` function above, can be used with *string*. Negative *x* or *y* position the cursor relative to the left or bottom edges of the screen, respectively. The function returns the number of lines written.

Binary Input/Output

The facility for binary file input and output allows users and sites to create arbitrary binary file formats for writing and reading `spec` data arrays. C source code for a number of formats is included in the `spec` distribution.

`fmt_read(file, fmt, arr [, header [, flags]])` –

`fmt_write(file, fmt, arr [, header [, flags]])` –

`fmt_close(file, fmt)` –

In these functions *file* is the name of the data file, *fmt* selects which format to use and *arr* is the data array. The optional *header* argument is an associative array that may contain identifying information to be saved with the binary values in the data array. The optional *flags* argument is reserved for future enhancements.

Commands For Variables

`global name ...` – Declares *name* to be a global symbol. A global symbol retains its value after each parsed program is executed. If *name* is used as an array name, each element of the array is global. By appending empty square brackets to *name* the type of the symbol can be forced to be an associative array, which may be useful if *name* is to be used as an argument to a macro function before its type has been established by usage.

`unglobal name ...` – Makes the global or constant symbol *name* no longer global.

`constant name expression` – Declares *name* to be a constant, global symbol having the value given by *expression*. A constant symbol cannot be changed by assignment.

`local name ...` – Allows reuse of a preexisting name and gives the new instance of that name scope only within the statement block in which it is defined. The name may be that of a macro, in which case the macro definition is unavailable within the statement block. By appending empty square brackets to *name* the type of the symbol can be forced to be an associative array, which may be useful if *name* is to be used as an argument to a macro function before its type has been established by usage.

`delete assoc-array[elem] ...` – Removes the element *elem* of the associative array *assoc-array*.

`[[extern] shared] [type] array var[cols]` – Declares a one-dimensional data array.

`[[extern] shared] [type] array var[rows][cols]` – Declares a two-dimensional data array.

(Text forthcoming. See the on-line *arrays* help file in the meantime.)

`syms [-v] [+|-BGLADNSIC] [pattern ...]` – Lists *spec*'s current variables. Without arguments, all the variables are listed, along with their memory consumption and type. With the `-v` flag, the variables are listed along with their values in a format that can be saved to a file and read back as commands. If arguments are given as *pattern*, only symbols matching the arguments are printed. Such arguments may contain the `?` and `*` metacharacters.

In addition, the type of symbols listed can be controlled using the flags in the following table where a `-` flag prevents symbols with the given attribute from being listed and a `+` flag includes symbols with the given attribute in the list.

B	Built-In
G	Global
L	Local
A	Associative array
D	Data array
N	Number type
S	String type
I	Immutable attribute
C	Constant attribute

Commands For Macros

Built-In Commands

`def name string` – Defines a macro named *name* to be *string*. Each time *name* occurs on input, it is replaced with *string*. The definition is made immediately, so the macro can be used later in the same statement block in which it is defined and can be redefined within the same statement block.

Note that the macro definition is made regardless of any surrounding flow control statements, since the enclosing mini-program is not yet completely parsed and is not executing.

`rdef name expression` – Defines a macro named *name* to be *expression*, which is almost always a string constant. Each time *name* occurs on input, the value *expression* is substituted. Unlike `def`, described above, the macro definition is not made until all the encompassing statement blocks are parsed and the resulting mini-program is executed. Consider the following example.

```

if (flag == 1)
    rdef plot "onp;offt;lp_plot;ont;plot_res;offp"
else if (flag == 2)
    rdef plot "splot;onp;offt;lp_plot;plot_res;ont;offp"
else if (flag == 3)
    rdef plot "onp;plot_res;offp"
else
    rdef plot ""

```

Clearly, it is necessary for the mini-program to be parsed and executed to decide which is the appropriate definition to assign to the `plot` macro.

`prdef` – Displays all macro definitions. The displayed definitions are prepended with `def name `` and terminated with ``` so if saved to a file, the definitions can be read back. (See the standard macro `savmac` on page 132.)

`prdef pattern ...` – As above, except only macro names matching *pattern* are listed, where *pattern* may contain the metacharacters `?` or `*`.

`lsdef` – Lists the name and the number of characters in each macro definition.

`lsdef pattern ...` – As above, except only macro names matching *pattern* are listed, where *pattern* may contain the metacharacters `?` or `*`.

`undef name ...` – Removes the named macros, which can be ordinary macros, macro functions or `cdef()` macros.

`cdef("name", string ["key" [, flags]])` – Defines parts of *chained* macros. A chained macro definition is maintained in pieces that can be selectively included to form the complete macro definition. The argument *name* is the name of the macro. The argument *string* contains a piece to add to the macro.

The chained macro can have three parts: a beginning, a middle and an end. Pieces included in each of the parts of the macros are sorted lexicographically by the keys when putting together the macro definition. Pieces without a key are placed in the middle in the order in which they were added, but after any middle pieces that include a key.

The *key* argument allows a piece to be selectively replaced or deleted, and also controls the order in which the piece is placed into the macro definition. The *flags* argument controls whether the pieces are added to the beginning or to the end of the macro, and also whether the pieces should be selectively included in the definition depending on whether *key* is the mnemonic of a configured motor or counter.

The bit meanings for *flags* are as follows:

- 0x01 – only include if key is a motor mnemonic
- 0x02 – only include if key is a counter mnemonic
- 0x10 – place in the beginning part of the macro
- 0x20 – place in the end part of the macro

If *flag* is the string "delete", the piece associated with *key* is deleted from the table. If the name is the null string, the piece associated with *key* is deleted from all the chained macros. If *key* is the null string, the *flags* have no effect.

The `cdef()` function will remove any existing macro defined using `def` or `rdef`. However, the commands `lsdef`, `prdef` and `undef` will function with chained macros.

When `spec` starts and when the `reconfig` command is run (or the `config` macro is invoked), all the chained macros are adjusted for the currently

configured motors and counters.

`cdef("?")` – Lists all the pieces of all the chained macros.

`cdef("name", "", "?")` – Lists the pieces of the macro named *name*.

Built-In Macro Names

The following macro names are built-in to `spec`. They are run at the specified times only if they have been given a definition.

`begin_mac` – If a macro by this name exists, it will be run after reading the hardware configuration file and all the start-up command files, but before reading commands from the keyboard.

`end_mac` – If a macro by this name exists, it will be run when `spec` exits from either a `^D` or a `quit` command.

`config_mac` – If a macro by this name exists, it will be run after reading the configuration file at start up and after the `reconfig` command is executed.

`prompt_mac` – If a macro by this name exists, it will always be run just before `spec` issues the main, level-zero prompt. If an error occurs while running `prompt_mac`, it will be automatically undefined.

`cleanup`, `cleanup1` – If either or both exists, they will be run whenever an error is encountered, the `exit` command is encountered, or a user types `^C`. The `cleanup` macros is run first. After running the clean-up macros, `spec` gives the standard prompt and waits for the next command from the keyboard.

`cleanup_once` – A clean-up macro (added in release 4.03.13) that is always deleted before a new `spec` main prompt is issued. If defined, its definition will be pushed on to the input stream whenever an error is encountered, the `exit` command is encountered, or a user types `^C`.

`cleanup_always` – Like `cleanup_once`, but its definition is not removed except by an explicit `undef` command.

Macro Arguments

Within ordinary macros (not macro functions), the symbols `$1`, `$2`, . . . , are replaced by the arguments with which the macro is invoked. Arguments are defined as strings of characters separated by spaces. Also, quotes can be used to include more than one space-separated string in a single argument. Use `\"` or `\'` to pass literal quotes. Arguments can be continued over more than one line by putting a backslash

at the end of the line.

Within a macro the following substitutions are also made:

- \$0 is replaced with the macro name,
- \$* is replaced with all the arguments,
- \$# is replaced with the number of arguments,
- \$\$ is a literal \$.

When a macro defined without arguments is invoked, only the macro name is replaced with the definition. On the other hand, when a macro defined with arguments is invoked, all characters on the input line up to a ;, a { or the end of the line are eaten up, whether or not the macro uses them as arguments.

When numbered arguments are referred to in the macro definition, but are missing when the macro is invoked, they are replaced with zeros. If \$* is used in the definition and there are no arguments, no characters are substituted.

As of spec release 3.03.10, there is no limit on the length of macro definitions. Macros can have a maximum of 25 arguments. The total combined size of all macro definitions is currently limited to 491,520 characters.

Beware of unwanted side affects when referencing the same argument more than once. For example, if i++ is an argument referenced more than once in a macro definition, the result is probably not what was desired.

Macro Functions

Macro functions (introduced in release 4) are a type of macro that can return values and can be used as an operand in expressions. The macro definition can include function arguments, which then become available to the body of the macro function. For example,

```
def factorial(n) '{
    if (n <= 1)
        return(1);
    return(n * factorial(n-1))
}'
```

The syntax of macro functions requires the macro name followed by a set of parenthesis which can contain a comma-separated list of argument names. The arguments names become local variables within the macro definition. The definition must be a statement block, that is, the statements must be enclosed in curly brackets.

Hardware Functions and Commands

MCA and Other Data Acquisition

`mca_sel(n)` – Selects which MCA-type device to use with subsequent `mca_get()`, `mca_put()` and `mca_par()` commands. The numbering of MCA-type devices is set in the *config* file. Returns `-1` if not configured for device *n*, otherwise returns zero. It is not necessary to use `mca_sel()` if only one MCA-type device is configured. The selected MCA-type device does not change when reading the *config* file with the `reconfig` command.

`mca_sel("?")` – Lists the configured MCA devices, with the currently selected device marked with an asterisk.

`mca_par(par [, val])` – A device-dependent function to access various features and parameters of the currently selected MCA-type device. The string *par* selects an option. The argument *val* contains an optional numeric value. See the help file for the particular device for implemented options and return values.

`mca_get(grp, elem)` or `mca_get(array)` – Gets data from the currently selected MCA-type device, and transfers it to element *elem* of data group *grp* or to the elements of the data array *array*. Generally returns the number of points read or `-1` for failure.

`mca_put(grp, elem)` or `mca_put(array)` – Sends data from data group *grp*, element *elem* or from the data array *array* to the currently selected MCA-type device. Generally returns the number of points written or `-1` for failure.

`mca_spar(sel par [, val])` – As `mca_par()` above, but selects which MCA device with the *sel* argument.

`mca_sget(sel grp, elem)` or `mca_sget(sel array)` – As `mca_get()` above, but selects which MCA device with the *sel* argument.

`mca_sput(sel grp, elem)` or `mca_sput(sel array)` – As `mca_put()` above, but selects which MCA device with the *sel* argument.

CAMAC (IEEE-583) Hardware Functions

CAMAC functions are available if the appropriate hardware devices and software drivers have been installed on the computer. The *config* file describes which CAMAC hardware is installed. Refer to the *Administrator's Guide* for information on the supported CAMAC controllers and how to install the corresponding drivers.

CAMAC modules are programmed with *FNA* codes where *F* is a function code, *N* is the slot number and *A* is a subaddress number. Slot numbers are assigned in the

config file. Built-in code for the specialized CAMAC devices used for controlling motors, clocks and scalers is accessed through commands such as `mcount()`, `tcoun`t() and `move_all`. However, simple devices, such as input or output registers, can be accessed directly by the user. These modules are also assigned slot numbers in the *config* file. They are also given device numbers, starting at 0, that are used in the following functions.

`ca_get(device, address)` – The CAMAC module having device number *device*, as set in the *config* file, is read using $F = 0$ and $A = address$ with the 24-bit value so obtained returned. Resets to command level if not configured for *device*.

`ca_put(x, device, address)` – This function is similar to `ca_get()` above, except the 24-bit value *x* is written using $F = 16$. The actual number written is returned, which is the 24-bit integer representation of *x*. Resets to command level if not configured for *device*.

`ca_fna(f, n, a [, v])` – Sends the arbitrary FNA command to the module in slot *n*. If the dataway command given by *f* is a write function, the 24-bit value to be written is contained in *v*. If the dataway command given by *f* is a read command, the function returns the 24-bit value obtained from the module. The user should avoid issuing commands that would cause a LAM and should certainly avoid issuing commands to slots that are being used for motor or counter control by *spec*'s internal hardware code.

`ca_cntl(cmd, [, arg])` – Performs the selected CAMAC crate command according to the parameter *cmd*, as follows:

"Z" or "init" – performs a crate initialize (reset).

"C" or "clear" – performs a crate clear.

"inhibit" – set crate inhibit if *arg* is 1 and clears crate inhibit if *arg* is 0.

During normal operation, you should not need to issue these commands. You should probably issue a `reconfig` after sending a crate initialize or clear.

GPIB (IEEE-488) Hardware Functions

GPIB functions are available if the appropriate hardware and software drivers have been installed on the computer. Information in the *config* file describes the particular GPIB hardware in use. Refer to the *Administrator's Guide* for information on the supported GPIB controllers and how to install the corresponding drivers.

spec allows up to four GPIB controllers to be in use at once. The controller unit numbers are set in the *config* file. For the functions below, there are two methods by which the unit number can be specified. If no unit number is specified, the default, unit 0, is used. The first method of addressing is of the form "unit:addr" where the

quotes are required. Alternatively, the unit number can be coded in the GPIB address as $unit \times 100 + addr$.

`gpib_get(addr)` – Returns a string from the GPIB device with address `addr`. The device must terminate the string with either a newline (`\n`) or a carriage return and a newline (`\r\n`). In either case, the terminator is removed before the string is returned. At most, 8,192 characters can be read at a time. (That limit is as of `spec` release 4.05.09 – the limit was 512 from release 4.05.02 and 255 since release 2.14.)

`gpib_get(addr, n)` – As above, but reads `n` bytes and does not look for or remove the terminator.

`gpib_get(addr, s)` – As above, but tries to read up to the terminator given by the first character of the string `s`, except for the special cases described below. The terminator is removed.

`gpib_get(addr, mode)` – If `mode` is the string "byte", reads and returns one unsigned binary byte. The following modes read short or long integers and work the same on both *big-endian* and *little-endian* platforms. If `mode` is the string "int2" reads two binary bytes and returns the short integer so formed. If `mode` is the string "int4" reads four binary bytes and returns the long integer so formed. By default, the incoming data is treated as *big endian*. If the incoming data is *little endian*, use "int2_swap" or "int4_swap".

`gpib_put(addr, string)` – Writes the string `string` to the device with GPIB address `addr`. Returns the number of bytes written. The length of the string is not limited, but null bytes cannot be sent.

`gpib_poll(addr)` – Returns the serial-poll status from the device with address `addr`.

`gpib_cntl(addr, cmd)` – Performs the selected GPIB command on the device with address `addr`. The string `cmd` is one of the following:

"gtl" – Go to local.

"llo" – Local lockout.

"sdc" – Selected device clear.

"dcl" – Device clear (sent to all devices).

"get" – Group execute trigger (sent to addressed device).

"responsive" – This special command will return a nonzero value if the GPIB controller (not device) associated with `addr` is present in `spec`'s configuration file and responds to `spec`'s presence test.

Returns zero if command is sent successfully, Otherwise returns `-1`.

RS-232 Serial Interfaces

Serial functions are available on all the systems `spec` supports. The device nodes are selected in the `config` file, where they are assigned `spec` device numbers starting from zero that are used in the functions below.

`ser_get(dev_num)` – Reads a string of characters, up to and including a newline, from serial device `dev_num` and returns the string so read. The routine will not return until the read is satisfied. Use `ser_get(dev_num, 0)` to read up to a newline with a timeout.

`ser_get(dev_num, n)` – If `n` is greater than 0, this routine reads `n` bytes from the serial device `dev_num` and returns the string so obtained. If `n` is zero, the routine reads characters up to and including a newline. If the read is not satisfied within two seconds, the routine returns the null string.

`ser_get(dev_num, mode)` – If `mode` is the string "byte", reads and returns one unsigned binary byte. The following modes read short or long integers and work the same on both *big-endian* and *little-endian* platforms. If `mode` is the string "int2" reads two binary bytes and returns the short integer so formed. If `mode` is the string "int4" reads four binary bytes and returns the long integer so formed. By default, the incoming data is treated as *big endian*. If the incoming data is *little endian*, use "int2_swap" or "int4_swap". If the device returns less than the required number of bytes, the characters so read are discarded, and `ser_get()` returns the null string. The tty mode for the device needs to be selected as `raw` in the `config` file for this command to work properly.

`ser_put(dev_num, string)` – Writes the string `string` to the serial device having device number `dev_num`, as set in the `config` file. Returns the number of bytes written.

`spec` uses the standard `tty(4)` driver interface on the UNIX systems. Serial line characteristics are set according to the device modes set in the `config` file. Available modes are either *raw* or *cooked*, with *cooked* mode also having *noflow*, *igncr* (no-op on non-System V) and *evenp* or *oddp* options. The baud rate for the serial line is also taken from the `config` file.

`ser_par(dev_num, par [, val])` – Sets parameters for the serial device `dev_num`. Values for `par` include:

"timeout" – Sets or returns the current value for the timeout for reads from *cooked* serial lines. The units are seconds. If `val` is zero, timeouts are disabled. If `val` is less than zero, the default timeout of two seconds is set.

PC Port I/O

The port I/O functions allow arbitrary access to I/O ports on a PC computer, and therefore, should be used with caution. To lessen the chance of writing data to a port that might damage the computer or corrupt data, valid addresses for the following functions must be explicitly assigned in the *config* file.

`port_get(addr)` – Reads one byte from the PC I/O port with the address *addr*. Resets to command level if the port has not been selected in the *config* file. Otherwise, returns the byte read.

`port_getw(addr)` – As above, but reads and returns a two-byte value.

`port_put(addr, byte)` – Writes the byte *byte* to the PC I/O port with the address *addr*. Resets to command level if the port has not been selected for writing in the *config* file. Otherwise, returns zero.

`port_putw(addr, word)` – As above, but writes the two-byte value *word* to the I/O port.

VME Hardware Functions

The type of data access and/or VME address modifier for the following functions can be selected with the optional argument *dmode* as follows (if more than one option is needed, make a comma-separated list in the single string argument):

"D8" – byte access

"D16" – short-word access

"D32" – long-word access, but only available with `vme_get32()` and `vme_put32()`.

"DPRT" – use the address modifier appropriate for dual-port memory access, on adapters that support it.

"amod=0xXX" – specify the hexadecimal value for the address modifier.

The default mode for the A16 access functions `vme_get()` and `vme_put()` is D8. The default mode for the A32 access functions `vme_get32()` and `vme_put32()` is D32. Not all VME adapters supported by *spec* support A32 access. There are currently no functions for A24 access.

`vme_get(addr [, dmode])` – Returns the data at *addr* in the 64K A16 address space.

`vme_put(addr, data [, dmode])` – Writes *data* to *addr* in the 64K A16 address space.

`vme_get32(addr [, dmode])` – Returns the data at *addr* in the A32 address space.

`vme_put32(addr, data [, dmode])` – Writes *data* to *addr* in the A32 address space.

`vme_move(from, to [cnt [, dmode]])` – Copies data between a *spec* data array and VME A32 address space. One of the *from* and *to* arguments must be the name of a *spec* data array while the other must be a VME address. If the optional argument *cnt* is present, it designates how many data items (not bytes) to copy. If missing or zero, the number of elements in the array is copied.

Miscellaneous

`reconfig` – Reconfigures the hardware. This command obtains any modified configuration information, including hardware devices and types, CAMAC slot assignments and motor parameters and settings, from the *config* and *settings* files. The sequence of events is as follows:

First, *spec* waits for all asynchronous activity (moving and counting) to finish. It then does a *sync* of the motor controller registers, comparing them with the internal program positions. Next, all open devices are closed. The *config* file is then read to obtain the configuration information, and the program opens and possibly does hardware presence tests on the selected devices. Finally, the internal program motor positions are updated from the *settings* file and then resynchronized with the motor hardware.

`set_sim(how)` – If *how* is 0, simulate mode is turned off. If *how* is 1 (or positive), simulate mode is turned on. In either case the program waits for moving and counting to finish before changing the mode, and the function returns the previous mode (0 or 1). If *how* is -1 (or negative) the function returns the current value of simulation mode without changing it. Whenever simulation mode is turned off, the motor settings file is reread to restore the motor positions. Simulation mode cannot be turned off if *spec* was invoked with the `-s` flag.

`wait()` – Waits for all asynchronous activity to complete. Returns **true**.

`wait(flag)` – If *flag* has bit one set (1), waits for moving to finish. If *flag* has bit two set (2), waits for counting to finish. If *flag* has bit three set (4), waits for other data acquisition (multi-channel scaling, for example) to finish. If *flag* has bit five set (32), the function returns **true** if any of the activities flagged by the first three bits are active, otherwise it returns **false**. Thus `wait(35)` (or `wait(0x23)`) returns **true** if either counting or moving is active.

`stop()` – Stops all asynchronous activity. Returns **true**.

`stop(flag)` – If *flag* has bit one set (1), stops all motors that are moving. If *flag* has bit two set (2), stops the timer, counters and any other data acquisition (multi-channel scaling, for example).

Counting

`mcount(counts)` – Starts the timer/clock counting for *counts* monitor counts. Returns zero.

`tcount(t)` – Starts the timer/clock counting for *t* seconds, where *t* may be nonintegral. Returns zero.

`getcounts` – Loads the built-in array `S[]` with the contents of the scalers.

`cnt_mne(counter)` – Returns the string mnemonic of counter number *counter* as given in the configuration file. (Mnemonics are, at most, 7 characters long.) Resets to command level if not configured for *counter*.

`cnt_name(counter)` – Returns the string name of counter number *counter* as given in the configuration file. (Names are, at most, 15 characters long.) Returns "?" if not configured for *counter*.

`cnt_num(mne)` – Returns the counter number corresponding to the counter mnemonic *mne*, or -1 if there is no such counter configured.

`counter_par(counter, par [, val])` – Returns or sets parameters associated with counter number *counter* as given in the configuration file. The following *par* arguments are supported for all counters:

"unit" – returns the unit number of the indicated counter.

"channel" – returns the channel number of the indicated counter.

"scale" – returns the value of the scale factor set in the *config* file for the indicated counter.

"responsive" – returns nonzero if the hardware appears to be working for the indicated counter.

"controller" – returns a string that indicates the controller type of the indicated counter.

"disable" – returns a nonzero value if the counter has been disabled by software. If *val* is given and is nonzero, then the counter is disabled. If *val* is given and is zero, the counter becomes no longer disabled. A disabled counter channel will not be accessed by any of `spec`'s counting commands. Any `cdef()`-defined macros will automatically exclude the portions of the macro keyed to the particular counter when the counter is software disabled.

In addition, device-dependent values for *par* are available for specific counter models. See the *Hardware Reference* for values for specific controllers.

The counting functions `mcount()` and `tclock()` both program the timer/clock for a specified count time. Before the count period begins, both functions clear and enable all configured scalars and MCAs. The routines return immediately, as counting is asynchronous. Use `wait()`, described below, to determine if counting has been completed. A `^c` will halt the timer/clock.

When the count time has expired, or counting is aborted by a `^c`, the scalars and MCAs are disabled. Normally, enable signals from the timer/clock are used to directly gate the scalars or MCAs. Software gating takes place whether or not hardware gating is in place and can be used in lieu of hardware gating, although the interval over which the gating occurs will not be as precisely controlled.

Controlling Motors

`move_all` — This command sets motors in motion. The sequence of events is as follows. For some motor controllers, `spec` first examines the controller registers of all nonbusy motors and makes sure the contents agree with the current positions in program memory. If there is a discrepancy, the user is asked to choose the correct position. Next, `spec` prepares to move all motors that are not already at the positions specified in the built-in `A[]` array, interpreted in user units. A motor will not be moved if it is currently moving or if it is marked as protected (via the configuration file) or unusable (due to failure of a hardware presence test). If the target position of any of the motors is outside the software limits, the entire move is canceled, and the program resets to command level. Otherwise, the motors are started, and the command returns.

The sequence of commands when using `move_all` should almost always be,

```
wait(1)           # Wait for moving to finish
read_motors(0)    # Put current positions of all motors in A[]
(Assign new values to elements of A[] to be moved)
move_all          # Move to those positions
```

If `read_motors()` is called before the motors have stopped, the values in `A[]` will reflect the motor positions before they stopped. If `read_motors()` is not called at all, or if you do not explicitly assign a value to each element of `A[]`, then you will not know for sure where some motors will be going when `move_all` is called.

A `^c` halts moving, as does the `sync` command.

`move_cnt` — This command is similar to `move_all`, described above, but with the following differences. Just before the motors are started, the clock/timer is

enabled and programmed to gate the scalers with a longer-than-necessary count time. The motors are then started at the base rate set in the *config* file, but are not accelerated to the steady-state rate. No backlash correction is done at the end of the move. When the move is completed, the clock/timer is stopped. The `move_cnt` command is used in powder-averaging scans. (See the powder-mode macros on page 148.)

`sync` – If any motors are moving, they are halted. The motor positions maintained by the motor controller are then compared with the motor positions currently set in the program. If there is a discrepancy, the user is asked which should be changed. The `sync` command is used to place the motor hardware in a known state and is supposed to fix any problems in communicating with the controllers.

`motor_mne(motor)` – Returns the string mnemonic of motor number *motor* as given in the configuration file. (Mnemonics are, at most, 7 characters long.) Resets to command level if not configured for *motor*.

`motor_name(motor)` – Returns the string name of motor number *motor* as given in the configuration file. (Names are, at most, 15 characters long.) Returns "?" if not configured for *motor*.

`motor_num(mne)` – Returns the motor number corresponding to the motor mnemonic *mne*, or -1 if there is no such motor configured.

`motor_par(motor, par [, val])` – Returns or sets configuration parameters for motor *motor*. Recognized values for the string *par* follow. Note, not all parameters are meaningful for all motor controllers.

"step_size" – returns the current step-size parameter. The units are generally in steps per degree or steps per millimeter. If *val* is given, then the parameter is set to that value, but only if changes to the step-size parameter have been enabled using `spec_par("modify_step_size", "yes")`.

"acceleration" – returns the value of the current acceleration parameter. The units of acceleration are the time in milliseconds for the motor to accelerate to full speed. If *val* is given, then the acceleration is set to that value.

"base_rate" – returns the current base-rate parameter. The units are steps per second. If *val* is given, then the base rate is set to that value.

"velocity" – returns the current steady-state velocity parameter. The units are steps per second. If *val* is given, then the steady-state velocity is set to that value.

"backlash" – returns the value of the backlash parameter. Its sign and magnitude determine the direction and extent of the motor's backlash

- correction. If *val* is given, then the backlash is set to that value. Setting the backlash to zero disables the backlash correction.
- "config_step_size" – returns the step-size parameter contained in the *config* file.
- "config_acceleration" – returns the acceleration parameter contained in the *config* file.
- "config_velocity" – returns the steady-state velocity parameter contained in the *config* file.
- "config_base_rate" – returns the base-rate parameter contained in the *config* file.
- "config_backlash" – returns the backlash parameter contained in the *config* file.
- "controller" – returns a string containing the controller name of the specified motor. The controller names are those used in *spec*'s *config* files.
- "unit" – returns the unit number of the specified motor. Each motor controller unit may contain more than one motor channel.
- "channel" – returns the channel number of the specified motor.
- "responsive" – returns a nonzero value if the motor responded to an initial presence test or appears otherwise to be working.
- "active" – returns a nonzero value if the motor is currently moving.
- "disable" – returns a nonzero value if the motor has been disabled by software. If *val* is given and is nonzero, then the motor is disabled. If *val* is given and is zero, the motor becomes no longer disabled. A disabled motor channel will not be accessed by any of *spec*'s commands, and, of course, cannot be moved. Any *cdef()*-defined macros will automatically exclude the portions of the macro keyed to the particular motor when the motor is software disabled.
- "slop" – returns the value of the slop parameter. If *val* is given, sets the slop parameter. When this parameter is present, discrepancies between hardware and software motors positions are silently resolved in favor of the the hardware when the number of steps in the discrepancy is less than the value of the slop parameter. (Not yet implemented for all motor controllers.)
- "home_slew_rate" – returns the value of the home-slew-rate parameter. If *val* is given, sets the parameter. This parameter is the steady-state velocity used during a home search. (Only available for selected controllers.)
- "home_base_rate" – returns the value of the home-base-rate parameter. If *val* is given, sets the parameter. This parameter is the base-rate velocity used during a home search. (Only available for selected controllers.)

"home_acceleration" – returns the value of the home-acceleration parameter. If *val* is given, sets the parameter. This parameter is the acceleration used during a home search. (Only available for selected controllers.)

"dc_dead_band" – returns the value of the dead-band parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_settle_time" – returns the value of the settle-time parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_gain" – returns the value of the gain parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_dynamic_gain" – returns the value of the dynamic-gain parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_damping_constant" – returns the value of the damping-constant parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_integration_constant" – returns the value of the integration-constant parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_integration_limit" – returns the value of the integration-limit parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_following_error" – returns the value of the dc-following parameter for certain DC motors. Sets the parameter if *val* is given.

"dc_sampling_interval" – returns the value of the sampling-interval parameter for certain DC motors. Sets the parameter if *val* is given.

"encoder_step_size" – returns the value of the encoder step size parameter. Sets the parameter if *val* is given.

"step_mode" – returns the value of the step-mode parameter. Sets the parameter if *val* is given. A zero indicates full-step mode, while a one indicates half-step mode.

"deceleration" – returns the value of the deceleration parameter. Sets the parameter if *val* is given.

"torque" – returns the value of the torque parameter. Sets the parameter if *val* is given.

Rereading the *config* file resets the values of all the motor parameters to the values in the *config* file. Little consistency checking is done by `spec` on the values programmed with `motor_par()`. You must be sure to use values meaningful to your particular motor controller.

In addition, device-dependent values for *par* are available for specific motor controllers. See the *Hardware Reference* for values for specific controllers.

`dial(motor, user_angle)` – Returns the dial angle for motor *motor* corresponding to user angle *user_angle* using the current *offset* between user and dial angles for *motor*. The value returned is $(user_angle - offset) / sign$, where *sign* is ± 1

and is set in the *config* file. The value is rounded to the motor resolution as set by the step-size parameter in the *config* file. Resets to command level if not configured for motor *motor*.

`read_motors(how)` – Reads the current motor positions from the motor controllers and places the values in the `A[]` array, depending on the value of the argument *how*. If bit 1 is set, the function returns dial values, otherwise user values are returned. If bit 2 is set, a forced read of all hardware takes place. (For efficiency, normally most motor controllers are not read if the position hasn't been changed by a move.) If bit 3 is set and if there is a discrepancy between the software and hardware, the software will be silently corrected to match the hardware. Note, the forced-read and “silent-sync” features are not yet implemented for all motor controllers. Check the *Hardware Reference* or contact CSS for hardware-specific information.

`chg_dial(motor, dial_angle)` – Sets the dial position of motor *motor* to *dial_angle*. Returns nonzero if not configured for *motor* or if the protection flags prevent the user from changing the limits on this motor. Resets to command level if any motors are moving.

`chg_dial(motor, cmd)` – Starts motor *motor* on a home or limit search, according to the value of *cmd* as follows:

"home+" – Move to home switch in positive direction.

"home-" – Move to home switch in negative direction.

"home" – Move to home switch in positive direction if current dial position is less than zero, otherwise move to home switch in negative direction.

"lim+" – Move to limit switch in positive direction.

"lim-" – Move to limit switch in negative direction.

Positive and negative direction are with respect to the dial position of the motor. (At present, most motor controllers do not implement the home or limit search feature.) Returns `-1` if not configured for motor *motor* or if the motor is protected, unusable or moving, else returns zero.

`get_lim(motor, flag)` – Returns the dial limit of motor number *motor*. If *flag* > 0 , returns the high limit. If *flag* < 0 , returns the low limit. Resets to command level if not configured for *motor*.

`user(motor, dial_angle)` – Returns the user angle for *motor* corresponding to dial angle *dial_angle* using the current *offset* between user and dial angles for *motor*. The value returned is $sign \times dial_angle + offset$, where *sign* is ± 1 and is set in the *config* file. The value is rounded to the motor resolution as set by the step-size parameter in the *config* file. Resets to command level if not configured for *motor*.

`chg_offset(motor, user_angle)` – Sets the offset between the dial angle and the user angle, using the current dial position and the argument *user_angle* for motor *motor* according to the relation $user_angle = offset + sign \times dial_angle$ where *sign* is ± 1 and is set in the *config* file. Returns nonzero if not configured for *motor*. Resets to command level if any motors are moving.

`set_lim(motor, low, high)` – Sets the low and high limits of motor number *motor*. *low* and *high* are in dial units. It does not actually matter in which order the limits are given. Returns nonzero if not configured for *motor* or if the protection flags prevent the user from changing the limits on this motor. Resets to command level if any motors are moving.

Data Handling and Analysis Functions

`spec` can store data in up to 256 independent data arrays called groups. Each group is configured (see below) to have a particular number of data elements per point. For example, each point in a group could have elements for H, K, L, and detector counts. Alternatively, each point could have just one element and be used to hold data obtained from an MCA.

Groups are configured using the `data_grp()` function. A group can have up to 2048 elements per point. The maximum number of points in a group is determined by the product of the number of elements per point and the number of points. That product can be no more than 65,536, and may be slightly less depending on how the number of elements divides into 2048. The maximum number of points for all groups is 262,144. (These limits are arbitrary and are set to control the size of static data arrays and auxiliary files. If requested, CSS can make the limits larger.)

When starting `spec` for the first time or with the `-f` (fresh) flag, one data group (group 0) is configured for 4096 points, with each point consisting of two elements.

When leaving `spec`, the current data group configuration and data points are saved.

`spec` has several functions to manipulate the internal data. These functions allow unary and binary arithmetic operations, math functions and analysis operations to be performed on all the elements of a group or among elements in different groups.

In the functions described below, if an element number is negative, the element number is obtained by adding the number of elements per point in the group to the negative element number. For example, `element-1` is the last element, `element-2` is the second to last, etc.

All functions reset to command level if an invalid group, point or element is given as an argument. Functions that don't need to return anything in particular return zero.

`data_grp(grp, npts, wid)` – Configures data group *grp*. The group will have *npts* points, each having *wid* elements. If *npts* and *wid* match the previous values for the group, the data in the group is unchanged. Otherwise, the data values of the reconfigured group are set to zero. If *wid* is zero, the group is eliminated. If *npts* is zero, as many points as possible are configured. If *npts* is negative, as many points as possible, but not more than $-npts$ are configured. If *grp* is -1 , the current group configuration is displayed.

`data_info(grp, what)` – Returns data group configuration information for group *grp*, according to the the string *what*. Values for *what* are:

"npts" – the number of configured points.

"elem" – the number of configured elements.

"last" – the number of the last point added to the group.

"precision" – the number of bytes per element, either 4 or 8.

If the group number is invalid, or if the string *what* is none of the above, returns -1 .

`data_get(grp, npt, elem)` – Returns the value of element *elem* of point *npt* in group *grp*.

`data_put(grp, npt, elem, val)` – Assigns the value *val* to element *elem* of point *npt* in group *grp*.

`data_nput(grp, npt, val0 [, val1 ...])` – Assigns values to point *npt* of group *grp*. Element 0 is assigned *val0*, element 1 is assigned *val1*, etc. Not all elements need be given, although elements are assigned successively, starting at element 0.

`data_uop(g_src, e_src, g_dst, e_dst, uop [, val])` – Performs the unary operation specified by the string *uop* on element *e_src* for all points in group *g_src*. The results are put in element *e_dst* of the corresponding points in group *g_dst*. The source and destination groups and/or elements may be the same. If the number of points in the groups differ, the operation is carried out on up to the smallest number of points among the groups. Values for *uop* are:

"clr" – clear to zero.

"fill" – each element is set to point number, starting at 0.

"neg" – Negative of source.

"abs" – Absolute value of source.

"inv" – Inverse of source.

"sin" – Sine of source.

"cos" – Cosine of source.

"tan" – Tangent of source.

- "asin" — Arcsine of source.
- "acos" — Arccosine of source.
- "atan" — Arctangent of source.
- "log" — Natural logarithm of source.
- "exp" — Exponential of source.
- "log10" — Log base 10 of source.
- "pow" — The *val* power of source.
- "copy" — Value of source.
- "rev" — Reversed copy of source.
- "sqrt" — Square root of source.
- "set" — All elements set to the value of *val*.
- "contract" — Every *val* points are averaged to make a new point.
- "add" — Source plus *val*.
- "sub" — Source minus *val*.
- "mul" — Source times *val*.
- "div" — Source divided by *val*.

If any of the operations would result in an exception (divide by zero, log or square root of a negative number, etc), the operation is not performed and a count of the operations skipped is printed as an error message.

`data_bop(g0_src, e0_src, g1_src, e1_src, g_dst, e_dst, bop)` — Performs the binary operation specified by the string *bop* on elements *e0_src* and *e1_src* for all points in the groups *g0_src* and *g1_src*. The results are put in element *e_dst* for the corresponding points of group *g_dst*. The source and destination groups and/or elements may be the same. If the number of points in the groups differ, the operation is carried out on up to the smallest number of points among the groups. Values for *bop* are:

- "add" — The sum of the source elements.
- "sub" — Source 0 minus source 1.
- "mul" — The product of the source elements.
- "div" — Source 0 divided by source 1.

If the divide would result in an exception, the operation is not performed and a count of the operations skipped is printed as an error message.

`data_anal(grp, start, npts, e1_0, e1_1, op [, val])` — Performs the operations indicated by *op* on *npts* points in group *grp*, starting at point *start*. The operations use the values in element *e1_0* (if applicable) and *e1_1*. If *npts* is zero, the operations are performed on points from *start* to the last point added using `data_nput()` or `data_put()`. The values for *op* are:

"min" — Returns the minimum value of e_{l_1} . (e_{l_0} is unused.)

"max" — Returns the maximum value of e_{l_1} . (e_{l_0} is unused.)

"i_at_min" — Returns the point number of the data point with the minimum value of e_{l_1} . (e_{l_0} is unused.)

"i_at_max" — Returns the point number of the data point with the maximum value of e_{l_1} . (e_{l_0} is unused.)

"i_<=_value" — Returns the point number of the nearest data point in e_{l_1} at or below val , starting from the first point. (e_{l_0} is unused.)

"i_>=_value" — Returns the point number of the nearest data point in e_{l_1} at or above val , starting at the last point. (e_{l_0} is unused.)

"uhmx" — Returns the value in e_{l_0} corresponding to half the maximum value in e_{l_1} and at a higher index.

"lhmx" — Returns the value in e_{l_0} corresponding to half the maximum value in e_{l_1} and at a lower index.

"sum" — Returns the sum of the values in e_{l_1} . (e_{l_0} is unused.)

"fwhm" — Returns the full-width in e_{l_0} at half the maximum value of e_{l_1} .

"cfwhm" — Returns the center of the full-width in e_{l_0} at half the maximum value of e_{l_1} .

"com" — Returns the center of mass in e_{l_0} with respect to e_{l_1} . The value is the sum of the products of each e_{l_0} and e_{l_1} divided by the number of points.

"x_at_min" — Returns the value of e_{l_0} at the minimum in e_{l_1} .

"x_at_max" — Returns the value of e_{l_0} at the maximum in e_{l_1} .

"sumsq" — Returns the sum of the squares in e_{l_1} . (e_{l_0} is unused.)

The following operations treat a data group as a two dimensional data array with rows indexed by the point number and the columns indexed by the element number. The operations work on the portion of the group determined by the starting row $start$ the number of rows $npts$, the starting column e_{l_0} and the end row e_{l_1} . As usual, if $npts$ is zero, all points (rows) from $start$ to the last are considered. A negative element (column) number is added to the group width to obtain the element (column) to use.

"gmin" — Returns the minimum value.

"gmax" — Returns the maximum value.

"gsum" — Returns the sum of all values.

"i_at_gmin" — Returns the index number of the minimum value. The index number is the row number times the group width plus the element number.

"i_at_gmax" — Returns the index number, as defined above, of the maximum value.

`data_read(file_name, grp, start, npts)` – Reads data from the ASCII file `file_name`, and stuffs the data into group `grp` starting at point `start`, reading up to `npts` points. If `npts` is zero, all the points in the file are read. The values on each line of the file are assigned into successive elements for each point in the group. If there are more elements on a line in the file than fit in the group, or if there are more points in the file than in the group, the extra values are ignored. Lines beginning with the # character are ignored. Returns -1 if the file can't be opened, otherwise returns the number of points read.

`data_fit(pars, grp, start, npts, el_data, el_data [, ...])` – Performs a linear fit of the data in element `el_data` to the terms in the elements specified by `el_pars`. The fitted parameters are returned in the array `pars` supplied by the user. The function returns the *chi-squared* value of the fit, if the fit was successful. A -1 is returned if there are insufficient arguments or the covariance matrix is singular. The fit algorithm is along the same lines as the `lfit()` routine in *Numerical Recipes* (W.H. Press, et al., Cambridge University Press, 1986, page 512).

`data_plot(grp, start, npts, el_0, el_1 [, el_2 ...])` – Plots the current data in group `grp` starting at point `start` and plotting `npts` points. Element `el_0` is used for x . Elements given by the subsequent arguments (up to a maximum of 64) are plotted along the y axis. The element arguments can be combined in a single space- or comma-delimited string, which can make creation of macros to plot a variable numbers of curves in the same plot window easier.

If `npts` is zero, only the points from `start` to the last point added using `data_nput()` or `data_put()` are plotted.

If preceded by a call of `plot_cntl("addpoint")` and the ranges have not changed, only point `start + npts - 1` is drawn. If preceded by a call of `plot_cntl("addline")` the current plot will not be erased, and the plot ranges will not be changed.

The plotting area is not automatically erased by a call of `data_plot()` – use `plot_cntl("erase")` for that. The axis ranges are set using the `plot_range()` function. See `plot_cntl()` for other options that affect drawing the plot.

`plot_cntl(cmd)` – Selects built-in plotting features. The argument `cmd` is a string of comma- or space-delimited options. The following options may be preceded by a minus sign to turn the associated feature off, nothing (or an optional plus sign) to turn the feature on or a question mark to return a value of one or a zero that indicates whether the associated feature is currently on or off:

- "xlog" — Use a logarithmic *x* axis.
- "ylog" — Use a logarithmic *y* axis.
- "xexact" — Force *x*-axis minimum and maximum to be set to the scan endpoints (as opposed to being rounded).
- "colors" — Enable the use of colors.
- "dots" — Draw graphics-mode points with large dots.
- "lines" — Connect graphics-mode points with lines.
- "ebars" — Draw vertical lines through each point of length equal to the twice the square root of the *y* value.
- "persist" — Keep graphics mode on after ordinary user input.
- "perpetual" — Keep graphics mode on continuously (appropriate if using X windows, for example).

Other `plot_cntl()` options are:

- "colors=*bgnd:win:text:axis:symb:...*" — Assigns colors for drawing the various graphics-mode elements. The values for *bgnd* (the background color of the area outside the axis), *win* (the background color of the area inside the axis), *text* (the color of the text), *axis* (the color of the axis) and *symb . . .* (the color of the plotting symbols) are integers.

The first 10 colors are standardized according to the following table:

0	background	(normally white or black)
1	foreground	(normally black or white)
2	blue	3 red
4	green	5 yellow
6	cyan	7 magenta
8	white	9 black

Other colors may be available depending on the particular device. You don't have to assign values to all colors.

- "colors[*numb*]" — Returns the current color assignments, where *numb* is a number between zero and 67. Numbers zero through three return the colors assigned to the *bgnd*, *win*, *text* and *axis* elements respectively. Numbers from 4 through 67 return the colors assigned to the symbols for data elements zero through 63.
- "filter*numb*" — Selects filter number *numb*, where *numb* can be any of the numerals from 1 through 5. All plotting commands are directed to this filter. The default filter is filter 1. Each filter is associated with a separate process. On an X windows display, each filter is associated with a separate window.

- "title=*string*" — On an X windows display, the title given by *string* is used in the `XSetWMName()` and `XSetWMIconName()` calls to set the window and icon labels. With most X11 window managers, that means the title will appear in the window's title bar.
- "geometry=*widthxheight+xoff+yoff*" — With the *x11* high-resolution plot windows, sets the size and position of the window. As with the conventional X11 syntax for specifying window geometry, not all parts of the geometry string are required.
- "open" — Turn on graphics mode. If there is no graphics filter program currently active for the current filter number (see above), the filter program associated with the current `GTERM` variable is started. Recognized `GTERM` values are `vga`, `ega`, `cga`, `herc`, `x11`, and `sun`.
- "close" — Turn off graphics mode, unless the perpetual or persistent mode has been selected.
- "kill" — Turn off graphics mode and terminate graphics process.
- "erase" — Clear the graphics-mode screen (or the text screen if graphics mode is off).
- "addpoint" — Before a call to `data_plot()` will cause the plot to be made with minimal redrawing. Used for updated plotting during scans.
- "addline" — Before a call to `data_plot()` will prevent the current data from being erased and the new data from changing the ranges when the new data points are drawn. Used for plotting several data sets from different data groups on top of each other.
- "mca" — Before a call to `data_plot()` will cause the data points to be displayed using a minimal redrawing algorithm, appropriate for displaying data actively being accumulated by an MCA-type device. The "dots" and "ebars" modes must be turned off for the algorithm to work effectively.
- "lp" — Before a call to `data_plot()` will generate printing instructions appropriate for plotting on a 132-column printer.

`data_dump(grp, start, npts, el_0, [, el_1 ...] [, fmt1] [, fmt2])` — Efficiently writes elements from group *grp* to turned on output devices. The starting point is *start* and the number of points is *npts*. The elements specified by *el_0*, *el_1*, etc., are printed. If *el_0* is the string "all", all the elements for each point are printed. If *npts* is zero, only the points from *start* to the last point added using `data_nput()` or `data_put()` are printed. The element arguments can be combined in a single space- or comma-delimited string.

The optional argument *fmt1* is a string, having the format "%#", that specifies how many data points (specified by the number #) are to be printed on each line. If the number # is followed by the letter C, a backslash is added to each

continued line, appropriate for saving MCA data in manageable length lines. New versions (since May 1, 1995) of the C-PLOT *scans.4* user function interpret the continued lines correctly for MCA data. The optional argument *fmt2* is a string that specifies an alternate *printf()*-style format for the values. Only *e*, *g* and *f* formats are recognized. For example, "%15.8f" uses fixed-point format with eight digits after the decimal point and a fifteen-character-wide field. The default output format is "%g". See *printf()* in a C manual for more information. Note that in the default installation, the internal data arrays use single-precision floating values, which contain only about 8 decimal digits of significance.

`plot_move(x, y [, string [, color]])` — Moves the current position to column *x* and row *y*, where column 0, row 0 is the upper left corner of the screen. If the third argument *string* is present, it is written as a label at the given position. If using color high-resolution graphics, the fourth argument, if present, is the color to use to draw the label. The background color for the entire label will be the background color at the starting position. If graphics mode is not on, `plot_move()` works just as `tty_move()`. Returns **true**.

`plot_range(xmin, xmax, ymin, ymax)` — Sets the ranges of the internally generated plots. If any of the arguments is the string "auto", the corresponding range limit is determined automatically from the data at the time the plot is drawn. If any of the arguments is the string "extend", the corresponding range limit is only changed if the current data decrease the minimum or increase the maximum. Returns **true**.

The Data-Pipe Facility

`spec`'s `data_pipe()` function allows integration of external code with `spec`. With the data-pipe facility, `spec` sends information to the external program, allows the external program to execute for a time, and then receives information back from the external program. The information can be in the form of a string or a number, and can also include the contents of a `spec` data group or data array. The handshaking and data transfer between `spec` and the data-pipe program is done in an overhead module included in the `spec` distribution that is linked with the external code.

From `spec`, access to the data-pipe facility is through the `data_pipe()` function called from the user level. Usage is as follows.

`data_pipe("?")` — Lists the currently running data-pipe processes with name and process id.

`data_pipe(program, "kill")` — Kills the process associated with *program*.

`data_pipe(program [, args [, grp_out|arr_out [, grp_in|arr_in]]])` – Initiates or resumes synchronous execution of the special process named *program*. If *program* contains a / character, then it contains the complete absolute or relative path name of the program to run. Otherwise the program must be in the `SPECD/data_pipe` directory, where `SPECD` is the built-in `spec` variable containing the path name of `spec`'s auxiliary file directory. You can use the string "." for *program* as an abbreviation for the same program name as used in the last call to `data_pipe()`.

The string value of *args* is made available to the user code in the program as described in the next section.

The optional arguments *grp_out* and *grp_in* are data group numbers. If *grp_out* is present, the contents of that group are sent to the data-pipe program. If *grp_in* is present, it is the number of the data group that will receive values from the data-pipe program. The data-pipe program configures the size of *grp_in* for an implicit call to `data_grp()` within `data_pipe()`. If the *grp_in* argument is absent, `spec` will not receive data-group data from the data-pipe program. If *grp_out* is also absent, group data won't be sent to the data-pipe program. Even without group arguments, the data-pipe program can still return values to `spec` in the form of assigning a number or string return value to `data_pipe()`.

Either or both of the data group arguments can be replaced with the array arguments *arr_out* and *arr_in*. The arrays referred to by these arguments must be the data arrays declared explicitly with the `array` keyword. When sending array data to the data-pipe program, the array data is first converted to double precision floating point format. The received data is always double, but is converted to fit the declared data type of *arr_in*. Only as much data as will fit into the array will be assigned. The number of columns in *arr_in* should match the width of the data sent over by the data-pipe program. If not, the data will still be assigned to the array, but will be misaligned.

Prior to `spec` release 4.03.13, only one `data_pipe()` function could be active at a time.

The user C code can be compiled and linked using the command

```
dpmake program [UOBJ=...] [LIBS=...] [optional_make_args]
```

The command `dpmake` is a short shell script which invokes the `make` utility using the makefile `data_pipe.mak` in the `SPECD/data_pipe` directory. The file `program.c` will be compiled and linked with the data-pipe overhead module, and the result placed in an executable file named *program*. If additional object modules or libraries need to be linked, they can be specified with the `UOBJ=` or

LIBS= parameters. If the tools provided are not sufficient, you can create custom makefiles based on the distributed *data_pipe.mak*.

After linking *program*, move it to the *SPECD/data_pipe* directory for easy access by all users.

The subroutines available from the user C code portion of the data-pipe program are described below.

The skeleton user C-code part of the data-pipe program should contain the following:

```
#include <user_pipe.h>

user_code(argc, argv)
char      **argv;
{
    ...
}
```

The include file *user_pipe.h* contains declarations of the subroutines available in the C code. The file resides in the *SPECD/data_pipe* directory.

The subroutine *user_code()* is called by the overhead part of the data-pipe program each time *data_pipe()* is invoked in *spec*. The parameter *argc* is set to the number of space-delimited words present in the string value of the *args* parameter to *data_pipe()*. The parameter *argv* is an array of character pointers that point to each of the *argc* space-delimited words in the *args* string. Alternatively, the *get_input_string()* function (see below) returns the *args* string in its entirety.

The *user_code()* routine will be called every time the *data_pipe()* function is called from *spec*. The data-pipe program does not exit between calls of *user_code()*, so you should be careful about allocating memory or opening files each time *user_code()* is called without freeing the memory or closing the files each time *user_code()* returns. Alternatively, you can make sure such things are only done the first time *user_code()* is called.

Besides the *argc, argv* technique for accessing the *args* typed in the *data_pipe()* call, the following function is available (as of *spec* release 4.03.13):

`char * get_input_string()` — Returns a pointer to memory holding a copy of the second argument *args* entered with the call to *data_pipe()*.

If *data_pipe()* is sending a data group or array to the user code, the following subroutines provide access to the data parameters and values.

`int get_group_number()` — Returns the group number specified as the *data_pipe()* *grp_out* argument. A -2 is returned if an array was

specified. A -1 is returned if neither data group or array was specified.

`int get_group_npts()` – Returns the number of points in the `data_pipe()` `grp_out` or the number of rows in `arr_out`.

`int get_group_width()` – Returns the number of elements per point in the `data_pipe()` `grp_out` or the number of columns in `arr_out`.

`int get_input_data(double *x, int pts, int wid)` – Transfers data from the `grp_out` or `arr_out` specified in the call to `data_pipe()` to the memory area indicated by the pointer `x`. The pointer `x` is treated as an array dimensioned as `x[pts][wid]`. If the data group/array has more points/rows than `pts` or more elements/columns than `wid`, only as many points/rows or elements/columns as are available in the data group/array are copied. Data from only a single element/column may be retrieved using one or more calls of `get_input_elem()` below. If the data in the data group from `spec` is float rather than double (which depends on `spec`'s installation configuration), float-to-double conversion is done within the call to `get_input_data()`. The return value is the number of points/rows copied.

`int get_input_elem(double *x, int pts, int el)` – Transfers one element of the data from the `grp_out` or `arr_out` specified in the call to `data_pipe()` to the memory area indicated by the pointer `x`. No more than `pts` points are copied from element/column `el` of the the data group/array. If the data in the data group from `spec` is float rather than double (which depends on `spec`'s installation configuration), float-to-double conversion is done within the call to `get_input_data()`. The return value is the number of points/rows copied.

The following subroutines allow you to send group/array data back to `spec` when `data_pipe()` is invoked with a `grp_in` or `arr_in` argument. For a data group, the call to `data_pipe()` will implicitly call `data_grp()` to configure the size of the return group according to the parameters set in the following subroutines. For an array, the array must already be declared and dimensioned.

There are two ways to send group/array data back to `spec`. The subroutine `set_return_data()` allows you to send the entire data group in one call that passes both a pointer to the data and the data group size to the data-pipe program overhead code. Alternatively, you can use the `set_return_group()` subroutine to configure the data group/array size, followed by one or more calls to `set_return_elem()` to set one element/column of the data group/array at a time.

`int get_return_group_number()` – Returns the group number specified as the `data_pipe()` `grp_in` argument. A -2 is returned if an array was specified. A -1 is returned if neither data group or array was specified.

`void set_return_data(double *x, int pts, int wid, int last)` – Configures the return data group and copies the data at the same time. The pointer `x` is considered as an array of dimension `x[pts][wid]` for the purpose of transferring data to the data group. The argument `last` sets the index number of the last element added to the group, which is used by the various data analysis and plotting functions available in `spec`.

`void set_return_group(int pts, int wid)` – Configures the size of the return data group without copying data. This subroutine must be called once before calling `set_return_elem()` below.

`void set_return_elem(double *x, int pts, int el, int last)` – Copies one element to the return data group, which must have been previously configured by a call of `set_return_group()`, above. If the parameters `pts` or `el` exceed the values configured, or if the return group hasn't been configured, the subroutine returns `-1`. Otherwise zero is returned.

You can set the value that the `data_pipe()` function returns in `spec` from the user C code in your *data-pipe* process. You can have `data_pipe()` return a number or a string or, if necessary, reset to command level. If no explicit return value is assigned in the user C code, `data_pipe()` returns zero.

`int set_return_string(char *s)` – Sets the return value of `data_pipe()` to the string `s`. This subroutine returns `-1` if memory could not be obtained for the string `s`, otherwise it returns zero.

`void set_return_value(double v)` – Sets the return value of `data_pipe()` to the value `v`.

`void do_error_return()` – Calling this subroutine from the user C code causes control to pass back to `spec` without returning data group or array values, if they have been set. The return value of `data_pipe()` will be the value set by `set_return_value()` above, if such a value has been set, otherwise the return value of `data_pipe()` will be `-1`. This subroutine does not return.

`void do_abort_return()` – Calling this subroutine from the user C code causes control to pass back to `spec` without returning data group or array values, if they has been set. In `spec`, there is no return from `data_pipe()`, rather `spec` resets to command level. This subroutine does not return.

`void do_quit_return()` – Calling this subroutine from the user C code causes control to pass back to `spec` normally as if `user_code()` returned normally, but the *data-pipe* program will then exit. This subroutine does not return.

String and Number Functions

Math Functions

`exp(x)`, `exp10(x)` – Returns e^x and 10^x respectively.

`log(x)`, `log10(x)` – Returns the natural logarithm and the base 10 logarithm of x respectively.

`sqrt(x)` – Returns the square root of x .

`pow(x, y)` – Returns x^y .

`fabs(x)` – Returns the absolute value of x .

`int(x)` – Returns the integer part of x . The integer part is formed by truncation towards zero.

`rand()` – Returns a random integer between 0 and 32767.

`rand(x)` – If x is positive, returns a random integer between 0 and x , inclusive. If x is negative, returns a random integer between $-x$ and x , inclusive. Values of x greater than 32767 or less than -16383 are set to those limits. If x is zero, zero is returned. The C-library `rand()` function is used to obtain the values. The seed is set to the time of day on the first call. The randomness (or lack thereof) of the numbers obtained is due to the C library implementation.

`srand(seed)` – Sets the seed value for the random number generator used by the `rand()` function to the integer value `seed`. This function allows the same sequence of random numbers to be generated reproducibly by resetting the seed to the same value.

`sin(x)`, `cos(x)`, `tan(x)` – Returns the sine, cosine and tangent, respectively, of the argument x , which must be in radians.

`asin(x)`, `acos(x)`, `atan(x)` – Returns the arc sine, arc cosine and arc tangent, respectively, of the argument x . The return value is in radians. `asin()` and `atan()` return values in the range $-\pi/2$ to $\pi/2$, while `acos()` returns values from 0 to π .

`atan2(y, x)` – Returns the arc tangent of y/x using the signs of the arguments to determine the quadrant of the return value. The return value is in the range $-\pi$ to π . Having both y and x zero is an error.

String Functions

`index(s1, s2)` – Returns an integer indicating the position of the first occurrence of string `s2` in string `s1` or zero if `s1` does not contain `s2`.

`split(string, array)` – Splits the string *string* at space characters and assigns the resulting substrings to successive elements of the associative array *array*, starting with element 0. The space characters are eliminated. The function returns the number of elements assigned.

`split(string, array, delimiter)` – As above, but splits the string into elements that are delimited by the string *delimiter*. The delimiting characters are eliminated.

`substr(string, m)` – Returns the portion of string *string* that begins at position *m*.

`substr(string, m, n)` – As above, but the returned string is no longer than *n*.

`length(string)` – Returns the length of the string *string*.

`sprintf(format [, a, b, ...])` – Returns a string containing the formatted print. See *sprintf()* in a C-language reference manual.

`sscanf(string, format, a [, b, ...])` – Scans the literal string or string variable *string* for data, where *format* contains a format specification in the same style as the C language *scanf()* function. Each subsequent argument is a variable name or array element that will be assigned the values scanned for. The function returns the number of items found in the string.

Conversion Functions

`asc(s)` – Returns the ASCII value of the first character of the string value of the argument.

`bcd(x)` – Returns a 24-bit integer that is the binary-coded decimal representation of the nonnegative integer *x*.

`dcb(x)` – Returns the nonnegative integer corresponding to the 24-bit binary-coded decimal representation *x*.

`deg(x)` – Returns the argument converted from radians to degrees.

`rad(x)` – Returns the argument converted from degrees to radians.

STANDARD MACRO GUIDE

Introduction

The standard macros included with the `spec` distribution are an integral part of the `spec` package. For some sites, the standard macros are sufficient for performing experiments. For others, the standard macros provide a starting point for custom development. This reference presents some of the standard macros, grouped by functionality. For many macros, the complete definition is printed. At the end of this reference, the construction of the scan macros is discussed in depth.

The following files, found in the `macros` subdirectory of the distribution, contain the definitions for all the macros in the standard library. If it is not obvious in which file a particular macro is stored, you can always change to the `macros` directory and type `grep macro_name *`, where `grep` is the standard UNIX file searching utility.

File	Contents
<i>count.mac</i>	Counting macros (<code>ct</code> , <code>uct</code> , <code>count</code> , <code>show_cnts</code> , ...).
<i>cplot.mac</i>	The <code>cplot_plot</code> macro.
<i>energy.mac</i>	For an energy-selecting monochromator (<code>Escan</code> , <code>set_E</code> , ...).
<i>file.mac</i>	The <code>newfile</code> macro.
<i>getscan.mac</i>	The <code>getscan</code> macro.
<i>hkl.mac</i>	General reciprocal space macros (<code>br</code> , <code>mk</code> , <code>ca</code> , <code>wh</code> , ...).
<i>motor.mac</i>	Motor moving and status (<code>mv</code> , <code>umv</code> , <code>wa</code> , <code>set</code> , <code>set_lm</code> , ...).
<i>plot.mac</i>	Data plotting (<code>plot</code> , <code>plot_res</code> , <code>rplot</code> , <code>splot</code> , <code>ansiplot</code> , ...).
<i>powder.mac</i>	Powder-mode macros (<code>setpowder</code> , <code>_pmove</code> and <code>_pcount</code>).
<i>region.mac</i>	Macros to define a series of scans(<code>setreg</code> and <code>doreg</code>).
<i>scans.mac</i>	Basic scan macros (<code>ascan</code> , <code>a2scan</code> , <code>hklscan</code> , ...).
<i>scans1.mac</i>	Scan helper macros (<code>_head</code> , <code>_loop</code> , <code>setscans</code> , ...).
<i>start.mac</i>	The <code>startup</code> macro.
<i>temper.mac</i>	Temperature control (<code>settemp</code> , <code>measuretemp</code> , <code>te</code> , ...).
<i>util.mac</i>	Misc. utility macros (<code>do</code> , <code>qdo</code> , <code>savmac</code> , <code>comment</code> , <code>u</code> , <code>help</code> , ...).
<i>fivec.src</i>	Five-circle geometry macros.
<i>fourc.src</i>	Four-circle geometry macros.
<i>sixc.src</i>	Six-circle geometry macros.
<i>spec.src</i>	Version for no diffractometer.
<i>surf.src</i>	Special liquid surface diffractometer macros.
<i>twoc.src</i>	Two-circle geometry macros.
<i>zaxis.src</i>	Z-axis geometry macros.

When installed, the `.mac` files above are combined into one file and placed (assuming the default auxiliary file directory) in `/usr/local/lib/spec.d/standard.mac`. A file formed from the first four letters of the geometry configuration contains the geometry

macros from one of the `.src` files above. For example, `/usr/lib/spec.d/four.mac` is created for the four-circle geometry and contains the macros from `fourc.src`.

The macro definitions are the least stable part of the `spec` package. The macros are easy to change — no recompilation of C code is necessary — and the intent of the designers of the `spec` package was to put its flexibility in the macros. Thus, the definitions presented below may differ with the macro definitions in your current version of `spec`.

You may find the existing macros do not accomplish what you want. A simple procedure for modifying a standard macro is to use the macro `savmac` to copy the definition of the existing macro into a file. You then edit the macro definition in that file and read it back in using the `do` macro. You can gather your customized macros into a file named `spec.mac` in your current directory or into the file `/usr/lib/spec.d/site.mac`. Both of these files are read every time you start `spec`, whether or not you are starting fresh.

If you have made generally useful modifications to the standard macros, or if you have developed your own macros, please send copies to Certified Scientific Software. We include many user-contributed macros in each new release of the `spec` package.

Some Tips

The syntax rules for defining macros are given in the *Reference Manual* on page 90. The suggestions that follow offer some additional guidance for writing macros that will fit in well with the standard library.

When a macro requires arguments, it is a good idea to check that the right number of arguments have been given, and if not, print a usage message and exit to command level. The symbol `$#` will be set to the number of arguments when the macro is run. For example,

```
def ascan '
    if ($# != 5) {
        print "Usage:  ascan  motor start finish  intervals time"
        exit
    }
    ...
,
```

If an argument is supposed to be a motor number or mnemonic, use the `_check0` macro before operating on the motor. The `_check0` macro exits to command level if

the argument is not a valid mnemonic. For instance, to check the first argument of a macro, use

```
_check0 "$1"
```

A mistyped mnemonic might otherwise become a variable with an arbitrary value (zero for a new variable) resulting in an operation on the wrong motor (usually motor zero).

It is good practice to refer to arguments just once when writing macros to avoid side effects that occur, for example, if the macro is invoked as `my_mac i++`. Here the variable `i` would be incremented each time `$1` is used in the macro. In the scan macros, the arguments are assigned to global variables just after the usage check:

```
def ascan '  
    ...  
    { _m1 = $1; _s1 = $2; _f1 = $3; _n1 = int($4); _ctime = $5 }  
    ...  
,
```

When a macro changes a parameter or mode that affects later data, it is a good idea to note that change in the data file and on the printer. Macros such as `comment`, `qcomment` and `gpset` are available for that purpose.

If possible, declare local variables `local` to avoid conflicts with other variables, especially when macros are nested or parsed together.

Watch out for name conflicts when naming new macros and using variables. You can prevent most conflicts by using the `local` keyword to explicitly declare local names within a statement block. Names declared that way can be used as symbols within the statement block even if they are already in use as macros. Otherwise, if you construct commands using a variable name that is really a macro name, when that intended variable is encountered, it will be replaced by the macro, making a mess of things.

Note that several one-letter names such as `d`, `h`, `p` and `l` are already in use as macro names. Don't use these names as variables, unless they are declared `local` inside a statement block. Typing `lsdef ?` will list all one letter macro names. Typing `lsdef _?` will list all two letter macro names that begin with an underscore.

Command files that define macros often assign default values to related global variables. You should always check if these global variables have already had a value assigned before assigning default values. If the user had assigned a new value to a variable, you do not want that assignment undone if the macro file is reread. The built-in `whatis()` function can be used to see if a variable has been assigned a value (see page 75 for an explanation of the `whatis()` return values),

```

if ((whatis("DATAFILE")>>16)&0x0800) {
    print "Warning: No open data file. Using \"/dev/null\".\n"
    open(DATAFILE = "/dev/null")
}

```

When writing macros that move motors, be careful with the `move_all` command. When moving motors, always do a `waitmove` and `getangles` first. Then assign new values to `A[]`, and finally call `move_all` (or `move_em`).

When obtaining input from the user, the functions `getval()` and `yesno()` are useful. For example,

```

_update = yesno("Show updated moving and counting", _update)
g_mode = getval("Geometry mode", g_mode)

```

results in the following screen output:

```

Show updated moving and counting (NO)?
Geometry mode (3)?

```

You can also use the `input()` built-in function to obtain user input. Remember, though, that `input()` returns a string. If the string contains a valid number, the automatic string-to-number conversion will take place, if context requires it. However, no expression simplification is done on the string, so a response of `2+2` will not have a number value of 4 when returned by `input()`.

When using `on()` and `off()` to control output, do the operations on `"tty"` last. Since `"tty"` is always turned back on if everything else is turned off, the commands

```

off("tty");on(PRINTER);print "hello world";on("tty");off(PRINTER)

```

will not have the desired effect. The first `off()` turns off everything, so `"tty"` is automatically turned back on, and the message goes to both `PRINTER` and `"tty"`.

Use existing UNIX utilities if they can be of help. For example, if you manipulate UNIX file names in your macros you can use the return value of the `test` utility to check for existence of a file. For example, the function `unix("test -r $1")` will return zero if the file specified by the argument exists and is readable.

Utility Macros

UNIX Commands

These simple macros are for commonly used UNIX commands.

```
def cd      'chdir("$*")'           # Change working directory
def pwd     'print CWD'           # Print working directory
def u       'unix("$*")'          # Execute arbitrary shell commands
def ls      'unix("ls $*")'       # List files
def l       'unix("ls -l $*")'    # Long listing of files
def cat     'unix("cat $*")'      # Show file contents
def less    'unix("less $*")'     # Peruse files with handy utility
def mail    'unix(sprintf("%s $*", MAIL))' # Send mail
def ed      'unix("ed $*")'       # Invoke an editor
def ned     'unix("ned $*")'      # Invoke another editor
def vi      'unix("vi $*")'      # Invoke another editor
```

The `u` macro, without arguments, spawns an interactive subshell, using your *SHELL* environment variable.

Note how the above macros supply parentheses and quotation marks around the arguments, as required by the parser's grammar rules.

Basic Aliases

The main purpose of these macros is to provide a shorthand way of typing some useful commands.

```
def d       'print date()'        # Print current time and date
def p       'print'               # Shorthand for print
def h       'help'                # Shorthand for help, below
def hi      'history'             # Shorthand for history
def beep    'printf("\a")'        # Sound the bell
def cl      'tty_cntl("cl")'      # Clear the screen
def com     'comment "$*"'        # Shorthand for comment, below
def ond     'if (DATAFILE)        # Send output to data file
{ on(DATAFILE) }'
def offd    'if (DATAFILE)        # Stop sending
{ off(DATAFILE) }'
def onp     'if (PRINTER)         # Send output to the printer
{ on(PRINTER) }'
def offp    'if (PRINTER != "")   # Stop sending
{ off(PRINTER) }'
def ont     'on("tty")'          # Send output to the terminal
def offt    'off("tty")'         # Stop sending
def fon     'if ($# == 1) on("$1")
else { print "Usage: fon filename"; on(); }'
def foff    'if ($# == 1) off("$1")'
```

```

        else { print "Usage: foff filename";on(); }'
def waitall  '{ user_waitall; wait(0) }'    # Wait for async activity
def waitmove '{ user_waitmove; wait(1) }'   # Wait for moving to end
def waitcount '{ user_waitcount; wait(2) }' # Wait for counting to end
def chk_move  '(wait(0x21) || USER_CHK_MOVE)'
def chk_count '(wait(0x22) || USER_CHK_COUNT)'
def chk_acq   '(wait(0x24) || USER_CHK_ACQ)'
def w         '{ waitall; beep }'          # Wait, and be audible when done

```

Basic Utility Macros

These straightforward macros combine a number of built-in functions and commands to provide a higher level of functionality with minimal input. First, here is their usage:

```

help [topic]           # Display help files
config                 # Edit hardware configuration
onsim                  # Turn on simulate mode
offsim                 # Turn off simulate mode
debug [value]         # Select debugging categories
bug                    # Mail a bug report
whats object          # Identify the object
gpset variable value  # Comment if a variable has changed

```

Here are the definitions for some:

```

# Examine help file, use default if no argument.
def help '
    if ($#)
        gethelp("$1");
    else {
        local t
        for (t="help";; )
            if (gethelp(t) || (t = input("\nSubject? ")) == "")
                break
    }
'

```

```

# View (and modify), then reread configuration file.
# Use -s flag if in simulate mode. Re-order motor numbers
# with _assign. Check for monochromator mnemonics.
def config '
    wait(-1)
    user_waitall
    sync
    unix(sprintf("%s/edconf %s %s/%s",\
                SPECD, set_sim(-1)? "-s:""", SPECD, SPEC))
    reconfig
    user_config
    _assign
    _assign_mono
,
def user_config ''

# Turn simulate mode on. Comment on printer and file if changed.
def onsim '{
    local t

    if (!(t = set_sim(1))) { qcomment "Simulate mode ON" }
    printf("Simulate was %s, is now %s.\n", t? "on":"off",\
          set_sim(-1)? "ON":"OFF")
}'

# Turn simulate mode off.
def offsim '{
    local t

    if (t = set_sim(0)) { qcomment "Simulate mode OFF" }
    printf("Simulate was %s, is now %s.\n", t? "on":"off",\
          set_sim(-1)? "ON":"OFF")
}'

# Easy way to set the debug level.
# +arg adds bits to DEBUG. -arg removes them.
def debug '{
    local t

    if ($# == 0) {
        gethelp("debug")
        t = input(sprintf("\nDebug value (%d)? ", DEBUG))
    } else
        t = "$*"
    if (index(t, "+"))      DEBUG |= 0+t
    else if (index(t, "-")) DEBUG &= ~(0-t)
    else                    DEBUG = 0+t
}'

```

```

# Send a bug report to the administrator.
def bug '
    print "The mail utility will be run for you.  Describe your"
    print "problem to the administrator.  When you are done, type ^D."
    {
        local s
        s = unix(sprintf("%s -s \"Bug from %s\" %s", MAIL, USER, ADMIN))
        printf("Bug report %ssent to %s.", s? "not ":"", ADMIN)
    }
,

# Set something and comment if it has changed.
def gpset '
    if ($1 != $2) {
        comment "$2 reset from %g to %g" "$2,$1"
        $2 = $1
    }
,

```

Reading From Command Files

```

do command_file           # Run a command file
qdo command_file        # Run a command file without echo
newmac                     # Reread standard command files

```

Since the `do` and `qdo` macros have nearly identical functionality, the commands for both are placed in a single macro named `_do`. This macro implements special functions, such as letting a dot stand for the previous command file or searching for a command file first in the current directory and then in a special command file directory.

```

# "do" a command file.
def do '_do $* do'

# Quietly "do" a command file.
def qdo '_do $* qdo'

# Run a command file.  If not in current directory, look for
# file in DO_DIR.  Save file name in DOFILE so "do ." works.
def _do '
    if ($# != 2 || ("${2}" != "do" && "${2}" != "qdo")) {
        print "Usage:  do file"
        print "          qdo file"
        exit
    }
    if ("${1}" == "." && DOFILE == "") {
        print "No previous do file."
        exit
    }

```

```

if ("$1" != ".") {
    DOFILE = "$1"
    if (DO_DIR != "." && unix(sprintf("test -r %s", DOFILE))) {
        local t

        t=sprintf("%s/%s", DO_DIR, DOFILE)
        if (!unix(sprintf("test -r %s", t)))
            DOFILE = t
    }
}
if (!unix(sprintf("test -r %s", DOFILE))) {
    qcomment "do %s" DOFILE
    $2file(DOFILE)
} else {
    printf("Can\'t read command file \"%s\".\n", DOFILE)
    exit
}
,

```

The `newmac` macro rereads the standard macro files that reside in `SPECD` the auxiliary file directory (usually `/usr/lib/spec.d`). Invoking `newmac` is useful if a new version of the standard macros has been installed but you do not want to start fresh, or if you have somehow corrupted the definition of a standard macro and want to get back the original definition.

Saving To Output Devices

```
comment format args                # Send a comment to output devices
qcomment format args              # Send comment to file and printer
prcmd command                     # Print the output of a command
savcmd command file               # Save a command to a file
savmac macro_name file           # Save a macro to a file

# Print a comment on the screen, printer and data file.
def comment '
    printf("\n%s. $1.\n", date(), $2)
    qcomment "$1" "$2"
,

# Print a comment on the printer and data file.
def qcomment '
    if (PRINTER != "")
        fprintf(PRINTER, "\n%s. $1.\n", date(), $2)
    if (DATAFILE != "")
        fprintf(DATAFILE, "#C %s. $1.\n", date(), $2)
,

# Have output of any command sent to the printer.
# Commands are all on one line to avoid outputting prompts.
def prcmd 'onp;offt;printf("\n%s\n", "$*");ont;$*
offp'

# Have output of any command sent to a file.
def savcmd '
    if ($# != 2) {
        print "Usage: savcmd command filename"
        exit
    }
    on("$2");offt;printf("\n%s\n", "$1");ont;$1
    close("$2")
,

# Save a macro definition to a file.
def savmac '
    if ($# != 2) {
        print "Usage: savmac macro_name filename"
        exit
    }
    on("$2"); offt
    prdef $1
    ont; close("$2")
,
```

Start-up Macros

These macros ask for all the initialization information that may be needed by the other macros.

```
startup           # Sets things up to start experiment
newsample         # Gets title and data file for experiment
newfile [name [num]] # Sets up a data file
starttemp         # Asks for temperature control parameters
setscans          # Asks for scan options
setplot           # Asks for plotting options
startgeo          # Queries user for all geometry parameters
save [file]       # Saves important variables to a file
savegeo           # Used by "save", saves geometry parameters
saveusr           # Used by "save", user defined
```

The `startup` macro prompts for values for a number of global variables and also invokes all the other initialization macros, leading to more than a screenful of questions. You can, however, always jump back to command level by typing the interrupt character (`^c`) if you do not need to change items at the bottom of the list. Be careful, though, since some of the initialization macros, (`setplot`, for example) don't save the entered information until all their questions are answered.

```
def startup '
  printf("\n(newsample)")
  newsample
  {
    local t
    t = PRINTER != "" && PRINTER != "/dev/null"
    if (yesno("\nUse a printer for scan output", t)) {
      PRINTER = getval("Printer device",PRINTER)
      if (index(PRINTER,"/") == 0)
        PRINTER = sprintf("/dev/%s",PRINTER)
      if (open(PRINTER))
        PRINTER = "/dev/null"
    } else
      PRINTER = "/dev/null"
  }
  if (substr(PRINTER,1,5) != "/dev/")
    PRINTER = sprintf("/dev/%s",PRINTER)
  if (open(PRINTER))
    PRINTER = "/dev/null"
  newfile
  DO_DIR = getval("\nCommand file directory",DO_DIR)
  COUNT = getval("Default count time for ct and uct",COUNT)
  UPDATE = getval("Update interval for umv, uct, etc. in seconds",UPDATE)
  if (whatis("starttemp") >> 16)
    printf("0arttemp")
  starttemp
```

```

setscans
setplot
startgeo

```

In the standard distribution, `starttemp` has a null definition.

The `newfile` macro creates, opens or reopens standard `spec` data files. The filename and scan number may be given as arguments. Otherwise, you are asked for the information. If you have a directory named *data* in your current directory, and there are no / characters in the file name you give, the data file will be placed in the *data* directory. If the file already exists, new scans will be appended to the file. The existing file is not removed.

The `startgeo` macro is defined differently for various geometries, but should query the user for values for all the relevant geometry parameters.

The `save` macro is not really an initialization macro, but it creates a file that can be used for initialization. The purpose of the macro is to save all the important global variables in a file that can be run as a command file at a later time to restore the values of those variables. For example, if the user anticipates starting fresh with a new version of the software, having a file created by the `save` macro will simplify creating a new program state.

```

# Save current globals to a save file
def save '{
    local f

    if ($# == 0)
        f = getval("File for saving globals", "saved")
    else if ($# == 1)
        f = "$1"
    else {
        print "Usage:  save [filename]"
        exit
    }
    unix(sprintf("if test -s %s ; then mv %s %s.bak ; fi", f, f, f))
    on(f); offt
    printf("PRINTER=\ "%s"\ "\n", PRINTER)
    savegeo
    saveusr
    ont; close(f)
    qcomment "Globals saved in \ "%s\ "" "f"
    printf("Type \ "do %s\ " to recover.\n", f)
}'

```

The macro `savegeo` saves all the geometry parameters for the particular configuration. You can define the macro `saveusr` to save whatever else is desired.

Motor Macros

```
mv motor pos           # Move a motor
mvr motor pos          # Move a motor, relatively
mvd motor dial_pos    # Move a motor to a dial position
tw motor inc           # Tweak a motor, interactively
umv motor pos          # Move while updating screen
umvr motor pos         # Move while updating screen
wa                     # Show positions of all motors
lm                     # Show limits of all motors
wm m1 m2 ...           # Show positions and limits of motors
uwm m1 m2 ...          # Show positions while motors are moving
set motor pos          # Set user angle for a motor
set_dial motor pos     # Set dial angle for a motor
set_lm motor low high # Set user limits for a motor
an tth_pos th_pos      # Move two theta and theta
pl chi_pos phi_pos     # Move chi and phi (four-circle)
uan tth_pos th_pos     # Move while updating screen
upl chi_pos phi_pos    # Move while updating screen
```

The following macro moves a single motor, adding a comment to the printer that the motor was moved:

```
# Move a single motor
def mv '_mv $*; move_poll'
def umv '_mv $*; _update1 $1' # "update" version of mv
def _mv '
    if ($# != 2) {
        print "Usage: mv motor position"
        exit
    }
    _check0 "$1"
    waitmove; getangles; A[$1]=$2
    if (PRINTER != "")
        fprintf(PRINTER, "\nmv $1 %g\n", A[$1])
    move_em
,
```

In `mv`, as in all the macros that move motors, the `move_em` macro is invoked, rather than the `move_all` command. Normally, `move_em` is defined as

```
def move_em '
    user_premove
    move_all
    user_postmove
```

One can define the `user_premove` and/or `user_postmove` macros to take into account special conditions. For example, to check for limits that depend on the relative position of motors, one could define `user_premove` as

```

def user_remove '
    if (fabs(A[tth] - A[th]) > 10) {
        print "Move exceeds Theta - Two Theta relative limit."
        exit
    }
    move_all
,

```

The `set` macro changes the offset between user and dial units.

```

# Define a new motor position
def set '
    if ($# != 2) {
        print "Usage:  set motor new_user_value"
        exit
    }
    {
        local old
        _check0 "$1"
        waitmove; getangles
        old = A[$1]
        if (chg_offset($1, $2))
            exit
        getangles
        if (old != A[$1]) {
            comment "%s reset from %g to %g" "motor_name($1), old, A[$1]"
        } else
            print "No change."
    }
,

```

The `set_dial` macro changes the dial position of the motor, which means a change to the contents of the motor controller register. `set_dial` refuses to set the dial beyond the current software limits for the motor. `set_dial` also changes the offset to maintain the prior value of the user angle. These two macros document the change in the data file and on the printer.

The `set_lm` macro converts the user-unit arguments to dial units for the call to `set_lim()`.

```

Change a motor limit
def set_lm '
    if ($# != 3) {
        print "Usage:  set_lm motor low high"
        exit
    }
    {
        _check0 "$1"
        if (!set_lim($1, dial($1, $2), dial($1, $3))) {
            onp
            printf("\n%s limits set to %g %g (dial units).\n",\
                motor_name($1), get_lim($1, -1), get_lim($1, +1))
            offp
        }
    }
,

```

The macros in the above list that begin with a `u` continuously read motor positions from the controller and show the positions on the screen. The frequency of screen updates is set by the global variable `UPDATE`, which is used as an argument to the `sleep()` function. Setting `UPDATE=.25` places a 1/4 second pause between updates. The `umv` macro first calls `_mv` and then calls the internal `_update1` macro. The other updated-move macros are defined similarly.

```

def umv _'mv $*; _update1 $1 ' # "update" version of mv

# Displays updated position of 1 motor while it is moving
def _update1 '
    if (chk_move) {
        printf("\n%10.9s\n", motor_name($1))
        while (wait(0x22)) {
            getangles
            printf("%10.4f\r", A[$1])
            sleep(UPDATE)
        }
        getangles
        printf("%10.4f\n", A[$1])
    }
,

```

The technique for displaying status information about all the motors is a little complicated. `spec` places no restriction on what order the motors are assigned to the controller, but does recognize that there is a preferred order for displaying motor information. To this end, the macros use an array `mA[]` which contains reordered motor numbers. The four-circle macro source file contains the following code, which is executed when the command file is read and when the `config` macro is run.

```

# Conventionally, the first four motors are tth, th, chi, phi.
# The following code guarantees this.
def _assign '{
    local i j
    mA[0]=tth
    mA[1]=th
    mA[2]=chi
    mA[3]=phi
    for (i = 4, j = 0; i < MOTORS; j++) {
        if (j == tth || j == th || j == chi || j == phi)
            continue
        mA[i++] = j
    }
}'

```

Similar code is contained in the macro source files for the other geometries.

An internal macro named `_mo_loop` exists to loop through all the motors printing selected fields. Its use is best illustrated by example. First here is its definition:

```

# Looping routine used in many macros.
# Normally k is set to MOTORS, but can be set to something else, e.g., 4
# (Kludge with printf(" ") avoids auto linefeed on 80th column.)
def _mo_loop '{
    local s
    for (j = i; j < i + 8 && j < k; j++)
        if (motor_name(mA[j]) != "unused") {
            s = s sprintf("%$1", $2)
            if (j < i + 7)
                s = s " "
        }
    print s
}'

```

It is within this macro that motors named `unused` are not used in printing motor information.

The `wa` macro that displays information for all motors is typical of a macro that calls the `_mo_loop` macro.

```

# Where - all motors
def wa '
    waitmove; get_angles
    onp
    printf("\nCurrent Positions (user, dial)\n")
    {
        local i j k
        for (i = 0, k = MOTORS; i < k; i += 8) {
            _mo_loop 9.9s "motor_name(mA[j])"
            _mo_loop 9.9s "motor_mne(mA[j])"
            _mo_loop 9.4f "A[mA[j]]"
            _mo_loop 9.4f "dial(mA[j], A[mA[j]])"
        }
    }
    offp
,

```

The first argument for `_mo_loop` is a `printf()` field specification, the second argument is the field value. The field values use the `mA[]` array to reorder the motor numbers.

Counting Macros

```
ct [time]           # Count, then show_cnts
count [time]       # Count for time
show_cnts          # Count, then display results
uct [time]         # Updated counting
```

When *time* is positive, counting is to seconds. When *time* is negative, counting is to monitor counts. If the counting macros are invoked without an argument, the count time used is that contained in the global variable `COUNT`.

Counting in `spec` combines timing generators and scalers. Three scaler channels are normally used in the standard macros. The first channel takes an accurate clock input (normally at 1 KHz), the second takes the monitor input, and the third takes the detector.

The scaler channels identifying the various inputs are determined by the values of the global constants, `sec`, `mon` and `det`. Their default values are 0, 1 and 2, respectively, matching the recommended hardware cabling. If you do not connect your counting sources to the default scaler channels, you must explicitly assign new values to `sec`, `mon` and `det` for the standard counting macros to work properly.

The additional global variables `MON` and `DET` are then set to the channels that are to be treated as the monitor and detector for particular scans, normally `mon` and `det`. If, while doing a line-up scan of a motor on which the monitor is mounted, you would want the monitor counts to be plotted as a function of the motor position, enter `DET=mon` before doing the scan. Do not forget to reassign `DET=det` at the end of the scan.

The global variable `COUNT` is set to a default count time (normally 1 second). If the macro `ct` is invoked without arguments, counting will last for the duration given by `COUNT`.

```
# A user calls "ct" to count for some interval and display results
def ct '{
    rdef cleanup \'
        undef cleanup
        onp; show_cnts; offp
    \'
    waitmove
    count_em $*
    waitcount
    undef cleanup
    onp; show_cnts; offp
}'
```

```

# "count" is the basic macro to count to monitor or to time.
# It runs the clock and reads the scalers.
def count '{
    waitmove
    if ($1) for (;;) {
        count_em $1
        waitcount
        get_counts
        chk_beam
    }
    if (S[sec] && MON >= 0)
        MON_RATE = S[MON]/S[sec]
}'
# The macro "show_cnts" reads the scalers and displays the results.
def show_cnts '{
    local i

    get_counts
    printf("\n%s\n\n", date())
    for (i=0;i<COUNTERS;i++)
        if (cnt_name(i) != "unused")
            printf("%12s = %g%s\n", cnt_name(i), S[i], \
                i != sec && S[sec]? sprintf(" (%g/s)", S[i] / S[sec]):"")
}'

```

Updated counting is done with the uct macro,

```

def uct '{
    waitmove
    count_em $*
    if (chk_count) {
        local i
        printf("\n")
        for (i=0;i<COUNTERS;i++)
            if (cnt_name(i) != "unused")
                printf("%12.12s ", cnt_name(i))
        printf("\n")
        while (chk_count) {
            get_counts
            for (i=0;i<COUNTERS;i++)
                if (cnt_name(i) != "unused")
                    printf("%12g ", S[i])
            printf("\r")
            sleep(UPDATE)
        }
        get_counts
        for (i=0;i<COUNTERS;i++)
            if (cnt_name(i) != "unused")
                printf("%12g ", S[i])
        printf("\n")
    }
}'

```

Plotting Macros

At present, `spec`'s plotting is done entirely at the macro level and does only character plots on the screen and printer. Several screen plotting macros are defined for various types of terminals.

```
setplot [mode]      # Select plotting options
plot                # Plot data on printer
rplot              # Plot updated data at each point of scan
splot              # Plot data on screen
pts                # List current data on the screen
lp_plot            # Primitive 132-column wide plot for printers
plot_res           # Show results after scans
splot_res          # Show results on screen plot
rplot_res          # Show results on updated plot during scans
```

The `setplot` macro defines the `plot` macro, depending on your choices of plot modes.

The `scan_plot` macro is invoked within the looping portion of all the scans.

The `setplot` macro assigns values to the global variable, `PLOT_MODE`, according to the values defined in the following table:

Bit Value	Description
1	Do updated plotting during scans.
2	Do screen plot after scan.
4	Do printer plot after scan.
8	Scale x -axis of screen plots to fit width of scan.
16	Force y -axis minimum to be zero.
32	Use logarithmic y -axis.
64	Do simple background subtraction
128	Use high-resolution plotting device
256	With high-res, don't use large dots
512	With high-res, don't connect points with lines
1024	With high-res, don't draw error bars

The `scan_plot` macro is called for each point of a scan, while `plot` is called at the end of each scan.

The `splot` macro draws a screen plot. The `rplot` macro is called to redraw the plot with minimal updating during data accumulation.

Reciprocal Space Macros

The following macros are general and applicable to most geometry configurations. Macros special to the four-circle geometry configuration are described in the *Four-Circle Reference* that follows this guide.

```
ca H K L           # Display calculated positions for H K L
cal H K L          # As above, but don't reset positions
ci tth th chi phi  # Display calculated H K L for angles
br H K L           # Move to H K L
mk H K L           # Move to H K L
ubr H K L          # Move to H K L while updating screen
umk H K L          # Move to H K L while updating screen
mi ALPHA BETA      # Move to ALPHA BETA
wh                 # Display H, K, L, tth, th, chi, phi, etc.
pa                 # Display geometry parameters
```

The difference between `ca` and `cal` is that the first macro restores the A[] angles and H, K and L to the current diffractometer position, while the second macro leaves them at the calculated values.

There is no difference between the `br` and `mk` macros, except their names.

The `ubr` and `umk` macros continuously read the motor positions from the controller and show the positions on the screen. The frequency of updates is set by the global variable `UPDATE`.

```
# Go to a Bragg position
def br '_br $*; move_poll'
def _br '
    if ($# != 3) {
        print "Usage:  br H K L"
        exit
    }
    waitmove; { H=$1; K=$2; L=$3 } getangles; calcA
    onp; offt; printf("\nbr %g %g %g\n", H, K, L); offp; ont
    move_em
,
```

```

# Calculate motor positions for a given H, K, and L
def cal '
    if ($# != 3) {
        print "Usage:  cal H K L"
        exit
    } ;
    {H = $1; K = $2; L = $3 } calcA; calcHKL
    onp
    printf("\nCalculated Positions:\n")
    _var
    offp
,

# As above but reset positions to diffractometer positions
def ca '
    if ($# != 3) {
        print "Usage:  ca H K L"
        exit
    } ;
    {H = $1; K = $2; L = $3 } calcA; calcHKL
    onp
    printf("\nCalculated Positions:\n")
    _var
    offp
    waitmove; getangles; calcHKL
,

# Where - reciprocal and real space
def wh '
    waitmove; getangles; calcHKL
    onp
    _var
    offp
,

# A macro called by "wh", "ca" and "ci" to display important
# geometry quantities.  (Four-circle version.)
def _var '
    printf("\nH K L = %.5g %.5g %.5g\n", H, K, L)
    printf("ALPHA = %.5g BETA = %.5g", ALPHA, BETA)
    printf("  AZIMUTH = %.5g LAMBDA = %g\n\n", AZIMUTH, LAMBDA)
    _mot 4
,

```

Scan Macros

The following sections summarize the usage of the standard scans in `spec`. This discussion is followed by a description of macros to customize the scan output sent to the printer and data file. Refer to page 154 for a detailed discussion of a single-motor scan.

All the scans use the same basic invocation syntax. For example, the single motor, absolute-position scan is invoked as

```
ascan motor start finish intervals time
```

The range of a scan is specified by the starting position `start`, the final position `finish` and the number of intervals `intervals` of the scanned parameters. Thus

```
ascan tth 26 28 20 60
```

would start with the `tth` motor at 26° , and move the motor to 28° , using 20 intervals of 0.1° . The total number of points scanned is one more than the number of intervals, in this case, 21 points. The time per point `time`, if a positive number, indicates counting to seconds. In the above example, each point takes 60 seconds. Using a negative `time` indicates counting to monitor counts.

Scan Miscellany

```
resume                # Restart an aborted scan
setscans              # Set scan-mode options
```

If a scan is halted by typing the interrupt character (`^C`), or because of some other error such as hitting a motor limit, you can normally restart the scan by typing `resume`. You must not have changed the value of any of the internal scan variables in the meantime. If you use `resume` to continue a relative position scan, such as `lup` or `dscan`, the scanned motors will not be returned to the center point when the scan ends, as they otherwise would be.

You also have the option to control how the scan data is displayed on the screen as it is taken. An option to the `setscans` macro allows the motor positions and scalers contents to be displayed while they are changing during a scan. The rate of updates is set by the `UPDATE` global variable, just as with the `umv` and `uct` macros. An option to the `setplot` macro, presented earlier, allows real-time plots of the data points to be displayed as they are measured.

Another option selected in the `setscans` macro lets you choose whether to do prescan motor limit checks with reciprocal space scans. The purpose of these checks is to avoid running into a software motor limit in the middle of a scan. For regular motor

scans, the limit checks are done only at the scan endpoints. For reciprocal space scans, the motor positions do not necessarily change monotonically and so the motor limits must be checked at each scan point. Since this checking requires a time-consuming loop at the macro level, you may choose to disable the feature if you are confident your reciprocal space scans will not send a motor outside the ranges defined by software limits.

Motor Scans

```

ascan motor start finish intervals time
a2scan m1 s1 f1 m2 s2 f2 intervals time
a3scan m1 s1 f1 m2 s2 f2 m3 s3 f3 intervals time
mesh m1 s1 f1 intervals1 m2 s2 f2 intervals2 time
lup motor start finish intervals time
dscan motor start finish intervals time
d2scan m1 s1 f1 m2 s2 f2 intervals time
d3scan m1 s1 f1 m2 s2 f2 m3 s3 f3 intervals time
th2th tth_start_rel tth_finish_rel intervals time

```

The argument *motor* (and *m1*, *m2* and *m3*) is a motor number or mnemonic, such as *th*, *tth*, *chi*, or *phi*. *ascan*, *a2scan* and *a3scan* are single-, two- and three-motor absolute-position scans. *mesh* is a nested two-motor scan, where the first motor scans through its range at each point of the second motor's scan. *lup* (or equivalently *dscan*), *d2scan* and *d3scan* are single- and two- and three-motor relative position scans. The starting and finishing positions are given relative to the current position, and the motors are returned to their starting position at the end of the scan. These relative position scans are defined in terms of the absolute-position scans.

The *th2th* macro is a special case of the *d2scan* that will scan the *tth* and *th* motors, with *th* moving half the range of *tth*. Its definition is,

```

def th2th '
  if ($# != 4) {
    print "Usage: th2th tth_start_rel tth_finish_rel intervals time"
    exit
  }
  d2scan tth $1 $2 th ($1)/2 ($2)/2 $3 $4
,

```

Basic Reciprocal Space Scans

```
hscan start_H finish_H intervals time
kscan start_K finish_K intervals time
lscan start_L finish_L intervals time
hklskan s_H f_H s_K f_K s_L f_L intervals time
hklnesh Q1 s_Q1 f_Q1 intervals1 Q2 s_Q2 f_Q2 intervals2 time
```

The first three scans are special cases of `hklskan`, as in

```
def hscan '
    if ($# != 4) {
        print "Usage: hscan start finish intervals time"
        exit
    }
    waitall; getangles; calcHKL
    hklskan $1 $2 K K L L $3 $4
,
```

`hklnesh` does a grid scan of two reciprocal coordinates, where $Q1$ and $Q2$ are literally H, K or L, and $Q1 \neq Q2$.

A `waitmove`, `getangles` and `calcHKL` are done at the start of the scan to obtain the current diffractometer position to determine the values of the unspecified coordinates in `hscan`, `kscan`, `lscan` and `hklnesh`.

Special Reciprocal Space Scans

```
klradial angle start_radius finish_radius intervals time [H=expr]
hlradial angle start_radius finish_radius intervals time [K=expr]
hkradial angle start_radius finish_radius intervals time [L=expr]
klcircle radius start_angle finish_angle intervals time [H=expr]
hlcircle radius start_angle finish_angle intervals time [K=expr]
hkcircle radius start_angle finish_angle intervals time [L=expr]
```

The first three scans trace a path that would form a radial cut if projected on to the K - L , H - L or H - K planes at the specified angle in degrees from the positive K , H or H axes, respectively. `start_radius` and `finish_radius` specify the radial distance from the origin. The unscanned coordinate will be reevaluated at each point according to the optional expression in the last argument, which can be a function of the other coordinates, for example, $H=L/300$. Otherwise, the unscanned coordinate will remain constant.

The second three scans differ only in that they trace out a circular arc in the projected plane at the radius given by the first argument. `start_angle` and `finish_angle` are the endpoints in degrees of the arc of the scan.

Temperature Scans

```
# Temperature scan
tscan start finish intervals time [sleep]
# Delta temperature scan
dtscan start finish intervals time [sleep]
```

These two macros scan the temperature setpoint. The macro `settemp` (see the temperature control macros, described below) is called to change the setpoint. A call of the `sleep()` function is done after calling `settemp`, but before counting, if the optional argument `sleep` is greater than zero.

Powder Mode

```
setpowder                # Select "powder" mode
setpowder off            # Turn "powder" mode off
setpowder motor full_width # Specify powder motor and width
```

Powder mode enables you to measure intensities while a motor is scanned through some range. When turned on, powder mode affects all scans where motors are moved. If invoked without arguments `setpowder` prompts for the powder motor and for the full width of the rocking movement to take place at each point of the scan. Invoked with the `off` argument, scans return to their normal mode.

Customizing Scan Output

To allow you to customize the scan headers and the information saved with each data point, several macros are available for redefinition. Their default definitions are:

```
def Pheader ''          # Printer header
def Fheader ''          # File header
def Plabel  '""'        # Printer/Video column labels
def Pout    '""'        # Printer/Video output at each point
def Flabel  '""'        # File column labels
def Fout    '""'        # File output at each point
```

Four of these must be defined as strings—in the default case, null strings. Here are examples of how these macros might save temperature set point and measurement information on the output devices.

```

def Pheader 'printf("\n Setpoint = %g (%g C)\n", TEMP_SP,DEGC_SP)'
def Fheader '_cols++;printf("#X %gKohm (%gC)\n",TEMP_SP,DEGC_SP)'
def Plabel 'sprintf("%7.7s %7.7s ", "T-set", "T-degC")'
def Pout 'sprintf("%7.5g %7.5g ",TEMP_SP,DEGC)'
def Flabel '"DegC "'
def Fout 'sprintf("%g ",DEGC)'

```

The `Pheader` and `Fheader` macros must print newline-terminated lines. More than one line is permitted, however. Besides adding scan header information to the data file, `Fheader` also adjusts the value of the global variable `_cols`, which is used to indicate the number of data columns in the data file. In the example shown, the `Flabel` and `Fout` definitions add one column to the data file, so `_cols` is incremented by one. The `Plabel` and `Pout` macros add columns to the printer (and screen) output. The columns in `Flabel` should be separated by double spaces (the data file convention). The columns in the other headers should be separated by single spaces. In each case, the spaces come after the label.

The `Ftail` macro is available for adding scan results to the data file at the conclusion of a scan. By default `Ftail` is defined as nothing,

```
def Ftail '' # File tail
```

You might define it to be

```
def Ftail 'printf("#R %d %g %g %g %g %g %g\n", \
                SCAN_N, pl_xMAX, pl_MAX, pl_FWHM, pl_CWHM, pl_COM, pl_SUM);'
```

where the values being printed are from the `pl_anal()` function described on page ##. The `#R` characters begin the data file control line for scan results.

Temperature Control Macros

```

te                # Read or set the temperature
settemp           # Set the temperature
measuretemp       # Measure the temperature
showtemp          # Show temperature parameters
teramp            # Ramp the temperature

```

Methods for handling temperature control and other experimental parameters are likely to vary greatly from lab to lab and experiment to experiment. You may be able to modify these standard macros to suit your specific needs.

The temperature control model assumed by these macros uses two independent instruments: one instrument to control the temperature and one instrument to measure the temperature. The following global variables are used by the macros:

TEMP_SP	The set point of the controller in ohms, volts, etc.
T_LO_SP	The lower limit for the controller set point.
T_HI_SP	The upper limit for the controller set point.
DEGC_SP	The temperature from which the set point is derived.
TEMP_CS	The value of the temperature sensor in ohms, volts, etc.
DEGC	The measured temperature.

The macro below displays the current set point and measured temperature.

```
# Display temperature parameters
def showtemp '
    measuretemp
    printf("Temperature Setpoint = %g (%gC)\n",TEMP_SP,DEGC_SP)
    printf("                Measured = %g (%gC)\n",TEMP_CS,DEGC)
,
```

You must supply the macro `measuretemp`. It should read `TEMP_CS` from the temperature sensor and convert it to `DEGC`. Sample `measuretemp` macros are given below.

The `te` macro is the one you would use most often to display or set the temperature set point.

```
# Simple read or set temperature
def te '
    if ($# == 1) {
        settemp $1
        qcomment "Temperature Setpoint at %g" "TEMP_SP"
    }
    onp; showtemp; offp
,
```

If invoked without arguments, it simply displays the current temperature parameters. Otherwise it invokes the `settemp` macro. The `settemp` macro checks its argument against the set point limits and then calls the `_settemp` macro, which you must supply.

```
# Assign the temperature setpoint
def settemp '
    if ($# != 1) {
        print "Usage:  settemp set_point"
        exit
    } else {
        local _1
        _1 = $1
        if (_1 < T_LO_SP || _1 > T_HI_SP) {
            printf("Temp limits are %g to %g.\n",T_LO_SP,T_HI_SP)
            exit
        }
        TEMP_SP = _1
        _settemp
    }

```



```

    }
,

```

Here are examples of `_settemp` macros from several installations (the symbol `_1` is defined in `settemp`):

```

# Write setpoint to a Lakeshore 82C Controller on GPIB bus
def _settemp '
    gpib_put(12, sprintf("S%6.4f", _1))
,
# Write setpoint to home-made GPIB device used at MIT
def _settemp '{
    local _s
    _s = int(32767*_1/10)
    gpib_put(4, sprintf("%c%c%c%c\160\200",\
        0x80|(_s &0xF), 0x90|((_s>>4)&0xF),\
        0xA0|((_s>>8)&0xF), 0xB0|((_s>>12)&0xF)))
}'
# Write setpoint to a home-made device used with CAMAC at Harvard
def _settemp '
    ca_put(bcd(10000*_1), 0, 0)
,

```

Here are examples of different `measuretemp` macros:

```

# Read parameters from a Lakeshore 82C Controller on GPIB bus
def measuretemp '{
    local _s
    gpib_put(12, "W0")
    _s=gpib_get(12)
    TEMP_SS=substr(_s,1,6)
    TEMP_CS=substr(_s,9,6)*100
    RtoT_0 DEGC TEMP_CS
    TEMP_SP=substr(_s,17,6)
    RtoT_0 DEGC_SP TEMP_SP
}'
# Read setpoint from CAMAC and temperature from GPIB device
def measuretemp '
    TEMP_CS = gpib_get(1)/1000
    RtoT_0 DEGC TEMP_CS
    TEMP_SP = dcb(ca_get(0, 0))/10000
    RtoT_0 DEGC_SP TEMP_SP
}'

```

Keep in mind that `measuretemp` is also called at each iteration of the standard scan macros.

The macro `RtoT_0`, used above, is one of several in the standard package that convert between degrees C and kilohms for common thermistors:

```
# Temperature to kohms
def TtoR_0 '
    local _k          # YSI 44011 (100kohm @ 25C)  20 to 120 C
    $1 = exp(-11.2942      +5.3483e3      /(_k = ($2) + 273.15)\
              -1.42016e5  /(_k*_k) -1.172e7      /(_k*_k*_k))
    ,
# Kohms to temperature
def RtoT_0 '
    local _l          # YSI 44011 (100kohm @ 25C)  20 to 120 C
    $1 = (1/(+2.2764e-3      +2.20116e-4 *(_l = log($2))\
              +2.61027e-6 *_l*_l      +9.02451e-8 *_l*_l*_l) - 273.15)
    ,
```

(The four parameters in each equation were obtained by fitting a table of values supplied by the manufacturer of the thermistors. No guarantees are made about the accuracy of the fitted parameters.)

The following macro will gradually change (or ramp) a temperature controller to a new set point. If the ramp time is greater than 500 seconds, the temperature is changed every 10 seconds, otherwise the temperature is changed every 2 seconds.

```
# Read or set or ramp the temperature
def teramp '{
    if ($# == 1) {
        te $1
    } else if ($# == 2) {
        local _i _s1 _f1 _d1 _rtime _stime
        _f1 = $1
        _rtime = $2
        _stime = _rtime < 500? 2:10
        _s1 = TEMP_SP
        _d1 = (_f1 - _s1) / _rtime * _stime
        qcomment "Ramp Temp Setpoint from %g to %g" "_s1,_f1"
        for (_i=0; _i<=_rtime; _i+=_stime, _s1 += _d1) {
            settemp _s1
            measuretemp
            printf("Set=%7.4f Meas=%7.4fC\r",TEMP_SP,DEGC)
            sleep(_stime)
        }
        showtemp
    } else {
        print "Usage:  teramp set_point  or  teramp set_point time"
        exit
    }
}'
```

Printer Initialization Macros

These macros send out the particular character sequences that put various printers into compressed mode to fit 132 columns of text on 8½" wide paper.

```
# Put DecWriter into compressed mode
def initdw 'onp; offt; printf("\033[4w"); offp; ont'

# Put Epson printer into compressed mode
def initfx 'onp; offt; printf("\017"); offp; ont'

# Put Okidata printer into compressed mode
def initoki 'onp; offt; printf("\035"); offp; ont'

# NEC P6/P7 printer, put into compressed mode
def initnec 'onp; offt; printf("\033!\004"); offp; ont'
```

The Scan Macros In Detail

All the scan macros in the standard package share a similar structure. To keep the format of the output sent to the data file, printer and screen consistent, common parts of each scan are defined as macros that are called by all the scans. For example, the `scan_head` macro is called by each scan to write scan headers on all the output files and devices. Certain macros are shared by all the scans for another reason. Special operating modes or options are implemented by redefining shared macros. For example, the `scan_move` macro, called within the looping portion of the scans, is normally defined as `_move`, which is:

```
def _move 'move_em; waitmove; getangles; calcHKL'
```

In powder mode, `scan_move` is defined as `_pmove`, a slightly more complicated macro, designed to move the designated powder averaging motor some width on alternating sides of the center trajectory of the scan,

```
def _pmove '  
  if (_stype&2)  
    _cp = A[_pmot]  
  A[_pmot] = _cp + _pwid/2  
  _pwid = -_pwid  
  move_em; waitmove; getangles; A[_pmot] = _cp; calcHKL  
,
```

The following paragraphs explain in detail the construction of the scan macros, using the single-motor scan, `ascan`, as an example. Here is its definition:

```
def ascan '  
  if ($# != 5) {  
    print "Usage:  ascan motor start finish intervals time"  
    exit  
  }  
  _check0 "$1"  
  { _m1 = $1; _s1 = $2; _f1 = $3; _n1 = int($4); _ctime = $5 }  
  
  if (_n1 <= 0) {  
    print "Intervals <= 0"  
    exit  
  }  
  
  _bad_lim = 0  
  _chk_lim _m1 _s1  
  _chk_lim _m1 _f1  
  if (_bad_lim) exit  
  
  HEADING = sprintf("ascan %s %g %g %g %g", "$1", $2, $3, $4, $5)  
  _d1 = (_f1 - _s1) / _n1++  
  _cols=4
```

```

X_L = motor_name(_m1)
_sx = _s1 ; _fx = _f1
_stype = 1|(1<<8)
FPRNT=sprintf("%s H K L", motor_name(_m1))
PPRNT=sprintf("%8.8s", motor_name(_m1))
VPRNT=sprintf("%9.9s", motor_name(_m1))
scan_head
def _scan_on \'
  for (; NPTS < _n1; NPTS++) {
    A[_m1] = _s1 + NPTS * _d1
    scan_move
    FPRNT=sprintf("%g %g %g %g",A[_m1],H,K,L)
    PPRNT=sprintf("%8.4f",A[_m1])
    VPRNT=sprintf("%9.4f",A[_m1])
    scan_loop
    pl_put(NPTS, A[_m1], S[DET])
    scan_plot
  }
  scan_tail
\'
_scan_on

```

In `ascan`, as in all scans, the first thing to do is to check the number of arguments, `$#`, and if incorrect, print a usage message:

```

if ($# != 5) {
  print "Usage:  ascan motor start finish intervals time"
  exit
}

```

Next, the `_check0` macro is called,

```
_check0 "$1"
```

as it is whenever a motor mnemonic is used as an argument in the standard macros. The macro checks its argument against all valid motor mnemonics and motor numbers. The purpose is to prevent unintentionally sending motors into motion if the user mistypes a mnemonic. The definition of `_check0` is

```

def _check0 '{
    local _i

    for (_i = 0; _i <= MOTORS; _i++)
        if (_i == MOTORS) {
            print "Invalid motor name: $1"
            exit
        } else if ($1 == _i) {
            if ("$1" != motor_mne(_i) && "$1" != _i) {
                print "Invalid motor name: $1"
                exit
            } else
                break
        }
    }
,

```

Next in `ascan`, the global variables used in the scan are initialized from the arguments.

```
{ _m1 = $1; _s1 = $2; _f1 = $3; _n1 = int($4); _ctime = $5 }
```

The global variables being assigned are shared by all the scans.

Next in `ascan`, a check is made to ensure the number of intervals is positive.

```

if (_n1 <= 0) {
    print "Intervals <= 0"
    exit
}

```

The next four lines do a motor limit check before the start of the scan.

```

_bad_lim = 0
_chk_lim _m1 _s1
_chk_lim _m1 _f1
if (_bad_lim) exit

```

The `_chk_lim` macro sets the flag `_bad_lim` if the position given by the second argument is outside the limits of the motor given by the first argument.

```

def _chk_lim '{
    local _u _t

    if ((_u = dial($1, $2)) < (_t = get_lim($1, -1))) {
        printf("%s will hit low limit at %g.\n", motor_name($1), _t)
        _bad_lim++
    } else if (_u > (_t = get_lim($1, 1))) {
        printf("%s will hit high limit at %g.\n", motor_name($1), _t)
        _bad_lim++
    }
}'

```

The prescan limit check is straightforward for simple motor scans. For reciprocal

space scans, the limit check must loop through all the points of the scan since the motor positions are not necessarily monotonic functions of the scan variables.

Next in `ascan`, the global variable `HEADING` is initialized.

```
HEADING = sprintf("ascan %s %g %g %g %g", "$1", $2, $3, $4, $5)
```

It is used in the scan headers written to the file, screen and printer, and records the arguments with which the scan was invoked.

Next, some global scan variables are initialized.

```
_d1 = (_f1 - _s1) / _n1++
_cols=4
X_L = motor_name(_m1)
_sx = _s1 ; _fx = _f1
_stype = 1|(1<<8)
```

The `_d1` variable is set to the step size for the scan. The number of intervals in `_n1` is incremented so its value will be the actual number of points. The `_cols` global variable is set to the number of extra columns this scan will use in the data file. Here it is four, for the motor position and values of H , K and L at each point.

`X_L` is set to the x -axis label to use on the plot of the scan. The globals `_sx` and `_fx` are set to the endpoints of the x axis to be used in plotting the data on the screen during the scan.

The variable `_stype` is treated as a two byte integer and holds a code representing the current scan type. The low-order byte is a bit flag, while the high order byte contains a number value. The expression `1|(1<<8)` use the bitwise-or and the bitwise-shift operators to put values in each byte. Currently, the following codes are used:

Code	Type Of Scan	High-Order Byte
1	motor	number of motors
2	HKL	nothing
4	temperature	nothing

Next in `ascan`, the global variables `FPRNT`, `PPRNT` and `VPRNT` are given string values to be used for file, printer and video-screen column labels particular to this scan.

```
FPRNT=sprintf("%s H K L", motor_name(_m1))
PPRNT=sprintf("%8.8s", motor_name(_m1))
VPRNT=sprintf("%9.9s", motor_name(_m1))
```

Each label contains the name of the motor being scanned, although printed with a different field width. Different widths are used to fit the widths and number of fields on the target devices. A challenge in constructing the scan macros is to fit all the desired columns of information within a single line. All the scans must limit the line

length to 132 columns for output sent to the printer. (80-column printers must be operated in compressed mode to make their carriages effectively 132 columns wide.) The video screen is 80 columns wide. For the data file, there is no restriction on width. Also for the data file, no attempt is made to line up items in columns.

Next in `ascan` is the `scan_head` macro, called to do the general initialization. All scan macros call `scan_head`. The default definition of `scan_head` is

```
def scan_head '_head'
```

where `_head` is defined as,

```
def _head '
  _scan_time
  waitall; get_angles; calcHKL
  NPTS = T_AV = MT_AV = 0
  DATE = date()
  TIME = TIME_END = time()
  _cp = A[_pmot]
  rdef cleanup "_scanabort"

  # DATA FILE HEADER
  ond; offt
  printf("\n#S %d %s\n#D %s\n",++SCAN_N,HEADING,DATE)
  if (_ctime < 0)
    printf("#M %g (%s)\n", -_ctime, S_NA[MON])
  else
    printf("#T %g (%s)\n", _ctime, S_NA[sec])
  printf("#G")
  for (_i=0; _i<NPARAM; _i++) printf(" %g", G[_i])
  printf("\n")
  printf("#Q %g %g %g\n", H, K, L)
  {
    local _i _j _k
    for (_i = 0, _k = MOTORS; _i < _k; _i += 8) {
      printf("#P%d ", _i/8)
      _mo_loop .6g "A[mA[_j]]"
    }
  }
  Fheader
  printf("#N %d\n", _cols + 3)
  printf("#L %s%s Epoch %s %s\n",FPRNT,Flabel,\
    S_NA[_ctime < 0? sec:MON],S_NA[DET])
  offd

  # PRINTER HEADER
  onp; offt
  printf("\n\nScan %3d %s file = %s %s user = %s\n%s\n\n",\
    SCAN_N,DATE,DATAFILE,TITLE,USER,HEADING)
  {
    local _i _j _k
    for (_i = 0, _k = MOTORS; _i < _k; _i += 8) {
```



```

        printf(" ")
        _mo_loop 9.9s "motor_name(mA[_j])"
        printf(" ")
        _mo_loop 9.6g "A[mA[_j]]"
    }
}
Pheader
printf("\n # %11.9s %11.9s %11.9s %8.8s %8.8s %8.8s %s%s\n",\
      "H", "K", "L", S_NA[sec], S_NA[MON], S_NA[DET], PPRNT, Plabel)
offp

# TTY HEADER
ont
printf("\nScan %3d  %s  file = %s  %s  user = %s\n%s\n\n",\
      SCAN_N, DATE, DATAFILE, TITLE, USER, HEADING)
printf(" # %s %8.8s %8.8s %10.10s%s\n",\
      VPRNT, S_NA[DET], S_NA[MON], S_NA[sec], Plabel)
,

```

The commands at the beginning of `_head`,

```
waitall; get_angles; calcHKL
```

insure the motors are stopped and positions current before proceeding. Next, `_head` initializes some variables. `NPTS` is the loop variable in the scans that will run from 0 to `_n1`. `T_AV` and `MT_AV` maintain the average temperature (from the global variable `DEGC`) and the average monitor counts or time per point during the scan. `DATE` and `TIME` are set to the current date and time. `TIME_END` is updated at each point with the current time. The `_cp` variable is used in powder mode and is set to the center position of the powder-average motor.

Next in the header macro, the real space motor positions and the reciprocal-space position are made current with `getangles` and `calcHKL`. The cleanup macro is defined to be the standard macro `_scanabort`. The macro named `cleanup` is special as `spec` automatically invokes that macro when a user types `^C` or on any other error, such as hitting motor limits, trying to go to an unreachable position or encountering a syntax error in a macro. The definition of `_scanabort` is,

```

def _scanabort '
    _cleanup2
    _cleanup3
    comment "Scan aborted after %g points" NPTS
    sync
    undef cleanup
,

```

The `_cleanup2` macro is defined for delta scans to move motors back to their starting positions. The `_cleanup3` macro is available to users for defining some kind of private clean up actions.

Finally, the headers are written to the file, printer and screen in turn. Included in the headers are the user-defined `Fheader`, `Flabel`, `Pheader` and `Plabel`.

Returning back to `ascan`, the next part of the macro is the loop:

```
def _scan_on \'
  for (; NPTS < _n1; NPTS++) {
    A[_m1] = _s1 + NPTS * _d1
    scan_move
    FPRNT=sprintf("%g %g %g %g",A[_m1],H,K,L)
    PPRNT=sprintf("%8.4f",A[_m1])
    VPRNT=sprintf("%9.4f",A[_m1])
    scan_loop
    pl_put(NPTS, A[_m1], S[DET])
    scan_plot
  }
  scan_tail
\'
_scan_on
```

The loop is implemented as a macro to enable the scan to be continued with the `resume` macro. The relevant global variables are initialized outside the loop, so that invoking `_scan_on` continues the scan where it had left off when interrupted. Here is the `resume` macro.

```
def resume '
  if (NPTS >= (index(HEADING, "mesh")? _n1*_n2 : _n1)) {
    print "Last scan appears to be finished."
    exit
  }
  def cleanup "_scanabort"
  comment "Scan continued"
  _scan_on
,
```

The `scan_move`, `scan_loop` and `scan_plot` macros are invoked by all the scans. In the loop, the motor array `A[]` is set to the target position for the scanned motor and the motor is moved using the `scan_move` macro, normally defined as `_move`:

```
def _move 'move_em; waitmove; getangles; calcHKL'
```

String variables are then assigned to values that will be written to the output devices using the `scan_loop` macro. The `scan_loop` macro is generally defined as `_loop` which has the definition,

```

# The loop macro, called by all the scans at each iteration
def _loop '
    scan_count _ctime
    measuretemp
    calcHKL
    z = _ctime < 0? S[sec]/1000:S[MON]
    T_AV += DEGC; MT_AV += z
    printf("%3d %s %8.0f %8.0f %10.6g%s\n",\
        NPTS,VPRNT,S[DET],S[MON],S[sec]/1000,Pout)
    onp; offt
    printf("%3d %11.5g %11.5g %11.5g %8.6g %8.0f %8.0f %s%s\n",\
        NPTS,H,K,L,S[sec]/1000,S[MON],S[DET],PPRNT,Pout)
    offp; ond; offt
    printf("%s%s %d %g %g\n",FPRNT,Fout,(TIME_END=time())-EPOCH,z,S[DET])
    offd; ont
,

```

This macro first counts by calling the `scan_count` macro, normally defined as `_count`, which is, in turn, defined as `count`. (In powder mode, or when using updated counting during scans, `_count` is defined differently.) The `_loop` macro then calls `measuretemp`. With this macro, you can have any per-point action done, not limited to, nor necessarily even including, measuring the temperature of the sample. Next in `_loop` the sums for computing the average temperature and monitor count rate are adjusted. Finally the video screen, printer and data file are updated with the results of the current iteration.

The last thing in `_scan_on` is a call to `scan_tail`, normally defined as `_tail`:

```

# The tail macro, called by all the scans when they complete
def _tail '
    undef cleanup
    TIME_END = time()
    if (!(_stype&8)) {
        ond; offt
        Ftail
        offd; ont
        plot
    }
,

```

This macro removes the definition of `cleanup`, since it is no longer needed, and if not a *mesh* scan, adds the user defined results to the file and calls the `plot` macro.

Standard Data-File Format

The data files created by the macros have a simple format. The files are ASCII. Control lines in the file begin with a # followed by an upper-case letter. Other lines are blank or contain scan data.

The control conventions are:

Code	Parameters	Description
#C	<i>comment</i>	Comments inserted by many of the standard macros.
#D	<i>date</i>	A string representing the current date, in the format Wed May 4 23:59:49 1988.
#E	<i>seconds</i>	The UNIX epoch at the time the file was opened.
#F	<i>filename</i>	The name by which the file was opened.
#G	<i>parameters</i>	The contents of the geometry parameter array $G[]$.
#L	<i>labels</i>	Labels for the data columns.
#M	<i>counts</i>	If counting to monitor counts, the number of counts.
#N	<i>columns</i>	The number of columns of data that follow.
#O0	<i>motor names</i>	The names of the first 8 motors.
#O1	<i>motor names</i>	The names of the second 8 motors, etc.
#P0	<i>motor positions</i>	The starting positions of the first 8 motors.
#P1	<i>motor positions</i>	The starting positions of the second 8 motors, etc.
#Q	<i>H K L</i>	The starting values of a scan in reciprocal space.
#R	<i>results</i>	Results from various <code>pl_anal()</code> operations optionally printed at the conclusion of a scan.
#S	<i>number</i>	A new scan having scan number <i>number</i> follows, normally preceded by a blank line.
#T	<i>seconds</i>	If counting to time, the time used.
#U		Reserved for user.
#X	<i>setpoint</i>	The temperature setpoint.

FOUR-CIRCLE REFERENCE

Introduction

When invoked by the name *fourc*, *spec* runs with code appropriate for a four-circle diffractometer. This section of the *Reference Manual* focuses on the features of *spec* unique to the *fourc* version.

The four circles of the standard four-circle diffractometer are: 2θ , the angle through which the beam is scattered, and θ , χ , and ϕ , the three Euler angles, which orient the sample. Of these three, θ is the outermost circle with its axis of rotation coincident with that of 2θ . The χ circle is mounted on the θ circle, with its axis of rotation perpendicular to the θ axis. The ϕ circle is mounted on the χ circle such that its axis of rotation lies in the plane of the χ circle.

From the keyboard and on the screen, the angles are named `tth`, `th`, `chi` and `phi`, respectively, and conventionally referred to in that order. For *fourc* to work properly, angles with these names must be configured.

In describing the operation of a four-circle diffractometer, it is convenient to consider three coordinate systems: 1) a frame fixed in the laboratory, 2) a frame fixed on the spectrometer and 3) the natural axes of the sample. Note that *fourc* uses right-handed coordinate systems. All rotations are right-handed except for the χ rotation.

- (1) The x - y plane of the laboratory coordinate system is called the scattering plane and contains the sample and the points reached by the detector as it rotates on the 2θ arm. A counter-clockwise rotation of the 2θ axis corresponds to increasing 2θ , with the 2θ rotation axis defining the positive z direction in the laboratory. The zero of 2θ is defined as the setting at which the undeflected X-ray beam hits the detector.

The positive y axis is along the line from the sample to the X-ray source. The position at which θ rotates the χ circle to put the χ rotation axis along the y axis defines the zero of θ . A clockwise rotation of χ corresponds to increasing χ .

The zero of χ is the position which puts the ϕ rotation axis along the positive z axis. The positive x axis direction is determined by the cross product of the y and z axes ($\hat{\mathbf{x}} = \hat{\mathbf{y}} \times \hat{\mathbf{z}}$), which completes the definition of the right-handed coordinate system.

It is important to note that the zeroes of 2θ , θ and χ and the direction of positive rotation of all the circles must be set as described above and cannot be freely redefined.

- (2) The spectrometer coordinate system is defined as a right-handed system fixed on the ϕ rotation stage at the sample position such the coordinate system is aligned with the laboratory coordinate system when all four spectrometer angles are zero. This definition determines the zero of ϕ .
- (3) The third coordinate system is aligned with specific directions in the sample. A common and useful example are coordinates defined as the lattice vectors of a crystalline sample. When placing a sample in the spectrometer, it is unlikely that its axes will line up with the spectrometer axes. Nevertheless, *fourc* allows the sample orientation to be fully specified by finding the angles at which two Bragg peaks are detected and giving the corresponding reciprocal lattice indices. This process is described fully in the section on the *Orientation Matrix*.

To orient a sample so as to measure the intensity at a particular reciprocal lattice position requires that the reciprocal lattice vector of interest is aligned with the scattering vector of the spectrometer. Since any rotation about the scattering vector does not change the diffraction condition, there is a high degree of degeneracy that must be resolved in order for *fourc* to determine unique angle settings. How the degeneracy is lifted described in the section on *Four-Circle Modes*.

Diffraction Alignment

This section presents a guide on how to set up a four-circle spectrometer. This summary applies whether the scattering plane is horizontal or vertical.

The first step, which should need only be done the first time *fourc* is used with the diffractometer, is to ensure each diffractometer motor is set up with the correct name, mnemonic, rotation sense, steps-per-degree, etc. The `config` macro is generally used for this purpose. For stepping motors, the rotation sense of an axis depends on the details of the motor controller and cable connections. If the rotation sense isn't as described in the *Introduction*, change the *sign_of_user_x_dial* parameter in the motor configuration file.

For each motor, *fourc* keeps track of both a dial and a user position. The dial position is meant to agree with the readout of the physical dial on the spectrometer. The value and the sign of the *steps_per_unit* parameter should be chosen so that the dial position and its direction in the computer agree with the physical dial reading. Use the `set_dial` macro to set the dial positions. The user positions should correspond to the underlying “true” orientation angles of the spectrometer that satisfy the constraints given above. Use the `set` macro to set the user positions.

Once properly configured, diffractometer alignment proceeds as follows.

- (1) Arrange for the X-ray beam to go through the center of rotation. Generally, the center of rotation is found with a pin and a telescope.
- (2) Arrange for the X-ray beam to be perpendicular to the 2θ axis. This condition is typically verified by comparing X-ray burns made on X-ray sensitive paper with 2θ near the undeflected beam direction and with 2θ offset by 180° .
- (3) Set 2θ so that the undeflected X-ray beam direction corresponds to the zero of 2θ .
- (4) Align the χ rotation axis with the laboratory y axis to set the zero of θ .
- (5) Align the ϕ rotation axis with the θ rotation axis to set the zero of χ .

One way to do (4) and (5) is as follows:

- (i) Mount a Si(111) wafer so that the (111) direction is (approximately) along the ϕ axis.
- (ii) Find the (111) Bragg reflection. Note the values of θ and χ . Call them θ_1 and χ_1 .
- (iii) Rotate ϕ by 180° .
- (iv) Find the Bragg reflection again. Note the values of θ and χ . Call them θ_2 and χ_2 .
- (v) $\frac{1}{2}(\chi_1 + \chi_2)$ corresponds to $\chi = 90^\circ$ in a correctly aligned spectrometer; $\frac{1}{2}(\theta_1 + \theta_2)$ corresponds to $\theta = \frac{1}{2}2\theta$.

The Huber four-circle diffractometer is an example of an instrument that works with *fourc*.¹ Another common spectrometer configuration has two crossed $\pm 20^\circ$ tilt stages on top of full θ and 2θ circles. This configuration is compatible with the four-circle code if the tilt stage immediately adjacent to the θ circle is χ and the other is ϕ . When the first tilt stage is zero, χ is at 90° .

¹The four-circle Huber has dial readings with all right-handed rotations, so the χ circle should have dial readings and user readings that are in opposite senses. If the θ circle is offset by 180° then the dial readings and the user readings of all angles can have the same sense.

Orientation Matrix

Angle calculations for the four-circle diffractometer are described in detail in Busing and Levy.² You may also refer to that paper to learn how to calculate the *orientation matrix*. The orientation matrix, \mathbf{UB} , describes the sample orientation with respect to the diffractometer angles. Given \mathbf{UB} , it is possible to calculate the diffractometer angles $(2\theta, \theta, \chi, \phi)$ necessary to rotate a particular scattering vector \mathbf{Q} , indexed by (H, K, L) , into the diffraction position. The matrix \mathbf{B} transforms the given (H, K, L) into an orthonormal coordinate system fixed in the crystal. The matrix \mathbf{U} is the rotation matrix that rotates the crystal's reference frame into the spectrometer's.

The first step in constructing an appropriate orientation matrix is to enter the sample crystal lattice parameters a, b, c, α, β and γ .³ These are real-space parameters, as might be found in Wychoff⁴ or Pearson.⁵ Use the macro `setlat` to assign values:

```
1.FOURC> setlat 3.61 3.61 3.61 90 90 90
```

```
2.FOURC>
```

Next, you must specify the sets of values of $(2\theta, \theta, \chi, \phi)$ at which two Bragg reflections are in the diffracting position. One of these is called the primary reflection. *Fourc* ensures that the values of (H, K, L) reported for the primary reflection agree (to within a scale factor) with the values entered. However, because of experimental errors and/or uncertainties in the unit cell parameters, the values of (H, K, L) reported for the other Bragg reflection, called the secondary reflection, may not agree perfectly with the entered values (although they should be close).

You can use the `or0` and `or1` macros to enter the parameters for the primary and secondary reflections, respectively. However, the `or0` and `or1` macros require that the diffractometer be moved to the associated reflections, as these macros use the current angles and the entered (H, K, L) in the calculation of the orientation matrix. Alternatively, you can use the `setor0` and `setor1` macros, which prompt for both (H, K, L) and the angles that define the orientation matrix, without moving the spectrometer to the given settings.

²W. R. Busing and H. A. Levy, *Acta Cryst.* **22**, 457 (1967). Note however, that this paper uses right-handed coordinates systems and left-handed rotations for all rotations except for χ , which is right-handed.

³The conventional symbols for the crystal lattice angles include α and β . These angles are unrelated to the orientation angles α and β defined in the introduction. The different meanings should be clear from context.

⁴R. W. G. Wychoff, *Crystal Structures* (Wiley, New York, 1964).

⁵P. Villars and L. D. Calvert, *Pearson's Handbook of Crystallographic Data for Intermetallic Phases* (American Society for Metals, Metals Park, Ohio, 1985).

Four-Circle Modes

As noted above, because there are three Euler angles (θ, χ, ϕ) while the direction of the scattering vector \mathbf{Q} is specified by only two angles, there is a degeneracy associated with the transformation from (H, K, L) to $(2\theta, \theta, \chi, \phi)$. The degeneracy is resolved in *fourc* by providing a constraint. In *fourc*, the different constraints are called *modes*. The value of the `g_mode` geometry parameter determines the prevailing mode.

Fourc defines several angles in order to specify certain of the modes.

The angle ω is defined as $\omega = \theta - (2\theta)/2$, and is referred to as `OMEGA`.

The angle ψ (referred to as `AZIMUTH`) specifies a counter-clockwise rotation about the diffraction vector. The zero of ψ is determined by a reference vector, different from the diffraction vector. The azimuthal angle ψ is defined to be zero when this vector is in the diffraction plane.

The angles α (`ALPHA`) and β (`BETA`) are defined such that the angles between the azimuthal reference vector and the incident and scattered X-rays are $90^\circ - \alpha$ and $90^\circ - \beta$, respectively. One commonly used azimuthal reference vector is the sample's surface normal, which then makes α and β correspond to the incident and exit angles of the X-rays on the surface.

Omega Equals Zero (`g_mode = 0`)

The simplest constraint is $\omega = 0$. Any (H, K, L) with a small enough 2θ can be reached in this mode.

Omega Fixed (`g_mode = 1`)

In this mode ω is fixed to a finite value. Suppose you want to go to the $(0, 0, 2)$ Bragg reflection but with $\omega = 10$. You would then type

```
2.FOURC> setmode 1
```

```
Now using Omega-Fixed mode.
```

```
3.FOURC> OMEGA=10
```

```
4.FOURC> br 0 0 2
```

```
5.FOURC>
```

Zone or χ and ϕ Fixed ($g_mode = 2$)

A zone of reciprocal space is a plane passing through the origin. A zone axis is the direction of the normal to this plane. The zone axis can be specified by the vector product of any two non-colinear points in the zone. **Zone** mode forces the specified zone axis to be normal to the scattering plane. In other words, a zone of the crystal is leveled into the scattering plane of the spectrometer. Any point in this plane can then be reached with θ and 2θ (i.e., χ and ϕ are fixed). The geometry code provides for the calculation of the χ and ϕ necessary to put any two reciprocal space positions in the scattering plane via the `cz`, `sz` or `mz` macros. These macros are explained later.

Phi Fixed or Three Circle ($g_mode = 3$)

The angle ϕ is fixed at some arbitrary value.

Azimuth Fixed ($g_mode = 4$)

This mode fixes the value of the rotation angle ψ of a reference vector about the scattering vector **Q**.

Azimuth-fixed mode provides a degree of control that is particularly useful in a surface-diffraction experiment. If the reference vector is chosen to be the surface normal, setting ψ to 90° rotates the surface normal into the plane defined by the scattering vector and the diffractometer θ - 2θ axis. This means that the incidence angle α will equal the exit angle β .

The reference vector is defined by the geometry parameters `g_haz`, `g_kaz`, and `g_laz`. Use the macro `setaz` to specify the (H, K, L) of the reference vector. For example, to set the reference vector to $(0, 0, 2)$, use

```
5.FOURC> setaz 0 0 2
```

```
6.FOURC>
```

Azimuth-fixed mode will fail if you try to make measurements with **Q** parallel to the reference vector, since there is then no way to define a rotation about **Q**. The remedy is either to switch to another mode (usually the best choice) or to switch to another reference vector.

Alpha Fixed ($g_mode = 5$)

This mode allows you to hold the value of α constant while moving to various values of (H, K, L) . This is useful in experiments in which it is necessary to control the X-ray penetration depth into a sample. More generally, if looking at a weak signal from

a surface, keeping α small will keep the background small with no restriction (in principle) on the momentum transfer normal to the surface.

Suppose you are studying a sample of copper and you want to fix the incidence angle to be equal to the critical angle for total external reflection. Set the reference vector to the surface normal of the crystal and then type the following commands:

```
6.FOURC> setmode 5  
Now using Alpha-fixed mode.
```

```
7.FOURC> ALPHA=0.4126
```

```
8.FOURC> br 0 2 2
```

```
9.FOURC>
```

To implement **alpha-fixed** mode, *fourc* calculates the value of ψ needed to fix the angle between the incident wave vector and the reference vector.⁶

Beta Fixed (`g_mode = 6`)

This is the same as **alpha-fixed** mode except that β , rather than α , may be fixed.

⁶S. G. J. Mochrie, *J. Appl. Cryst.* **21**, 1-3 (1988).

Freezing Angles

For the **omega-fixed**, **phi-fixed**, **zone**, **azimuth-fixed**, **alpha-fixed** and **beta-fixed** modes you may *freeze* the value of the associated angle (or angles), so when calculating motor positions corresponding to an arbitrary (H, K, L) using `calcA` (within `br`, for example), the angle (or angles) will be reset to the frozen value before the calculation no matter what the current position of the diffractometer.

The macro `freeze` waits until all motors have stopped moving, then sets a variable (`g_frz`) indicating frozen mode is on and saves the current position of the frozen angle in another variable. Usage might be:

```
9.FOURC> setmode 1
Now using Omega-Fixed mode.

10.FOURC> freeze 5

Sun Jan 28 12:16:23 1990. Freezing Omega at 5.

11.FOURC>
```

If the value to freeze for the current mode is not given as an argument to the macro, the current value of the related angle or parameter is used to set the frozen value. In **zone** mode, both the χ and ϕ values need to be given as arguments.

The macro `unfreeze` sets `g_frz` to zero. Subsequent angle calculations will use whatever the current value of the associated constrained angle or angles for the current mode.

Sectors

Sectors correspond to different symmetry transformations of $(2\theta, \omega, \chi, \phi)$ and may be of help in avoiding blind spots. Sectors can also be useful for samples in cryostats or ovens. All modes can have the positions chosen for the motors further influenced by the choice of sector. However, the modes **azimuth-fixed**, **alpha-fixed** and **beta-fixed** only allow sectors numbered zero through three, below.

For a given pair of incident and scattered X-ray beams, \mathbf{k}_i and \mathbf{k}_f , there are eight orientations of the crystal in the spectrometer that give the same scattering since they present equivalent projections to the incident and scattered beams. The eight orientations are labeled as sectors 0 through 7.

Four of the orientations come from the symmetries of a pair of vectors. They correspond to the identity operation (i.e., the current diffraction angles), a rotation of 180° about \mathbf{k}_i , a rotation of 180° about the bisector of \mathbf{k}_i and \mathbf{k}_f , and a rotation about z by $180^\circ - 2\theta$. The last two symmetries are based on interchanging the role of the entrance angle and the exit angle of the X-rays on the sample. These four symmetries give sectors 0, 2, 4 and 6 in *fourc*. For each of these four positions another position can be obtained by increasing θ by 180° , decreasing ϕ by -180° and changing the sign of χ . Since a rotation of 180 and of -180 are the same, both θ and ϕ can be increased by 180 . These orientations give sectors 1, 3, 5 and 7. Studying these operations shows that sectors 2, 3, 6 and 7 have opposite signs of 2θ from the current position. Sectors 2, 3, 4 and 6 have flipped the *up* direction of the sample (normal to the scattering plane) to the *down* direction.

The value of `g_sect` determines in which sector of reciprocal space the diffractometer operates.

The actual transformations of the angles are:

	0	1	2	3	4	5	6	7
$2\theta \rightarrow$	2θ	2θ	-2θ	-2θ	2θ	2θ	-2θ	-2θ
$\omega \rightarrow$	ω	$\omega - 180^\circ$	$-\omega$	$180^\circ - \omega$	$-\omega$	$180^\circ - \omega$	ω	$\omega - 180^\circ$
$\chi \rightarrow$	χ	$-\chi$	$\chi - 180^\circ$	$180^\circ - \chi$	$180^\circ - \chi$	$\chi - 180^\circ$	$-\chi$	χ
$\phi \rightarrow$	ϕ	$\phi - 180^\circ$	ϕ	$\phi - 180^\circ$	$\phi - 180^\circ$	ϕ	$\phi - 180^\circ$	ϕ

In addition, a sector 8 is defined that can be used when the χ and ϕ circles of the diffractometer are arc segments in contrast to the complete circle. Sector 8 minimizes $|\chi - 90^\circ|$ and $|\phi|$.

Cut Points

The angles -90° and $+270^\circ$ are at the same position on a circle, but if a motor in `spec` is sitting at 0° , it will move counter-clockwise by 90° to get to the -90° position and clockwise by 270° to get to the $+270^\circ$ position.

By setting a lower cut point for a particular motor, you can choose what value on the circle the four-circle angle calculations will produce. That is, the calculations will place the angle between the lower cut point and that value plus 360° .

Cut points can be set for the θ , χ and ϕ four-circle angles. For the 2θ angle, the lower cut point is fixed at -180° . For the **azimuth-fixed** modes, there is pseudo-cut point for the azimuthal angle. If the cut point is less than zero, the calculated azimuth will be between -180° and 0° , otherwise the azimuth will be between 0° and $+180^\circ$.

The default cut points are shown below.

Angle	Default Cut Point
θ	-180
χ	-180
ϕ	-180
ψ	0

The macro `cuts` can be used to set or display the cut points. With no arguments, it displays the current cut points. Two arguments are used to set a single cut point. Four arguments set all cut points.

```
11.FOURC> cuts
```

```
Cut Points:
```

```
   th   chi   phi  azimuth
-180 -180 -180         0
```

```
12.FOURC> cuts phi 90
```

or

```
13.FOURC> cuts -180 -180 90 0
```

```
14.FOURC>
```


Four-Circle Files

Most questions regarding the behavior of `spec` when used with a four-circle diffractometer can ultimately be resolved by consulting the appropriate source code. The file `geo_fourc.c` in the standard distribution of `spec` contains the code for all the four-circle calculations. The file `u_hook.c` contains a few lines of code that connect the code in `geo_fourc.c` with the rest of the program. Finally, the files `macros/fourc.src` and `macros/ub.mac` in the standard distribution contains the definitions for all the four-circle and orientation matrix macros.

Four-Circle Variables

The four-circle coordinate variables (H , K and L , in particular) are stored in a built-in array named `Q[]`. The four-circle geometry calculations either use the motor positions contained in the `A[]` array to calculate values for the `Q[]` parameters or place motor positions in `A[]` based on the current values in `Q[]`. Each four-circle variable has a descriptive macro definition as an alias, such as `def OMEGA 'Q[6]'`.

Variable	Alias	Description
<code>Q[0]</code>	<code>H</code>	x component of the scattering vector.
<code>Q[1]</code>	<code>K</code>	y component of the scattering vector.
<code>Q[2]</code>	<code>L</code>	z component of the scattering vector.
<code>Q[3]</code>	<code>LAMBDA</code>	Incident X-ray wavelength λ .
<code>Q[4]</code>	<code>ALPHA</code>	Incident angle α .
<code>Q[5]</code>	<code>BETA</code>	Exiting angle β .
<code>Q[6]</code>	<code>OMEGA</code>	$\omega = \theta - (2\theta)/2$.
<code>Q[7]</code>	<code>AZIMUTH</code>	Azimuthal angle.
<code>Q[8]</code>	<code>F_ALPHA</code>	Frozen value of α for alpha-fixed mode.
<code>Q[9]</code>	<code>F_BETA</code>	Frozen value of β for beta-fixed mode.
<code>Q[10]</code>	<code>F_OMEGA</code>	Frozen value of ω for omega-fixed mode.
<code>Q[11]</code>	<code>F_AZIMUTH</code>	Frozen value of ψ for azimuth-fixed mode.
<code>Q[12]</code>	<code>F_PHI</code>	Frozen value of ϕ for phi-fixed mode.
<code>Q[13]</code>	<code>F_CHI_Z</code>	Frozen value of χ for zone mode.
<code>Q[14]</code>	<code>F_PHI_Z</code>	Frozen value of ϕ for zone mode.

The geometry parameters in the table below affect the geometry calculations in various ways. Although the parameters can be changed by assignment, the preferred method is to use the indicated macro for setting the parameters.

Variable	Alias	Related Macro	Description
G[0]	g_mode	setmode	Specifies the four-circle <i>mode</i> .
G[1]	g_sect	setsector	Specifies the <i>sector</i> .
G[2]	g_frz	freeze	Nonzero when an angle is <i>frozen</i> .
G[3]	g_haz	setaz	<i>H</i> of the azimuthal reference vector.
G[4]	g_kaz	setaz	<i>K</i> of the azimuthal reference vector.
G[5]	g_laz	setaz	<i>L</i> of the azimuthal reference vector.
G[6]	g_zh0	mz, sz	<i>H</i> of first zone -mode vector.
G[7]	g_zk0	mz, sz	<i>K</i> of first zone -mode vector.
G[8]	g_zl0	mz, sz	<i>L</i> of first zone -mode vector.
G[9]	g_zh1	mz, sz	<i>H</i> of second zone -mode vector.
G[10]	g_zk1	mz, sz	<i>K</i> of second zone -mode vector.
G[11]	g_zl1	mz, sz	<i>L</i> of second zone -mode vector.
U[0]	g_aa	setlat	<i>a</i> lattice constant in Angstroms.
U[1]	g_bb	setlat	<i>b</i> lattice constant.
U[2]	g_cc	setlat	<i>c</i> lattice constant.
U[3]	g_al	setlat	α lattice angle.
U[4]	g_be	setlat	β lattice angle.
U[5]	g_ga	setlat	γ lattice angle.
U[6]	g_aa_s	setrlat	a^* reciprocal lattice constant.
U[7]	g_bb_s	setrlat	b^* reciprocal lattice constant.
U[8]	g_cc_s	setrlat	c^* reciprocal lattice constant.
U[9]	g_al_s	setrlat	α^* reciprocal lattice angle.
U[10]	g_be_s	setrlat	β^* reciprocal lattice angle.
U[11]	g_ga_s	setrlat	γ^* reciprocal lattice angle.
U[18]	g_u00	or0, setor0	Observed 2θ of primary reflection.
U[19]	g_u01	or0, setor0	Observed θ of primary reflection.
U[20]	g_u02	or0, setor0	Observed χ of primary reflection.
U[21]	g_u03	or0, setor0	Observed ϕ of primary reflection.
U[24]	g_u10	or1, setor1	Observed 2θ of secondary reflection.
U[25]	g_u11	or1, setor1	Observed θ of secondary reflection.
U[26]	g_u12	or1, setor1	Observed χ of secondary reflection.
U[27]	g_u13	or1, setor1	Observed ϕ of secondary reflection.
U[12]	g_h0	or0, setor0	<i>H</i> of primary reflection.
U[13]	g_k0	or0, setor0	<i>K</i> of primary reflection.

U[14]	g_l0	or0, setor0	<i>L</i> of primary reflection.
U[15]	g_h1	or1, setor1	<i>H</i> of secondary reflection.
U[16]	g_k1	or1, setor1	<i>K</i> of secondary reflection.
U[17]	g_l1	or1, setor1	<i>L</i> of secondary reflection.

The first three parameters select modes and set flags. The next three parameters describe the components of the azimuthal reference vector. The six after that describe the *zone*-mode vectors.

The next sets of parameters describe the orientation matrix, including the lattice constants of the sample and the parameters of the primary and secondary orientation reflections. Remember that the `calcG` macro, described below, must be called to make sure the orientation matrix is recalculated after changing any of these related values above. The `or0`, `setor0`, `or1`, `setor1`, `or_swap` and `setlat` macros do just that.

Four-Circle Functions

You can access the four-circle calculations through `spec`'s user-hook routine `calc()`. The table below summarizes the calculations available.

Function	Alias	Description
<code>calc(1)</code>	<code>calcA</code>	Calculate motor positions for current <i>H K L</i> .
<code>calc(2)</code>	<code>calcHKL</code>	Calculate <i>H K L</i> for motor positions in <code>A[]</code> .
<code>calc(4)</code>	<code>calcG</code>	Recalculate orientation matrix.
<code>calc(5)</code>	<code>calcZ</code>	Calculate χ and ϕ for zone feature.
<code>calc(7,0)</code>	<code>calcD</code>	Calculate direct lattice from reciprocal parameters.
<code>calc(7,1)</code>	<code>calcR</code>	Calculate reciprocal lattice from direct parameters.
<code>calc(8)</code>	<code>calcE</code>	Calculate λ for current monochromator positions.
<code>calc(9)</code>	<code>calcM</code>	Calculate monochromator position for current λ .
<code>calc(10)</code>	<code>_begUB</code>	Initialize sums for fitting UB .
<code>calc(11)</code>	<code>_addUB</code>	Add a reflection to fitting sums.
<code>calc(12)</code>	<code>_fitUB</code>	Fit UB .
<code>calc(13)</code>	<code>calcL</code>	Calculate lattice parameters from UB .

Four-Circle Macros

The macros below are used in setting the parameters and selecting modes.

Name	Arguments	Description
setmode	1 optional	Choose geometry mode.
setsector	1 optional	Choose sector.
setlat	6 optional	Set lattice parameters.
setaz	3 optional	Set azimuthal reference vector.
setmono	2 optional	Set beam-line monochromator parameters.
or0	3 optional	Set primary orientation reflection.
or1	3 optional	Set secondary orientation reflection.
setor0	none	Alternative to set primary orientation reflection.
setor1	none	Alternative to set secondary orientation reflection.
or_swap	none	Swap values for primary and secondary vectors.
freeze	none	Turn on <i>freeze</i> mode.
unfreeze	none	Turn <i>freeze</i> mode off.
cz	6	Calculate zone.
sz	6	Set zone parameters.
mz	6	Move to zone.
cuts	2/6 optional	Show or set cut points.

Most of the macros with optional arguments will prompt for the required values if invoked without arguments. Using any of these macros to change a parameter will produce a comment on the printer and in the data file.

The `or0` macro is a typical example of these macros.

```
14.FOURC> prdef or0
def or0 '
  if ($# == 3) {
    _1 = $1; _2 = $2; _3 = $3
  } else if ($# == 0) {
    print "\nEnter primary-reflection HKL coordinates:"
    _1 = getval(" H", g_h0)
    _2 = getval(" K", g_k0)
    _3 = getval(" L", g_l0)
  } else {
    print "Usage:  or0  or  or0 H K L"
    exit
  }
  waitmove; get_angles
  gpset _1      g_h0      # gpset documents the change
  gpset _2      g_k0
  gpset _3      g_l0
  gpset A[mA[0]] g_u00
```

```

gpset A[mA[1]] g_u01
gpset A[mA[2]] g_u02
if (_numgeo > 3) { gpset A[mA[3]] g_u03 }
if (_numgeo > 4) { gpset A[mA[4]] g_u04 }
if (_numgeo > 5) { gpset A[mA[5]] g_u05 }
gpset LAMBDA g_lambda0
calcG

```

15.FOURC>

Zone Macros

Zone mode is controlled with the `cz` (calculate zone), `sz` (set zone) and `mz` (move zone) macros. Given two Bragg reflections, `cz` will calculate and display the values of χ and ϕ necessary to put both of these reflections in the scattering plane. To find the angles needed to put (0, 0, 2) and (0, 2, 2) in the scattering plane, type

```

15.FOURC> cz 0 0 2 0 2 2
Chi = 45 Phi = 90

```

16.FOURC>

Once appropriate values of χ and ϕ have been calculated, the scattering plane can be set using the `p1` (plane) macro, which moves the 2θ and θ motors together,

```

16.FOURC> p1 45 90

```

17.FOURC>

Alternatively, you can use the `mz` macro, which calculates the necessary χ and ϕ , moves there, sets **zone** mode, if not already in it, and saves the values of the zone vectors in the `G[]` geometry parameter array.

```

17.FOURC> mz 0 0 2 0 2 2

```

```

18.FOURC> p A[chi], A[phi]
45 90

```

19.FOURC>

The `sz` macro calculates and displays the χ and ϕ values, sets **zone** mode, if not already in it, saves the values of the zone vectors, sets the frozen values of **zone**-mode χ and ϕ , but does not move the diffractometer.

The `cz`, `sz` and `mz` macros make use of the `Z[]` array variables to pass the zone vectors to the geometry code.

```

19.FOURC> prdef cz
def cz '
    if ($# != 6) {
        print "Usage:  cz h0 k0 l0 h1 k1 l1"
        exit
    }
    Z[0]=$1; Z[1]=$2; Z[2]=$3; Z[3]=$4; Z[4]=$5; Z[5]=$6
    calcZ
    printf("Chi = %g  Phi = %g\n", A[chi], A[phi])
    waitmove; get_angles; calcHKL
'

20.FOURC>

```

Least-Squares Refinement of Lattice Parameters

In the previous sections, the procedure described for determining the orientation matrix required a knowledge of the lattice parameters of the crystal and the position of two reflections. When such information is unknown, the orientation matrix can be fit to an unlimited number of observed peak positions using a least-squares procedure.⁷ Lattice parameters derived from the fitted orientation matrix can then be calculated, although such lattice parameters are not constrained to exhibit any symmetry whatsoever.

In spec's implementation of least square refinement, three macros are used to create a file that contains the observed peak positions. That file is eventually run as a command file, and the least squares analysis is performed.

The `reflex_beg` macro initializes the reflections file:

⁷J. Matthews and R. L. Walker, *Mathematical Methods of Physics*, (Benjamin, Menlo Park, 1970), p. 391.

```

global REFLEX # Variable for file name
# Open the file, save old one as .bak and write header
def reflex_beg '{
    if ("$1" == "0") {
        if (REFLEX == "")
            REFLEX = "reflex"
        REFLEX = getval("Reflections file", REFLEX)
    } else
        REFLEX = "$1"

    if (open(REFLEX))
        exit
    close(REFLEX)
    if (file_info(REFLEX, "-s"))
        unix(sprintf("mv %s %s.bak", REFLEX, REFLEX))
    fprintf(REFLEX, "# %s\n\n_begUB\n\n", date())
}'

```

The reflex macro adds lines to the file that contain the (H, K, L) and $(2\theta, \theta, \chi, \phi)$ of each reflection:

```

# Add reflection to the file
def reflex '
    if ($# != 3) {
        print "Usage:  reflex H K L"
        exit
    }
    if (REFLEX == "") {
        REFLEX = getval("Reflections file", "reflex")
        if (REFLEX == "")
            exit
    }
    waitmove; get_angles; calcHKL
    fprintf(REFLEX, "H = %g;  K = %g;  L = %g\n", $1, $2, $3)
    {
        local i

        for (i=0; i<_numgeo; i++)
            fprintf(REFLEX, "A[%s]=%9.4f;  ", motor_mne(mA[i]), A[mA[i]])
        fprintf(REFLEX, "\n")
    }
    fprintf(REFLEX, "# counts = %g\n", S[DET])
    fprintf(REFLEX, "_addUB\n\n")
'

```

Finally, the `reflex_end` macro puts the proper trailer on the file:

```
# Add trailer to file
def reflex_end '
    fprintf(REFLEX, "_fitUB\n")
    printf("Type \"qdo %s\" to calculate new orientation matrix\n",\
          REFLEX)
',
```

When you are ready to calculate the orientation matrix, simply run the command file. There is no limit to the number of reflections contained in the file. You can also edit the file by hand to add or subtract reflections.

The calculated orientation matrix will remain valid until you type `calcG`, or invoke a macro that calls `calcG`. Those macros are `or0`, `or1`, `or_swap` and `setlat`.

The six original lattice parameters will remain unchanged when using the above macros to fit the orientation matrix to the reflections. The `calcL` macro can be invoked to calculate the lattice parameters derived from the fitted orientation matrix and place their values in the appropriate elements of the parameter array. The old lattice parameters will be lost.

Here is a sketch of the commands you use to perform the least squares refinement of the lattice parameters.

```
20.FOURC> reflex_beg
Reflections file (reflex)? <return>

21.FOURC> (find and move to a reflection ...)

22.FOURC> reflex 2 2 0

23.FOURC> (find and move to another reflection ...)

24.FOURC> reflex 2 0 2

25.FOURC> ...

26.FOURC> reflex_end
Type "qdo reflex" to recalculate orientation matrix.

27.FOURC> qdo reflex
Opened command file 'reflex' at level 1.

28.FOURC>
```

At least three reflections must be used for the least-squares fitting to work.

To calculate the new lattice parameters, use the `calcL` macro:

```
28.FOURC> calcL

29.FOURC> pa

4-Circle Geometry, Omega-Fixed (mode 1), frozen coordinate, sector 0

Primary Reflection (at lambda 1.54):
  Omega Chi Phi = 5.00025 90.102 -82.5455
  H K L = 2 2 0
Secondary (at lambda 1.54):
  Omega Chi Phi = 0 0 90
  H K L = 0 1 0
Azimuthal Reference:
  H K L = 2 2 0
Lattice Constants:
  a b c = 4.12644 4.12295 4.11078
  Alpha Beta Gamma = 90.0435 89.9763 90.1221

30.FOURC>
```

Here is a typical reflections file created by the above macros:

```
# Wed Jan 31 21:55:01 1990

_begUB

H = 2; K = 2; L = 0
A[tth]= 63.8185; A[th]= 36.9320; A[chi]= 89.8765; A[phi]= -80.0815
# counts = 3456
_addUB

H = 2; K = 0; L = 2
A[tth]= 63.8335; A[th]= 36.8920; A[chi]= 145.4185; A[phi]= 42.8145
# counts = 6345
_addUB

H = 0; K = 2; L = -1
A[tth]= 49.4100; A[th]= 29.6725; A[chi]= 35.8180; A[phi]= 45.9070
# counts = 5634
_addUB

H = 0; K = 2; L = -1
A[tth]= 49.3550; A[th]= 29.7225; A[chi]= 35.9380; A[phi]= 45.9070
# counts = 4563
_addUB

_fitUB
```

Within this file, the `calc()` function codes defined by the macros `_begUB`, `_addUB` and `_fitUB` are used to access the C code that performs the least squares operations.

ADMINISTRATOR'S GUIDE

Introduction

The first section of this guide outlines the procedure for installing `spec` on your computer. Later sections describe the format of some of the installed files.

Quick Install

For those who need little explanation, here are minimal installations instructions based on the customary configuration:

```
# First time, create a specadm user account, then
cd ~specadm
mkdir spec4.04.01          # choose name based on release
cd spec4.04.01
tar xvf ../spec_XXX.tar   # use pathname of distribution
# If an update, start with last version's parameters
cp ../specYYY/install_data . # use pathname of previous distribution
./Install                 # as root
```

The *Install* program will display the current installation parameters and prompt for changes. Once the software is installed, run the `spec` executable and type `config` to invoke the hardware configuration editor.

For those needing more detailed instructions, read on.

Steps For Installing spec

To install `spec` on your computer, be sure you have the software development tools available. `spec`'s requirements include the *make* utility, a compatible C compiler, and compatible runtime libraries for linking.

Before installing `spec` for the first time, you need to make several decisions:

Decide who will own the files. Most sites create a special user account, usually with the name *specadm*, that has a home directory used to hold the `spec` distribution files and is the user account assigned to own the installed `spec` files. Having a special *specadm* account allows you or other users to configure and edit `spec` files without invoking super-user powers, thus lowering the risk of making catastrophic errors.

Decide where to put the distribution. The distribution files are what come directly from CSS, usually as a *tar* file sent via internet, but sometimes on floppy disk, magnetic tape or CD-ROM. These distribution files need to be extracted into a distribution directory that serves as a “staging” area for the installation. The files needed

while running `spec` will be copied elsewhere during installation, so the distribution directory need not be always accessible to the `spec` users. However, it is usually advised to keep the distribution files available for reference. You will likely receive updated distributions from time to time, and it is a good idea to keep each distribution in a separate directory. The usual choice is to put the distribution in the `specadm` home directory in a subdirectory named after the `spec` version number, such as `~specadm/spec4.04.01`.

Decide where the installed files will go. Two directories are needed for the installed files. One is a directory for the executable programs. Most sites choose `/usr/local/bin`. This directory needs to be in each `spec` user's search path. The second directory is for `spec`'s auxiliary files. The usual choice is `/usr/local/lib/spec.d`.

Once these decisions are made and you have created a `spec` administrator user account (if used) and the directories mentioned above (if needed), you are ready to perform the installation. In brief, to install `spec` you:

- extract the distribution from the supplied media or tar file,
- run the installation program to install the `spec` files,
- install optional kernel drivers (seldom necessary),
- enter the hardware configuration for your particular experimental set-up.

In addition to the information presented below, be sure to consult any *README* files in the `spec` distribution directories for up-to-date information on installation procedures.

Extracting the Distribution

If you have made a `spec` administrator's account, you should become that user (or the root user), either by logging in or by using the `set-user id` command `su specadm`. You can then change to the `spec` administrator's home directory and use it as a place to hold the distribution files.

Make a subdirectory to receive the current version of `spec`. If the distribution is numbered release 4.04.01, you might make a directory called `spec4.04.01` using the command `mkdir spec4.04.01`. Change to the new directory with `cd spec4.04.01`.

The distribution will be in *tar* format. Usually the distribution is obtained as a tar file via internet. If the distribution arrives via magnetic media or a CD-ROM, the command to extract the tar file should be printed on the distribution label. For a tar file, the command to extract the files is

```
tar xvf specdist.tar
```

Installing the spec Program Files

To install the spec files, particularly for a first time installation, you may need to install as root. Use the `su` command to gain super-user privileges. If you are updating from previous distributions, you can copy the most recent *install_data* file containing your default installation parameters to the new distribution directory. Then, from the current spec distribution directory, type `./Install` to run the installation program.

The *Install* program will first indicate the current installation parameters. You may either accept those or enter new parameters. When entering new parameters, the default response to each question is given in parenthesis. Most questions present a number of choices. You can either type the number of the choice or you can type out the literal selection. If one of the options is the word “other”, such as for the name of a directory, you can directly type your selection when prompted. For most questions, the first choice listed is probably the best response. For example,

```
Choices for binaries directory are:
```

- 1) /usr/local/bin
- 2) /usr/local
- 3) /usr/local/spec/bin
- 4) /u/bin
- 5) /LocalLibrary/Spec/bin
- 6) other

```
Choose binaries directory (/usr/local/bin)?
```

Entering a single minus sign (-) will move back to the question for the previous parameter, allowing you to enter a different value.

The installation questions ask for the following parameters:

platform – Selects from the supported combinations of computer and operating system.

geometry – Selects from the supported diffractometer configurations.

installed name – Selects a name for the installed program. For configurations with special geometry code, the first four letters of the name must match the first four letters of the geometry configuration. Thus, *fourcL* and *fourcR* might be the names for the left- and right-hand sides of a rotating anode lab with two four-circle diffractometers.

additional geometries – Asks if you want to enter more combinations of the previous two items.

file ownership – Selects the name of the owner of the *spec* files.

binaries directory – Selects where the programs that users run directly from the shell will go. (It's better to put *spec* files in some place other than the standard */bin* or */usr/bin*, in order to be able to distinguish files that are standard UNIX from those that have been added locally.) This directory should be in each *spec* user's search path.

auxiliary directory – Selects where *spec* puts its auxiliary files. Among the auxiliary files are user's state files, which can occupy significant disk space if there are many users running *spec*. Thus, this directory should not be on a file system that is in danger of running out of space.

ESRF library directory – Gives the location of the ESRF device server libraries if installing on an ESRF device server platform.

EPICS library directory – Gives the location of the EPICS channel access libraries if installing on an EPICS platform.

config file permissions – Selects who can change the hardware configuration file. On a low-security site, select the first choice, which lets all *spec* users make changes as needed. See security notes on page 210 for additional considerations.

data points type – Selects the precision for *spec*'s internal data-points arrays. The choice affects both the virtual memory consumption of *spec* while it is running and the disk space consumed by user state files. Single precision data is suitable for most installations.

alternate history code – Selects an alternative library that can be linked with *spec* to provide a more powerful history recall syntax than the standard *spec* history mechanism. The alternative library also includes command line editing features. See the files in the *readline* subdirectory of the distribution for additional information.

CAMAC hardware support – Selects whether to include CAMAC hardware support in the installed programs.

KSC 3929 SCSI-CAMAC file location – Gives the location of Kinetic Systems supplied interface file for their model 3929 SCSI-to-CAMAC crate controller. If using the CSS-supplied software on the HP platform, enter `no`. If using the CSS-supplied driver on the SunOS 4.1 platform, enter `no`.

Nat Instr GPIB file location – Gives the path name of the *cib.o* file that may be needed to be linked with `spec` when using National Instruments GPIB boards. You may enter `no` for the PCII and AT-GPIB models of the board for PC computers, for the SB-GPIB board with versions 1.3 and 2.1 of the driver for Sun computers, for the 1014-1S board for Sun computers and for the boards for IBM PS/2 and RS/6000 AIX systems.

SICL GPIB library flags – Sets the flags required to load the *libsicl.a* file, necessary if you are using the SICL GPIB presently supported on HP 700 series platforms with either HP or IOtech devices. Note, CSS has an alternative to the IOtech version of *libsicl.a* for the SCIS488/H product. Contact CSS for more information.

VME hardware support – Selects whether to include direct support for VME hardware in the installed program. Note, this option does not refer to the VME hardware used in ESRF and EPICS environments.

Nat Instr VME file location – Gives the location of the *nivxi.o* file that will be needed to be linked in with `spec` if you are using a National Instruments VME controller.

Extra compiler flags – Allows you to add extra compiler flags for both compiling the site-dependent source files and linking.

Extra object files – Allows you to specify extra site-dependent object files to include when linking the `spec` executable.

Extra library flags – Allows you to specify extra site-dependent libraries to be searched during the link phase when producing the `spec` executable.

After answering the questions, the installation should then continue automatically, producing output similar to the following:

```
Checking if u_hook.c needs compiling ...
Compiling u_hook.c ...
Checking if u_hdw.c needs compiling ...
Compiling u_hdw.c ...
Checking if readline library is up to date (may take some time) ...
Checking if geo_fourc.c needs compiling ...
Checking if spec needs to be linked ...
Linking spec ...
Making directory /usr/local/lib/spec.d
Making directory /usr/local/lib/spec.d/help
Making directory /usr/local/lib/spec.d/hgr
Making directory /usr/local/lib/spec.d/fourc
Making directory /usr/local/lib/spec.d/fourc/userfiles
Installing config auxiliary files ...

Installing fourc ...
  ln /usr/local/bin/spec /usr/local/bin/fourc
Installing prototype config file ...
  Fixing "fourc" config and settings file permissions ...
  No settings file. Creating one.
Installing macros ...
  Installing /usr/local/lib/spec.d/standard.mac ...
  Installing /usr/local/lib/spec.d/four.mac ...

Installing help files ...
  Clearing out old help files ...
  Making the "help" help file ...
  Making the "changes" help file ...
  Copying help files ...
  Changing ownership of help files to specadm ...
  cp chk_file /usr/local/lib/spec.d
Installing the "chelp" program ...
Installing the "showscans" package ...
  cp show.awk /usr/local/lib/spec.d
  cp scans /usr/local/bin
Installing the "contents" program ...
  cp contents /usr/local/bin
Installing the "tidy_spec" program ...
  cp tidy_spec /usr/local/bin
```

If you change certain parameters that require relinking `spec` and don't see the `Linking spec ...` message when rerunning the *Install* program, simply remove the *spec* file and run *Install* again.

Installing Drivers

In the early history of `spec`, much of the hardware communication was done through drivers. A driver is a piece of software that becomes part of the operating system and allows `spec` to communicate with the particular hardware the driver supports or performs some other function not available with the standard operating system. Most current configurations of `spec` do not require any CSS-provided drivers.

Note that for many PC configurations where `spec` once required use of a driver, `spec` now supports driverless configurations. The hardware is controlled directly from user level using `inb()/outb()` type of calls. Frequent polling of status registers replaces the use of interrupts. The performance of `spec` using polled mode seems to be as effective as the interrupt-driven drivers. Those configurations that do require `spec` drivers are described below.

The files for `spec` drivers are in the *drivers* subdirectory of the standard distribution. If the file *drive.tar.Z* exists in the directory, run the *TYPEME* script to extract the driver files. A *README* file in that directory contains up-to-date information about the currently available drivers and instructions on installing them.

Some hardware devices on some platforms, such as the National Instruments GPIB controllers on most platforms, require a vendor-supplied driver.

`spec` drivers are available for the items in the following list.

- The DSP 6001/6002 CAMAC crate controller with the DSP PC004 IBM PC interface on System V release 3.2 and 4.x, SCO XENIX 386 3.2.3 and SCO UNIX platforms. Note, this controller can also be used in a driverless configuration, where all I/O to the CAMAC is through user-level `inb()/outb()` calls. Contact CSS to discuss performance considerations.
- The Kinetic Systems 3922 CAMAC controller and 2926 PC board on on System V release 3.2 and 4.x platforms. Note, this controller can also be used in a driverless configuration, where all I/O to the CAMAC is through user-level `inb()/outb()` calls. Contact CSS to discuss performance considerations.
- The Kinetic Systems 3912 CAMAC controller on BSD and ULTRIX platforms.
- The Kinetic Systems 3929 SCSI CAMAC crate controller on SunOS 4.1 and Solaris 2.x platforms. Note, this controller can be used on HP-UX platforms with special code built into `spec` – no driver is necessary.
- The Oregon Micro Systems PCX, PC38 or PC39 PC board motor controller on System V release 3.2 and 4.x platforms. It is also possible to use this board in a driverless, polled mode. Contact CSS to discuss performance considerations.

- The Oxford/Tennelec/Nucleus PCA II multichannel analyzer card on System V release 3.2, 4.x and *linux* platforms. It is also possible to use this board in a driverless, polled mode. However, dead-time correction will not be available in that case. The newer model PCA-3 board doesn't require a driver in any case.
- On PC platforms using System V Release 3 UNIX, a special "nap" driver must be installed to implement sub-second sleeps. This driver is not needed on System V Release 4 platforms.
- On SCO XENIX, SCO UNIX and AIX PS/2 platforms, a special driver must be installed to enable access to I/O ports. If any of the hardware devices you use require I/O port access, install this driver.
- On *linux* platforms, *spec* is installed as a set-user-id root program, in order that the *iopl()* system call can be executed. Note, though, that the root privilege is only enabled for the duration of that call.

For the hardware controllers that require an interrupt on ISA-bus platforms, you have to carefully select from the limited interrupts available. With 8-bit boards, such as the DSP PC004 and the OMS motor controllers, there are only four interrupts available: 3, 4, 5 and 7. The conventional PC device assignments are COM2, COM1, LP2 and LP1 respectively. (Interrupts 0, 1, 2 and 6 are used by the clock, keyboard, slave interrupt controller and floppy controller.) When you install a new board that requires an interrupt vector, you may have to disable one of the above devices. They are generally eliminated in the order LP2, COM2, COM1, LP1.

If you run out of 8-bit interrupts, it is possible to wire the interrupt signal from an 8-bit board to an unused interrupt line in a 16-bit slot. Contact CSS for advice.

Selecting the Hardware Configuration

The final step in the initial *spec* installation is to set the hardware configuration specific to your site. You can either do that by starting the *spec* program and typing *config*, which is macro that executes the configuration editor, or you can run the configuration editor directly. For the latter, first change to */usr/local/lib/spec.d* (or to the auxiliary file directory specified when you did the installation). If you are installing the normal four-circle version of the program, type the command *edconf fourc*. If you installed a different geometry, give the name of that geometry as the argument. This command starts a spread-sheet styled program that lets you select motor parameters, devices names, CAMAC slot assignments, etc.

Refer to the notes on the configuration editor that follow for instructions on using *edconf*.

Fine Tuning the Hardware Configuration

You may wish to fine tune the installation by avoiding linking in hardware support for devices you will never be using at your site. The advantage is only that the `spec` executable images will be smaller and will require a bit less memory when running. The savings will be only on the order of 100 kbytes or so, depending on how many hardware modules are eliminated. You select which hardware modules are included by editing the file `u_hdw.c` and undefining the appropriate hardware keywords. After doing so, you will need to redo the installation, although you may use the `-p` flag for a partial installation to save time.

Adding Site-Dependent Help Files

If a file named `.local` in the `help` subdirectory of the auxiliary file directory exists and contains a list of file names, those names will be added to the topics contained in the `help` help file when `spec` is installed.

When `spec` starts up, the help file `news`, which is provided by CSS, and the file `local`, if it exists, will be displayed. The help file format is described on page 70 in the *Reference Manual*.

Adding Site-Dependent C Code

This step applies only to sophisticated end users of `spec` who understand the C language and need to customize `spec` for specific, site-dependent uses. Most readers can skip to the next section. Note also, local code can be accessed using the *data-pipe* facility explained on page 113 in the *Reference Manual*.

`spec` has provisions for end users to add their own C code to the program. User-added code is accessed using the built-in `calc()` function. If you wish to incorporate non-standard calculations within the `spec` program, you can do so by adding hooks for the code in the `u_hook.c` source file. C code that you add should, in general, be limited to calculations. You should avoid I/O, signal catching, etc. Consult CSS for specific information about what is appropriate for including in user-added C code. The `geo_*.c` files in the standard `spec` distribution that contain the X-ray diffractometer geometry code are examples of site-dependent code.

Within *u_hook.c* there is a routine called `init_calc()`. This routine is called once when `spec` starts up. Within `init_calc()`, calls to the routine

```
ins_calc(num, func)
int      (*func)();
```

insert the C routine *func* in a table of functions. These functions are called when `calc(num)` or `calc(num, arg)` is typed as a command to `spec`. The routine *func()* should be specified as either

```
func(num)
```

or

```
func(num, arg)
double arg;
```

depending on whether `calc()` is to be invoked with one or two arguments.

Any return value from *func()* is ignored. However, you can have the `calc()` routine return a value by assigning a number to the variable

```
extern double calc_return;
```

in *func()*. If no explicit assignment is made to `calc_return`, `calc()` returns zero.

The argument *num* can be from 0 to 63, but must be chosen not to conflict with any of the other `ins_calc()` entries already existing in *u_hook.c*.

You can also create built-in arrays of double precision, floating point numbers that can be used to communicate values between your C code and the user of the program. The routine

```
ins_asym(x, n, s)
double **x;
char    *s;
```

inserts the array *x* consisting of *n* elements into the table of built-in symbols. The character pointer *s* points to a string containing the name used to refer to the array from `spec` command level. For example,

```
#define N_PARAM 28
double *gparam[N_PARAM];
init_calc() {
    ...
    ins_asym(gparam, N_PARAM, "G");
    ...
}
```

inserts the 28-element array referred to as `G[]` into the program. Since the array `gparam[]` is an array of pointers, you must use the indirection operator (*) when referring to the values of the floating point numbers in your C code, as in

```
...
*gparam[3] = 1.54;
...
if (*gparam[2] == 0)
    ...
```

If you make any changes to *u_hook.c*, you must relink and reinstall the `spec` binary.

Updating spec

`spec` updates are normally extracted into a directory named after the `spec` version number, as in `/usr/specadm/spec4.04.01` for version 4.04.01 of the software. The instructions above for extracting the distribution and installing the files also apply to updates.

Existing *settings* and *config* files from a previous installation will not be disturbed during an update. Each user's state file will also remain intact, although it is recommended that users start out fresh (by typing `fourc -f`) after new versions are installed to incorporate improvements to the standard macros in their state files. Occasionally, new versions will not be compatible with previous state files, and `spec` will automatically throw out the old state files, print a message and start fresh anyway. Use the *tidy_spec* utility to clean out old state files to free up disk space, especially if the old state files are obsolete with the newer version of `spec`.

Update notes that accompany each new distribution indicate whether new versions of any drivers need to be installed. Also, the help file *changes* will contain summaries of the significant bug fixes and improvements to `spec` included in the update.

Installed Files

File Hierarchy

After installation and site configuration on a computer running, for example, both a four-circle and a z-axis diffractometer, the spec file hierarchies would appear something like the following,

```

/usr/local/bin-----|
                    |-camac
                    |-chelp
                    |-contents
                    |-dpmake
                    |-four
                    |-scans
                    |-showscans
                    |-show_state
                    |-tidy_spec
                    |-wiz_passwd
                    |-zaxis
                    |
                    |-README
                    |-chk_file
                    |-edconf          |-config
                    |-four.mac       |-settings
                    |-fourc-----|-conf.mac
                    |               |-userfiles-----|
                    |
                    |-data_pipe-----|-data_pipe.mak
                    |                 |-data_pipe.o
                    |                 |-pipe_test.c
                    |                 |-user_pipe.h
                    |
                    |-help-----|-angles
                    |            |-ackno
                    |            |-changes
                    |            |-config
                    |            |- ...
                    |
                    |-hgr-----|-x11filt
                    |           |-pcfilt
                    |           |- ...
/usr/local/lib/spec.d-|-io_enable
                    |-show.awk
                    |-site.mac
                    |-site_f.mac
                    |-standard.mac  |-config
                    |-zaxi.mac       |-settings          |-user_ttyL
                    |-zaxis-----|-conf.mac           |-user_ttyM
                    |               |-userfiles-----|-user_ttyP
                    |               |                 |-user_ttyS
                    |               |                 |- ...

```

where */usr/local/bin* is the installation directory, configured as INSDIR in the

Makefile, and */usr/local/lib/spec.d* is *SPECD*, the auxiliary file directory. Of the programs installed in */usr/local/bin*, *contents* and *showscans* are described in the *User Manual*, while the *camac* utility program is described below. The *chelp* utility is a stand-alone help file viewer that allows browsing of the *spec* help files without having to run *spec*.

The subdirectory *fourc* contains the files specific to the four-circle diffractometer. The name *fourc* matches the name by which the program is invoked. The first four letters of the name determine the geometry configuration. If a single computer is to control two spectrometers, they could be called *fourcL* and *fourcR*, and the *Install* program would create separate subdirectories called by those names for each.

Within the diffractometer directory is the associated configuration file, *config*, which specifies the hardware and the motor parameters to be used. Also associated with each diffractometer is the *settings* file that tracks changes in the motor position and limit settings. The *edconf* program can be used to to modify the contents of these two files.

The subdirectory *userfiles* contains each user's state files on a per terminal basis. These files allow the user to exit *spec* and restart at a later time, retaining macro definitions, variable assignments, etc. The *spec* administrator may, from time to time, delete old state files for users not expected back again, especially if disk space is a problem. The *tidy_spec* program reports on the disk usage of all the *userfiles* directories and provides options for removing files by age, user, tty or spectrometer geometry. Type `tidy_spec -` from the shell for usage options.

The CAMAC Debugging Program

The utility *camac* is a program that allows you to debug some CAMAC modules, but only on systems that use a CSS-provided CAMAC kernel driver. Invoked as `camac -z`, it will initialize the crate by issuing a dataway Z. Without arguments, the program lets you send arbitrary commands to any module. See the source code in *drivers/camac.c* for details on how the program works. The *camac* program only works when the CAMAC access is through CSS supplied drivers. Within *spec*, the functions `ca_cntl()` and `ca_fna()` allow similar functionality for debugging CAMAC devices.

Enabling I/O Ports On PC Platforms (not *linux*)

The *spec* distribution includes a program called *io_enable* for use on systems derived from AT&T System V/386 UNIX. This program enables user programs (such as *spec*) to write to selected I/O ports on the adapter cards plugged into the motherboard. *spec* automatically invokes *io_enable* if the hardware configuration requires it. Ports only need to be enabled once after booting the computer. The total number of ports that can be enabled is quite limited, though. The *io_enable* program must be owned by root and have set-user-id execute mode.

You can invoke *io_enable* from the shell using the following syntax:

```
io_enable [-d] port_address [port_address ...]
```

The `-d` option causes the specified ports to be removed from the systems list of enabled ports. Port addresses must be given as hexadecimal numbers with a `0x` prefix. The addresses may be specified either individually, as a hyphen-separated range, or as a base address followed by a comma and a count. Usage might be,

```
io_enable 0x300-0x302 0x312,4
```

to add 0x300, 0x301, 0x302, 0x312, 0x313, 0x314 and 0x315 to the enabled-port list.

Enabling I/O Ports On PC Platforms Running *linux*

On PC systems running *linux*, I/O port access is achieved by changing the process's I/O privilege level using the system call *iopl()*. For this method to work, the *spec* executable must be set-user-id root. Note, though, that the root privilege is only enabled for the duration of the system call. The rest of the time, the process only has the access privileges of the real user. (Prior to *spec* release 4.03.01, I/O port access on *linux* required a kernel patch.)

Enabling E/ISA I/O Ports On HP 700 Systems

The process for making the E/ISA I/O ports on HP 700 platforms involves a few steps. A config file for the HP *eisa_config* utility is required, such as the following, which is appropriate for the Oregon Micro Systems PC boards:

```
BOARD
    ID = "OMS0001"
    NAME = "OMS PCX/PC38 Motor Controllers"
    MFR = "Oregon Micro Systems"
    CATEGORY = "OTH"
    SLOT = ISA8
    IOCHECK = INVALID
    DISABLE = UNSUPPORTED

FUNCTION = "I/O Registers"
    TYPE="OTH"
    CHOICE="default"
    LINK
        PORT = 330h
        SIZE = BYTE
        SHARE = NO
    LINK
        PORT = 331h
        SIZE = BYTE
        SHARE = NO
    LINK
        PORT = 332h
        SIZE = BYTE
        SHARE = NO
    LINK
        PORT = 333h
        SIZE = BYTE
        SHARE = NO
```

A similar file needs to be created for each ISA card you intend to use with *spec* and placed in the directory */etc/eisa* with a name that matches the board ID, as in *!OMS0001.CFG* for ID = "OMS0001". The IDs are arbitrary for ISA boards. Be sure to enter the PORT addresses for all the I/O ports to be used on the card.

After the configuration files are installed, you must run the HP *eisa_config* utility to save the configuration. You must also reboot the computer to make the boards available.

When *spec* starts up, it will use the CSS *hp_ports* utility to create a special file that will be used to access the ports. For example, the special file */dev/ioport/0x330* will be created to access the ports described in the example above. *spec* will create one special file for each eight consecutive I/O ports.

The Configuration Editor

The *edconf* program is the primary means for maintaining the hardware configuration. When running *spec*, *edconf* is usually run by invoking the `config` macro. Without arguments, *edconf* will use the *config* and *settings* files in the current directory. If given a directory name as an argument, it will use the files in that directory. If invoked with the `-s` flag, *edconf* will run in *simulate* mode, allowing you to view but not modify the files. If you do not have write permission for the *config* file, *edconf* will automatically run in *simulate* mode.

To get a list of the available commands while running *edconf*, type a question mark (?). The following commands are available:

<i>Arrow keys</i>	Move around.
h j k l	Move left, down, up, right.
<code><return></code>	Enter data or move down one row.
<code><space></code>	Move right.
^F	Scroll forward through motors on the Motor screen, counters on the Scalers screen, devices on the CAMAC screen, items on the drop-down menus, etc.
^B	Scroll backwards, as above.
'	Enter to change string-valued cell (e.g. motor names).
+ - > <	Step through list of choices (if <code><></code> appears in label).
r	Reread <i>settings</i> and <i>config</i> files.
w	Write <i>settings</i> and <i>config</i> files.
R	Read from backup <i>settings</i> and <i>config</i> files.
c	Change to next screen.
M	Change to Motor screen.
m	Step through motor parameter screens.
C	Change to CAMAC screen.
I	Change to Interfaces screen.
D	Change to Devices screen.
A	Change to MCA/CCD Acquisition screen
s	Change to Scalers (counters) screen.
i	Insert a motor at current position.
d	Delete the motor at current position.
G	Toggle all-motors mode with linked configurations.
^D	Blank out optional motor parameter fields.
^L	Refresh the screen.

- `^w` Gain wizard access to set protections.
- `H` Print help information for current screen.
- `?` Print command help window.
- `^c` Exit.

There are five types of data cells in the configuration spread sheet. For number-valued cells, simply enter the number value when the cell is highlighted. For string-valued cells, such as motor names, you must first type a single quote (') before entering the string. For YES/NO cells, type a `y` or an `n`. The fourth type of cell has the characters `<>` before the cell label. For these cells, use the `+`, `-`, `<`, or `>` keys to step through the possible choices. For all types of cells, the `<return>` key enters the values. The fifth type of cell is used for entering motor unit/channel values. For these cells, type just the channel number for unit zero, or type both numbers separated by a literal `/`.

Use the `w` command to write out the configuration, and use `^c` to exit the program.

The Settings File

The *settings* file is a binary file that contains consecutive data for each motor according to the following structure:

```

struct sav_mot {
    long    sm_pos;           /* Current dial position */
    float   sm_off;          /* Current user/dial offset */
    double  sm_low;          /* Software low limit */
    double  sm_high;         /* Software high limit */
};

```

The *settings* file must have write permission for everybody who runs `spec`, as it is updated every time someone moves a motor or changes an offset or limit. When `spec` starts out, it checks to see if there is `At` present, there is no fail-safe file locking mechanism to prevent more than one program from writing to the *settings* file at the same time. It is up to the users to make sure they are not running two versions of `spec` on the same computer using the same diffractometer or that they are not running *edconf* and moving motors with `spec` at the same time.

The Config File

The *config* file is an ASCII file that describes the diffractometer hardware configuration. Although the *config* file can be edited by hand, you will be safer using the *edconf* program to make modifications as *edconf* insures the *config* file obeys the structuring rules required by `spec`.

Comment lines in the *config* file begin with a #. Other lines contain key words specifying devices, CAMAC slots or motor and counter parameters. Key words are followed by a space-delimited equals sign and one or more parameters.

The *config_adm* help file contains up-to-date information about currently recognized key words and supported hardware devices. The *Hardware Reference* section of the manual describes the specific hardware devices recognized by *spec* and indicates the *config* file syntax required to specify each device.

CAMAC Slots

CAMAC slot assignments in the *config* file consist of a module code on the left and a slot number on the right. For example,

```
CA_KS3610 = 2
```

tells the program a Kinetic Systems 3610 hex scaler is in slot 2.

The following modules names are recognized by *spec*. More than one of the modules marked with an asterisk are allowed. Append *_#* to number modules consecutively, where *#* is 0, 1, 2, etc.

CA_BR5302*	BiRa 5302 ADC as counters
CA_DSP2190	DSP Technology MCS Averager
CA_DXP*	XIA DXP MCA
CA_E250*	DSP E250 12-Bit D/A as Motor Controller
CA_E500*	DSP Technology E500A Stepper Motor Controller
CA_IOM1	BiRa 2601 I/O For E500 Multiplexing
CA_IOM2	F16,A0 I/O For E500 Multiplexing
CA_IOM3	F16,A1 I/O For E500 Multiplexing
CA_IO*	Any module to be accessed with F codes of 0 or 16
CA_KS3112*	Kinetic Systems 3112 D to A (as motor controller)
CA_KS3116*	Kinetic Systems 3116 16-Bit D/A as Motor Controller
CA_KS3195*	Kinetic Systems 3195 16-Bit D/A as Motor Controller
CA_KS3388	Kinetic Systems 3388 GPIB interface
CA_KS3512*	Kinetic Systems 3512/14 ADC as counters
CA_KS3610	Kinetic Systems 3610 6-Channel, 50 MHz Counter
CA_KS3640C*	Kinetic Systems 3640 Up/Down Counter as Counter
CA_KS3640M*	Kinetic Systems 3640 Up/Down Counter (for SMC's)
CA_KS3640T	Kinetic Systems 3640 Up/Down Counter as Timer
CA_KS3655	Kinetic Systems 3655 8-Channel Timing Pulse Generator
CA_KS3929	Kinetic Systems 3929 SCSI Crate Controller
CA_KS3929_HP	Kinetic Systems 3929 SCSI to CAMAC on HP

CA_LC2301	LeCroy 2301 interface for QVT MCA
CA_LC3512	LeCroy 3512 Spectroscopy ADC
CA_LC3521	LeCroy 3521A Multichannel Scaling
CA_LC3588	LeCroy 3588 Multichannel Scaler
CA_LC4434*	LeCroy 4434 32-Channel Scaler
CA_LC8206	LeCroy MM8206A Histogramming Memory
CA_QS450*	DSP Technology QS-450 4-Channel Counter
CA_RTC018	DSP Technology RTC-018 Real Time Clock
CA_RTC018M	DSP RTC-018 2nd Unit For Monitor
CA_SMC*	Joerger Stepper Motor Controller SMC-L or SMC-24
CA_TS201	DSP Technology TS-201 Dual Timer/Scaler

Motor Parameters

Motor parameter assignment consists of key words of the form MOT000, MOT001, ..., followed by 11 values. The MOT key words must be numbered consecutively starting at zero. The values are:

- 1 Controller type (E500, SMC, OMS, ...)
- 2 Steps per unit (degrees, mm, ...) (sign changes direction of motion)
- 3 Sign be twe e nse ra rd d à lu n ist(+l o r- l)
- 4 Steady state rate (Hz) (must be positive)
- 5 Base rate (Hz) (must be positive) (also used as backlash rate)
- 6 Steps for backlash (sign changes direction of motion)
- 7 Acceleration time (msec)
- 8 Not used
- 9 Motor flags in hexadecimal (protection, units, etc.)
- 10 Motor mnemonic (th, phi, sl1, ...)
- 11 Motor name (Theta, Phi, Slit 1, unused, ...)

An example is:

```
# Motor  cntrl  steps  sign  slew  base  backl  accel  nada  flags  mne  name
MOT000 = E500 -2000   1 2000   200    50   125    0 0x003 tth  Two Theta
```

Valid controller types currently include:

18011	Oriel Encoder Mike Controller 18011
18092	Oriel Encoder Mike 18092
ANORAD	Anorad I-Series Controller
ANORAD_E	As above, but with encoder
CM3000	Compumotor 3000
CM4000	Compumotor 4000
CMSX	Compumoter SX
CMSX_E	As above, but with encoder
DAC_B12	PC DAC 12-Bit D/A (binary output)
DAC_B16	PC DAC 16-Bit D/A (binary output)
DAC_T12	PC DAC 12-Bit D/A (two's complement)
DAC_T16	PC DAC 16-Bit D/A (two's complement)
E250	DSP E250 12-Bit D/A as Motor Controller
E500	DSP Technology E500A
E500_M	As above, but with multiplexor
ECB_M	RISO ECB Motors
EPICS_M1	EPICS using spec's <i>config</i> motor parameters
EPICS_M2	EPICS using EPICS' database motor parameters
ES_OMS	ESRF using VME OMS
ES_VPAP	ESRF using VME Vpap
GALIL	Galil DMC-1000 PC Board
HLV544	Highland Technology VME V544
HUB9000	Huber 9000
IP28	Micro-Controle IP28
ITL09	Micro-Controle ITL09
ITL09_E	As above, but with encoder
IXE	Phytron IXE
KS3112	Kinetic Systems 3112 12-Bit D/A
KS3116	Kinetic Systems 3116 16-Bit D/A
KS3195	Kinetic SystemsS 3195 16-Bit D/A
MAXE	ESRF VME MAXE Motor Controller
MAXE_E	As above, but with encoder
MAXE_S	As above, but with servo
MAXE_DC	ESRF VME DC Motor Controller
MC4	Klinger MC4
MCB	Advanced Control Systems MCB
MCU	Advanced Control Systems MCU
MCU_E	As above, but with encoder
MCU_H	As above, with with Heidenhain encoder

MCU_O	As above, but with old PROMs
MM2000	Newport MM2000/3000
MM2000_E	As above, but with encoder
MM2500	Newport MM2500
MM2500_E	As above, but with encoder
MM4000	Newport MM4000/4005
MM4000_E	As above, but with encoder
MMC32	NSLS homemade
MURR	Missouri Research Reactor Motors
MURR_E	As above, but with encoder
NF8732	New Focus Picomotor 8732
NONE	Pseudo controller
NSK	NSK Motor Controller
NT2400	Laboratory Equipment Corporation Model
OMS	Oregon Micro Systems PCX/34/38/39/48
OMS_E	As above, but with encoder
PI	PI DC Motor Controller
PM500	Newport PM500
RIGAKU	Rigaku RINT-2000 Motor Controller
SCICE_A	SCICE Actuator Device
SIX19	Micro-Controle SIX19
SMC	Joerger Single Motor Controller
SPI8	Advanced Control System SPI-8
TSUJI	Tsuji PM16C-02N
XIAHSC	XIA HSC-1
XRGCI_M	Inel XRGCI motor controller/timer

Field 2, the steps per unit, may be non-integral, and the units can be in degrees, millimeters or whatever. The rest of the numeric fields must be integral. The motor names should be kept to nine characters or less, as the standard macros truncate them to fit a nine-character field when printing them out.

Field 8 is reserved.

Field 9, the flags field, contains several kinds of information. The lowest order two bits are used to enable particular operations on the selected motor. If bit 0 is set, the user can move the motor. If bit 1 is set, the user can change the software limits of the motor. Bits 2 and 3 are used by the *edconf* program to prevent users from changing certain configuration information. Bits 8 through 12 are used with the shared *config* file feature described below.

Optional motor parameters appear on lines following the MOT keywords and are of the form

```
MOTPAR:dcgain = 1500
```

Each MOTPAR refers to the immediately preceding motor. Possible parameters are:

```
encoder_step_size
step_mode
slop
home_slew_rate
home_base_rate
home_acceleration
dc_dead_band
dc_settle_time
dc_gain
dc_dynamic_gain
dc_damping_constant
dc_integration_constant
dc_integration_limit
dc_following_error
dc_sampling_interval
deceleration
read_mode
torque
misc_par_1
misc_par_2
misc_par_3
misc_par_4
misc_par_5
misc_par_6
```

Most parameters are not used by most motor controllers.

Linked Configurations

An installation such as a synchrotron beamline uses many motors with most associated with beamline control. Spectrometers used for particular experiments have motors that aren't used in other experiments. To avoid having to merge the motor configurations and settings from one set of files to another when the spectrometer is changed, you can set things up so that a single version of the *config* and *settings* files will describe a number of different spectrometers. Here is how to set up the files:

- (1) If you already have several geometry configurations installed, you should make backup copies of the *config* and *settings* files from the current geometries.
- (2) If you already have several geometry configurations installed, remove the *config* and *settings* files from all but one of the geometry directories. Save the

config file that has the most motors, as you will have to add motors from the other geometries to the remaining *config* file.

- (3) Set up hard links in all the geometry directories so that the *config* and *settings* in all the geometry directories refer to the same file. For example, if the files already exist in the *fourc* directory, use the commands

```
ln fourc/config surf/config
ln fourc/settings surf/settings
```

to create hard links in the *surf* directory. Don't use symbolic links.

- (4) Edit the *config* file by hand to add new control lines that assign numbers to the different geometries. These control lines must be before the lines that assign motor information. The format of the geometry control lines is as follows:

```
GEO0 = common
GEO1 = fourc
GEO2 = surf
GEO3 = fivec
etc.
```

The parameter `GEO0` always refers to the motors that are common to all the geometries. Subsequent lines assign consecutive numbers to the other geometries.

- (5) Now run *edconf*. The motor screen will have a new field that lets you assign a spectrometer geometry to each motor or to make the motor in common with all the spectrometers. You can do the same for each scaler on the scaler screen (as of release 4.03.12).

The hard links must be maintained for the shared *config* and *settings* file scheme to work. You can safely use *vi* and *cp* to manipulate the files. However, using *mv* will destroy the links. Also, the editor *ned* will destroy the links.

When running *edconf* with a geometry directory as an argument or when invoking the *config* macro from *spec*, use the `G` command to toggle between displaying all the motors and scalers in the *config* file and just those motors and scalers used by the given geometry.

Security Issues

At some installations, you may wish to prevent ordinary users from accessing selected motors. `spec` offers several levels of security. The security works with the UNIX file ownership and protection mechanisms, so it is important that user and group ownership of the configuration files and file write permission be properly set.

To restrict configuration modification to a single user or group, you must set the write permission of the diffractometer's associated *config* file accordingly. Type `chmod 644 /usr/lib/spec.d/fourc/config` to allow only the owner of the file to modify it. Setting the mode to 664 allows users in the owner's group to also modify the file.

There are three basic levels of security for each motor. The first level is the most restrictive, as it prevents the motor from being moved and prevents changes to the position being made in the *settings* file. The motor position can still be read from the motor controller, though, and the user angle can still be changed using the `chg_offset()` function (invoked by the `set` macro). If there is ever a conflict between the current position and the position in the *settings* file, such as might happen if the power was turned off to the motor controller, the controller registers are automatically adjusted to match the position in the *settings* file without moving the motor.

The second level of security allows a motor to be moved, but prevents the software limits from being changed. Not only is the `set_lim()` command restricted, but also the `chg_dial()` command, as a change in the dial position would effectively change the position of the limits.

The third level offers no security and allows any operation on a motor.

When these motor restrictions are set in the *config* file, the restrictions apply to everyone, even the owner of the *config* file. To move a restricted motor, you must first change the *config* file.

Extra Protection

At some `spec` installations, the administrators need to prevent users from accessing or modifying the configuration of certain motors. The *edconf* program supports a wizard mode that allows such protection. You will need to create a user in the system password file named *specwiz*. If you type `^w` while running *edconf* you will be prompted for the wizard's password. If you enter it properly, you will be able to select additional levels of configuration protection.

To prevent users from disabling the wizard protections by editing the *config* file by hand, you can use file protection features built in to UNIX. One possibility is to make

the *edconf* program set-user id *specwiz* or *specadm*, change the ownership of the *config* files to *specwiz* or *specadm*, and change the modes of the *config* files to *rw-r--r--*. Do that using commands (as super user) like

```
chown specadm edconf fourc/config surf/config ...
chmod u+s edconf
chmod 644 fourc/config surf/config ...
```


HARDWARE REFERENCE

Introduction

spec includes built-in support for a wide variety of motor controllers, counters, timers and other data acquisition devices, allowing great flexibility in a site's hardware configuration. Information on currently supported devices follow. Support for additional devices is continually being added.

Interface Controllers and General Input/Output

The interface screen of the configuration editor is selected using the I command. Before any interfaces have been configured the screen looks like:

Interface Configuration

```

CAMAC          DEVICE  ADDR  <>MODE          <>TYPE
  NO
GPIB           DEVICE  ADDR  <>MODE          <>TYPE
0   NO
VME            DEVICE  ADDR          <>TYPE
  NO
SERIAL         DEVICE <>TYPE  <>BAUD          <>MODE
0   NO
1   NO
2   NO
3   NO
IO PORT                ADDR  <>MODE  NUM
  NO
  NO
  NO
  NO
```

Type ? or H for help, ^C to quit

The following sections explain the choices for each type of interface. To select a particular interface, use the arrow keys to move to the cell containing the word NO and type y for yes and then <return>. For CAMAC, GPIB and VME devices, move the cursor to the last column and use the < or > keys to select the correct device and then enter <return>.

CAMAC Controllers

`spec` supports only one CAMAC controller at a time. The following CAMAC crate controllers are available:

<i>config</i> file	Description
PC_DSP6001	DSP 6001 with PC004 (no driver)
GP_CC488	DSP CC-488 GPIB Crate Controller
CDEV	DSP-6001/DCC-11/KS-3912 Boards
CA_JOR73A	Jorway 73A SCSI to CAMAC
CA_KS3929_HP	KS 3929 SCSI to CAMAC on HP
CA_KS3929	KS 3929 SCSI to CAMAC on Sun
GP_KS3988	KS 3988 GPIB Crate Controller
PC_KSC2926	KSC 2926 with 3922 (no driver)
CA_KSC	Kinetic Systems CAMAC driver

To select one, move the cursor to the NO box under CAMAC and type `y` and `<return>`. Then move the cursor to the rightmost column and type `<` or `>` until the appropriate controller appears and then enter `<return>`. Finally, select appropriate parameters from the other columns.

CAMAC Controllers That Use spec Drivers

config file:

```
CDEV = device_name INTR|POLL
```

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES	/dev/ca00		INTR	DSP-6001/DCC-11/KS-3912 Boards

`spec` drivers are available for the DSP 6001/6002 with PC004 for PC platforms, the Kinetic Systems 3922 with 2926 for PC platforms, or the Kinetic Systems 3912 CAMAC controller for BSD and Ultrix platforms. The appropriate CSS driver must be installed in each case. The PC platform controllers may be used in a polled or interrupt-driven mode. In interrupt-driven mode, a CAMAC look-at-me (LAM) will generate a call to a `spec` interrupt service routine.

To use the DSP 6001/6002 controllers in interrupt-driven mode, the boards must be modified to give the module a software programmable interrupt-enable capability. The modifications involve cutting four traces and soldering four jumper wires on one of the boards in the 6001/2 module. Contact CSS to obtain the explicit instructions. No modifications are required to operate the controller in polled mode.

DSP 6001/6002 CAMAC With No Driver

config file:

```
PC_DSP6001 = base_address
```

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES		0x240		DSP 6001 with PC004 (no driver)

Use this entry to select the DSP 6001/6002 CAMAC controller on a PC, if you aren't using a driver. No interrupt is used in this configuration.

Kinetic Systems 3988 GPIB To CAMAC

config file:

```
GP_KS3988 = gpib_address INTR|POLL
```

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES		7	POLL	KS 3988 GPIB Crate Controller

Use this configuration for the Kinetic Systems 3988 GPIB-to-CAMAC controller. This controller may be used in a polled or interrupt-driven mode. In interrupt-driven mode, a CAMAC look-at-me (LAM) generates a GPIB service request (SRQ), which in turn, will call a *spec* interrupt service routine. Interrupt-driven mode is currently only available with National Instruments GPIB controllers, and only when not using the *cib.o* GPIB configuration (see below). If multiple versions of *spec* are sharing the controller on the same computer, the controller must be operated in polled mode.

DSP CC-488 GPIB To CAMAC

config file:

```
GP_CC488 = gpib_address INTR|POLL
```

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES		7	POLL	DSP CC-488 GPIB Crate Controller

This configuration is for the DSP CC-488 GPIB-to-CAMAC controller. This controller may be used in a polled or interrupt-driven mode. In interrupt-driven mode, a CAMAC look-at-me (LAM) generates a GPIB service request (SRQ), which in turn, will call a *spec* interrupt service routine. Interrupt-driven mode is currently only available with National Instruments GPIB controllers. If multiple versions of *spec*

are sharing the controller on the same computer, the controller must be operated in polled mode.

Jorway 73A SCSI To CAMAC

config file:

CA_JOR73A = *device_name*

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES	/dev/sga			Jorway 73A SCSI to CAMAC
YES	/dev/scsi/1			Jorway 73A SCSI to CAMAC

The Jorway 73A SCSI-to-CAMAC controller is supported on HP 700 Series and *linux* platforms. Note, the four-position “piano” switch should be left in the factory configuration, with positions 1, 2 and 3 off, and position 4 on.

Kinetic Systems 3929 SCSI To CAMAC

config file:

CA_KS3929 = *device_name*

edconf interfaces screen:

CAMAC	DEVICE	ADDR	<>MODE	<>TYPE
YES	/dev/sga			KS-3929 SCSI to CAMAC
YES	/dev/ksc0			KS-3929 SCSI to CAMAC
YES	/dev/scsi/1			KS-3929 SCSI to CAMAC

This Kinetic Systems SCSI-to-CAMAC controller is available on SunOS 4.x SBus platforms, where a CSS provided driver can be installed. On *linux* and HP 700 Series workstations, *spec* provides direct software support. In all cases, this controller only operates with *spec* in a polled mode. In addition the software interface available from Kinetic Systems is supported (see below).

Kinetic Systems CAMAC Software

config file:

```
CA_KSC = device_name
```

edconf interfaces screen:

```

CAMAC      DEVICE  ADDR  <>MODE                                <>TYPE
  YES    /dev/rcamac                                Kinetic Systems CAMAC Driver
```

Kinetic Systems sells software interfaces for some of their CAMAC controllers on some UNIX platforms. Presently, only the package for the 3929 SCSI-to-CAMAC controller on the HP700 platform has been used with *spec*. To use the Kinetic Systems software, the location of the their object module must be given in response to the "KSC 3929 SCSI-CAMAC file location" query when running the *Install* program. (Note, however, there is no reason to use the rather expensive KSC software for the 3929 on the HP platform as CSS provides bundled support.)

GPIB Controllers

spec works with a variety of GPIB controllers as described below. (In addition, the Kinetic Systems 3388 CAMAC-to-GPIB controller is available.) Up to four GPIB controllers can be configured simultaneously.

On platforms that implement the System V interprocess communications (IPC) semaphore and shared-memory system calls, more than one *spec* process can share a single GPIB controller. For those systems, each *spec* must have the *shared* version of the GPIB controller selected. In the *config* file, a *_L* is appended to the module key word to indicate the shared version. If multiple GPIB controllers are configured, the controller unit number of the shared controller must be the same in each version of *spec*.

spec uses the GPIB controllers at *board* (as opposed to *device*) level, which makes it unlikely that other programs can use a GPIB controller while *spec* is using it.

National Instruments GPIB with National Instruments Drivers

The National Instruments GPIB boards and drivers should be installed according to National Instruments instructions. `spec` communicates with the boards using only the device node `/dev/gpib0`.

(On System V PC Platforms, where the GPIB driver is linked into the kernel, and the kernel is patched using the `ibconf` program, you should run `ibconf` directly on the driver file, so that each time you rebuild a kernel, you won't need to rerun `ibconf`. Thus, you might run

```
ibconf /etc/conf/pack.d/ib2/Driver.o
```

after the driver has been installed.)

When you do configure `/dev/gpib0` with the `ibconf` program, set the controller primary address to 0, the secondary address to none, board-is-system-controller mode to yes and disable-auto-serial-polling mode to yes.

Other board configuration parameters are programmed by `spec` each time it is run, thus overriding any values you may set using the `ibconf` program. Those parameters are: timeout setting, EOS byte, terminate-read-on-EOS mode, type of compare on EOS, set-EOI-w/last-byte-of-write mode and UNIX signal. The special nodes that may be created for each individual device by the National Instruments installation program are not used at all.

National Instruments GPIB with `cib.o`

config file:

```
PC_GPIBPC4 = device_name
PC_GPIBPC4_L = device_name
```

edconf interfaces screen:

GPIB	DEVICE	ADDR	<>TYPE
YES	/dev/gpib0		Nat Inst with cib.o
YES	/dev/gpib0		Nat Inst with cib.o (shared)

National Instruments provides a C language interface to its driver on most platforms in a file called `cib.c`. For some older versions of the National Instruments driver, `spec` has built-in C code that can be used instead of the National Instruments C interface. When using the National Instruments `cib.o` file, it isn't possible for GPIB devices that generate service requests (SRQ) to generate interrupts. Such devices must be used in polled mode.

In addition to specifying this choice in the *config* file, the location of the `cib.o` file must be entered when `spec` is installed. (See Page 191 in the *Administrator's Guide*.)

Any of the National Instruments boards (except the DEC MicroVax board) can be used with the *cib.o* file.

National Instruments GPIB on linux

config file:

```
PC_GPIBPC = device_name
PC_GPIBPC_L = device_name
```

edconf interfaces screen:

```

GPIB          DEVICE  ADDR                                     <>TYPE
YES /dev/gpib0/master                                     National Instruments GPIB
YES /dev/gpib0/master                                     Nat Inst GPIB (shared)
```

On *linux* platforms, there is a freely available GPIB driver for National Instruments boards available by anonymous ftp from the site `koala.chemie.fu-berlin.de` in the directory `/pub/linux/LINUX-LAB/IEEE488`. That driver should be installed and configured according to the documentation in the driver package. Note, however, that the DMA option in the driver installation should *not* be selected, as the implementation of DMA in the driver is notoriously flakey. The only part of the driver package needed by *spec* is the *driver/gpib0.o* module. *spec* does not use any of the application library included with the driver package. The file `/etc/gpib.conf` associated with the application library, so also is not used by and does not influence *spec*. Also, *spec* communicates with the driver directly through the `/dev/gpib0/master` device node. No other device nodes are used.

National Instruments GPIB-ENET

config file:

```
PC_GPIBPC5 = hostname
PC_GPIBPC5_L = hostname
```

edconf interfaces screen:

```

GPIB          DEVICE  ADDR                                     <>TYPE
YES           gpib0                                     Nat Inst GPIB-ENET
YES           gpib0                                     Nat Inst GPIB-ENET (shared)
```

spec must be linked with the National Instruments *cib.o* module that comes with the GPIB ethernet device. On the *edconf* interfaces screen, enter the name of the board as indicated by the National Instruments *ibconf* utility.

National Instruments PCII GPIB on PC UNIX System V Platforms

config file:

```
PC_GPIBPC = device_name
PC_GPIBPC_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                     <>TYPE
YES       /dev/gpib0                                     Nat Inst GPIB PCII
YES       /dev/gpib0                                     Nat Inst GPIB PCII (shared)
```

The PCII board is an old model, but is still supported by `spec`. This selection uses `spec`'s built-in code. The *cib.o* file configuration described earlier may also be used.

National Instruments AT-GPIB on PC UNIX System V Platforms

config file:

```
PC_GPIBPC2 = device_name
PC_GPIBPC2_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                     <>TYPE
YES       /dev/gpib0                                     Nat Inst AT-GPIB
YES       /dev/gpib0                                     Nat Inst AT-GPIB (shared)
```

The AT-GPIB board is a current model. This selection uses `spec`'s built-in code. The *cib.o* file configuration described earlier may also be used.

National Instruments GPIB on SCO UNIX and IBM AIX Platforms

config file:

```
PC_GPIBPC = device_name
PC_GPIBPC_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                     <>TYPE
YES       /dev/gpib0                                     National Instruments GPIB
YES       /dev/gpib0                                     Nat Inst GPIB (shared)
```

`spec`'s built-in code should work with SCO UNIX, IBM PS/2 and IBM RS/6000 platforms with this configuration choice. If new versions of the National Instrument drivers don't work, switch to the *cib.o* file configuration, described previously.

National Instruments SB-GPIB Ver 1.3 on SunOS 4.x Platforms

config file:

```
PC_GPIBPC = device_name
PC_GPIBPC_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       /dev/gpib0                                Nat Inst SB-GPIB Ver 1.3
YES       /dev/gpib0                                Nat Inst SB-GPIB Ver 1.3 (shared)
```

This old version of the driver is supported with built-in code.

National Instruments GPIB 1024-1S on SunOS 4.x Platforms

config file:

```
PC_GPIBPC2 = device_name
PC_GPIBPC2_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       /dev/gpib0                                Nat Inst GPIB 1014-1S
YES       /dev/gpib0                                Nat Inst GPIB 1024-1S (shared)
```

An old version of the driver for this board is supported with built-in code.

National Instruments SB-GPIB Ver 2.1 on SunOS 4.x Platforms

config file:

```
PC_GPIBPC3 = device_name
PC_GPIBPC3_L = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       /dev/gpib0                                Nat Inst SB-GPIB Ver 2.1
YES       /dev/gpib0                                Nat Inst SB-GPIB Ver 2.1 (shared)
```

This old version of the driver for this board is supported with built-in code.

National Instruments GPIB on DEC MicroVax

config file:

```
PC_GPIB11 = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       /dev/ib                                DEC GPIB11V-2 For Q-Bus
```

National Instruments distributes the driver source code for this module and platform. The driver is unlikely to be updated, so *spec* should continue to work with this module and platform indefinitely. No *cib.o* file is available for this platform.

HP SICL GPIB On HP Platforms

config file:

```
PC_SICL_H = sicl_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       hpib                                HP SICL GPIB
```

This configuration choice supports HP's GPIB using HP's SICL interface library. When the *spec Install* script is run, the question regarding GPIB SICL must be answered with "yes", and a *libsicl.a* or *libsicl.sl* must be available on the system for the SICL GPIB board to be available. If using the HP E2050 LAN/HP-IB Gateway, the device name is of the form *lan[hostname]:interface* where *lan* is the symbolic name for the device set in the */usr/pil/etc/hwconfig.hw* file, *hostname* is the host name or IP address of the gateway and *interface* is the interface as set for the host-name parameter in the gateway on-board configuration.

IOtech SCSI To GPIB On HP Platforms

config file:

```
PC_SICL_H = device_name
```

edconf interfaces screen:

```

GPIB      DEVICE  ADDR                                <>TYPE
YES       /dev/IOtech1                       IOtech SCSI488/H SICL GPIB
```

This configuration uses IOtech's SCSI-to-GPIB interface on an HP 700 Series workstation with the IOtech SICL software. When the *spec Install* script is run, the flags required to load the GPIB SICL library must be entered. Note, CSS has available a

modified version of the standard *libsicl.a* that doesn't use the *ieee488* daemon program, but does allow multiple processes to access the GPIB controller. Contact CSS to obtain the modified version of the library.

IOtech SCSI To GPIB on Sun Platforms

config file:

```
PC_IOTECH = device_name
```

edconf interfaces screen:

GPIB	DEVICE	ADDR	<>TYPE
YES	<eee488/ieee		IOtech SCSI488/S GPIB

Use this configuration for the IOtech SCSI-to-GPIB on a SunOS 4.x platform.

The full name of the default device is */dev/ieee488/ieee*. (The *edconf* program only displays the last twelve characters of long device names unless the cursor is on the device name cell.)

Scientific Solutions IEEE-488 on PC Platforms

config file:

```
PC_TEC488 = base_address  
PC_TEC488_L = base_address
```

edconf interfaces screen:

GPIB	DEVICE	ADDR	<>TYPE
YES		0x300	Scientific Solutions IEEE-488
YES		0x300	Scien Solut IEEE-488 (shared)

Only the Scientific Solutions (Tecmar) GPIB board (old style) is currently supported by *spec* (not the GPIB-LM model). Very old models of the board do not work. The card is accessed completely through user-level I/O. No kernel driver is needed.

Kinetic Systems 3388 CAMAC-To-GPIB Module

config file:

```
CA_KS3388 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	KS3388		KS 3388 GPIB Interface

When using the Kinetic Systems 3388 GPIB-to-CAMAC module, you must set the talk/listen address switch inside the module to correspond to address zero.

VME Controllers

National Instruments VME with National Instruments Drivers

config file: *edconf*

```
PC_NIVME = /dev/null
```

interfaces screen:

VME	DEVICE	ADDR	<>TYPE
YES			National Instruments VME

spec supports the National Instruments VXI-SB2020 on SunOS 4.1 and the VXI-AT2023 on System V PC platforms. The NI drivers must be installed, and the location of the NI *cvxi.o* must be specified when *spec* is installed.

The National Instruments software needs several patches when used on SVR3 PC platforms. Contact CSS for details.

Bit 3 Model 403 ISA-VME

Bit 3 Model 616/617 PCI-VME

Bit 3 Model 487-1 with Model 933 Driver Software

Bit 3 Model 466-1/467-1 with Model 944 Driver Software

Serial (RS-232C) Ports

config file:

```
SDEV_# = device_name baud_rate tty_modes
```

edconf interfaces screen:

SERIAL		DEVICE	<>TYPE	<>BAUD	<>MODE
0	YES	/dev/ttya1	<>	9600	cooked igncr
1	YES	/dev/ttya2	<>	2400	raw
2	NO				
3	NO				

Serial ports for use with the user-level `ser_get()` and `ser_put()` built-in functions are also selected on the interfaces screen. The device name, baud rate and serial line modes are selected for up to four serial devices. The number of the device is used as the first argument to the `ser_put()` and `ser_get()` functions. Available tty modes are either *raw* or *cooked*, with *cooked* mode also having *noflow*, *igncr* (a no-op on non-System V systems) and *evenp* or *oddp* options.

The <>TYPE field allows serial devices connected through special software servers used at ESRF or with EPICS to be selected. For normal serial devices, the field should contain the characters <>. See the *Reference Manual* for a description of the `spec` functions that access the serial ports.

Generalized CAMAC I/O

config file:

```
CA_IO = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	IO	0	Generalized CAMAC I/O

CAMAC modules configured for Generalized CAMAC I/O can be accessed using `spec`'s `ca_get()` and `ca_put()` functions. (See page 95 in the *Reference Manual*.) The module address for those functions is the unit number on the CAMAC configuration screen. Note that even more generalized access is available using the `ca_fna()` function. (See page 95 in the *Reference Manual*.) Arbitrary commands can be sent to any CAMAC module, whether or not the module is listed in the *config* file.

PC Port Input/Output

config file:

```
PC_PORT_# = base_address number_of_ports read_write_flag
```

edconf interfaces screen:

IO	PORT	ADDR	<>MODE	NUM
	YES	0x300	Read	1
	YES	0x310	R/W	4
	NO			
	NO			

On ISA bus systems on 80x86-compatible systems and on HP 700 platforms with E/ISA bus support, the ports available for the built-in `port_get()`, `port_getw()`, `port_put()` and `port_putw()` functions are configured on the interfaces screen. The board's hexadecimal base address is given, along with the number of contiguous 8-bit ports (maximum of 16) that can be accessed. The ports can be configured for read-only access or for read-write access. Be careful not to select port addresses associated with standard PC hardware such as the video board or the hard disk! Also be sure to include enough 8-bit ports to handle 16-bit word access, if that is how you will be using the ports. On the HP platforms, a config file must also be set up in the */etc/eisa* directory using the HP *eisa_config* utility.

Motor Controllers

Advanced Control System MCB (GPIB and Serial)

config file:

```
RS_MCB = device_name baud_rate number_of_motors
GP_MCB = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		6		4	Advanced Control System MCB (GPIB)
YES	/dev/tty00	<>	9600	4	Advanced Control System MCB (Serial)

When used on the serial interface, the MCB appears to output characters using even parity in spite of the board's jumpers being set to generate no parity. To accommodate that hardware idiosyncrasy, *spec* opens the port using even parity. If future versions of the MCB generate parity properly according to the jumper settings, the jumpers should be set for even parity to accommodate *spec*.

Advanced Control System MCU-2 (Serial)

config file:

```
RS_MCU = device_name baud_rate number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/ttyS1	<>	9600	12	Advanced Control System MCU

On the *Motors* screen of the configuration editor, one can select among *MCU*, *MCU_E*, *MCU_H* and *MCU_O* for the controller type.

The *MCU_O* controller type indicates that the controllers have the old-style firmware (which doesn't implement the # start character for message sent to the controllers), so that *spec* will not test first for the new-style firmware, thus avoiding timeouts.

The *MCU_H* controller type is for a special version of the controller which includes a hardware tie-in for a Heidenhain encoder. For such controllers, the `motor_par("encoder_step_size")` parameter is relevant in order to set the ratio between the Heidenhain encoder readings and the *MCU-2* step size. The default value is 131072 / 360.

The *MCU_E* type indicates an encoder is used. However, currently, there is no difference in the software whether *MCU_E* or *MCU* is selected.

Please note, many users have had problems establishing initial communication between the MCU-2 controllers and their computers. The problem is almost always related to the cable. Standard RS-232C cables do not appear to work. A custom cable wired according to the diagram in the MCU-2 manual seems to be needed. Note also, that the connections for pins 2 and 3 may need to be swapped from what is shown in the manual. Be prepared to try the cable both ways before adding the final touches.

Command pass through for the MCU-2 controllers is available using the following:

`motor_par(motor, "send", cmd)` – Sends the string `cmd` to the MCU channel associated with `motor`. For example, set `cmd` to "J500" to set the jog rate for `motor` to 500 steps per second.

`motor_par(motor, "read", cmd)` – Sends the string `cmd` to the MCU channel associated with `motor`, as above, and returns a string containing the response.

Compumotor 3000 (GPIB and Serial)

config file:

```
RS_CM3000 = device_name baud_rate number_of_motors
GP_CM3000 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00	<>	9600	4	Compumotor 3000 (Serial)
YES		6		4	Compumotor 3000 (GPIB)

Use of this old motor controller is not recommended!

Compumotor 4000 (GPIB and Serial)

config file:

```
RS1_CM4000 = device_name baud_rate number_of_motors
RS2_CM4000 = device_name baud_rate number_of_motors
GP_CM4000 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00	<>	9600	4	Compumotor 4000 (Serial port 1)
YES	/dev/tty00	<>	9600	4	Compumotor 4000 (Serial port 2)
YES		3		4	Compumotor 4000 (GPIB)

The Compumotor 4000 can be used on either an RS-232 or GPIB interface. You must program the RS-232 baud rate or the GPIB address using the Compumotor front

panel controls. You should consult the Compumotor manual for details, but in brief, the procedure is as follows. You must first enter the *ACCESS* code (the factory default is *4000*). You then choose the *IMMED* function, and then the *DEFINE GPIB ADDR* statement to select the the GPIB address. Alternatively, choose the *IMMED* function, and the the *RS232 PORT1* or *RS232 PORT2* statement to select the port and configure the baud rate for the RS-232 interface.

There are many other configuration options with this controller. Other than the GPIB address and the baud rate, you should probably not change any of these others. You can reestablish the factory defaults using the *RESET* function from the main menu.

Output pins 46 or 47 on the programmable output connector can be used to gate a counter during powder-mode scans. While the powder-mode motor is moving during these scans, *SPEC* sets pin 46 high and pin 47 low. Use the one appropriate for your particular counter. To gate the Ortec 994 counter/timer, for example, pin 46 and an even-numbered pin (all are logic ground) are connected to the front panel *enable* BNC connector on the Ortec module.

Compumotor AX (Serial)

config file:

```
RS_CMAX = device_name baud_rate number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00	<>	9600	4	Compumotor AX Motor Controller (Serial)

Command Pass Through

Command pass through is available using the following functions. Command pass through should be used with caution to avoid interfering with the built-in programming commands *SPEC* sends to the controllers.

`motor_par(motor, "send", cmd)` — Sends the string *cmd* to the motor channel associated with *motor*.

`motor_par(motor, "read", cmd)` — Sends the string *cmd* to the motor channel associated with *motor*, as above, and returns a string containing the response.

Compumotor SX (Serial)

DSP E250 12-Bit DAC as Motor Controller (CAMAC)

config file:

```
CA_E250 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	E250	0	DSP E250 12-Bit D/A as Motor Control

Some **spec** users use this DAC to control piezo-electric motion devices. Commanding such a device to move from **spec** results in an instantaneous change in the output voltage of the DAC.

DSP E500 Stepper Motor Controller (CAMAC)

config file:

```
CA_E500 = slot_number  
CA_E500M = slot_number  
CA_IOM1 = slot_number  
CA_IOM2 = slot_number  
CA_IOM3 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	E500	0	DSP E500 Stepper Motor Controller
2	IOM1		BiRa 2601 I/O For E500 Multiplexing
3	IOM2		F16,A0 I/O For E500 Multiplexing
4	IOM3		F16,A1 I/O For E500 Multiplexing

Selecting one of the three I/O module configurations above allows **spec** to multiplex an E500 motor channel. Currently only one channel of one E500 can be multiplexed. However, up to sixteen motors can be multiplexed on that channel. Contact CSS for additional information on the multiplexing circuitry required.

Huber SMC 9000 (GPIB)

config file:

```
GP_HUB9000 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		6		8	Huber 9000 Motor Controller (GPIB)

Note, *spec* does not support the RS-232C interface for this controller, as CSS was unable to make it work reliably with *spec*. Note also, a feature of this motor controller is that if a limit switch is hit, communication with the remote computer is shutdown so that one must manually move the motor off the limit and reset the controller in order to reestablish remote communication.

Inel XRGCI as Motor Controller (Serial)

config file:

```
RS_XRGCI_M = device_name baud_rate number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/ttyh3	<>	4800	3	Inel XRGCI as Motor Controller

Joerger SMC Stepper Motor Controllers (CAMAC)

config file:

```
CA_SMC = slot_number
```

```
CA_KS3640M = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	SMC	0	Joerger Single Motor Controller
2	KS3640M	0	KS 3640 Counter with Joerger SMC

Joerger controller models SMC-L, SMC-24 and SMC-LP can be used with *spec*. None of these models contain absolute motor position registers, so they are generally supplemented with additional CAMAC counter modules to keep track of motor position. *spec* supports the Kinetic Systems 3640 Up/Down Counter for this purpose. This module has four 16-bit counters that must be connected in series to form two 32-bit counters. One 3640 module is thus needed for each two SMC modules.

The SMC controllers can be used without a supplemental counter, but `spec` will not do as well in keeping track of absolute motor positions.

The Joerger controllers have no provision for a soft abort. If you type `^C` to abort moving, the controller will simply stop sending pulses to the motors. Inertia may cause the motors to continue to turn a bit, and absolute positions will be lost. The only work around for this problem is to keep motor velocities low.

The model SMC-LP allows programmable motor speed and acceleration. According to the module's documentation, if the internal frequency adjustments are at the factory settings, motor speed can be programmed between 50 and 2000 steps per second, while acceleration time can be varied from 20 to 2000 msec. For the other models, these parameters are set manually using potentiometers on the module —values entered in the `config` file are ignored.

Klinger MC-4 Stepping Motor Controller

config file:

```
RS_MC4 = device_name baud_rate number_of_motors
GP_MC4 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00	<>	4800	4	Klinger MC4 (Serial)
YES		3		4	Klinger MC4 (GPIB)

The Klinger MC-4 can be used on either an RS-232 or GPIB interface.

Back panel switches SW, SX, SY and SZ should each have locations 1 and 2 up and locations 3 and 4 down (`spec` does not currently support an origin switch, top zero or encoders). If you use the MC-4 with a GPIB interface, back panel switch S3 location 4 should be down, indicating a line feed terminator. Use switch S1 locations 1 to 5 select the GPIB address.

If you use the MC-4 with an RS-232 interface, set back panel switch S2 locations 1 to 3 to select the baud rate. Location 4 should be up, indicating software handshake. Switch S3 location 1 should be down and locations 2 and 3 up, to set 8 bit data words with 2 stop bits. Switch S3 location 4 should be down, indicating a line feed terminator.

You must connect certain pins together on the *general purpose I/O* connector on the back panel of the MC-4 to use the MC-4 with `spec`. Wire one end of a 1K to 10K ohm resistor to pin 35 (+5V). The other end should be connected to both pins 26 (B8) and 29 (E3). Pin 26 is an programmable output pin, and pin 29 is an external input. `spec` sets the state of pin 26 at the beginning of each move (which may involve all the

motors on the board) and clears the state of the pin at the end of the move. The status of pin 29 is read to determine when all activity, including backlash, is completed. The connection to pin 35 (+5V) is necessary to pull up the output.

Output pins 24 or 25 on the *general purpose I/O* connector can be used to gate a counter during powder-mode scans. While the powder-mode motor is moving during these scans, *spec* sets pin 24 high and pin 25 low. Use the one appropriate for your particular counter. To gate the Ortec 994 counter/timer, for example, pin 24 and pin 36 (logic ground) are connected to the front panel *enable* BNC connector on the Ortec module.

Kinetic Systems 3112 12-Bit DAC as Motor Controller (CAMAC)

config file:

```
CA_KS3112 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	KS3112	0	KS 3112 12-Bit D/A as Motor Control

Some *spec* users use this DAC to control piezo-electric motion devices. Commanding such a device to move from *spec* results in an instantaneous change in the output voltage of the DAC.

Micro-Controle IP28 (GPIB and Serial)

config file:

```
RS_IP28 = device_name baud_rate number_of_motors
GP_IP28 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00	<>	9600	4	Micro-Controle IP28 (Serial)
YES		6		4	Micro-Controle IP28 (GPIB)

Code for this motor controller was developed by a *spec* user. If new users plan on using this controller, please contact CSS first.

MicroControle SIX19 (Serial)

config file:

```
RS_SIX19 = device_name baud_rate number_of_motors
```

edconf devices screen:

```
MOTORS          DEVICE  ADDR  <>MODE  NUM                                <>TYPE
   YES                                     13                                MicroControle SIX19
```

Missouri University Research Reactor Motor Controller (GPIB)

config file:

```
HW_MURR = number_of_motors
```

edconf devices screen:

```
MOTORS          DEVICE  ADDR  <>MODE  NUM                                <>TYPE
   YES                                     13                                Missouri Research Reactor Motors
```

edconf motor screen:

```
Number: <>Controller      0:  MURR_E   1:  MURR_E   2:  MURR    3:  MURR
Unit/Channel                0/1        0/2        0/3        0/4
Name                        Two Theta   Theta      Chi        Phi
Mnemonic                     tth        th         chi        phi
```

The University of Missouri Research Reactor uses custom motor controllers. Each channel of the motor controller requires one GPIB address. The GPIB address is set using the channel number in the required unit/channel configuration on the motor screen of the configuration editor. The unit number is not relevant for these motors.

The controller type MURR_E or MURR is selected depending on whether or not the controller channel uses an encoder.

There are several unique parameters associated with each motor channel. The parameters assume default values when the controllers are powered up. Alternate values can be set in the *config* file that will be programmed by *spec*. The values are

Name	Parameter Name	Power-On Value
Modulo	Generic Parameter 1	360000
Grain	Generic Parameter 2	5
Direction	Generic Parameter 3	1
Drive Mode	Generic Parameter 4	1
Cut Point	Generic Parameter 5	0

For the *modulo* parameter, the value 360000 is appropriate for rotation stages. For

translation stages, the maximum value (999999?) would be appropriate. The *grain* parameter is a multiplier for the steps sent to the motor by *spec* and should be simply related to a gear reducer value. The *direction* parameter should be set to zero to reverse the direction of the motor so that its position agrees with the controller display. The *drive-mode* parameter changes the meaning of the output signals. For a value of one, the output signals are count-up/count-down. For a value of zero, the output signals are step/direction.

The *cut-point* parameter allows negative positions to be reported by *spec*, even though the controller only reports positive positions. Positions reported by *spec* will be between the cut point and the cut point plus the modulo parameter multiplied by the step size parameter. For a modulo of 360000 and a cut point of -180, positions will be between -180 and +180, for example.

New Focus Model 8732 Picomotor Controller (GPIB and Serial)

config file:

```
RS_NF8732 = device_name baud_rate number_of_motors
GP_NF8732 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00		9600	4	New Focus Picomotor 8732 (Serial)
YES		6		4	New Focus Picomotor 8732 (GPIB)

Each New Focus Picomotor controller unit has five card slots available. Each card has four connectors, and each connector can control up to three channels. The slot, connector and channel numbers need to be encoded in the configuration file by entering unit/channel information with the configuration editor. The unit number selects which Picomotor controller. The channel number selects the slot/connector/channel number of the controller encoded as XYZ, where the slot number is $1 \leq X \leq 5$, the connector is $1 \leq Y \leq 4$ and the channel is $1 \leq Z \leq 3$.

The only parameter from the *config* file used to program these controllers is the steady-state rate, which is sent to the controller as the pulse frequency.

The New Focus Picomotors are unlike most motor controllers in that there is no way to read the motor positions, and no way to know how far the motors move when commanded. *spec* attempts to guess how far the motor has moved if the positions are read while the motor is active, or if the move was aborted or stopped, based on the elapsed time of the move and the pulse frequency programmed into the controller. However, the positions reported by *spec* for the New Focus controllers should not be taken too seriously.

The motor controller can only move 65,535 steps at a time. For larger moves, `spec` will automatically send 65,535 steps at a time until the complete move is performed.

Only the up-to-three motors on the same connector can be moved simultaneously. `spec` prints an error message if you try to move a motor on a one connector while motors on another connector are active. (The moves could be automatically queued in software if users think such additional code development is warranted.)

Note also, the commands `chg_dial(mne, "lim+")` and `chg_dial(mne, "lim-")` can be used to start continuous moves, which can be stopped either with `^C` or the `stop()` function. Please beware also, there is no limit switch capability with this controller, so moves must be stopped by user intervention.

Command pass through is available using the following functions.

`motor_par(motor, "send", cmd)` – Sends the string `cmd` to the New Focus channel associated with `motor`.

`motor_par(motor, "read", cmd)` – Sends the string `cmd` to the New Focus channel associated with `motor`, as above, and returns a string containing the response.

The following special commands are also available:

`motor_par(motor, "always_address", mode)` – By default, `spec` currently prefixes each command sent to the controller with an eleven-character channel-select instruction. Communication can be made more efficient by turning off the "always_address" feature by calling this function with `mode` set to 0. With the mode off, the channel-select instruction is only sent when needed. The reason for always sending the channel-select information is to avoid `spec` losing track of the current channel in the event a user manually switches the controller to local mode.

`motor_par(motor, "was_local")` – This command will force `spec` to issue a channel-select instruction for the next command, allowing recovery after manually switching the controller to local mode when the "always_address" mode is off.

Newport (Klinger) Motion Master 2000/3000 (GPIB, Serial and PC Board)

config file:

```
RS_MM2000 = device_name baud_rate number_of_motors
RS2_MM2000 = device_name baud_rate number_of_motors
GP_MM2000 = gpib_address number_of_motors
PC_MM2000 = base_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00		9600	4	Newport MM2000/3000 (Serial)
YES	/dev/tty00		9600	4	Newport MM2000/3000 (Daisy Chain)
YES		6		4	Newport MM2000/3000 (GPIB)
YES		0x280		4	Newport MM2000 (AT bus)

Optional parameters:

```

MOTPAR:dc_proportional_gain
MOTPAR:dc_derivative_gain
MOTPAR:dc_integral_gain
MOTPAR:dc_integration_limit
MOTPAR:dc_sampling_interval
MOTPAR:dc_following_error
MOTPAR:home_base_rate
MOTPAR:home_slew_rate
MOTPAR:home_acceleration
MOTPAR:slop

```

The Newport (formerly Klinger) MM2000 and MM3000 motor controllers are supported by `spec` on both RS-232C and GPIB interfaces. The MM2000 is also supported on the ISA bus interface. On the serial interface, `spec` supports the daisy chaining available on the MM2000 and MM3000 controllers. All these controllers can be used both with DC motors (with encoders) and with the 1.5M-type stepper motors.

Newport Motion Master 4000/4005 (GPIB and Serial)

config file:

```

RS_MM4000 = device_name baud_rate number_of_motors
GP_MM4000 = gpib_address number_of_motors

```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00		9600	4	Newport MM4000/4005 (Serial)
YES		6		4	Newport MM4000/4005 (GPIB)

Optional parameters:

```

MOTPAR:dc_gain
MOTPAR:dc_damping_constant
MOTPAR:dc_integration_constant
MOTPAR:dc_following_error
MOTPAR:home_slew_rate
MOTPAR:home_acceleration
MOTPAR:slop

```

Before using the MM4000/4005 with `spec`, you need to set the communication parameters using the front panel buttons and display. The default communication

timeout of 0.5 seconds should be fine. Choose CR as the communication terminator for both GPIB and RS-232C interfaces. The SRQ feature of the GPIB interface is not used by spec, so the IEEE SRQ setting must be set to NO. For the RS-232C interface, the factory defaults for a parity setting of none, a word length of 8 bits and a stop-bits setting of 1 bit should be appropriate.

NSLS Brand MMC32 Controller (GPIB)

config file:

```
GP_MMC32 = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		5		32	NSLS Brand MMC32 Controller (GPIB)

Oregon Micro Systems (PC Board and VME)

config file:

```
PC_OMS = device_name number_of_motors INTR|POLL
PC_OMSP = base_address number_of_motors POLL
PC_OMSP58 = base_address memory_address number_of_motors POLL
PC_OMSV = VME_address number_of_motors IRQ_number|POLL
PC_OMSV58 = VME_address number_of_motors IRQ_number|POLL
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/oms00		INTR	4	Oregon Micro Systems PCX/38/39
YES		0x330	POLL	4	Oregon Micro Systems PCX/38/39 polled
YES	0xe000	0x300	POLL	4	Oregon Micro Systems PC58 polled
YES		0xfc00	IRQ5	8	Oregon Micro Systems VME8
YES		0xf000	POLL	8	Oregon Micro Systems VME58

spec currently supports PC board and VME module Oregon Micro Systems motor controllers.

For the PC versions of spec, the PCX, PC38 or PC39 models may be used in two, four, six or eight motor configurations. (Note, the newer PC34 and PC48 models should be used with the PCX/38/39 configuration.) spec can operate with a CSS-supplied driver (on certain platforms) or completely from user level using I/O port polling. The driver does require a dedicated PC interrupt, and at present, is limited to support of only one board. If I/O port polling is used, spec allows use of multiple boards.

The driver is contained in the file *oms.c* in the *drivers* subdirectory of the *spec* distribution. See the *README* file in the that directory for instructions on installing the driver into the UNIX kernel.

There is no driver for the PC58 board. It is only supported in polled mode.

For VME, the VME8, VME44 and VME58 models are supported. The VME8 operates eight motors, while the VME44 operates four motors with encoders, although software options for encoders are not currently implemented in *spec*. Multiple OMS VME motor controllers can be used simultaneously and are generally operated in polled mode. Interrupt-driven mode is currently only supported with the National Instruments MXI-VXI controllers.

OMS motor controllers can have from two to eight motors. *spec* numbers the motors the OMS manuals designate X, Y, Z, T, U, V, R and S as 0 through 7, in that order.

The first example above selects the PC board with the driver node */dev/oms00*. The driver may be used in either interrupt or polled mode. Interrupt mode means the *spec* program will be interrupted when motors complete their motions or hit a limit. In polled mode, the *wait()* function must be called repeatedly to check the status of the motor. Interrupt mode generally gives better performance, although in earlier versions of *spec*, software problems could be overcome by using polled mode. A PC interrupt is always required when the driver is used, even when polled mode is selected.

The second example selects the PC board with I/O port polling, with the board's base address at *0x330*, and with four motors on the board.

The third example selects the PC58 board at I/O port *0x300*. The PC58 also require 4,096 bytes of low memory. The example configuration with the address entered as *0xE000* in the DEVICE column selects a real memory address of *0xE0000*, as the value in the configuration is multiplied by 16.

The fourth example selects the VME8 and VME44 modules, with the board's A16 base address jumpered at *0xFC00* and with the VME interrupt request jumpered for *IRQ5*. Any of the VME IRQ vectors may be selected as can be polled mode. If more than one VME OMS controller is being used, all must be in polled mode or all must use interrupts. Different boards may use the same interrupt, though.

The last example selects the VME58 model. Note, this model requires 4096 bytes of A16 address space, so valid addresses have one hexadecimal digit followed by three zeroes.

On the *motor* screen (M) of the configuration editor, all of the OMS controllers use the symbol *OMS* or *OMS_E* in the controller field of the screen. The latter indicates the motor is being used with an encoder.

Special Commands

The following special commands are available through the `motor_par()` function. The two letter commands are direct implementations of commands described in the OMS manual. Refer to that manual for more information. Not all commands are available on all versions of the OMS controllers or on all firmware versions for a particular controller.

`motor_par(motor, "PA", mode)` – If `mode` is 1, the controller turns motor power on before each move and off after the move (assuming motor power is controlled by the auxiliary output pins). If `mode` is 0, motor power stays on.

`motor_par(motor, "SE", msec)` – Sets the settling time in milliseconds to be used before the power is reduced in PA mode.

`motor_par(motor, "AF")` – Turns auxiliary power off.

`motor_par(motor, "AN")` – Turns auxiliary power on.

`motor_par(motor, "BH", mask)` – Sets general purpose output pins high, according to which of bits 0-13 in `mask` are set.

`motor_par(motor, "BL", mask)` – Sets general purpose output pins low, according to which of bits 0-13 in `mask` are set.

`motor_par(motor, "BX")` – Returns the state of the general purpose input pins. A one in any binary position in the value returned indicates that the corresponding pin is low.

`motor_par(motor, "RB")` – Returns the direction of the general purpose I/O lines. Output bits return a one, while input bits return a zero.

Command Pass Through

Command pass through is available using the following functions. Command pass through should be used with caution to avoid interfering with the built-in programming commands `spec` sends to the OMS controllers.

`motor_par(motor, "send", cmd)` – Sends the string `cmd` to the OMS channel associated with `motor`. For example, set `cmd` to "LF" to disable hardware limits on the associated motor.

`motor_par(motor, "read", cmd)` – Sends the string `cmd` to the OMS channel associated with `motor`, as above, and returns a string containing the response. For example,


```
30.FOURC> print motor_par(tth, "read", "RP")
240000
```

```
31.FOURC>
```

results in the string "AX RP\n" being sent to the controller.

Asynchronous Surface Scanning

The following commands implement a special asynchronous, two-dimensional scanning mode available with newer versions of the OMS firmware. The scan is in the form of a repeating square wave, as illustrated below.



The scan starts at the point o , as specified with the commands below, and continues in the x and y directions in the range as specified with the commands below. At the end of the range, the motors are returned to the starting position and the scan is repeated.

Two motors must be configured with the mnemonics dx and dy in order for the asynchronous scanning mode to be available. When not in scanning mode, these motors may be moved normally.

Once started, the scanning will continue until explicitly stopped either with the `stop_scan` command (shown below), with a $\wedge C$ typed at the keyboard or with a `sync` command (which aborts the motors, but doesn't update `spec`'s positions). While scanning, the `wait()` function will not indicate these motors are moving. The `getangles` command will, however, return the current positions of these motors.

When the dx and dy motors are scanning, the remaining motors may be moved independently.

`motor_par(motor, "x_start", value)` – Sets the starting position for the dx motor.

`motor_par(motor, "x_range", value)` – Sets the extent of the motion in the x direction.

`motor_par(motor, "x_stepsize", value)` – Sets the size of each step in x . The number of steps is determined by dividing this number into the range for x .

`motor_par(motor, "y_start", value)` – Sets the starting position for the dy motor.

`motor_par(motor, "y_range", value)` – Sets the extent of the motion in the *y* direction.

`motor_par(motor, "start_scan")` – Starts the asynchronous scan.

`motor_par(motor, "stop_scan")` – Stops the asynchronous scan.

Oriel Encoder Mike Controller 18011 (Serial)

config file:

```
RS_18011 = device_name baud_rate number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/ttyh3	<>	4800	3	Oriel Encoder Mike Controller 18011

Oriel Encoder Mike Controller 18092 (Serial)

config file:

```
RS_18092 = device_name baud_rate number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/ttyh3	<>	4800	3	Oriel Encoder Mike Controller 18092

Phytron IXE α -C (GPIB and Serial)

config file:

```
RS_IXE = device_name baud_rate number_of_motors
```

```
GP_IXE = gpib_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00		9600	2	Phytron IXE (Serial)
YES		6		4	Phytron IXE (GPIB)

Finding Reference Positions

The Phytron's axis initialization command "0", which searches for the *negative* limit switch as described in the Phytron manual, is sent to a specific motor when the `chg_dial()` function is invoked with either the "home+" or "home-" arguments. The

free-axis-displacement commands "L+" and "L-" are sent with the "lim+" and "lim-" arguments, respectively.

There is no fixed hardware signal for a home switch on the Phytron controller, but there are twelve digital inputs available on the input connector. There is also a command which will perform a relative move at the base rate until one of the inputs goes high or low or the magnitude of the move is reached. There are four parameters in this command: the direction of the move, the magnitude of the move, the binary input number and the sense of the input switch. All four parameters are set by entering a string as *generic parameter 1* on the optional motor parameter screen of the configuration editor. (Get there by typing `m` twice from the standard motor parameter screen.) The string you enter will be sent to the motor when the "home" argument is used with the `chg_dial()` function. The string is the actual command sent to the Phytron, and is of the form

*S*magvE *nnD*

where *S* is a + or a - for the sign of the move, *mag* is the magnitude of the move (maximum of 65535 steps), *nn* is the input number (01 through 12) and *D* is the sense of the input where 0 means the motor stops if the input goes off and 1 means the motor stops if the input goes on. For example,

+200vE071

would command the motor to move no more than 200 steps in the plus direction, or until input 7 goes ON.

Use the `m` command twice from the motor screen of the configuration editor to reach the screen where you can enter *generic parameter 1*. Type an initial single quote to enter a string.

Special Commands

On faster computers, the Phytron apparently cannot keep up with commands sent by the computer at full speed. You can slow down the communication between `spec` and the Phytron controllers with the following commands:

`motor_par(motor, "rdelay" [, value])` – If *value* is given, sets the delay before reading a response from the Phytron to *value* seconds, otherwise returns the current value. The default value is 0.015.

`motor_par(motor, "wdelay" [, value])` – If *value* is given, sets the delay before sending a command to the Phytron to *value* seconds, otherwise returns the current value. The default value is 0.015.

Only one copy of the *rdelay* and *wdelay* parameters is kept for all the Phytron controllers. The motor mnemonic *motor* can be associated with any of the Phytron controllers. The values for the parameters are saved in the state file, so should only need to be reset after starting fresh. (See page 51 in the *Reference Manual*.)

Command pass through is available using the following functions.

`motor_par(motor, "send", cmd)` — Sends the string *cmd* to the Phytron channel associated with *motor*.

`motor_par(motor, "read", cmd)` — Sends the string *cmd* to the Phytron channel associated with *motor*, as above, and returns a string containing the response.

`motor_par(motor, "usend", cmd)` — Sends the string *cmd* to the Phytron controller associated with *motor*.

`motor_par(motor, "uread", cmd)` — Sends the string *cmd* to the Phytron controller associated with *motor*, as above, and returns a string containing the response.

For example,

```
31.FOURC> print motor_par(tth, "read", "P20R")
240000
```

```
32.FOURC>
```

results in the string "\002XP20R\003\r\n" being sent to the controller. Command pass through should be used with caution to avoid interfering with the built-in programming commands *spec* sends to the Phytron controllers.

The following command is also available to help with debugging:

`motor_par(motor, "dump")` — Displays the values of Phytron parameters P01 through P10 and P12 through P17 for the channel associated with *motor*.

PC DAC as Motor Controller

config file:

```
PC_DAC_B12 = base_address number_of_motors
PC_DAC_T12 = base_address number_of_motors
PC_DAC_B16 = base_address number_of_motors
PC_DAC_T16 = base_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		0x300		1	PC DAC 12-Bit D/A (binary output)
YES		0x310		1	PC DAC 12-Bit D/A (two's complement)
YES		0x320		1	PC DAC 16-Bit D/A (binary output)
YES		0x330		1	PC DAC 16-Bit D/A (two's complement)

PMC Corporation DCX-100 (Serial and PC Board)

config file:

```
RS_DCX = device_name baud_rate number_of_motors
PC_DCX = base_address number_of_motors
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty00		9600	4	PMC Corp DCX-100 (Serial)
YES		0xD000		4	PMC Corp DCX-100 (PC Board)

Optional parameters:

```
MOTPAR:dc_proportional_gain
MOTPAR:dc_derivative_gain
MOTPAR:dc_integral_gain
MOTPAR:dc_integration_limit
MOTPAR:dc_sampling_interval
MOTPAR:dc_following_error
MOTPAR:slop
```

XIA HSC (Huber Slit Controller)

config file:

```
RS_XIAHSC = device_name baud_rate number_of_slits
```

edconf devices screen:

MOTORS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/ttyS0	<>	9600	2	XIA HSC-1 (Serial)

sample *edconf* motor screen:

Number: <>Controller	0:	XIAHSC	1:	XIAHSC	2:	XIAHSC	3:	XIAHSC
Unit/Channel		0/0		0/1		0/2		0/3
Name		Slit1 A		Slit1 B		Slit1 Gap		Slit1 Cen
Mnemonic		s1a		s1b		s1g		s1c
Steps per degree/mm		400		400		400		400

The X-Ray Instrumentation Associates Huber Slit Controller Model HSC-1 is a specialized device only used to control Huber slits. The HSC-1 communicates through a serial port, and several HSC-1 modules can be daisy-chained together and run off a single serial port. On the device screen of the configuration editor, the NUM field is set to the number of HSC-1 modules on the serial port.

The HSC-1 module needs to be sent a calibration command before the HSC-1 motors can be moved with the normal move commands. The HSC-1 manual describes a manual calibration procedure. It is also possible to set the controller to the calibrated state by sending the "calibrate" command with `motor_par()` as described

below.

Each HSC-1 contains two motors that control the slit blades. Each blade can be moved independently. The HSC-1 also implements commands to move both blades simultaneously to change either the gap or the center-of-the-gap position. `spec` can be configured to control just the two blades, just the gap and the center position, or all four motions. When all four motions are configured, moving either blade changes the positions reported for gap and center, and moving either the gap or the center causes the positions reported for each blade to change.

Configuration for the HSC-1 requires the unit/channel field on the second line of the motor screen of the configuration editor to be filled in according to the following special format. The unit number corresponds to successive entries on the devices screen – each unit is associated with a different serial port. The channel number combines two values. Each HSC-1 module requires an arbitrary module number N (see below). This number is multiplied by 10 and added to the channel number that identifies the motion, as follows: For motor controller N , channel $N \times 10 + 0$ corresponds to blade A, channel $N \times 10 + 1$ corresponds to blade B, channel $N \times 10 + 2$ corresponds to the gap and channel $N \times 10 + 3$ corresponds to the center of the gap. Channel numbers ending in 4 through 9 are invalid.

Entering the Serial Number

The module number N (see above) is used only for internal bookkeeping and does not designate a particular HSC-1 module. Each module is identified by a thirteen-character unique serial number of the form `XIAHSC-B-0014`. The serial number needs to be entered as a string in the field *generic parameter 1* on the optional motor parameter screen of the configuration editor. Use the `m` command twice from the motor screen of the configuration editor to reach the screen where you can enter *generic parameter 1*. Type an initial single quote to enter a string. Note, the serial number can also be entered as `B-0014`, `0014` or `14` if such a string is sufficient to distinguish among modules. Also note, the serial number should be entered for just one of the motors associated with module N .

If the alias feature of the HSC-1 is used, and bit 6 of the *control word* (see below) is set for “use alias as ID”, `spec` requires there be no space characters in the alias. Of course, the alias, rather than the serial numbers would need to be entered as *generic parameter 1*. Setting and changing aliases requires establishing serial communication with the modules, which may be difficult for novice `spec` administrators, so CSS recommends simply using the serial numbers as the modules come from the factory.

Motor Parameters

The *steps per deg/mm* parameter should be set to 400 for the HSC-1 modules.

Normally, the positions for each blade become more positive as the blade is opened. However, if the *sign of user * dial* parameter is negative for either blade (or both), the motor position will become more negative as the blade opens. The sense of the center-of-the-gap motion can also be changed by changing the sign of the *sign of user * dial* parameter. The gap motion is always positive as the gap increases, though.

The backlash, speed and acceleration parameters in the *config* file are ignored.

CSS recommends using the calibration feature of the HSC modules to set the zero positions, rather than using the `set` macro to set the user *offset* parameter. That is, it is best to keep the user and dial positions the same. The `chg_dial()` function will, in fact, send the “immediate calibration” command to the controller, but only when setting the position to zero. Note, the gap should be physically at zero before using `set_dial`. Use of the standard `set_dial` macro should be followed by the `set` macro to set the user-dial offset back to zero.

The HSC-1 controller stores a number of parameters in nonvolatile memory. `spec` will read and display them with the command `motor_par(motor, "dump")`, where *motor* is the mnemonic for any of the motions on the particular HSC-1 module. The display format is as follows:

```
1   Outer motion limit (rw) = 4400 (11 mm)
2   Origin position (rw) = 400 (1 mm)
3   Motor A position (ro) = 900 (2.25 mm)
4   Motor B position (ro) = 900 (2.25 mm)
5   Motor step delay (rw) = 200 (roughly 0.272 mm/sec)
6   Gear backlash (rw) = 10 (0.25 mm)
7   Control word (rw) = 142 (0x8e)
8   Escape character (rw) = 33
9   Arbitration priority (rw) = 8
10  Motor A phase (ro) = 0
11  Motor B phase (ro) = 0
12  Calibration complete (ro) = 150
13  EEPROM signature (ro) = 23205
14  EEPROM version (ro) = 4
```

These parameters can be modified using a command such as

```
motor_par(motor, "send", "W 6 20")
```

which changes memory location 6 (gear backlash) to 20.

Special Commands

The `motor_par()` options implemented for the HSC-1 are as follows:

- `motor_par(motor, "calibrate")` – Sends the “immediate calibration” command to the unit. The effect is to set the current position of each blade in the controller to the *origin* parameter. It also sets *spec*’s positions for the gap, center and blades to zero. Thus the gap should physically be at zero before sending this command.
- `motor_par(motor, "origin", value)` – Sets the controller’s *origin* parameter (parameter 2 of the controller’s memory map) to *value*. The units of *value* are steps, where 400 steps corresponds to 1 mm. The origin parameter determines how far beyond the zero position each slit blade can be moved. Note, changing the origin parameter will change the setting of the current position of the blades. The blades should thus be both at zero before sending the "origin" command, and a "calibrate" command should be sent immediately afterwards. The factory default value for *origin* is 400.
- `motor_par(motor, "range" [, value])` – Sets the controller’s “outer limit” parameter (parameter 1 of the controller’s memory map) to *value*, where the units of *value* are steps. This parameter controls how far each blade can be moved. The factory default value for this parameter is 4400.
- `motor_par(motor, "step+")` or `motor_par(motor, "step-")` – Moves blade A or blade B one step in the specified direction. This command can be used to position the slits whether or not they have been calibrated.
- `motor_par(motor, "send", cmd)` – Sends the string *cmd* to the HSC-1 unit associated with *motor*. The module serial number will be included automatically.
- `motor_par(motor, "read", cmd)` – Sends the string *cmd* to the HSC-1 unit associated with *motor*, as above, and returns a string containing the response.
- `motor_par(motor, "usend", cmd)` – Sends the string *cmd* to the serial port connected to the HSC-1 unit associated with *motor*. The *cmd* must include the full HSC-1 command syntax.
- `motor_par(motor, "uread", cmd)` – Sends the string *cmd* to the serial port connected to the HSC-1 unit associated with *motor*, as above, and returns a string containing the response.

Timers and Counters

Am9513-based Counter/Timer PC Boards

config file:

```
PC_AM9513 = base_address number_of_counters
```

edconf devices screen:

SCALERS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		0x340		3	Am9513 Counter/Timer PC Boards

edconf scalers screen:

NUMBER	NAME	MNEMONIC	<>DEVICE	UNIT	CHAN	<>USE AS	SCALE	FACTOR
0	Seconds	sec	AM9513	0	0	timebase		1000
1	Monitor	mon	AM9513	0	1	monitor		1
2	Detector	det	AM9513	0	2	counter		1

The ComputerBoards CIO-CTR05/10/20 cards, the Keithley-Metrabyte Model CTM-05/10 cards, the Scientific Solutions Labmaster series cards, and similar models from other manufactures all use the Advanced Micro Devices Am9513 System Timing Controller chip. The chip contains five 16-bit counters that can be programmed in a wide range of configurations. *spec*'s programming uses two of the counters for a 32-bit detector counter, two for a 32-bit monitor counter and one for a 16-bit elapsed time counter. On boards with two or four chips, the additional chips are each programmed for two more 32-bit detector counters. You can program the chip through *spec* to count to either a time preset using the `tcount()` function or a monitor-count preset using `mcount()`.

You must connect the detector to the input connector pin labeled *source 3*. Counts received from the monitor go to the pin labeled *source 5*. In addition, you must wire the connector pin labeled *output 1* to the pins *gate 2*, *gate 4* and *gate 5*. (In the new Keithley-Metrabyte CTM-05A manual, the *source* pins are now labeled *ACLKIN*, the *output* pins are now labeled *ATIMEROUT*, and the *gate* pins are now labeled *AGATE*.)

If it is a two- or four-chip board, the additional detectors are connected to the *source 3* and *source 5* pins of the chips. In addition, the *output 1* from the first chip must be also connected to *gate 2* and *gate 4* of the additional chips.

The counter boards are accessed from user level and are polled to determine when count intervals have elapsed. Thus, interrupts should be disabled on the boards.

You will need to enter the base address of the counter chip in the *config* file. Note that for the Labmaster board, the base address of the counter chip is eight plus the base address of the board itself.

When counting to time, the resolution of the clock depends on the length of the count interval. The maximum count time is 71.5 minutes. The time base resolution (in seconds) is set according to the following table:

0.00001	for $t < 0.6$ sec
0.0001	for $t < 6$ sec
0.001	for $t < 60$ sec
0.01	for $t < 655.35$ sec (10.9 min)
0.0655	for $t < 71.5$ min

When counting to monitor counts, the 0.01 second time base is used, and the value returned for the time channel will be corrected to account for the rollovers that occur every 655.36 seconds.

Bi Ra 5302 64-Channel ADC (CAMAC)

config file:

```
CA_BR5302 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	BR5302	0	BiRa 5302 ADC as counters

scalars screen:

NUMBER	NAME	MNEMONIC	<>DEVICE	UNIT	CHAN	<>USE AS	SCALE FACTOR
0	Sensor 1	sen1	BR5302	0	0	counter	2

Up to 64 counters may be configured per ADC module. Each channel is 12-bits. Currently *spec* assumes a ± 10 Volt range on each channel and scales the readings to that range. The scale factor from the *config* file is used to program the gain on the corresponding channel. The values returned by *getcounts* in the *S[]* array are scaled by the gain value. Allowed values for the gain are from 1 to 1024 in powers of two. If an illegal value is entered, *spec* uses the next lower legal value.

DSP RTC-018 Real Time Clock (CAMAC)

config file:

CA_RTC018 = *slot_number*

edconf CAMAC screen:

Slot	Module	Unit	Description
1	RTC018		DSP RTC-018 Real Time Clock

The Standard Engineering DSP RTC018 Real-Time Clock is wired as follows:

- (1) If counting to time, the crystal oscillator output (2^{18} Hz) is connected to input A. If counting to the monitor, the monitor cable from one of the scaler inputs (usually channel 1) is connected to input A.
- (2) The *preset out* output is connected to *start*.
- (3) If using the DSP QS-450 scaler, connect the *busy* output to the *gate* input on the scaler module. If using the Kinetic Systems 3610 hex scaler, connect the *end* output to the *inhibit* input of the scaler module.
- (4) A 1 KHz signal should be fed into a scaler input (normally channel 0).
- (5) The detector signal should be fed into a third scaler input (normally channel 2).

DSP QS-450 4-Channel Counter (CAMAC)

config file:

CA_QS450 = *slot_number*

edconf CAMAC screen:

Slot	Module	Unit	Description
1	QS450	0	DSP QS-450 4-Channel Counter

DSP TS-201 Dual Timer/Scaler (CAMAC)

config file:

CA_TS201 = *slot_number*

edconf CAMAC screen:

Slot	Module	Unit	Description
1	TS201		DSP TS-201 Dual Timer/Scaler

Inel 715 Dual Scaler

config file:

```
RS_INEL = device_name baud_rate number_of_counters
```

Inel XRGCI as Timer/Counter

config file:

```
RS_XRGCI_T = device_name baud_rate number_of_counters
```

Joerger VSC16/8 Timer/Counter (VME)

config file:

```
PC_VSC16T = base_address number_of_counters INTR|POLL
```

edconf devices screen:

SCALERS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES		0x1000	POLL	8	Joerger VSC16/8 as Timer/Counter
YES		0x1100	POLL	8	Joerger VSC16/8 as Counters

Note, you need to add two zeros to the value of the settings of the six hexadecimal digits on the module's address switches to form the A32 address entered in spec's configuration editor.

Note also, the *ARM IN* connector needs to be jumpered to the *ARM OUT* on the Joerger front panel. If more than one Joerger module is used, the one module designated as Timer/Counter is the master, and the *ARM OUT* from that module needs to be connected to the *ARM IN* of all the modules.

CSS recommends users order the Joerger scaler with a 1 MHz crystal oscillator rather than the 10 MHz oscillator normally provided. The 10 MHz oscillator only allows preset counting times of a bit more than seven minutes before the counter overflows. The oscillator can also easily be changed in the field. Its only purpose is to provide the front panel time-base output. In either case, a value corresponding to the oscillator rate must be entered into the configuration editor on the scalers screen for the scale factor for the channel corresponding to seconds.

Kinetic Systems 3610 6-Channel 50 MHz Counter (CAMAC)

config file:

CA_KS3610 = *slot_number*

edconf CAMAC screen:

Slot	Module	Unit	Description
1	KS3610	0	KS 3610 6-Channel 50 MHz Counter

Kinetic Systems 3640 Used as Counter or Timer (CAMAC)

config file:

CA_KS3640T = *slot_number*

CA_KS3640C = *slot_number*

edconf CAMAC screen:

Slot	Module	Unit	Description
1	KS3640T		KS 3640 Counter used as Timer
2	KS3640C	0	KS 3640 Counter used as Counter

In order to use a 3640 as a timer to gate other 3640 modules, you need to make modifications. The modifications will leave one of the front panel inhibit inputs alone and convert the other to an inhibit output. One way to do this is to add an LM311 comparator IC to the circuit. The negative input (pin 2) of the 311 is connected to the LAM signal from pin 9 of IC 30 of the 3640 module. The positive input (pin 3) is held high at about 3.2V through a 1200 ohm over 2200 ohm voltage divider between +5V and ground. Pins 1 and 4 of the 311 are connected to ground. Pin 8 is connected to +5V. The output of the 311 (pin ?) is connected to the front panel LIMO connector, which must have the factory connection cut. In addition, a 2K Ohm resistor is connected between the output and +5V to pull up the output.

In operation, the inhibit inputs of all the 3640s are connected to the inhibit outputs of all the 3640s. Thus any module can be used as the gate. One module should be fed a fixed time base, say 1KHz or 10KHz which must come from some external source.

Kinetic Systems 3655 Timing Generator (CAMAC)

config file:

```
CA_KS3655 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	KS3655		KS 3655 8-Channel Timing Generator

Current users of `spec` have made the following modifications to the Kinetic Systems Model 3655 Timing Generators. This timing generator is customarily used with the Kinetic Systems Model 3610 Hex Scaler. All modifications are made on the component side of the board.

- (1) Bring the internal 1 KHz timing signal out through the channel 7 front-panel lemo connector. Do this by first unsoldering the end of the wire that connects the center pin of the channel 7 lemo connector from the feed through on the circuit board. Do not unsolder the wire from the connector, as it will be difficult to solder on a new wire. Instead splice a longer wire to the one already attached to the connector and solder the other end of that wire to pin 11 of chip BJ. (Pin 11 is the center pin on the front-panel side of the chip.)
- (2) Disconnect the internal inhibit signal from the CAMAC dataway, and bring it out through the channel 8 lemo connector. Do this by folding up and/or snipping pin 6 of the socketed 7407 chip in position BX, near the front panel. Next either fold up pin 3 of chip AU or cut the long trace that leads from that pin to the dataway Inhibit connector. Then solder a wire from pin 3 of chip AU to the bottom lead of resistor R45. The resistor is located near the top-right corner of the circuit board, and the bottom lead is the one nearest the letters R45.

Channel 7 is then connected to channel 0 of the scaler module, while Channel 8 is connected to the *inhibit* input. The signal from the source of monitor counts is connected, using a tee, to the *clock* input of the 3655 and the monitor scaler channel (usually channel 1; the detector is usually channel 2).

Ortec 974/994/995/997 NIM Timers and Counters

config file:

```
RS_OR9XT = device_name baud_rate number_of_counters
GP_OR9XT = gpib_address number_of_counters
RS_OR9XC = device_name baud_rate number_of_counters
GP_OR9XC = gpib_address number_of_counters
RS_OR9XB = device_name baud_rate number_of_counters
GP_OR9XB = gpib_address number_of_counters
```

edconf devices screen:

SCALERS	DEVICE	ADDR	<>MODE	NUM	<>TYPE
YES	/dev/tty1		9600	4	Ortec 974/994 Counter/Timer (Serial)
YES		3		4	Ortec 974/994 Counter/Timer (GPIO)
YES	/dev/tty2		9600	2	Ortec 974/994/995/997 Counter (Serial)
YES		3		2	Ortec 974/994/995/997 Counter (GPIO)
YES	/dev/tty3		9600	3	Ortec 994 Blind Timer/Counter (Serial)
YES		3		3	Ortec 994 Blind Timer/Counter (GPIO)

edconf scalers screen:

NUMBER	NAME	MNEMONIC	<>DEVICE	UNIT	CHAN	<>USE AS	SCALE	FACTOR
0	Seconds	sec	OR9XX	0	0	timebase		1000
1	Monitor	mon	OR9XX	0	1	monitor		1
2	Detector	det	OR9XX	0	2	counter		1

spec supports the Ortec 974, 994, 995 and 997 counter and counter-timer NIM modules over both GPIO and RS-232 interfaces. When running the configuration editor, select from the above descriptions on the device configuration screen to specify which Ortec modules you are using and how you are using them.

Only one module can be selected as a counter/timer. The 974 module can be assigned a maximum of four channels. The 994 should be assigned two channels normally and three channels when used as a blind timer. The 995 has two channels and the 997 has one.

On the scaler configuration screen, choose OR9XX as the controller for all channels associated with an Ortec module. The unit numbers selected for each channel correspond to the order the Ortec modules appear on the device configuration screen. When using the 994 as a blind timer, you must select channel number 2 for the timebase.

Using the 974

The 974 is a four-channel counter/timer having a minimum 0.1 second time base. You should connect the monitor counts through a tee to the EXT IN connector on the back of module and to the COUNTER INPUT 2 connector on the front of the module. Use the COUNTER INPUT 3 and 4 connectors for one or two detector input

channels. Also, make sure that the internal dip switch S-1 has position 6 set to one-cycle.

Using the 994 as a Normal Timer

The 994 is a two-channel counter/timer with a minimum 0.01 second time base. In order to obtain accurate elapsed time readings, one counter channel is used to count time and the other is used to count monitor counts. An additional counter, such as the 995 or 997 is normally used to accumulate detector counts and is gated by the 994. The monitor count source should be connected to both the IN A and IN B front panel connectors of the 994 using a tee. The internal jumpers W3 and W4 must both be set to the TIME position. Jumper W1 must be set to the NORMAL position. Also, make sure the internal dip switch S-1 has position 6 set to one-cycle and position 7 set to COUNTER/TIMER. Finally, make sure the front panel DWELL switch is turned all the way off.

Using the 994 as a Blind Timer

In the blind timer mode, the 994 has the monitor counts connected to IN A and detector counts connected to IN B. The internal jumpers W3 and W4 must both be set to the COUNTS position. Jumper W1 must be set to the NORMAL position. Also, make sure the internal dip switch S-1 has position 6 set to one-cycle and position 7 set to COUNTER/TIMER. Finally, make sure the front panel DWELL switch is turned all the way off.

When operated as a blind timer, `spec` cannot read back the elapsed time from the module. Instead, when counting to monitor counts, when counting in powder mode, when reading the counters during updated counting and when counting is aborted with a `^C`, the elapsed count time is estimated from the software clock.

Gating

An external enable signal from certain motor controllers may be fed into the rear-panel gate BNC input on the 974 or the front panel enable BNC input on the 994 for precise counter gating in powder-mode scans.

If using a second Ortec module as a counter, you must connect the INTERVAL BNC connector (rear panel on 974, front panel on 994) to the master GATE on the 974 rear panel or to the ENABLE or individual GATE inputs on the 994, 995 or 997 modules.

Setting Operational Parameters

The `counter_par()` function can be used to set various parameters associated with the Ortec module code in `spec`. The first argument to `counter_par()` is a channel number, although all the commands affect all channels of the associated module, or all of the Ortec modules, if appropriate.

`counter_par(counter, "alarm", mode)` – If `mode` is zero, turns off the more efficient **ALARM** mode of operation of the timer, and turns on a slower polled mode. If `mode` is one, **ALARM** mode is turned on. The default operation is for **ALARM** mode to be turned on, and there is generally no reason to turn it off.

`counter_par(counter, "alarm")` – Returns one if **ALARM** mode is on. Otherwise returns zero.

`counter_par(counter, "display", channel)` – Sets the counter channel that will be displayed on the associated module. For the 974 modules, valid values for `channel` are 1 to 4. For the the 994 and 995 modules, valid values for `channel` are 0 and 1.

`counter_par(counter, "display")` – Returns the channel number currently being displayed.

`counter_par(counter, "local", mode)` – If `mode` is nonzero, will force the associated module to go into local mode to allow front panel operation. In addition, the module will be placed in local mode after each count interval. If `mode` is zero, the module will be set to remote mode at the start of the next count interval, and will not be set back to local mode after counting. When `spec` starts up, not switching to local mode is the default behavior to minimize overhead.

`counter_par(counter, "local")` – Returns zero if the associated module is to be kept in remote mode. Otherwise, returns one.

Software Timer

config file:

```
SW_SFTWARE = 1
```

If no hardware timer is available, the system clock can be used as a timer. Only counting to time is allowed, as counting to monitor makes no sense. The nominal resolution depends on the underlying operating system, although 10 msec is typical. The accuracy, though, is certainly less than that.

Multichannel Data Acquisition Devices

DSP 2190 MCS Averager

config file:

```
CA_DSP2190 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	DSP2190		DSP 2190 MCS Averager

The DSP Technology 2190 Multichannel Scaling Averager consists of a pair of CAMAC modules: the 2090S Multichannel Scaling module and the 4101 Averaging Memory module. These two modules must occupy consecutive slots in the CAMAC crate, with the 2090S in the lower-numbered slot. There is no entry for the 4101 module in the *config* file.

Functions

The `mca_par()` function controls the module's behavior as follows:

`mca_par("run")` – programs the MCS for the number of bins and sweeps set with the functions described below, then enables any other counters and starts the averaging process. When the programmed number of sweeps is completed, the MCS will generate a CAMAC LAM, which will cause the other counters to be disabled. Use the `wait()` function to determine when the programmed number of sweeps are complete.

`mca_par("halt")` – halts the MCS and disables the other counters.

`mca_par("bins")` – returns the number of bins in each sweep. (Referred to in the module documentation as “record length”).

`mca_par("bins", value)` – sets the number of bins in each sweep to *value*. The number of bins can range from 8 to 32,767.

`mca_par("sweeps")` – returns the number of sweeps to be summed.

`mca_par("sweeps", value)` – sets the number of sweeps to be summed in the next scan to *value*. The number of sweeps can range from 1 to 65,536.

`mca_par("sweeps_comp")` – returns the number of sweeps completed in the previous scan. An error message is printed if this function is called while a scan is in progress.

`mca_par("first_ch")` – returns the first channel to be read out using `mca_get()`.

`mca_par("first_ch", value)` – sets the first channel to be read out using `mca_get()` to *value*.

`mca_par("npts")` – returns the number of channels to be read out using `mca_get()`.

`mca_par("npts", value)` – sets the number of channels to be read out using `mca_get()` to *value*.

Note that the `mca_get()` function cannot be used while the MCS module is taking data.

Note that the 4101 doesn't actually average the sweeps, but only accumulates sums in each channel. To obtain an average, you must divide the data in each channel by the number of sweeps. The averaging scan will halt before the programmed number of sweeps is completed if any of the channels overflow.

The module expects an external trigger and the external trigger is required to begin each sweep.

LeCroy 2301 interface for qVT MCA

config file:

```
CA_LC2301 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	LC2301		LeCroy 2301 interface for qVT MCA

Functions

The `mca_par()` function controls the MCA module's behavior as follows:

`mca_par("clear")` – clears the MCA. `spec` inserts a 1.5 second delay to give the device time to clear.

`mca_par("run")` – starts the MCA.

`mca_par("halt")` – stops the MCA.

`mca_par("first_ch")` – returns the first channel to be read out.

`mca_par("first_ch", value)` – sets the first channel to be read out to *value*.

`mca_par("npts")` – returns the number of channels to be read out.

`mca_par("npts", value)` – sets the number of channels to be read out to *value*. The maximum number of channels is 1,024.

`mca_par("delay")` – returns the delay time in seconds that `spec` sleeps after the MCA is cleared.

`mca_par("delay", value)` – sets the time for `spec` to delay after sending the clear command. The hardware does require some delay. Some users have reported 1.5 seconds are needed, others report 0.1 seconds is adequate. The default value is 0.1 seconds.

LeCroy 3512 Spectroscopy ADC

config file:

```
CA_LC3512 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	LC3512		LeCroy 2301 interface for qVT MCA

LeCroy 3588 Fast Histogram Memory

config file:

```
CA_LC3588 = slot_number
```

edconf CAMAC screen:

Slot	Module	Unit	Description
1	LC3588		LeCroy 3588 Fast Histogram Memory

Keithley 2001 Multimeter (GPIB)

config file:

```
GP_K2001
```


Oxford/Tennelec/Nucleus PCA Mutliport, PCA II, PCA-3

config file:

```
GP_PCA_M = gpib_address
PC_PCA3  = base_address
PC_PCAII = device_name base_address INTR|POLL
```

edconf devices screen:

MCAs	DEVICE	ADDR	<>MODE	<>TYPE
YES		7		The Nucleus PCA Multport (GPIB)
YES		0x210		The Nucleus PCA-3 MCA Board
YES		0x1e0	POLL	The Nucleus PCA II MCA Board
YES	/dev/pca	0x1e0	INTR	The Nucleus PCA II MCA Board

The PCA II MCA can be used in either a user-level I/O mode or in an interrupt-driven mode with the CSS provided driver. The interrupt-driven mode allows dead-time corrections and more accurate counting times. See the *drivers/README* file in the spec distribution for information on installing the driver.

If using the interrupt-driven mode, note the following: Apparently, the PCA II doesn't trigger an interrupt on some PC mother boards. This problem can be fixed by changing the value of the resistor labeled R12 on the "PCA2 Memory Card" circuit diagram. This resistor is located near the lower left corner of the main board when viewed from the component side with the connector fingers pointing down and the input BNC to the right. R12 is about a centimeter down and to the left of the U26 IC. The circuit diagram indicates the resistor's value is 2K, however the boards seem to be shipped with a 1K resistor (brown-black-red stripes). Soldering a second 1K resistor alongside R12 and in parallel electrically will lower the resistance to 0.5K, which seems to work. (This modification was suggested by the manufacturer.)

Functions

The `mca_par()` function controls the board's behavior as follows:

`mca_par("clear")` — clears the channels of the current group.

`mca_par("run")` — programs the board with the current parameters and starts acquisition. Note that the `tcount()` and `mcount()` functions, as used in the various counting macros will also start PCA II acquisition.

`mca_par("halt")` — stops acquisition. Note that the PCA II will also be halted when the `tcount()` and `mcount()` functions, as used in the various counting macros, complete their count intervals or are aborted.

`mca_par("group_size")` — returns the current group size.

`mca_par("group_size", size)` – sets the group size to *size*. Legal values are 256, 512, 1024, 2048, 4096 and 8192. Values above 1024 may not be legal if insufficient memory is installed on the board.

`mca_par("select_group")` – returns the currently active group. Groups are numbered starting at zero.

`mca_par("select_group", group)` – set the active group to *group*. The number of possible groups is given by the total number of channels on the board divided by the group size. If the *group* passed to the function is greater than the maximum number of groups (based on the current group size and total number of channels), the current group selected is *group* modulus the maximum number of groups.

`mca_par("pha")` – selects pulse-height analysis mode on the board.

`mca_par("gain")` – returns the current gain value used in pulse-height analysis mode.

`mca_par("gain", value)` – sets the pulse-height analysis gain to *value*. Legal values are 256, 512, 1024, 2048, 4096 and 8192.

`mca_par("offset")` – returns the current channel offset used in pulse-height analysis mode.

`mca_par("offset", value)` – sets the pulse-height analysis offset to *value*. Legal values are multiples of 256 from 0 to 7936.

`mca_par("mcs")` – selects multichannel scaling mode on the board.

`mca_par("dwell")` – returns the current multichannel scaling dwell time.

`mca_par("dwell", value)` – set the multichannel scaling dwell time. Allowed values are numbers between 1e-6 and 60 seconds with mantissa of 1, 2, 4 or 8. A value of -1 selects external dwell. If *value* isn't an allowed value, it is rounded to the nearest allowed value.

`mca_par("mode")` – returns two if the board is in PHA live-time mode, one if the board is in PHA real-time mode and zero if the board is in MCS mode.

`mca_par("readone", channel)` – returns the contents of channel number *channel*.

`mca_par("chan#")` – returns the contents of channel number #. The channel number is with respect to the current group.

`mca_par("chan#", value)` – sets channel # to *value*. The channel number is with respect to the current group.

The following `mca_par()` functions are only valid when the board is used with the interrupt-driven driver.

`mca_par("preset")` – in PHA mode, returns the current live-time or real-time preset value in seconds.

`mca_par("preset", value)` – in PHA mode, sets the current live-time or real-time preset value to *value* seconds.

`mca_par("passes")` – in MCS mode, returns the number of preset passes.

`mca_par("passes", value)` – in MCS mode, sets the number of passes to *value*.

`mca_par("live")` – in PHA mode, selects live-time counting.

`mca_par("real")` – in PHA mode, selects real-time counting.

`mca_par("dead")` – in PHA mode, returns the percent dead time, if accumulating in live-time mode.

`mca_par("elapsed_live")` – in PHA live-time mode, returns the elapsed live time in seconds.

`mca_par("elapsed_real")` – in PHA mode, returns the elapsed real time in seconds.

`mca_par("elapsed_passes")` – in MCS mode, returns the elapsed number of passes.

Silena CATO MCA (Serial)

config file:

RS_CATO

Nicomp TC-100 Autocorrelator (Serial)

config file:

RS_TC100 = *device_name* *baud_rate*

The Nicomp TC-100 Autocorrelator is selected in the *config* file in the MCA section of the devices screen.

Functions

The `mca_par()` function controls the correlator behavior as follows:

`mca_par("clock")` – returns the value of the current clock time parameter in microseconds.

`mca_par("clock", value)` – sets the clock time parameter. The units for *value* are microseconds. Valid clock times are of the form *X.XeY* where *X.X* ranges from

0.1 to 1.6 and Y ranges from 0 to 5. Values outside these bounds will be rounded to the closest allowed value. The new value takes effect on the next run command.

`mca_par("prescale")` – returns the value of the prescale factor.

`mca_par("prescale", value)` – sets the value of the prescale factor. Valid prescale values are from 1 to 99. The new value takes effect on the next run command.

`mca_par("dbase_mode")` – returns the state of the baseline mode. A return value of 1 means delayed baseline mode is in effect. A return value of 0 means delayed baseline mode is off.

`mca_par("dbase_mode", 1|0)` – sets the state of the baseline mode. A value of 1 turns on delayed-baseline mode. A value of 0 turns it off. The new mode takes effect on the next run command.

`mca_par("dbase")` – returns the value of the delayed baseline from the last data obtained using `mca_get()`.

`mca_par("cbase")` – returns the value of the calculated baseline from the last data obtained using `mca_get()`.

`mca_par("tcnts")` – returns the value of the total-counts monitor channel from the last data obtained using `mca_get()`.

`mca_par("pcnts")` – returns the value of the total-prescaled-counts monitor channel from the last data obtained using `mca_get()`.

`mca_par("rtime")` – returns the value of the run-time monitor channel from the last data obtained using `mca_get()` in seconds.

`mca_par("clear")` – clears the correlator.

`mca_par("run")` – sends the current clock-time, prescale and delayed-baseline parameters to the correlator and starts the correlator. The `tcount()` and `mcount()` functions also start the correlator.

`mca_par("halt")` – stops the correlator. The correlator is also halted when count intervals specified by `tcount()` or `mcount()` have elapsed, or when counting is aborted using a ^C .

`mca_par("plot")` – reads off the real-time data plot from the running correlator. The data obtained is a very low resolution version of the correlation function.

`mca_get(grp, el)` – reads the current data from the correlator, and stuffs the data into the data group `grp` element `el`.

REFERENCES

Journal articles describing the various supported X-ray diffractometers are available in the following references.

The four-circle diffractometer is discussed in

W. R. Busing and H. A. Levy, *Acta Cryst.* **22**, 457 (1967).

There is an error in Equation (48) of the above paper. The last line of the equation should be $\omega = \text{atan}(-R_{23}, R_{13})$.

Surface diffraction using a four-circle diffractometer is discussed in

S. G. J. Mochrie, *J. Appl. Cryst.* **21**, 1-3 (1988).

The z-axis diffractometer is described in

J. M. Bloch, *J. Appl. Cryst* **18**, 33 (1985).

The liquid-surface diffractometer supported by spec is described in

J. Als-Nielsen and P. S. Pershan, *Nucl. Instrum. Methods* **208**, 545 (1983).

A. H. Weiss, M. Deutsch, A. Braslau, B. M. Ocko, and P. S. Pershan, *Rev. Sci. Instrum.* **57**(10), 2554 (1986).

Angle calculations and operating modes for a six-circle diffractometer are presented in

M. Lohmeier and E. Vlieg, *J. Appl. Cryst.* **26**, 706 (1993).

A description of the CAMAC driver can be found in

G. Swislow, A. Braslau and S. G. J. Mochrie, *Interrupt-Driven CAMAC Software for UNIX-Based Computers*, AT&T Bell TM# 11115-870817-39.

INDEX

special characters

- !! or !-1 to recall previous command, 16, 51
- # to begin a comment line, 26, 29, 46, 204
- * metacharacter, 26, 50, 83, 89, 91
- ?
 - as metacharacter, 26, 50, 83, 89, 91
 - to list *edconf* commands, 202
- [and] to form arrays, 47
- \
 - to continue a line, 50
 - to introduce special characters, 49
- ^ to substitute in most recent command, 16, 51
- _check0 macro, 124, 154
- _chk_lim macro, 154
- _cleanup2 macro, 159
- _cleanup3 macro, 159
- _do macro, 130
- _loop macro, 160
- _mo_loop macro, 138
- _move macro, 154, 160
- _pmove macro, 154
- _scan_on macro, 160
- _scanabort macro, 159
- { and }
 - to delimit block, 18
 - to group lines as a parse tree, 45
- 80386 systems
 - fractional-second sleeps on, 76

A

- A[]
 - as built-in variable, 56
 - motor positions in, 19, 29–29, 101, 126, 175
 - placing values in, 105
- a2scan macro, 12, 146
- a3scan macro, 12, 146
- acos() function, 68, 118
- Administrator, *spec*
 - dial and user settings set by, 30
 - manual for, 187–111
 - news file updated by, 5
- AIX, use of *spec* with, 4
- Alpha-Fixed geometry mode, 170–72
- Am9513 counter chip, boards that use, 251
- an (angle) macro, 135
- Angles
 - dial. *See* Dial positions (angles)
 - freezing, 172
 - user. *See* User positions (angles)

Arithmetic operators, 59–62

array command, 68, 89

Array. *See also* A[], G[], mA[], Q[], S[], S_NA[] and Z[]

- adding built-in to *u_hook.c*, 196
- syntax of, 47

array_dump() function, 68

array_fit() function, 68

array_op() function, 68

array_pipe() function, 68

array_plot() function, 68

array_read() function, 68

asc() function, 68, 119

ascan macro, 12, 123, 124, 145, 146, 154–58

asin() function, 68, 118

atan() function, 68, 118

atan2() function, 68, 118

awk UNIX utility

- script, *show.awk*, 41–41

- to manipulate *spec* data files, 35

Azimuth-Fixed geometry mode, 170–72

B

Background subtraction in *scans.4*, 38

bcd() function, 68, 119

beep macro, 128

Beta-Fixed geometry mode, 171–72

Bitwise operators, 14, 59–62

Boolean operators, 14, 59–61

br (Bragg) macro, 9, 123, 143–44

break statement, 22, 63

bug macro, 128

C

C code, adding site-dependent, 195–97

C math functions, 14

ca (calculate) macro, 10, 123, 143–44

ca_cntl() function, 68, 95

ca_fna() function, 68, 95

ca_get() function, 34, 68, 95

ca_put() function, 34, 68, 95

cal (calculate) macro, 143–44

calc() function, 31, 56, 68, 79–79, 177, 195

calcG macro, 177, 182

calcHKL macro, 32, 159

calcL macro, 182–83

CAMAC (IEEE-583) interface, 3, 34

- hardware functions, 94–95

- installing driver for, 193

- slot assignments, 194, 204–105

camac program, 199
 cat macro, 127
 cd (change directory) macro, 28, 127
 cdef() function, 54, 68, 91, 100, 103
 chdir() function, 28, 56, 68, 69
 chg_dial() function, 30, 68, 105, 210
 chg_offset() function, 30, 68, 106, 210–110
 ci (calculate inverse) macro, 10, 143
 cleanup macro, 25, 54, 159–59
 Clock. *See* counting
 close() function, 17, 68, 76, 82, 86
 cnt_mne() function, 68, 100
 cnt_name() function, 68, 100
 cnt_num() function, 68, 100
 COLS built-in variable, 56
 com (comment) macro, 128
 Command files. *See* File (command)
 Command recall feature (history), 15–16, 79
 syntax for, 50
 Commands, spec
 listing, 26, 48, 80
 types of
 built-in, 68–119
 diagnostic, 26–27
 hardware, 68, 94–97
 macro, 68, 90–93
 program state, 68
 utility, 69–81
 comment macro, 7, 123, 125, 132
 Comments
 in a command file, 26
 pound sign to begin, 26, 29, 46, 204
 syntax for, 46
 config file. *See* File (configuration)
 config macro, 128
 Configuration editor. *See* edconf
 constant command, 55, 68, 89
 Constants
 numeric, 45
 decimal, 49
 floating point, 49
 hexadecimal, 49
 integer, 49
 octal, 49
 syntax of, 49
 string, 45
 escape sequences for, 49
 single or double quotation marks as delimiters
 for, 49
 syntax of, 49
 contents program, 40
 continue statement, 22, 63
 Control key actions
 ^ \ to quit, 5

^C
 cleanup macro run automatically after, 25,
 159–60
 control to command level after, 63
 files turned off after, 17
 to exit *edconf* program, 203
 to halt timer/clock, 12, 101
 to interrupt or abort, 10, 12, 13, 54, 133, 145
 to reset spec, 46
 to stop motors, 10
 ^D to terminate spec session, 55
 ^V to quit, 5
 Conversion
 between degrees C and kilohms, 152
 functions, 68, 119
 cos() function, 68, 118
 count macro, 140–41
 count.mac file, 123
 counter_par() function, 68, 100
 COUNTERS built-in variable, 56
 Counting, 11, 32–33, 100, 140–42
 C-PLOT package
 spec used with, 35–41
 csh UNIX history mechanism, 15
 ct (count) macro, 11, 33, 123, 140
 Cut points for a four-circle diffractometer, 174
 cuts macro, 174, 178
 CWD built-in variable, 56
 cz (calculate zone) macro, 178–79
D
 d (date) macro, 125, 128
 d2scan macro, 12, 146
 d3scan macro, 12, 146
 Data analysis features, 106–17
 Data file. *See* File (data)
 data_anal() function, 68, 108
 data_bop() function, 68, 108
 data_dump() function, 68, 112
 data_fit() function, 68, 110
 data_get() function, 68, 107
 data_grp() function, 68, 106, 107, 114
 data_info() function, 68, 107
 data_nput() function, 68, 107, 108, 110, 112
 data_pipe() function, 68, 113, 113–17
 data_plot() function, 68, 110, 112
 data_put() function, 68, 107, 108, 110, 112
 data_read() function, 68, 110
 data_uop() function, 68, 107
 date() function, 14–15, 68, 69, 128
 Date, returning the current, 69
 dcb() function, 68, 119
 DEBUG built-in variable, 19, 57

debug macro, 128
 def command, 45, 68, 90–93
 Default count time, 140
 deg() function, 68, 119
 delete command, 68, 89
 Device names, specifying–204. *See also* Hardware configuration
 Diagnostic commands, 26
 Dial positions (angles), 8, 28
 listed in degrees, 30
 returning, 104
 setting, 105
 dial() function, 68, 104
 Diffractometer
 angle settings, 5
 configuration, maintaining, 202
 four-circle
 alignment, 166
 cut points, 174
 functions, 177
 geometry for operating, 31
 macros, 178–79
 modes–71
 orientation matrix for, 168–68
 reference manual for, 165–83
 sectors, 173
 spec support of, 22
 variables, 175–77
 geometry, 31–32
 liquid surface
 geometry for operating, 31
 spec support of, 22
 operation, beginner’s guide to, 4–13
 two-circle
 operated by angles alone, 31
 spec support of, 22
 z-axis
 geometry for operating, 31
 spec support of, 22
 Directory, spec, 198–99
 /usr/lib/spec.d for auxiliary files, 58, 199
 /usr/local/lib/spec.d for auxiliary files, 4
 changing, 28, 69
 data, 134
 distribution, 188, 197
 drivers for driver files, 193
 help, 27
 macros for macro source files, 123, 175
 DISPLAY built-in variable, 57
 do macro, 26, 123, 130
 dofile() function, 25, 68, 83
 Drivers, installing, 193
 dscan macro, 12, 36, 146

DSP 6001/6002 CAMAC crate controller, 193
 DSP PC004 IBM PC interface card, 193
 DSP RTC018 Real-Time Clock, wiring, 253

E

ed macro, 127
 edconf program (configuration editor)
 to maintain diffractometer configuration, 202
 to set dial and user settings, 30, 194
 else statement, 21
 energy.mac file, 123
 eprint command, 68, 83, 86
 eprintf() function, 68, 83, 86
 Escape sequences for string constants, 49
 Executor, 45
 exit statement, 63
 exp() function, 68, 118
 exp10() function, 68, 118
 Experiments, automating, 3

F

fabs() function, 68, 118
 Fheader macro, 148, 149, 160
 File hierarchy, typical spec, 198–99
 File(s)
 ASCII
 command file as, 25
 configuration file as, 203
 data file as, 22, 35, 162
 auxiliary
 directory for, 4
 command, 25
 for least-squares refinement of lattice parameters, 180–83
 input, 45
 reading from, 83–84, 130
 startup, 5, 124
 configuration (*config*)
 baud rate for serial lines in, 97
 device numbers set in, 97
 installed hardware described in, 94
 modifying, 128, 203
 motor mnemonics in, 59
 permission levels for security of, 210
 purpose of, 199
 reading, 5–5, 99
 security through motor restrictions in, 210
 slot assignments in, 95, 204–105
 updating, 128, 197
 data
 adding scan results to, 149
 ASCII, standard format for, 22, 162
 controlling output to, 17–18
 inserting comments in, 7

- opening, 6
- selecting, 6
- standard format for, 35–36, 162
- summary scan information from, 40
- width for columns, 157
- distribution, 188
- functions for opening and closing, 81
- hardware configuration. *See* File (configuration)
- help, 27, 70
- index for *scans.4*, 40, 41
- log, 17, 82
- macro source, 123
- news, 5
- reflections, 180
- settings*, 5
 - preventing changes made in, 210
 - reading, 99
 - structure, 203
 - updating, 55, 197
- state, user's, 80, 197
- file.mac* file, 123
- `file_info()` function, 68, 74
- Filers)
 - help, 197
- Flabel macro, 148, 160
- Flow control, 20–22
 - with break statement, 63
 - with conditional statements, 62
 - with continue statement, 63
 - with exit statement, 63
 - with for statement, 62
 - with while statement, 62
- `fmt_close()` function, 68, 88
- `fmt_read()` function, 68, 88
- `fmt_write()` function, 68, 88
- for statement, 21, 62
- Four-circle diffractometer. *See* Diffractometer (four-circle)
- fourc.src* file, 123
- Fout macro, 148
- `fprintf()` function, 17, 68, 76, 77, 86
- freeze macro, 172, 178
- FRESH built-in variable, 57
- Ftail macro, 149
- Functions
 - calling user-added, 79–79, 195
 - types of
 - built-in, 68–119
 - CAMAC, 68, 94–95
 - command file, 68, 83–84
 - conversion, 68, 119
 - counting, 68, 101
 - four-circle, 177
 - GPIB, 68, 95

- hardware, 68, 94–97
- keyboard input and formatted output, 68, 81, 84
- miscellaneous, 68, 70–81
- number, 68, 118
- output control, 68, 81
- plotting and analysis, 68, 106
- serial, 68, 97–97
- string, 14–15, 68, 118
- system, 68, 69
- utility, 69–81

G

- G[]
 - built-in variable, 59
 - geometry parameters stored in, 162, 176, 196
- geo_fourc.c* file, 175
- Geometry
 - configurations, 143–44
 - diffractometer, 31–32
 - four circle, 162–83
- `get_lim()` function, 30, 68, 105
- getangles command, 29, 32–32, 126, 147, 159
- getcounts command, 58, 68, 100
- `getenv()` function, 68, 74
- `gethelp()` function, 27, 68, 70
- `getline()` function, 56, 68, 83, 84
- `getval` function, 126
- `getval()` function, 68, 78, 85
- global command, 18, 68, 89
- Global symbols, 18–20, 89
- GPIB (IEEE-488) interface, 3, 34
 - driver for, 193
 - hardware functions, 95
- `gpib_cntl()` function, 68, 96
- `gpib_get()` function, 56, 68, 96
- `gpib_poll()` function, 68, 96
- `gpib_put()` function, 34, 68, 96
- gpset macro, 125, 128
- Grammar rules
 - of keywords, operators, and commands, 64–68
 - of parser, 45
- grep* utility
 - for file searching, 123
 - to manipulate spec data files, 35
- GTERM built-in variable, 58

H

- h (help) macro, 27, 125, 128
- Hardware configuration, 3–5. *See also* File (configuration), *edconf*
 - reconfiguring, 99
 - selecting, 194
- Help facility, 27

help macro, 27, 123, 128
hi (history) macro, 128
history command, 15, 68, 79
History feature. *See* Command recall feature
hkcircle macro, 12, 147
hkl.mac file, 123
hklmesh macro, 12, 147
hklscan macro, 12, 36, 147
hkradial macro, 12, 147
hlcircle macro, 147
hlradical macro, 147
HOME built-in variable, 58
hscan macro, 12–12, 36, 147

I

I/O ports, enabling on PC platforms, 200
Identifiers (names)
 identifying with `whatis()`, 75–76
 syntax of, 46
if statement, 21, 62
image_get() function, 68
image_par() function, 68
image_put() function, 68
index() function, 15, 68, 118
init_calc() function, 196
initdw macro, 7, 153
initfx macro, 7, 153
initnec macro, 153
initoki macro, 153
Input preprocessor, 45
input() function, 15, 56, 68, 84, 126, 130
Input, translation of keyboard or command file, 45
int() function, 68, 118
Interfaces to user devices. *See* CAMAC interface,
 GPIB interface, and RS-232 interface
io_enable program, 200

J

Joerger SMC Stepper Motor Controller module, 233

K

Keyboard
 interrupts, 54–54
Keyboard
 reading input from, 84
Keywords
 as tokens, 45
 listing, 48, 80
Kinetic Systems Model 3655 Timing Generator modifi-
 cations, 256
klcircle macro, 147
klradical macro, 147

kscan macro, 12, 36, 147

L

l (list files) macro, 125, 127
Lattice parameters, calculating, 180–83
Least-squares refinement of lattice parameters,
 180–83
length() function, 15, 68, 119
less macro, 127
Lexical analyzer, 45
Limits. *See* Motor (limits)
lm macro, 9, 135
local command, 68, 89, 125
log() function, 68, 118
log10() function, 68, 118
Loop
 implemented as a macro in scans, 160–61
 while or for, 63
lp_plot macro, 142
ls (list files) macro, 127
lscan macro, 12, 147
lscmd command, 26, 48, 50, 68, 80
lsdef command, 68, 91
lsdef macro, 26, 50, 125
lup (lineup) macro, 12, 36, 146

M

mA[], motor numbers recorded in, 139
mA[], motor numbers reordered in, 137
Macro(s)
 arguments, style in manual for, 6
 cleanup, 54–54
 defining, 22–23, 90–90
 definition
 argument substitution in, 23, 93
 displaying, 90
 limits of, 93
 listing name and size of, 23–24, 91
 printing, 23, 90, 132
 removing, 24, 91
 library of predefined, 3, 23, 123
 listing all currently defined, 23–24, 91
 output devices used by, 17
 tips for writing, 124–26
 types of
 basic aliases, 127
 basic utility, 128–30
 command file, 130
 counting, 140–42
 four-circle, 178–79
 motor, 135–39
 plotting, 142–42
 printer initialization, 153
 reciprocal space, 143–44

- saving to output device, 132
 - scan, 145–49, 154–61
 - start-up, 133
 - temperature control, 149–52
 - utility, 127
 - zone, 179–80
- mail macro, 127
- Manual
 - administrator's, 187–111
 - conventions of type styles in, 6
 - four-circle reference, 165–83
 - reference, 45–119
 - standard macro reference, 123–62
 - user, 3–41
- Math functions, 14, 68, 118
- MCA. *See* Multichannel Analyzers
- mca_get() function, 68, 94
- mca_par() function, 68, 94
- mca_put() function, 68, 94
- mca_sel() function, 68, 94
- mca_sget() function, 68, 94
- mca_spar() function, 68, 94
- mca_sput() function, 68, 94
- mcount() function, 33, 68, 95, 100, 101
- measuretemp macro, 149–51, 161
- Memory
 - usage, showing, 80
- memstat command, 68, 80
- mesh macro, 12, 146
- Metacharacters, ? and *, 26, 50
- mi (move incident) macro, 143
- mk (move HKL) macro, 123, 143
- Motor(s)
 - controller registers, 29–30
 - controller types, 206
 - controlling, 101–6
 - limits
 - getting, 105, 135
 - setting, 30–31, 106, 135–37
 - software, 9–9
 - listing information for, 8
 - macros, 135–39
 - moving, 7–10, 28–31, 101, 135
 - parameter
 - assignment, 205
 - returning, 102
 - positions
 - (HKL) corresponding to set of, 10
 - displayed on screen, 10
 - reading, 137
 - setting, 7–10
 - storage of, 30
 - returning the mnemonic or name of, 102
 - securing from unauthorized use, 210–110

- stopping, 10, 101, 159, 234
 - unusable, 101
- motor.mac file, 123
- motor_mne() function, 68, 102
- motor_name() function, 30, 68, 102
- motor_num() function, 68, 102
- motor_par() function, 68, 78, 102
- MOTORS built-in variable, 58
- move_all command, 29, 31, 32, 68, 95, 101, 126, 135
- move_cnt command, 68, 101
- move_em macro, 126, 135, 160
- Multichannel analyzers (MCAs), 32, 101
- mv (move) macro, 10, 123, 135, 137
- mvd (move dial) macro, 135
- mvr (move relative) macro, 135
- mz (move zone) macro, 178–79

N

- National Instruments GPIB controller, 193
- ned macro, 127
- newfile macro, 6, 35, 123, 133
- newmac macro, 130
- newsample macro, 133
- NPTS loop variable, 159
- Number
 - functions, 68, 118
 - notation, 14, 49

O

- off() function, 17, 68, 76, 82, 126
- offd (off data file) macro, 17, 128
- offp (off printer) macro, 17, 128
- offsim (off simulate mode) macro, 128
- offt (off tty) macro, 17, 128
- Omega Equals Zero geometry mode, 169
- Omega-Fixed geometry mode, 169, 172
- on() function, 17, 68, 76, 77, 82, 126
- ond (on datafile) macro, 17, 128
- onp (on printer) macro, 17, 128
- onsim (on simulate mode) macro, 128
- ont (on tty) macro, 17, 128
- open() function, 17, 68, 76, 77, 81, 82
- Operators
 - tokens as, 45
 - types of
 - assignment, 61
 - binary, 60
 - ternary, 61
 - unary, 59
- or0 macro, 178–78
- or1 macro, 178
- Oregon Micro Systems PCX motor controller PC interface card, 193

Orientation matrix, 168–68
Output devices, commands for saving to, 132
Output files, controlling, 81

P

p (print) macro, 14, 125, 128
pa (parameters) macro, 143
PAGER built-in variable, 58
Parse tree, 45–46, 63
Parser, grammar rules of, 45
PC platforms
 enabling I/O ports on, 200
Pheader macro, 148, 149, 160
Phi-Fixed geometry mode, 170, 172
PI built-in variable, 18, 58
pl (plane) macro, 135
Plabel macro, 148, 149, 160
plot macro, 142, 142, 161
plot.mac file, 123
plot_cntl() function, 68, 110
plot_move() function, 68, 113
plot_range() function, 68, 110, 113
plot_res macro, 123, 142
Plotting
 functions, 106–17
 macros, 142–42
 scans, 13, 41
Points, maximum number of data, 106
port_get() function, 68, 98
port_getw() function, 68, 98
port_put() function, 68, 98
port_putw() function, 68, 98
Pout macro, 148, 149
pow() function, 68, 118
powder.mac file, 123
prcmd macro, 132
prdef command, 23, 26, 50, 68, 90
print command, 14, 55, 68, 82, 86
Printer
 controlling output to, 17–18
 initialization macros, 153
 selecting, 6
 setting top-of-form position on, 7
printf() function, 17–18, 22, 68, 82, 86, 139
Printing, formatted, 17–18, 86
Propagation of errors formalism in *scans.4*, 40
pts (points) macro, 13, 142
pwd (print working directory) macro, 127

Q

Q[]
 built-in variable, 59
 four-circle coordinate variables stored in, 31, 175

qcomment macro, 125, 132
qdo macro, 26, 123, 130
qdofile() function, 25, 68, 84

R

r2d2.src file, 123
rad() function, 68, 119
rand() function, 68, 118
rdef command, 45, 68, 90–93
read_motors() command, 68
read_motors() function, 56, 101, 105
README files for up-to-date information on
 devices supported in the *config* file, 204
 drivers currently available, 193
reconfig command, 30, 68, 99
reflex macro, 181
reflex_beg macro, 180
reflex_end macro, 182
Relational operators, 14, 59–62
resume macro, 145, 160
ROWS built-in variable, 58
rplot_res macro, 142
RS-232 (serial) interface, 3, 34
 hardware functions, 97–97
RTOT_0 macro, 152

S

S[]
 accessing contents of scalers through, 32
 as built-in variable, 58
 loading, 100
S_NA[], identifying scaler through, 32
savcmd macro, 132
save macro, 133–34
savegeo macro, 133, 134
saveusr macro, 133, 134
savmac macro, 124, 132
savstate command, 68, 80
Scaler channel assignments, 140
Scan header, 13, 35–36, 148, 154, 157
Scan types
 absolute-position motor, 12, 146
 powder-averaging, 102, 148
 reciprocal space, 12–12, 147
 relative-position motor, 12, 146
 temperature, 148
Scan(s)
 aborting, 13, 145, 159
 built of macros, 12
 grid, 147
 invocation syntax, 145
 macros, 145–49, 154–61
 merging in *scans.4*, 38
 motor, 146, 156

- number, 35, 37
- output, customizing, 148–49
- powder mode, 148
- reciprocal space, 146–46, 157
- restarting an aborted, 13, 145
- retrieving with *scans.4*, 37
- sample output, 12
- summary utilities, 40
- temperature, 148
- scan_count* macro, 161
- scan_head* macro, 154, 158
- scan_loop* macro, 160
- scan_move* macro, 154, 160
- scan_plot* macro—42, 160
- scan_tail* macro, 161
- scans.4* C-PLOT user function, 36–40
 - background subtraction with, 38
 - data columns used by, 39–39
 - error bars returned by, 40
 - file conventions, 38
 - file indexing by, 40
 - invoking, 36–37
 - memory for strings and scan numbers, 39
 - merging scans with, 38
 - options, 37
 - retrieving scans with, 37
- scans.mac* file, 123
- scans1.mac* file, 123
- Sectors for four-circle diffractometers, 173
- Security features of *spec*, 3, 210–110
- sed* utility to manipulate *spec* data files, 35
- ser_get()* function, 34, 56, 68, 97
- ser_par()* function, 68, 97
- ser_put()* function, 34, 68, 97
- set* macro, 8, 30, 123, 135, 136, 210
- set_dial* macro, 8, 135, 136
- set_lim()* function, 30, 68, 106, 136, 210
- set_lm* macro, 9, 31, 123, 135–36
- set_sim()* function, 68, 99, 129
- setaz* macro, 178
- setlat* macro, 178
- setmode* macro, 178
- setmono* macro, 178
- setplot* macro, 13, 133, 142–42, 145
- setpowder* macro, 148
- setscans* macro, 133, 145
- setsector* macro, 178
- setslits* macro, 133
- settemp* macro, 149, 149–52
- settings* file. *See* File (*settings*)
- shell escapes, *See* Subshells
- show_cnts* macro, 12, 123, 140–41
- showscans* program, 41–41
- showtemp* macro, 149–50
- Simulation mode, 99
- sin()* function, 68, 118
- Site-dependent C code, adding, 195–97
- site.mac* file, 124
- sleep()* function, 68, 76, 137
- slit.mac* file, 123
- sock_io()* function, 68, 80
- Software motor limits, 9–9. *See also* Motor (limits)
- spec*
 - as a calculator, 14–15
 - C-PLOT package used with, 35–41
 - customized with C code, 195–97
 - exiting, 5
 - features, 3
 - installation, 187–94
 - internal structure, 45–46
 - motor security of, 3, 210–110
 - purpose of, 3, 28
 - standard scans in, 145
 - start-up of four-circle version from a UNIX shell, 4
 - terminating, 5
 - UNIX utilities used with, 35–41
 - updating, 197–97
 - user interface, 14–35
 - welcome message, 4
- SPEC built-in variable, 58
- spec.mac* command file, 25, 124
- spec_par()* function, 52, 68, 76–79, 81, 83, 84, 102
- specadm* user account, 187, 188
- SPECD built-in variable, 58, 131
- Special characters in string constants, listing of, 49
- split()* function, 68, 119
- splot* macro, 13, 142–42
- splot_res* macro, 142
- sprintf()* function, 15, 68, 119
- sqrt()* function, 68, 118
- srnd()* function, 68, 118
- sscanf()* function, 68, 119
- start.mac* file, 123
- startgeo* macro, 134
- starttemp* macro, 133
- startup* macro, 5–6, 133
- stop()* function, 68, 99
- String
 - functions, 68, 118
 - patterns, 50
- stty* UNIX command, 54
- su* command, 188
- Subshells, spawning, 27–28, 69, 127
- substr()* function, 15, 68, 119
- Sun computers, use of *spec* with, 4
- surf.src* file, 123

syms command, 19, 26, 50, 68, 89
sync command, 30, 68, 102
Syntax conventions, 46–68
Syntax error, 45
System V/386 compatible UNIX
 serial line modes for, 97
 use of `spec` under, 200
sz (set zone) macro, 178–79

T

tan() function, 68, 118
tar command, 189
tcount() function, 32–33, 68, 95, 100, 101
te macro, 149–50
temper.mac file, 123
Temperature control, macros for, 149–52
teramp macro, 149, 152
TERM built-in variable, 58
Ternary operator for `spec` calculator, 14. *See also*
 Operators
test UNIX utility to check for file's existence, 126
th2th macro, 146
Three Circle geometry mode, 170
time() function, 15, 35, 68, 69
Timer/clock. *See also* Counting
 halted with ^C, 101
 starting, 100
Tokens, input text broken into, 45
tty_cntl() function, 50, 68, 87, 88
tty_fmt() function, 50, 68, 88
tty_move() function, 50, 68, 87, 113
tw (tweak) macro, 11, 135
twoc.mac file, 123

U

u macro, 28, 123
U[]
 built-in variable, 59
u_hook.c file, 175, 195
uan macro, 135
ubr macro, 143
uct macro, 12, 123, 140, 141, 145
umk macro, 143
umv (updated-move) macro, 10, 123, 135, 137, 145
umvr macro, 135
undef command, 68, 91
unfreeze macro, 172, 178
unglobal command, 68, 89
UNIX commands
 in macro definitions, 27–28
 macros for common, 127
UNIX epoch, 69
UNIX utilities, `spec` used with, 35–41

unix() function, 27–28, 53, 68, 69

Updated activities

 counting, 141
 moving, 10, 137, 143
 plotting, 142–42
 scans, 145
 setting UPDATE, 137, 143, 145

upl macro, 135

User account for administering `spec`, 187

USER built-in variable, 59

User positions (angles), 8, 28

 listed in degrees, 30
 listing, 29
 offset between dial angle and, 104
 returning, 104

user() function, 68, 105

util.mac file, 123

uwm macro, 135

V

Variable arguments, style in manual for, 6

Variables

 as tokens, 45
 attributes
 built-in, 55–59
 constant, 19, 55, 68, 89
 global, 18–19, 55, 68, 89, 133
 immutable, 55
 local, 55, 68, 89, 125
 changing, 19
 defined through usage, 18, 55
 four-circle, 175–77
 nonglobal, 68, 89
 limits of, 46
 symbols for, listing, 19–20
 syms, 68, 89
 types
 array, 55
 number, 55
 string, 55

VENIX

 quit control character on, 5
VERSION built-in variable, 59
vi (visual editor) macro, 127
vme_get() function, 68, 98
vme_get32() function, 68, 98
vme_move() function, 68, 99
vme_put() function, 68, 98
vme_put32() function, 68, 99

W

w (wait) macro, 10, 128

wa (where all) macro, 8, 123, 135, 138

`wait()` function, 33, 68, 77, 78, 99, 101
`waitall` macro, 128
`waitcount` macro, 128
`waitmove` macro, 32, 126, 128
Warning messages, 5
`wh` macro, 7, 123, 143–44
`whatis()` function, 21, 68, 75–76, 125
`whats` macro, 128
while statement, 21, 62
`wm` (where motors) macro, 9, 135

X

X rays, counting. *See* counting

Y

`yesno` macro, 126, 128
`yesno()` function, 68, 85

Z

`z[]`
built-in variable, 59
four-circle geometry zone mode, 177, 179
`zaxis.src` file, 123
Zone geometry mode, 170, 172, 179–80

NAME

spec – X-ray diffractometer control and general data acquisition package

SYNOPSIS

```
spec [ -ffsyq ] [ -g geometry ] [ -N my_name ] [ -u user ] [ -t tty ] [ -T fake_tty ] [ -d
debug ] [ -o option=value ] [ -p fd pid ] [ -D directory ]
```

DESCRIPTION

spec provides a software environment for the operation of an X-ray diffractometer and other data-acquisition instruments. **spec** contains a sophisticated command interpreter that uses a C-like grammar and is partially modeled on the standard UNIX utility *awk*. **spec** supports a variety of X-ray diffractometer configurations. The diffractometer geometry is chosen by the program name. Those currently supported include:

- spec** – Generic instrument control
- fourc** – Standard four-circle diffractometer
- twoc** – Standard two-circle diffractometer
- sixc** – Six-circle diffractometer ($\delta, \theta, \chi, \phi, \mu, \gamma$)
- psic** – An S4-D2 six-circle diffractometer
- kappa** – Kappa diffractometer
- surf** – Various liquid surface diffractometers
- zaxis** – Standard *z*-axis diffractometer

The following options are recognized:

- f** Fresh start. All symbols are set to their default values and the standard macros are read to establish the default state.
- F** Clean and fresh start. All symbols are set to their default values but no command files are read and no macros are defined. Only the built-in commands are available.
- s** Simulation mode. No hardware commands are issued. If started in simulation mode, simulation mode cannot be turned off without restarting the program.
- y** Yes, change motor controller registers initially if they disagree with the *settings* file. Normally, **spec** requires you to confirm such a change. This flag would be useful if you know controller power had been turned off.
- q** Indicates that **spec** should operate in quiet mode and allow output to all devices to be turned off. This option is only valid when used with the **-p** option.
- g** *geometry*
Force loading of macro files and activation of geometry calculations for the specified geometry, while using the configuration files taken from the name by which **spec** is invoked.
- N** *my_name*
Use *my_name* for setting the interactive prompt and the name of the directory containing the *config*, *settings* and state files. Normally the name by which **spec** is invoked is used.
- u** *user*
Use *user*'s last saved state as the current user's initial state.
- t** *tty*
Use the current user (or *user*'s) last saved state from the terminal specified by *tty*. The terminal can be specified as **-t** */dev/tty01* or **-t** *tty01*.

- T** *fake_tty*
This option creates a user state associated with *fake_tty*, which may be any name. This option allows you to bypass the locking feature that prevents multiple instances of **spec** to be started by the same user from the same terminal.
- d** *debug* Sets the initial value of the debugging variable `DEBUG` to *debug*. The available debugging categories are described on page 57 in the *Reference Manual*. A value of 192 is useful for debugging hardware problems.
- o** *option=value*
Initializes the `spec_par()` *option* to *value*. The available `spec_par()` options are described on page 76 in the *Reference Manual*.
- p** *fd pid* Indicates that **spec** input is coming from a pipe from another program. The argument *fd* is the file descriptor that **spec** should use for standard input. The argument *pid* is the process ID of the spawning process. If *fd* is zero, **spec** will not re-echo input from the file descriptor to **spec**'s standard output.
- D** *direc* Use *direc* instead of the compiled-in name (usually */usr/local/lib/spec.d*) or the **SPECD** environment name as the auxiliary file directory.

FILES

(**SPECD** is the auxiliary file directory, normally */usr/local/lib/spec.d*.)
(**geom** is the first four letters of the name by which **spec** was invoked.)
(**spec** is the complete name by which **spec** was invoked, as in *fourc*, *twoc*, etc.)

<i>./spec.mac</i>	Optional private command file always read at start-up.
SPECD / <i>site.mac</i>	Optional site command file always read at start-up.
SPECD / <i>site_f.mac</i>	Optional site command file only read when starting fresh.
SPECD / <i>standard.mac</i>	Standard macro definitions.
SPECD / <i>geom.mac</i>	Geometry macros.
SPECD / <i>spec/geom.mac</i>	More geometry macros.
SPECD / <i>spec/config</i>	Hardware configuration file.
SPECD / <i>spec/settings</i>	Motor settings file.
SPECD / <i>spec/conf.mac</i>	Optional configuration command file always read at start-up.
SPECD / <i>spec/userfiles/hdw_lock</i>	Spectrometer lock file.
SPECD / <i>spec/userfiles/user_tty.S</i>	User's state file. Uses only first 6 letters of <i>user</i> and <i>tty</i> .
SPECD / <i>spec/userfiles/user_tty.P</i>	User's data points file.
SPECD / <i>spec/userfiles/user_tty.L</i>	User's lock file.
SPECD / <i>help/*</i>	Help files.